



Global Measuring Standard

Guideline for measuring JOBST®
made-to-measure flat-knit garments



Compression to uplift your patients' life

At JOBST, we are passionate about providing compression therapy solutions that **empower** patients to **confidently** live the life they want, **without compromise** - despite their lymphatic or venous condition, or lipedema.

Everyone has different needs, so we design our products to support those differences. Our knowledge, insights and years of experience mean that we keep providing innovative products and services that suit your patients' lifestyle, giving them the comfort, support and peace of mind that they truly deserve while wearing compression.

Discover the world of JOBST and how we can not only help your patients to achieve their health goals, but more importantly, enable them to "**experience more of life's special moments**".

Overview of JOBST® made-to-measure flat-knit garments




Table of contents

	Introduction		
2	Compression to uplift your patients' life	20	JOBST Confidence Lower extremities - At a glance
3	Overview of JOBST made-to-measure flat-knit garments	51	JOBST Elvarex Upper extremities - At a glance
4	Table of contents	21	JOBST Relax Lower extremities - At a glance
5	Purpose of this document	52	JOBST Elvarex Soft Upper extremities - At a glance
6	When using this document, please keep in mind	53	JOBST Confidence Upper extremities - At a glance
7	How to use this document	54	JOBST Relax Upper extremities - At a glance
8	JOBST Measuring Standard		
10	Before you get started		Measuring instructions
13	JOBST eShop	22	Measuring for JOBST Lower extremities
		55	Measuring for JOBST Upper extremities
		48	JOBST Elvarex/Soft/Confidence Lower extremities - Overlapping garment instructions
	Tension Matrix		
15	Explanation of the Tension Matrix		Supporting video material
16	JOBST flat-knit Lower extremities Overview of tension per measuring point	49	JOBST Elvarex Lower extremities - Donning & Doffing videos
17	JOBST flat-knit Upper extremities Overview of tension per measuring point	50	JOBST Confidence Lower extremities - Donning & Doffing videos
		80	JOBST Confidence Upper extremities - Donning & Doffing videos
	JOBST flat-knit garments	81	JOBST Elvarex Plus and JOBST Elvarex Soft Seamless
14	How to choose between JOBST Elvarex and JOBST Confidence - Decision tree	82	Wash and care information
18	JOBST Elvarex Lower extremities - At a glance	84	Tips for better measurement - JOBST Confidence
19	JOBST Elvarex Soft Lower extremities - At a glance	86	Contact details

Introduction (1/3)

Purpose of this document



Measuring a flat-knit garment can be quite challenging, as there are many aspects that need to be considered during the process.



Being your trusted partner in compression therapy, we would like to ease the process of measuring and give as much support as possible.



The goal is to come close to a standardized measuring for JOBST flat-knit garments, in order to create compression hosiery that will have the best possible fit straight from the beginning.

Introduction (2/3)

When using this document, please keep in mind:



This booklet is a guide to measuring for JOBST made-to-measure flat-knit compression garments and must be used as a guide only. Please use alongside the JOBST order forms or the JOBST eShop.



Note: this Measuring Standard does not replace any formal training on treatment of lymphatic and venous diseases or lipedema.

Introduction (3/3)

How to use this document



Our JOBST Measuring Standards have been developed as an interactive PDF, which you can use on any digital device.



Using this document in full-screen mode enables you to interactively click through the slides, following the links while clicking on the respective buttons.



Clicking on the footer, you will always be able to go back to the table of contents or the tension matrix.



You also have the possibility to print the JOBST Measuring Standard completely or partly.

JOBST Measuring Standard

General information (1/2)



The limb should be decongested before measuring for a compression garment and edema reduced through multi-layer lymphedema bandaging or JOBST® FarrowWrap® compression system.



The best time to measure a patient for a compression garment is in the morning, or after a therapy session.

JOBST Measuring Standard

General information (2/2)



Measure hands, feet, fingers and toes accurately to the nearest millimeter.



Photo documentation is helpful in difficult cases.



Documentation of body measurements (no tension) and order measurements are important in order to track changes of the edema.



Do not measure into skin folds but place the measuring tape over the top of the skin fold.

Before you get started

This is what you need (1/3)

JOBST measuring tapes

Measuring tapes of different widths are recommended for taking measurements, depending on which areas of the extremities are being measured:

- A 0.5 cm wide measuring tape or appropriate finger / toe tape is used to ensure measurements to the nearest millimeter in the hand and foot areas.
- A regular JOBST measuring tape to measure the remaining circumference and length measurements (feet, legs, possibly abdomen, hand and arm).

JOBST skin pencil

A skin pencil is recommended for marking the measurement points.

Our JOBST skin pencils are particularly skin-compatible and are easy to remove from the skin due to their fat solubility.

Before you get started

This is what you need (2/3)

JOBST measuring card

The measuring card is a measurement aid developed exclusively for arm sleeves, gloves and toe caps. The card allows finger and toe lengths to be determined very easily and precisely. Moreover, the measurement aid is used for positioning measurement point G on the arm.

JOBST measuring board (Optional)

A measuring board is recommended for accurately measuring the lower extremities. It is particularly convenient to take combined length and circumference measurements without marking the patient's skin.

JOBST order form or JOBST eShop

In order to note down the measurements and place the order, please have either an order form or the JOBST eShop prepared.

Before you get started

This is what you need (3/3)

JOBST foam pad (Optional)

The foam pad is a measuring aid developed for measuring length K.

The pad allows accurate measurement to be taken for point K. It is placed in between the groin area.


JOBST elastic belt (Optional)

The elastic belt helps the patient find the best position for themselves and delivers a secured position for the measurement. It aids in the measurement of circumference and length at point H and T, especially for JOBST Confidence panty (AT) garment.

JOBST tape clipper (Optional)

The tape clipper is used on a measuring tape. This is needed for taking rise measurement for JOBST Confidence panty (AT) garment.

JOBST® eShop



Are you looking to make your work routine easier and digitize your processes?
JOBST eShop is here to help!




Intelligent patient management



Easy measurement documentation


Discover JOBST
eShop now!
eshop.jobst.us



Time efficient re-ordering

How to choose between JOBST Elvarex and JOBST Confidence in patients with lymphedema*

Decision tree



* Developed by Essity based on technical features of the products; individual patient needs and conditions must be taken into consideration

Tension Matrix

In order to measure for an optimal fit of our flat-knit garments, circumferential measurements should be carried out under different levels of tension.

As shown here, the tension can range from 0 (lay on measurement) to +++ (maximum tension).

Measuring for JOBST flat-knit garments involves tensions 0 to ++ pull.


The maximum tension of +++ is shown here to give you an indication and is actually not used while measuring patients.

-  **Lay on Measurement**
-  **Pull until first resistance**
-  **Pull until second resistance**
- Maximum tension**

JOBST flat-knit | Lower extremities

Overview of tension per measuring point


- For each of the **circumferential measurements** required for a **JOBST flat-knit garment**, a specific tension level will be used.
- The graphic shows you the overall picture and includes the tension matrix icons next to the measuring point.
- Additionally, you will find the **amount of tension** to be used for every **circumferential measuring point** in the step-by-step measuring instruction on the following slides.
- This tension matrix is valid for all JOBST flat-knit garments and shall be applied for Lymphedema and Lipedema patients.
- Please note that the **pain tolerance** of the patient must be taken into account.



JOBST flat-knit | Upper extremities

Overview of tension per measuring point

- For each of the **circumferential measurements** required for a **JOBST flat-knit garment**, a specific tension level will be used.
- The graphic shows you the overall picture and includes the tension matrix icons next to the measuring point.
- Additionally, you will find the **amount of tension** to be used for every **circumferential measuring point** in the step-by-step measuring instructions on the following slides.
- This tension matrix is valid for **all JOBST flat-knit garments** and shall be applied for Lymphedema and Lipedema patients.
- Please note that the **pain tolerance** of the patient must be taken into account.



* with silicone band & straight ending

** for JOBST Confidence

JOBST Elvarex* | Lower extremities

At a glance

Popular styles



AD
Knee High



AG
Thigh High



AT
Panty



BIT/BT
Leggings/Capri



CT
Capri



ET
Bermuda



FT
Biker Shorts



AA
Foot Cap

Compression classes

CCL 1

CCL 2

CCL 3 ▶ 3 Forte

CCL 4 ▶ 4 Super

Popular options



Functional Zone
T-Heel



Slant Foot



Functional Zone
Knee



Slip Form

JOBST Elvarex Soft | Lower extremities

At a glance

Popular styles



AD
Knee High



AG
Thigh High



AT
Panty



BIT/BT
Leggings



CT
Capri



ET
Bermuda

Compression classes



CCL 1



CCL 2



CCL 3

Popular options



Slant Foot



Slip Form

JOBST Confidence | Lower extremities

At a glance

Popular styles



AD
Knee High



AG
Thigh High



AT
Panty



BT/BT
Leggings



CT
Capri



ET
Bermuda

Compression classes



CCL 1



CCL 2



CCL 3

Popular options



Functional Zone
Ankle



Slant Foot




Functional Zone
Knee

JOBST Relax | Lower extremities

At a glance

Popular styles



AD
Knee High

AG
Thigh High

Compression classes

CCL 1

CCL 2


This garment has been developed to be soft and comfortable to the skin while sleeping or relaxing, forgoing the usage of a top band. It is not constructed for a lot of movement, especially strong bending of the knees. As a consequence the garment might move from its original position. If needed, adjust before and after sleeping. Furthermore, some patients might experience a slight volume reduction during the night and find the garment loose in the morning.

Measuring for JOBST | Lower extremities

General 1/2

Please keep in mind:

- Medial measurements are taken on the inside of the limb; lateral measurements on the outside.
- Proximal refers to points closer to the head; distal closer to the feet.



Measuring for JOBST | Lower extremities

General 2/2

When the patient is in a supine position:

- Position patient carefully on the measuring board.
- The foot must be flat on the board.
- The heel must be as close to the board as possible.
- If you have a measuring board or use JOBST Xpert Voice, it is faster to measure pairwise (length and circumference at the same time).
- If a measuring board is not available, make sure that the position of the foot is flexed and that the measuring tape is kept in a straight line from heel to pubic bone.

Measuring for JOBST | Lower extremities

Measuring the foot

- First, measure the foot lengths with the patient standing. Mark all length measuring points on the foot.
- The toe can be open or closed.





- Please note that a slant foot offers the best anatomical fit.
- Please keep in mind that you can also measure for a slant foot with closed toe.

Measuring for JOBST | Lower extremities

Lengths of foot (A, Z)

- **Open toe straight foot:** First measure the medial length from the back of the heel, to just under the base of the big toe. Continue with the lateral length from the back of the heel to just under the base of the small toe. Take the measurement in between to define length A.
- **Open toe slant foot:** First measure the medial length from the back of the heel, to just under the base of the big toe. Continue with the lateral length from the back of the heel to just under the base of the small toe. Note both measurements down.
- **Closed toe:** Also measure the total foot length from the back of the heel to the end of the longest digit on the medial side.
- Please note that for JOBST Relax, only medial length A is taken (IA-medial).




Measuring for JOBST | Lower extremities

A Circumference

- This measurement is taken with the patient standing.
- Measure circumference A at the base of the toes at the widest part of the foot.



- If measuring is done in lying position, please apply 0 tension.




Measuring for JOBST | Lower extremities

Y Circumference

- Position the patient in a supine position.
- Y circumference is taken around the heel to the front crease, at maximal dorsal flexion.



- If there are deep skin folds or irregular contours, measure cY more towards proximal.





Measuring for JOBST | Lower extremities

B Circumference and Length

- Circumference B is measured just above the malleoli, at the narrowest part of the ankle.
- Length B is measured where the B circumference is taken.
- For JOBST Relax, if cB is below 20cm, cB should be measured with 0 tension.




- If the limb definition is unclear, we recommend taking the measurement 10 cm from the heel.



Measuring for JOBST | Lower extremities

B1 Circumference and Length

- Circumference B1 is measured at the insertion of the Achilles tendon into the calf muscle.
- Length B1 is measured where the B1 circumference is taken.



- If you cannot feel the Achilles tendon, measure B1 at the midway point of length D.



Measuring for JOBST | Lower extremities

C Circumference and Length

- This is an anchor-point for an AD garment (knee-high).
- Circumference C is measured at the greatest calf circumference.
- Length C is measured where C circumference is taken.
- For CT (Capri Style) - if one wishes, the end point of the garment may be measured a little lower than the actual point C. But, the circumference has to be measured at exact point C.



- When calf circumference is unclear, C is located in-between B1 and D.



Measuring for JOBST | Lower extremities



D Circumference and Length

- Circumference D is measured two fingers below patella at pre-patellar bursa.
- Length D is measured where circumference D is taken.

0 For silicone band and straight ending



cD




- If the patella is hard to feel, measure two fingers below the knee crease on the medial side (length).

Measuring for JOBST | Lower extremities

E Length


- E length is measured across the middle of the patella.
- The leg must be extended and flat against the measuring board, with the foot in the correct flexed position (90 degree towards the board) and the heel as close to the board as possible.



Measuring for JOBST | Lower extremities

E Circumference


- Circumference E is measured over the middle of the patella.
- The leg must be flexed (30-45 degrees).
- This allows comfort and range of movement while wearing the garment.



Measuring for JOBST | Lower extremities

E Circumference


- Circumference E can be measured differently for severe cases and shall then be measured approx. 4 cm above the patella.



Measuring for JOBST | Lower extremities

E1 Circumference and E1 length for Bermuda


- Point E1 is required for Bermuda style.
- Circumference E1 measured at approx. 5 cm above point E.
- Length E1 is measured where circumference E1 is taken.



Measuring for JOBST | Lower extremities

F Circumference and Length


- Circumference F is located mid-thigh.
- This is the anchor point which holds the garment in place for AG (thigh-high) and AT (tights).
- Length F is measured where circumference F is taken.
- For AG-HT one legged panty, point F of the shorter leg (non-compression leg) should be measured with 0 tension. Also, length F of non-compression leg can differ from length F of compression leg.
- For JOBST Relax, if needed, slightly more tension can be used.



Measuring for JOBST | Lower extremities

G Circumference and Length


- Measure point G with the patient standing.
- Circumference G is measured at the widest part of the thigh.
- Length G is measured up to the gluteal fold.
- These measurements may not be taken at the same point.
- For AG-HT one legged panty, point G of the shorter leg (non-compression leg) should be measured with 0 tension.
- For JOBST Relax, if needed, slightly more tension can be used.



Measuring for JOBST | Lower extremities

K Length


- Measure point K with the patient standing.
- Place the JOBST foam pad in between the groin area.
- Now measure the length K from floor to the upper edge of the foam pad in a straight line.
- This measurement is essential for JOBST Confidence AG garments and all AT flat-knit garments.
- Length K is not necessary for JOBST Elvarex and JOBST Elvarex Soft AG garment, unless you want to order a to p comfort zone for JOBST Elvarex.



Measuring for JOBST | Lower extremities

H Circumference for JOBST Elvarex and Elvarex Soft


- H circumference is measured over the widest part of the buttocks area.



Measuring for JOBST | Lower extremities

H Circumference and Length for JOBST

- You can use an elastic belt to define and secure point H. Then, use a measuring tape and measure the circumference at the upper edge of the belt.
- Length H is measured at the side of the leg, from the feet up to the upper edge of the belt. Measure along the contour of the hip at the upper part.




Measuring for JOBST | Lower extremities

Body bandage for JOBST Elvarex and Elvarex Soft

How to measure GT and KT Length

- K / K1 to T: Measure length point T from eye-line of pubic bone to waist, pulling tape tension on any contours.
- G / K2 to T: Measure from gluteal crease to T, pulling tape tension on contours.



- If the belly of the patient is heavy, kindly ask your patient to flatten the belly with their hands when measuring length. This creates a 'lift'-effect.

Measuring for JOBST | Lower extremities

T Circumference for JOBST Elvarex and Elvarex Soft


- Point T is located at the waist and indicates the end of the garment.
- To measure the circumference at point T, place the elastic belt comfortably around the waist and measure the circumference along the upper edge of the belt.



Measuring for JOBST | Lower extremities

T Circumference and T Length for JOBST Confidence

- Point T is located at the waist and indicates the end of the garment.
- To measure the circumference at point T, place the elastic belt comfortably around the waist and measure the circumference along the upper edge of the belt.
- Length T is measured at the side of the leg from the feet to the upper edge of the belt, following the contour from G to T.



- For an optimum fit, find the natural waistline, consider support point such as the smallest circumference, skin crease, etc.

Measuring for JOBST | Lower extremities

Rise measurement TT (L) and K3 (L1) for JOBST Confidence

TT (L) and K3 (L1) need to be measured in one go using the JOBST tape clipper and the elastic belt.

Measuring length TT (L):


- Place the elastic belt comfortably around the waist.
- Place the tape with clipper at the middle of the belt at the front side of the body.
- Gently extend the tape between the legs, through the crotch to the middle of the belt at the back side of the body. This is length TT (L).

Measuring length K3 (L1):

- Slide the tape clipper from the top of the waistband frontside to the middle point between the legs. Take the measuring tape away. The number at the clipper is length K3 (L1).




- Avoid either securing the tape too tight or too loose.
- For men, the tape should be placed in the groin, beside the genitals and not on top of genitals.



JOBST Elvarex | Toe cap

A-A1 Length

- Measure the foot lengths with the patient standing.
Mark all length measuring points on the foot.
- The length is measured from A to A1.
- This length is taken for the lateral side and the medial side.
- These measurements will differ.



- It is recommended to have a length difference between medial and lateral of at least 2 cm.
- For a better anatomical fit, the recommended medial length is at least 8 cm.


JOBST Elvarex | Toe cap

A and A1 Circumference

- This measurement is taken with the patient standing.
- A is measured at the base of the toes or the widest part of the foot.
- A1 is measured at the end of the toe cap.



- When the toe cap and leg garment are worn together, it's recommended to order a slant foot to avoid pain on the bones and callus growth.
- When the compression class of the leg garment and the toecap are the same, add 1 cm to cA and cA1 to avoid double compression.




JOBST Elvarex | Toe cap

X and Z digit Circumferences

- Measure circumference X at the base of the toe.
- Measure circumference Z at the end of the toe.
- The X should not be less than the Z measurement.
- These measurements are taken without tension.
- All 5th toe circumferences are required for Elvarex® Plus, even if choosing open 5th toe option.
- 5th toe cut-back is only available in Elvarex® and Elvarex® Plus. No cut-back in Elvarex® Soft Seamless.

Digit Length


- Measure the toe length using the measuring card.
- Place the measuring card between the toes, measuring the side closest to the second toe.
- Measure at least up to the nail-bed to avoid swelling at the end of the toe.




- It is easier to measure this circumference using a small measuring tape.

JOBST Elvarex/Soft/Confidence | Lower extremities

Overlapping garment instructions



AD + CT Capri




AD + BT Leggings




AD + B1T Leggings

Add 2-3 cm to Capri/Leggings at the points of overlap.



AG + E(1)T Bermuda



AG + FT Biker Shorts

Add 1-2 cm to Bermuda/Biker at the points of overlap.

- In case of patients with severe oedema, you **may not** add extra measurements.


JOBST Elvarex | Lower extremities

Donning & Doffing videos

Tips for easy application




Tips for easy removal




JOBST Confidence | Lower extremities

Donning & Doffing videos

Tips for easy application



Tips for easy removal



JOBST Elvarex* | Upper extremities

At a glance

Popular styles



Combined
Armsleeve



AE
Long Gauntlet



AE
Long Glove



AC1
Glove



AC1
Gauntlet



CG1
Armsleeve

Compression classes

CCL 1

CCL 2

Popular option



Functional
Zone Elbow

*Caution: This product contains natural rubber latex which may cause allergic reactions.

JOBST Elvarex Soft | Upper extremities

At a glance

Style



CG1
Armsleeve

Compression classes

CCL 1

CCL 2

JOBST Confidence | Upper extremities

At a glance

Popular styles



CG1
Armsleeve



AG1
Combined
Armsleeve

Compression classes

CCL 1

CCL 2

Options



Functional Zone
Elbow

JOBST Relax | Upper extremities

At a glance

Popular styles



Combined
Armsleeve



CG1
Armsleeve

Compression classes

CCL 1



This garment has been developed to be soft and comfortable to the skin while sleeping or relaxing, forgoing the usage of a top band. It is not constructed for a lot of movement, especially strong bending of the knees. As a consequence the garment might move from its original position. If needed, adjust before and after sleeping. Furthermore, some patients might experience a slight volume reduction during the night and find the garment loose in the morning.

Measuring for JOBST | Upper extremities

General

Position the patient:

- Preferably on a chair, beside a table or bench. The arm lays on the surface and has a bent elbow as shown in the image below.
- First, mark ALL measuring points with a skin pencil.
- After that, proceed with measuring all circumferences. Then continue with measuring all lengths.



* with silicone band & straight ending

** for JOBST Confidence

Measuring for JOBST | Upper extremities

C Circumference

- Ask the patient to flex the hand towards inner wrist.
- At second crease, mark point C using a skin pencil.
- Ask the patient to lay the hand flat on the table (palm down) and measure circumference C around the wrist bone.

Note: the circumference is not taken at the exact position of the length mark.

When measuring for JOBST Elvarex and JOBST Elvarex Soft, please keep in mind:


- If proximal from C is bigger, always take the biggest circumference.
- If circumference is between 15 and 18 cm, add 1 cm to the measurement taken.
- If circumference is < 15 cm, add 2 cm to the measurement taken.



Measuring for JOBST | Upper extremities

C1 Circumference (Only for JOBST Confidence)


- Point C1 is located max. 6 cm proximal from C.
- Measure C1 circumference.



Measuring for JOBST | Upper extremities

D Circumference


- Point D is located at 2/3 length between point C and E.
- Measure D circumference here.



Measuring for JOBST | Upper extremities

E Circumference

- At the elbow crease, mark point E.
- Measure E circumference with the elbow flexed to ensure enough width for movement (45°).




- If the tissue from proximal E is bigger, take the bigger circumference.

Measuring for JOBST | Upper extremities

F Circumference


- At mid upper arm, mark point F.
This is 1/2 way between landmark E and G.
- This is the anchor point which holds the garment in place.
- Measure F circumference here.
- For JOBST Relax, if needed, slightly more tension can be used.



Measuring for JOBST | Upper extremities

G Circumference

- Use the measuring card to define point G.
Ask the patient to comfortably place it high in the axilla.
- Fold the measuring card over and mark point G.
This should be on the axillary line.
- Measure G circumference.
- For JOBST Relax, if needed, slightly more tension can be used.




* with silicone band & straight ending

Measuring for JOBST | Upper extremities


CC1, CD, CE Length

Now, continue with taking the length measurements as described below:


- Measure length C to C1 on inner aspect of the arm (only for JOBST Confidence).



- Measure length C to D on inner aspect of the arm.




- Measure length C to E on inner aspect of the arm and fix the measuring tape with your finger, before moving to the next measuring point (see following slide).




Measuring for JOBST | Upper extremities

CF, CG Length

- Measure length C to F on inner aspect of arm.




- Measure length from C to G on inner aspect of arm.



Measuring for JOBST | Upper extremities

Shoulder Cap for Bra Loop or Body Strap


- Measure the length from point G1 to outer aspect of bra strap point H.
- Apply the measuring tape firmly over the contour of the shoulder.



Measuring for JOBST | Upper extremities

Shoulder Cap for Bra Loop or Body Strap


- Measure the width of the bra strap for the VELCRO brand hook and loop fastener system.
- A minimum 2.5 cm strap width is recommended, to be effective.



Measuring for JOBST | Upper extremities

Body strap


- Measure from the shoulder across the body to where the strap should sit comfortably.
- Do not measure around the body, only on the frontside and in a straight line.
- Then, multiply the measured length by 2. This is the body strap length.



Measuring for JOBST | Gloves and Gauntlets

Defining measuring points A and B


- Place the hand flat on the table, palm up.
- Use the JOBST measuring card to identify point A.
- Place the measuring card in line with the base of the 5th and index finger. Just below the middle finger, mark point A
- After this, slide the card up, towards the thumb, and position the card in to the thumb web space.
- At middle of the hand mark point B.



Measuring for JOBST | Gloves and Gauntlets

Defining measuring points C and C1

- Ask the patient to flex the hand towards inner aspect of wrist.
- At second crease, mark point C.
- Point C1 is the ending point of the glove and should be located at least 5 cm from C.




- The maximum ending point is at 20.5 cm proximal from point A.

Measuring for JOBST | Gloves and Gauntlets

Defining measuring points D and E (long glove, long gauntlet)


- Point D is located at 2/3 length between landmark C and E.
- Point E is located 1 cm distal of the elbow crease.
- Point E is the ending point of the glove.



Measuring for JOBST | Gloves and Gauntlets

A Circumference


- Place the hand flat on a table (palm down).
- Ask the patient to apply pressure through to hand.
- Measure A circumference.



Measuring for JOBST | Gloves and Gauntlets

B Circumference


- Place the hand flat on a table (palm down).
- Ask the patient to apply pressure through to hand.
- Measure B circumference.



Measuring for JOBST | Gloves and Gauntlets

C Circumference

- Place the hand flat on a table (palm down).
- Measure circumference C around the wrist bone.



- If the glove is worn **together with an armsleeve**, you will have an **overlap** and should **add 1 cm to C**.


Measuring for JOBST | Gloves and Gauntlets

C1 Circumference

- Measure C1 circumference.




- If the glove is worn **together with an armsleeve**, you will have an **overlap** and should **add 1 cm to C1**.



Measuring for JOBST | Gloves and Gauntlets

D and E Circumference (long glove and long gauntlet)




- Measure D circumference.
- Measure E circumference 1 cm below the elbow crease.



Measuring for JOBST | Gloves and Gauntlets

AB, AC, AC1 Length

- Measure length AB on inner aspect of the palm.
- Measure length AC on inner aspect of the hand.
- Measure length AC1 on inner aspect of the hand and arm.




- For JOBST Relax gauntlet, add 1 cm measured lengths AB and AC.

Measuring for JOBST | Gloves and Gauntlets

AD, AE Length (long glove, long gauntlet)

- Measure length AD and length AE on inner aspect of the hand and arm.



Measuring for JOBST | Gloves and Gauntlets

Thumb Circumference

- Measure the thumb circumference using a relaxed figure-8 measurement, at the base of the joint of the thumb (MCP).



MCP - Metacarpophalangeal


Measuring for JOBST | Gloves and Gauntlets

X and Z digit Circumferences

- Measure X circumference at the base of the finger.
- Measure Z circumference at the point of the desired finger length.
- The X circumference should not be less than the Z circumference.



- When the desired finger length is less than 1.5 cm, use circumference X for Z.
- Please note that the minimum finger length is 1 cm.
- Please avoid an ending on the joint.




Measuring for JOBST | Gloves and Gauntlets

Digit Length

- Measure the finger length using the JOBST measuring card.
- Place the measuring card between the fingers and measure the length on the side which is closest to the middle finger.
- Measure up to the nailbed to avoid swelling at the end.



- Long fingers are advised if there is swelling in the fingers.
- Please note that the minimum finger length is 1 cm.



JOBST Confidence | Upper extremities

Donning & Doffing videos

Tips for easy application



Tips for easy removal



JOBST Elvarex Plus* and JOBST Elvarex Soft Seamless Upper & Lower extremities

Styles



AC1
Glove



AC1
Gauntlet



AE
Long Glove



AE
Long Gauntlet



AA
Toe Cap

Compression classes

CCL 1

CCL 2

Wash and care information (1/2)

Please remind the user to keep in mind the following:



Washing the garments every day, or at least every other day is not only important for hygienic reasons. It also keeps them elastic and maintains their function.



There are a few things to keep in mind when washing the garments:

- Wash the garments together with garments of the same color.
- A laundry net can protect the garments from damage.
- Wash the garments at the temperature indicated in the Instructions For Use, please refer to the washing label.
- Choose a gentle wash cycle.
- Do not use a brightener or fabric softener.



How to dry garments after washing:

Some products are suitable for tumble-drying.

Garments should not be dried in the sunlight as this harms the fabric.

Wash and care information (2/2)

Please remind the user to keep in mind the following:



What else should be kept in mind?

- Long fingernails and jewelry can damage the fabric when putting garments on or off.
- Do not use body lotions or greasy skin care ointments when wearing the garments.
JOBST offers skin care products that can be used without damaging the garments.
- When any loose threads or loops are noticed, they shall not be cut off - this can cause holes or ladders in the garments. Instead, the medical supply store shall be contacted to get a professional repair, when possible.




In addition, always **advise to follow the Instructions For Use for the respective garments** because they can differ from product to product.




JOBST Confidence is a flat-knit compression garment for patients with mild to moderate lymphedema and lipedema without pronounced skin folds or shape distortions. In order for the garment to exert its optimal compression pressure, precise measuring is necessary for an accurate fit. Please find below the measurement instructions to be used together with JOBST Global Measurement Guide.

Measure length IK




- Place the JOBST foam pad as far as comfortably possible into the crotch, between the legs.
- Measure perpendicularly from the floor to the top edge of the foam pad.

Measure length IH




- Put the belt around the **widest part** of the buttocks
- Measure **perpendicularly** from the bottom to the top edge of the belt

Measure circumference cT




- The measuring belt should be at the waist and around the hips in a straight line form
- Too high on the belly part or too low at the back results in an inaccurate fit

Measure length cT



- Place the belt at the desired waist place (approx. belly button)
- Length T should be measured **perpendicularly** from the floor to point G. And then, following the **contour** from G to T.
- In other words, measuring tape should be placed in a straight line from floor to the hips and shall follow the contours from hip to the upper edge of the belt.

Measure length TT/K3



- Place the measuring tape with tape slider on the upper edge of the belt.
- The measuring tape shall now pass through the gluteal cleft (and not on buttocks) - to the upper edge of the belt on back. This is length TT.
- The patient slides the tape slider to the center of the crotch.
- Remove the tape from patient's body. The measurement read at the tape slider position is length K3.

Check out the video on JOBST Confidence measurement



JOBST Confidence garments should not be pulled up too high! Please use the size of the woven-inn logo as reference.



JOBST Logo
approx 3cm long



JOBST Logo
> 3cm long



Which patients are suitable for JOBST Confidence?*

Due to the softness of the material, some patients are more suitable for JOBST Confidence products and others are less suitable.

JOBST Confidence is suitable for patients who...

- ... have mild to moderate lymphedema and lipedema without pronounced skin folds or shape distortions.
- ... have fragile or tender to touch skin.
- ... seek comfort and are new to flat-knit compression.
- ... are sporty and active.

JOBST Confidence is not suitable for patients who ...

- ... want to feel a strong pressure.
- ... hide their belly or have body shaping trousers.
- ... have considerably different circumferences in one limb, if you wish to have two legged panty.
- ... have belly overhangs.
- ... have large adipose tissue lobules.

Note: In these cases, possibly JOBST Elvarex or JOBST Elvarex Soft may be suitable. Please feel free to contact us.

Measuring points	min./max. Meas	Notes
circumferences		
cG	max. 90cm	We recommend using JOBST Elvarex or JOBST Elvarex Soft for larger circumferences
cH	up to 140cm	We recommend using JOBST Elvarex or JOBST Elvarex Soft for larger circumferences.
Length		
ℓB	min. 7cm	A shorter length is unrealistic.
TT	max. 100cm	We recommend using JOBST Elvarex or JOBST Elvarex Soft for larger length measurements.
K ³	max. 45cm	K3 must always be less than TT and must not exceed 45% of the total length (TT).
Length btw two points		
ℓG to ℓK	up to 4cm	Longer measurements than 6cm to 7cm are uncommon and may result in more material on the leg and therefore wrinkling.
ℓK to ℓH	approx 10cm, min 4cm	Usual length is between 4cm - 12cm.
ℓH to ℓT	min. 7cm	The waist band needs about 4cm to 5cm. At least 2cm of knitted fabric must also be added.

* This recommendation is to help healthcare professionals to navigate through JOBST flat-knit assortment and choose a suitable product for patients based on material properties as well as the indications and precautions in the Instructions For Use. For further information please read the Instructions For Use.



Discover more at www.JOBST.com

Available On-Line 24 hours per day, 365 Days per Year

eshop.jobst-usa



[/JOBSTUSA](https://www.facebook.com/JOBSTUSA)



[@JOBST_USA](https://www.instagram.com/JOBST_USA)



[@JOBSTforUSA](https://twitter.com/JOBSTforUSA)



jobst-usa.com

BSN Medical Inc., an Essity company
5825 Carnegie Blvd.
Charlotte, NC 28209-4633
Tel. 1-800-537-1063 Fax 1-800-835-4325



JOBST®,
an Essity brand