

# FarrowWrap® Leg piece Application Techniques

Applicable for the following fabrics

JOBST® FarrowWrap® Lite

JOBST® FarrowWrap® Strong

JOBST® FarrowWrap® Classic



SCAN ME  
for application video



1  
Apply a JOBST® Farrow liner and footpiece if applicable.  
Unwrap the legpiece and lay the wrap on the floor behind your leg with the label side furthest back and facing up. The shortest bands will be closer to the foot.



2  
Place the wrap with its spine positioned behind the leg starting from the ankle and finishing just below the knee.

Apply the top wrap as an anchor to hold the wrap in place. Then commence wrapping from the ankle up towards the knee.



3  
Grab each side of the band closest to your ankle and pull forward wrapping the band around the lowest part of the leg at or near end stretch and affix with the VELCRO®. Undo the top band before applying the second last band to ensure the final band sits over the top. If wearing a footpiece, the bottom band should overlap the footpiece.



4  
Repeat with each band moving up the leg, smoothing out any wrinkles as you go.

Try to ensure a 50% overlap between adjacent straps.



JOBST®,  
an Essity brand



Always follow the  
directions for use.

ESS0776 V2  
05/2024

**Essity Health & Medical**  
PO Box 337, Mount Waverley,  
Victoria, Australia 3149  
[medical.essity.com.au](http://medical.essity.com.au)  
[jobst.com.au](http://jobst.com.au)

**JOBST® Customer Service AU**  
T 1300 998 810 F 1300 998 820  
**JOBST® Customer Service NZ**  
T 0508 998 810 F 0508 998 820