EN Synthetic Splint System

Product Description

Dynacast® Prelude is a pre-padded synthetic splint system.

Intended use and indications

Dynacast® Prelude splints and pre-cut splints are intended for external rigid immobilization through splinting procedures. The purposes for external rigid immobilization include:

- Immobilization of bone fractures.
- Treatment of soft tissue and joint injuries.

Directions for Use

The products are intended to be used by qualified healthcare professionals.

1. In case of dispenser box cut the length needed for the treatment (Pictogram 1+2).
2. In case of dispenser box push back splint in pouch and tightly seal the remaining roll with the closing device (Pictogram 3+4). Be careful to not puncture the pouch as this will lead to premature hardening of the material.
3. Remove splint from pouch.
4. Open padding and cut the splint with a scissor to ensure complete overlap of padding. Padding is stretchable and can be pulled for better cover (Pictogram 5).
5. Apply water, passing under cool water outlet or using water or spray bottle (Pictogram 6).
6. Roll splint and squeeze out excess water.
7. Lay splint on towel, rolling in towel and squeezing to remove additional water, repeating if necessary (Pictogram 7).
8. Place splint on patient and wrap with a support bandage (Pictogram 8).
9. Mold splint to limb for 2–4 minutes until splint sets. Monitor patient according to standard procedure. Initial set: 3–5 minutes. Functional set: approximately 15 minutes (Pictogram 9).
10. Do not cover the splint/ support bandage during setting (e.g. do not place on a pillow or cover the splint with clothes).
11. Instruct patient to keep the splint dry.

Warnings and Precautions

- Any wound must be dressed prior to application of the splint.
- Skin contact with uncured polymer must be avoided. Wear gloves and make sure the padding completely covers the splint before applying on patient.
- Use water not warmer than room temperature to avoid burns during setting reaction.
- The setting reaction is exothermic, instruct the patient that the splint might feel warm during setting.
- Apply extra padding if required, e.g. on bony prominences.
- Do not use if the package is already open or damaged.
- Do not reuse the device.

Cleaning and care

A healthcare professional must confirm that the splint may be removed for cleaning. The cured splint can be cleaned by hand washing with mild detergent. Before re-application the splint must be completely dried to avoid skin macerations. Proper drying techniques are wiping with a towel first and then using a hairdryer on cold setting for complete drying.

Storage conditions

It is recommended to store the product at room temperature (15–25 °C/ 59–77 °F). Short term exposition to higher temperatures is tolerated.

Notice

Any serious incident that has occurred in relation to this device should be reported to the BSN medical Inc. and the competent authority of your state.

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Essity UK Ltd · Southfields Road
BSN medical GmbH · Schützenstraße 1–3
Essity Distribution B.V., Piet Heinkade 55