(ENG) Bacteria binding dressing with hydrogel

Device Description

Cutimed® Sorbact® Gel Dressing is a sterile, bacteria and fungi binding dressing covered with gel. It consists of a green Sorbact® wound contact layer with a water-based gel containing carbomer and propylene glycol (10%). Cutimed® Sorbact® Gel Dressing donates moisture, enables a moist wound environment and facilitates autolytic debridement.

Mode of Action

Cutimed® Sorbact® binds microorganisms, such as Staphylococcus aureus (including MRSA), Streptococcus species, Escherichia coli, Pseudomonas aeruginosa and Candida albicans, as shown in vitro. These microorganisms are removed from the wound each time the dressing is changed.

Intended Use

Cutimed® Sorbact® Gel Dressing is intended for use in the management of clean, contaminated, colonized or infected wounds with dry to low exudate levels, such as traumatic wounds, burns, cavity wounds, fistulas, pressure ulcers, diabetic foot ulcers and leg ulcers.

Contraindications

There are no known contraindications to the use of Cutimed® Sorbact® Gel Dressing. Attention is drawn to the warnings and precautions.

Warnings and Precautions

If the treated condition deteriorates, fails to improve or a side effect is observed, consult a physician or an appropriate health care provider.

Do not use on patients with a known allergy to propylene glycol or to any other components of the dressing. Caution should be taken when exposing children below 5 years and patients suffering from liver and kidney disease, due to limited capacity to metabolize propylene glycol. High exposure, e.g. in case of simultaneous use of other products containing propylene glycol, may in rare cases lead to adverse events, e.g. hyposmolality or metabolic acidosis.

The dressing should always be applied in direct contact with the wound area to allow microorganisms to bind to the dressing. Should the dressing adhere to the wound, moisten or soak the dressing to assist removal and avoid disruption of the healing wound.

Do not apply Cutimed® Sorbact® Gel Dressing on wound edges, to avoid maceration of the surrounding skin.

Do not use in combination with oily products, such as ointments, creams and solutions as they may decrease the binding of microorganisms.

Cutimed® Sorbact® Gel Dressing is for single use only and should not be re-used. Re-use may lead to increased risk of infection or cross-contamination.

Do not re-sterilize.

Do not use if the peel pouch is already open or damaged, as the sterility of the dressing is guaranteed only when the package is unopened and undamaged prior to use.

Prior to radiation therapy, remove Cutimed® Sorbact® Gel Dressing. A new dressing can be applied following treatment.

Cutimed® Sorbact® Gel is magnetic resonance imaging (MRI) safe.

Instructions for use for wound management

- Prepare the wound according to local clinical practice. If desired, gently press a Cutimed® Sorbact® round swab against the moist wound surface to remove superficial microbial contamination.
- Select an appropriate dressing size for the wound.
- Remove the dressing from the pouch using an aseptic technique.
- If cutting the dressing, use an aseptic technique. Discard any open and unused dressing.
- Apply the dressing. Ensure that the dressing comes into direct contact with the complete wound surface. The dressing should not overlap the surrounding skin. Fold back any potentially overlapping material towards the wound to avoid maceration of the surrounding skin. If desired, add remaining hydrogel from the peel pouch onto the applied dressing.
- Direct wound contact is a prerequisite for the effective binding of microorganisms. In deep or cavity wounds, Cutimed® Sorbact® gel dressings should be applied in so as to ensure direct contact with the whole wound surface. Remove excess gel with a sterile swab if required.
- In cavity wounds and fistulas, always leave a part of the dressing protruding so that it can be easily removed.
- Apply an appropriate secondary dressing. Depending on the level of exudate, cover Cutimed® Sorbact® gel either with an absorbent secondary dressing or with a transparent film dressing.
- The dressing change frequency depends on exudate levels, overall condition of the wound and surrounding skin. Should the clinical condition allow, the dressing can be left in place for up to 7 days. The final decision on the therapy regimen needs to be taken by health care professionals.

Storage and Disposal

Dispose of waste according to local environmental procedures.

Notice

Any serious incident that has occurred in relation to this device should be reported to BSN medical GmbH and the competent authority of your state.

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 essity		Artwork ID: JBR12201.02	Review Loop: 3	Printing Colours: C M Y K
Product Title: Cutimed® Sorbact® Gel			Date of Review Loop: 09.08.2019	
Contact: Artwork.Management@essity.com			Replaces Artwork ID: JBR12201.01	
Contact Address: BSN medical GmbH Quickbornstrasse24 20253 Hamburg Germany			Technical Drawing for print ID:	
Ref Number(s): 72611-00; 72611-01; 72611-16			Comment on Diecut: Die cut for control purpose only - not for production purpose. Please use current die cut, as agreed with BSN medical. Technical proof, not binding in colour.	
Selling Countries:				