Gelostretch®

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Datum: 05/2022

EN Zinc oxide gel bandages for semi-rigid support and compression dressings.

Indications:

For compression bandages for phlebological indications such as chronic venous insufficiency, including: venous leg ulcer and phleboedema, thrombophlebitis, phlebotrombosis. For support and relief bandages for diseases of the ligament and postural system such as distorsions, contusions, luxations, tendovaginitis, follow-up treatment after plaster cast removal, reduction of post-traumatic oedemas.

Contraindications:

Advanced peripheral arterial occlusive disease, decompensated cardiac insufficiency, septic phlebitis, phlegmasia coerulea dolens. Children under the age of 3 years.

Cautions:

- Incompatibility to fabric
- Weeping dermatosis
- Cutaneous infections
- Open wounds must be covered with an appropriate dressing before compression bandaging.
- Severely compromised skin sensibility and impaired sensitivity of the limb
- Advanced peripheral neuropathy
- Conditions in which increased venous and lymphatic return is not desired.
- The attending physician decides whether compression in general can be applied to patients suffering from the following diseases: rheumatoid arthritis, complex regional pain syndrome, CRPS (M. Sudeck), malignant lymphedema, gangrene.
- If pain or complications, e.g. sustained discolouration of the toes, occur, remove the bandage immediately and contact your healthcare professional.

Composition of substrate:

Gelostretch®:	76% viscose, 24% cotton
Gelocast®:	100% cotton
Gelocast® elastic:	71% viscose, 29% polyamide

Composition of the mass:

Zinc oxide, water, gel-former, iodine, parabens, glycerin.

Elasticity:

Gelostretch®:	stretchable in both longitudinal and transverse directions
Gelocast®:	inelastic
Gelocast® elastic:	elastic in the longitudinal direction

Special information:

Storage temperature: do not store above 25°C. Ready-to-use impregnated, quick-drying zinc oxide gel bandages.

Notice:

Any serious incident that has occurred in relation to this device should be reported to BSN medical Inc. (BSN medical GmbH within the European Union) and the competent authority of your state.

After opening the packaging, the bandage has to be applied immediately in a moist state.

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