

# Actimove®

Feel good. Move better.



*Actimove®*

*Orthopedic  
Supports  
Catalog*





# Orthopedic Supports Catalog

## Actimove® – for you and your patients' needs

Actimove® offers a broad portfolio of medically engineered braces and supports for most needs and patient groups, coming in five distinctive ranges. Actimove® Professional Line, Actimove® Sports Edition, Actimove® Everyday Supports, Actimove® Arthritis Care and Actimove® Kids give patients the confidence to manage their pain and to take on their challenges of daily life.

Actimove's packing design brings a fresh look with high impact into stores. The design also supports easy product choice – for both professionals and those self-selecting.

Drawing on over 80 years of experience in non-invasive orthopedics<sup>1</sup>, Actimove® prides itself on providing medical devices that are reliable and thoroughly tested. Actimove® braces and supports help to manage pain with excellent functional performance based on biomechanics. They are tested for biocompatibility, are highly comfortable and made from carefully chosen materials. Driven by constant innovation, Actimove® keeps working on improving its offers for you.

Actimove® is a brand of Essity, a leading hygiene and health company dedicated to improving people's well-being. Actimove® stands for trusted products and services to support sustainable business growth.

Whether users are running a marathon, recovering from an injury, overcoming daily arthritis, or recurring pain when moving, they can rely on Actimove® to help them with optimized support and pain-relief.

**Actimove®. Feel good. Move better.**

<sup>1</sup> In the market since 1940: the leading Plaster of Paris brand Gypsona®

# TABLE OF CONTENTS

## Actimove® SPORTS EDITION

### Wrist & Thumb

Actimove® Wrist Support, Adjustable	11
Actimove® Wrist Support, Elastic Wrap Around	12
Actimove® Wrist Stabilizer, Removable Metal Stay, Left or Right	13
Actimove® Thumb Stabilizer, Extra Stays	14

### Knee

Actimove® Patella Strap, Adjustable	15
Actimove® Dual Knee Strap, Adjustable	16
Actimove® Knee Support, Closed Patella	17
Actimove® Knee Support, Open Patella	18
Actimove® Knee Support, Open Patella Adjustable	19
Actimove® Knee Stabilizer, Adjustable Horseshoe and Stays	20
Actimove® Knee Brace, Wrap Around, Simple Hinges, Condyle Pads	21
Actimove® Knee Brace, Wrap Around, Polycentric Hinges, Condyle Pads	22
Actimove® Knee Brace, Adjustable Horseshoe, Simple Hinges, Condyle Pads	23
Actimove® PF Knee Brace – Lateral Support, Simple Hinges, Condyle Pads, J-shaped Buttress	24

### Ankle

Actimove® Ankle Support	25
Actimove® Ankle Support, Adjustable	26
Actimove® Ankle Support, Elastic Wrap Around	27
Actimove® Ankle Stabilizer, Criss-Cross Straps	29

### Back & Lumbal

Actimove® Back Support, 4 Stays, Adjustable Double Layer Compression	30
Actimove® Back Stabilizer, Rigid Panel with Pressure Pads SMART Easy-Closing-Pulley-System	31

### Elbow

Actimove® Elbow Strap, Hot/Cold Pack	32
Actimove® Elbow Support, Adjustable	33

### Shoulder

Actimove® Shoulder Support, Extra Pocket for Optional Hot/Cold Pack	34
---	----

### Thigh

Actimove® Thigh Support, Adjustable	35
-------------------------------------	----

## Actimove® EVERYDAY SUPPORTS

### Wrist & Thumb

Actimove® Wrist Stabilizer Carpal, Pre-shaped Metal Stay	37
--	----

### Knee

Actimove® Mild Knee Support	38
Actimove® Knee Support, Closed Patella	39



Actimove® Knee Support, Closed Patella, 2 Stays	40
Actimove® Knee Support, Open Patella, 4 Stays	41
<b>Ankle</b>	
Actimove® Mild Ankle Support	42
Actimove® Ankle Support	43
<b>Back &amp; Lumbal</b>	
Actimove® Lumbar Sacral Support	44
Actimove® Back Support, High-Density Foam Panel, Adjustable Double Layer Compression	46
<b>Elbow</b>	
Actimove® Mild Elbow Support	47
Actimove® Elbow Support	48
Actimove® Elbow Support, Pressure Pads and Strap	49
<b>Actimove® ARTHRITIS CARE</b>	
<b>Wrist</b>	
Actimove® Wrist Support	52
<b>Hands</b>	
Actimove® Arthritis Gloves	53
<b>Knee</b>	
Actimove® Knee Support	54
<b>Ankle</b>	
Actimove® Ankle Support	55
<b>Elbow</b>	
Actimove® Elbow Support	56
<b>Actimove® KIDS</b>	
<b>Wrist</b>	
Actimove® Wrist Stabilizer, Removable Metal Stay, Left or Right	58
<b>Knee</b>	
Actimove® Knee Support, Open Patella	59
<b>Ankle</b>	
Actimove® Ankle Support, Slim fit sleeve for shoes	60
<b>Shoulders</b>	
Actimove® Mitella, Arm Sling	61
<b>Cervical</b>	
Actimove® Cervical, Cervical Collar	62
<b>Elbow</b>	
Actimove® Elbow Support, Adjustable	63

# TABLE OF CONTENTS

## Actimove® PROFESSIONAL LINE

### Wrist & Thumb

Actimove® Rhizo Forte, Thumb Orthosis	66
Actimove® ManuMotion, Wrist Support	68
Actimove® Manus Forte, Wrist Brace	70
Actimove® Wrist Splint 8"	72
Actimove® Wrist Splint 6"	73
Actimove® Wrist Brace, Suede Finish, Lace-Closure	74
Actimove® Manus, Wrist Stabilizer	75
Actimove® Manus Forte Plus, Wrist and Thumb Brace	76
Actimove® Wrist Splint, Abducted Thumb, Compression-Formed Composite	78
Actimove® Gauntlet, Wrist and Thumb Stabilizer	80

### Knee

Actimove® Genu Tri-Panel, Knee Immobilizer	82
Actimove® Genu Clips, Knee Immobilizer	84
Actimove® Genu Eco, Knee Immobilizer	84
Actimove® GenuMotion, Knee Support	86
Actimove® Knee Brace, Composite Polycentric Hinges	88

### Ankle

Actimove® TaloMotion, Ankle Support	90
Actimove® AchilloMotion, Achilles Tendon Support	92
Actimove® TaloBrace, Comfort Ankle Brace	94
Actimove® TaloBrace Air, Comfort Ankle Brace	95
Actimove® TaloCast Air, Ankle Brace	96
Actimove® TaloCast AirGel, Ankle Brace	97
Actimove® TaloCast, Ankle Brace	98
Actimove® Standard Walker, Ankle Walker	100
Actimove® Standard Air Walker, Ankle Walker	101
Actimove® Closed Shell Air Walker, Ankle Walker/ Actimove® Heel Wedges, For Closed Shell Walker	102

### Back & Lumbal

Actimove® Lumbar Sacral Support Comfort, Additional Support Belt	106
Actimove® LumbaMotion, Back Support with Pressure Pad	108
Actimove® LumbaMotion Lady, Back Support with Pressure Pad	110
Actimove® Deluxe Lumbar Sacral Support, Rigid Panel, For Women	112
Actimove® LO•PLUS™ Lumbar Support, Stabilizing Posterior Panel	114

### Shoulder

Actimove® Sling Comfort, Arm Sling – Single Pack	117
Actimove® Mitella Comfort, Comfort Arm Sling	118

Actimove® Umerus Comfort, Comfort Shoulder Immobilizer	120
Actimove® Shoulder Abduction Sling, Shoulder Immobilizer	122

### Cervical

Actimove® Cervical, Cervical Collar Soft Density, Narrow Width	124
Actimove® Cervical, Cervical Collar Soft Density, Regular Width	124
Actimove® Cervical Comfort, Comfort Cervical Collar	126

### Abdominal & Ribs

Actimove® Abdominal Binder, Three-Panel	130
Actimove® Abdominal Binder, Four-Panel	130
Actimove® Abdominal Binder Comfort, Soft Pad	132
Actimove® Rib Belt, Two-Panel, For Men	134
Actimove® Rib Belt, Two-Panel, For Women	135

### Elbow

Actimove® EpiMotion, Elbow Support	138
Actimove® EpiSport, Elbow Clasp, Adjustable	140

### Thigh & Calf

Actimove® PowerMotion, Thigh Support	142
Actimove® PowerMotion, Calf Support	143

### Speciality Supports

Actimove® Athletic Supporter	145
Actimove® Hernia Support, Support Belt	146
Actimove® Hernia Support, Support Brief	147

# Choose the right product for your patients.

## Actimove® SPORTS EDITION

Keep cool and dry for great comfort



- Aids pain relief after acute sports injuries
- For the active younger people
- To help prevent re-injuries
- Made with COOLMAX® AIR Technology

## Actimove® EVERYDAY SUPPORT

Stability for daily mobility



- Aids pain relief in chronic degenerative joint diseases
- Designed for those looking for stability
- Medical compression to help reduce edema
- Motion-Comfort-Zone for sensitive areas
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex



## Actimove® ARTHRITIS CARE

Soothing warmth for pain relief



- Aids pain relief for chronic osteoarthritis and rheumatoid arthritis
- Designed to aid in helping to increase mobility
- Ceramic yarns retain and reflect body heat
- Validated compression for support and stability
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.

## Actimove® KIDS

Especially for kids



- Aids pain relief after acute injuries after play and sports
- Designed for children
- Includes products with COOLMAX® AIR Technology
- Neoprene-free, suitable for those sensitive to latex, as this product is not made with natural rubber latex.

## Actimove® PROFESSIONAL LINE

Excellent functionality for pain relief



- Physician recommended or prescribed
- For professional treatment of acute injuries, post-op use and rehabilitation of chronic orthopedic issues
- Pain relief from functional design based on biomechanics
- Control and guide joint movement
- Support - Stabilize - Immobilize


# Actimove® SPORTS EDITION

Younger, sporty people often wear neoprene braces & supports after sports injuries or to prevent re-injury. A major complaint about these is that they are too hot to wear and cause excessive sweating.<sup>1</sup>

Actimove® Sports Edition aids pain relief and healing through balanced warmth and medical compression. They offer great comfort from excellent moisture management and skin-friendly materials.

- Made with **COOLMAX® AIR<sup>2</sup>** technology, Actimove® Sports Edition is engineered for quick drying and exceptional breathability to keep people cool and dry.
- The perforated **performance material** of Actimove® Sports Edition is neoprene-free, friendly to the skin, and suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Furthermore, **validated compression** properties for support and stability help reduce swelling and manage pain relief.

**COOLMAX®**  
PRO



3D Diamond Lining with COOLMAX® AIR



Perforated Performance Material



<sup>1</sup> Essity End-User Validation Research 10/2012; n= 400 OSG shoppers in USA

<sup>2</sup> COOLMAX® AIR is a trademark of The LYCRA company

# Actimove® Wrist Support

## Adjustable

### Durable, adjustable support for weak or injured wrists

- Pain relief and healing are supported through balanced warmth and medical compression.
- Great comfort with exceptional breathability through quick drying perforated performance material.
- Neoprene-free, made with COOLMAX® AIR technology.
- Adjust the level of support for comfort no matter the activity.
- Allows for comfortable, dry hands based on the open palm design.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Universal sizing.
- Can be worn left or right.
- For stock keeping, long shelf life of five years.

This product is intended to provide support to the wrist.


Fields of application include:

- Weak or injured wrist
- Overuse syndromes
- Mild sprains and strains
- Arthritis
- Support during everyday activities



3D Diamond Lining with COOLMAX® AIR



	Code	Size	Measurements	Color
	7562610	Universal	5 <sup>1</sup> / <sub>2</sub> " – 8" / 14 – 20 cm	Black

# Actimove® Wrist Support\*

## Elastic Wrap Around

**Firm and adjustable support, padded for comfort for aching and fatigued wrists**

- Enjoy sports, as pain relief and healing are supported through balanced warmth and medical compression.
- Great comfort with exceptional breathability through quick drying perforated performance material.
- Neoprene-free, made with COOLMAX® AIR technology.
- Firm and adjustable support from powerful elastic strap no matter the activity.
- Comfortable to wear because of terrycloth padding.
- Universal sizing.
- Can be worn left or right.
- Contains latex.
- For stock keeping, long shelf life of five years.

This product is intended to provide support to the wrist.

Fields of application include:

- Weak, aching and fatigued wrists
- Support during sports

\*Caution: This product contains natural rubber latex which may cause allergic reactions.



3D Diamond Lining with COOLMAX® AIR



	Code	Size	Measurements	Color
	7341680	Universal	5 1/2" - 8" / 14 - 20 cm	Black





# Actimove® Wrist Stabilizer

Removable Metal Stay, Left or Right

Anatomically shaped stay for secure stabilization, managing sprains and strains

- Pain relief and healing are supported through balanced warmth and medical compression.
- Comfort and good breathability through neoprene-free performance material.
- Stabilizes wrist through an anatomically shaped removable stay for left and right hand application.
- Quick and easy to fit with slip-on design and individually adjustable straps.
- Special comfort in a premium soft material and seamless design.
- Can be easily worn under clothing due to low profile.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Available in three sizes.
- Can be worn left or right.
- For stock keeping, long shelf life of five years.

This product is intended to provide support to the wrist.


Fields of application include:

- Tendonitis
- Sprains, strains, soft tissue injuries
- Osteoarthritis
- Rheumatoid arthritis
- Post-operative care
- After cast removal



Slip-On Design



	Code	Size	Measurements	Color
	7572931	S	5" - 6 <sup>1</sup> / <sub>4</sub> " / 13 - 16 cm	Black
	7572932	M	6 <sup>1</sup> / <sub>4</sub> " - 7 <sup>1</sup> / <sub>2</sub> " / 16 - 19 cm	Black
	7572933	L	7 <sup>1</sup> / <sub>2</sub> " - 9" / 19 - 23 cm	Black



# Actimove® Thumb Stabilizer

## Extra Stays

### Anatomically shaped stays, managing thumb injuries, keeping finger dexterity

- Pain relief and healing are supported through balanced warmth and medical compression.
- Comfort and good breathability through neoprene-free performance material.
- Stabilization of the thumb through two integrated anatomically shaped stays.
- Comfortable to wear because of premium quality soft material.
- Firm and adjustable wrist support from powerful elastic strap.
- Full finger movement and hand dexterity from special open design.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Available in two sizes.
- Comes in black or beige.
- For stock keeping, long shelf life of five years.

This product is intended to provide support/stabilization to the thumb and wrist.


**Fields of application include:**

- Mild to moderate strains and sprains (e.g., skier's thumb)
- Inflammation
- Onset of osteoarthritis
- Thumb pain from repetitive motion activities, like texting or gaming



Full Finger Movement



	Code	Color	Code	Color	Size	Measurements
	7563630	Black	7563632	Beige	S/M	5 <sup>1</sup> / <sub>2</sub> " – 7 <sup>1</sup> / <sub>2</sub> " / 14 – 19 cm
	7563631	Black	7563633	Beige	L/XL	7 <sup>1</sup> / <sub>2</sub> " – 9" / 19 – 23 cm

# Actimove® Patella Strap

## Adjustable

### Targeted patella pain relief that stays in place

- Enjoy sports, as pain relief and healing are supported through balanced warmth and medical compression.
- Great comfort with exceptional breathability through quick drying perforated performance material.
- Neoprene-free, made with COOLMAX® AIR technology.
- Targeted pain relief by increased pressure on patella tendon from the high density EVA shaped insert.
- Stays in place under activity due to silicone dot lining.
- Adjust the level of support, for comfort no matter the activity.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Universal sizing.
- Can be worn left or right.
- For stock keeping, long shelf life of five years.

This product is intended to provide support/stabilization to the patella including lateral tracking.

**Fields of application include:**

- Patella tendonitis
- Patellofemoral pain syndrome
- Chondromalacia patella
- Osgood-Schlatter disease



Silicone Dot Lining



	Code	Size	Measurements	Color
	7558920	Universal	11 <sup>1/2</sup> " - 16 <sup>1/3</sup> " / 29.5 - 41 cm	Black



# Actimove® Dual Knee Strap

## Adjustable

### Dual strap design for additional stability and patella pain relief

- Enjoy sports, as pain relief and healing are supported through balanced warmth and medical compression.
- Great comfort with exceptional breathability through quick drying perforated performance material.
- Neoprene-free, made with COOLMAX® AIR technology.
- Additional stability from dual strap design supporting the painful patella.
- Targeted pain relief by increased pressure on patella tendon from the high density EVA shaped insert.
- Adjust the level of support, for comfort no matter the activity.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Available in four sizes.
- Can be worn left or right.
- For stock keeping, long shelf life of five years.

This product is intended to provide support/stabilization to the patella including lateral tracking.

**Fields of application include:**

- Patella and quadriceps tendonitis
- Patellofemoral pain syndrome
- Chondromalacia patella
- Osgood-Schlatter disease



3D Diamond Lining with COOLMAX® AIR



	Code	Size	Measurements	Color
	7559110	S	13" - 14 <sup>1/2</sup> " / 33 - 37 cm	Black
	7559111	M	14 <sup>1/2</sup> " - 16" / 37 - 41 cm	Black
	7559112	L	16" - 17 <sup>3/4</sup> " / 41 - 45 cm	Black
	7559113	XL	17 <sup>3/4</sup> " - 19 <sup>1/4</sup> " / 45 - 49 cm	Black

# Actimove® Knee Support

## Closed Patella

### Slip-on style for secure protection for minor knee injuries

- Enjoy sports, as pain relief and healing are supported through balanced warmth and medical compression.
- Great comfort with exceptional breathability through quick drying perforated performance material.
- Neoprene-free, made with COOLMAX® AIR technology.
- Feel secure as knee is protected by a snug fit that stays in place.
- Fast & easy slip-on style.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Available in seven sizes.
- Available in black or navy
- For stock keeping, long shelf life of five years.

This product is intended to provide support to the knee.

**Fields of application include:**

- Painful knee injuries and prevention
- Slight strains, sprains and contusions
- Minor meniscal damages
- Joint swelling



Perforated Performance Material

3D Diamond Lining with COOLMAX® AIR



	Code	Color	Code	Color	Size	Measurements
	7558620	Navy	7558627	Black	XS	12" - 14" / 30.5 - 35.5 cm
	7558621	Navy	7558628	Black	S	14" - 16" / 35.5 - 40.5 cm
	7558622	Navy	7558629	Black	M	16" - 18" / 40.5 - 46 cm
	7558623	Navy	7558630	Black	L	18" - 20" / 46 - 51 cm
	7558624	Navy	7558631	Black	XL	20" - 22" / 51 - 56 cm
	7558625	Navy	7558632	Black	XXL	22" - 24" / 56 - 61 cm
	7558626	Navy	7558633	Black	XXXL	24" - 26" / 61 - 66 cm

# Actimove® Knee Support

## Open Patella

### Slip-on style for secure protection for minor knee injuries

- Enjoy sports, as pain relief and healing are supported through balanced warmth and medical compression.
- Great comfort with exceptional breathability through quick drying perforated performance material.
- Neoprene-free, made with COOLMAX® AIR technology.
- Helps to relieve patella pain. Knee cap pressure is reduced by open design.
- Feel secure as knee is protected by a snug fit that stays in place.
- Fast & easy slip-on style.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Available in seven sizes.
- Available in black or navy.
- For stock keeping, long shelf life of five years.

This product is intended to provide support to the knee.

Fields of application include:

- Painful knee injuries and prevention
- Slight strains, sprains and contusions
- Minor meniscal damages
- Joint swelling



Oval-Shaped Patella Opening

3D Diamond Lining with COOLMAX® AIR



	Code	Color	Code	Color	Size	Measurements	
	7558512	Navy	7558519	Black	XS	12" - 14" / 30.5 - 35.5 cm	
	7558513	Navy	7558520	Black	S	14" - 16" / 35.5 - 40.5 cm	
	7558514	Navy	7558521	Black	M	16" - 18" / 40.5 - 46 cm	
	7558515	Navy	7558522	Black	L	18" - 20" / 46 - 51 cm	
	7558516	Navy	7558523	Black	XL	20" - 22" / 51 - 56 cm	
	7558517	Navy	7558524	Black	XXL	22" - 24" / 56 - 61 cm	
	7558518	Navy	7558525	Black	XXXL	24" - 26" / 61 - 66 cm	
				7558548	Black	XXXXL	26" - 28" / 66 - 71 cm



# Actimove® Knee Support

## Open Patella Adjustable

### Adjustable and secure protection, relieving patella pain in minor knee injuries

- Enjoy sports, as pain relief and healing are supported through balanced warmth and medical compression.
- Great comfort with exceptional breathability through quick drying perforated performance material.
- Neoprene-free, made with COOLMAX® AIR technology.
- Supports ease of application with wrap-around style for sore knees.
- Adjust the level of support for comfort no matter the activity.
- Helps to relieve patella pain. Knee cap pressure is reduced by open design.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Universal sizing.
- Can be worn left or right.
- For stock keeping, long shelf life of five years.

This product is intended to provide support to the knee.


Fields of application include:

- Painful knee injuries and prevention.
- Slight strains, sprains and contusions
- Minor meniscal damages
- Joint swelling



Oval-Shaped Patella Opening      3D Diamond Lining with COOLMAX® AIR



	Code	Size	Measurements	Color
	7559310	Universal	11 <sup>1</sup> / <sub>2</sub> " – 16 <sup>1</sup> / <sub>8</sub> " / 29.5 – 41 cm	Black
	7559313	Universal Plus	12 <sup>1</sup> / <sub>4</sub> " – 17" / 31 – 43 cm	Black

# Actimove® Knee Stabilizer

## Adjustable Horseshoe and Stays

### Sleeve design with adjustable horseshoe for improved patellar tracking

- Enjoy sports, as pain relief and healing are supported through balanced warmth and medical compression.
- Great comfort with exceptional breathability through quick drying perforated performance material.
- Neoprene-free, made with COOLMAX® AIR technology.
- Manage patellar pain syndrome with an adjustable horseshoe.
- Added knee support from medial and lateral spiral stays.
- Stays in place with improved fit through elastic hook and loop support strap.
- Easy application with convenient finger-loops.
- Improved patellar tracking through increased pressure by adjusting the horseshoe either inferiorly, superiorly, medially or laterally.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Available in seven sizes.
- Can be worn left or right.
- For stock keeping, long shelf life of five years.


This product is intended to provide support/stabilization to the patella including lateral tracking.

**Fields of application include:**

- Patellar pain syndrome
- Patellar tracking and subluxation
- Patellar instability
- Chondromalacia patella
- Patellar tendonitis



	Code	Size	Measurements	Color
	7245300	XS	12" - 14" / 30.5 - 35.5 cm	Black
	7245301	S	14" - 16" / 35.5 - 40.5 cm	Black
	7245302	M	16" - 18" / 40.5 - 46 cm	Black
	7245303	L	18" - 20" / 46 - 51 cm	Black
	7245304	XL	20" - 22" / 51 - 56 cm	Black
	7245305	XXL	22" - 24" / 56 - 61 cm	Black
	7245306	XXXL	24" - 26" / 61 - 66 cm	Black
	7245307	XXXXL	26" - 28" / 66 - 71 cm	Black



# Actimove® Knee Brace

Wrap Around, Simple Hinges, Condyle Pads

Front wrap closure system supports application, simple hinges for extra stability

- Pain relief and healing are supported through balanced warmth and medical compression.
- Great comfort with exceptional breathability through quick drying perforated performance material.
- Neoprene-free, made with COOLMAX® AIR technology.
- Manage knee pain from soft tissue, ligament injuries or osteoarthritis through extra stability from medial and lateral hinges.
- Supports application through front wrap closure system.
- Stays in place with improved fit through elastic hook and loop support strap.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Available in seven sizes.
- Can be worn left or right.
- For stock keeping, long shelf life of five years.

This product is intended to stabilize the knee.

Fields of application include:

- Light to moderate ACL, LCL, MCL, PCL sprains and tears
- Medial / lateral meniscal tears
- Osteoarthritis (OA)
- Knee joint instability



Removable Condyle Pads



	Code	Size	Measurements	Color
	7550010	XS	12" - 14" / 30.5 - 35.5 cm	Black
	7550011	S	14" - 16" / 35.5 - 40.5 cm	Black
	7550012	M	16" - 18" / 40.5 - 46 cm	Black
	7550013	L	18" - 20" / 46 - 51 cm	Black
	7550014	XL	20" - 22" / 51 - 56 cm	Black
	7550015	XXL	22" - 24" / 56 - 61 cm	Black
	7550016	XXXL	24" - 26" / 61 - 66 cm	Black

# Actimove® Knee Brace

Wrap Around, Polycentric Hinges, Condyle Pads

Front wrap closure system supports application, polycentric hinges for improved gait

- Pain relief and healing are supported through balanced warmth and medical compression.
- Great comfort with exceptional breathability through quick drying perforated performance material.
- Neoprene-free, made with COOLMAX® AIR technology.
- Manage knee pain from soft tissue, ligament injuries or osteoarthritis through extra stability from medial and lateral hinges.
- Supports application through front wrap closure system.
- Stays in place with improved fit through elastic hook and loop support straps.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Available in seven sizes.
- Can be worn left or right.
- For stock keeping, long shelf life of five years

This product is intended to stabilize the knee.


Fields of application include:

- Light to moderate ACL, LCL, MCL, PCL sprains and tears
- Medial / lateral meniscal tears
- Osteoarthritis (OA)
- Knee joint instability

Removable Condyle Pads



COOLMAX®  
PRO

	Code	Size	Measurements	Color
	7311800	XS	12" - 14" / 30.5 - 35.5 cm	Black
	7311801	S	14" - 16" / 35.5 - 40.5 cm	Black
	7311802	M	16" - 18" / 40.5 - 46 cm	Black
	7311803	L	18" - 20" / 46 - 51 cm	Black
	7311804	XL	20" - 22" / 51 - 56 cm	Black
	7311805	XXL	22" - 24" / 56 - 61 cm	Black
	7311806	XXXL	24" - 26" / 61 - 66 cm	Black
	7311828	XXXXL	26" - 28" / 66 - 71 cm	Black

# Actimove® Knee Brace

Adjustable Horseshoe, Simple Hinges, Condyle Pads

**Hinged knee sleeve for extra stability, adjustable horseshoe for improved patellar tracking**

- Pain relief and healing are supported through balanced warmth and medical compression.
- Great comfort with exceptional breathability through quick drying perforated performance material.
- Neoprene-free, made with COOLMAX® AIR technology.
- Manage knee pain from soft tissue and ligament injuries through extra stability with medial and lateral hinges.
- Optional patellar support with an adjustable horseshoe.
- Stays in place with improved fit through elastic hook and loop support strap.
- Easy application with convenient finger-loops.
- Improved patellar tracking through increased pressure by adjusting the horseshoe either inferiorly, superiorly, medially or laterally.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Available in seven sizes.
- Can be worn left or right.
- For stock keeping, long shelf life of five years.

**This product is intended to stabilize the knee.**

**Fields of application include:**

- ACL, LCL, MCL, PCL sprains and tears
- Medial/lateral meniscal tears
- Patellar instability
- Patellar tendonitis
- Chondromalacia patella
- Osgood-Schlatter disease



	Code	Size	Measurements	Color
	7568750	XS	12" - 14" / 30.5 - 35.5 cm	Black
	7568751	S	14" - 16" / 35.5 - 40.5 cm	Black
	7568752	M	16" - 18" / 40.5 - 46 cm	Black
	7568753	L	18" - 20" / 46 - 51 cm	Black
	7568754	XL	20" - 22" / 51 - 56 cm	Black
	7568755	XXL	22" - 24" / 56 - 61 cm	Black
	7568756	XXXL	24" - 26" / 61 - 66 cm	Black
	7568757	XXXXL	26" - 28" / 66 - 71 cm	Black



# Actimove® PF Knee Brace – Lateral Support

Simple Hinges, Condyle Pads, J-shaped Buttress

**Sleeve design with special lateral patellar femoral support, simple hinges for extra stability**

- Pain relief and healing are supported through balanced warmth and medical compression.
- Great comfort with exceptional breathability through quick drying perforated performance material.
- Neoprene-free, made with COOLMAX® AIR technology.
- Lateral patellar support and stabilization using special angled straps combined with the "J" shaped buttress.
- Supports side-to-side knee stability through medial and lateral hinges.
- Easy application with convenient finger-loops.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Available in seven sizes.
- For stock keeping, long shelf life of five years.

**This product is intended to stabilize the knee.**

**Fields of application include:**

- Chondromalacia patella
- Lateral pressure syndrome pain
- Patellar instability
- Lateral dislocation
- Lateral patellofemoral mal alignment

J-Shaped Buttress



	Code	Side	Code	Side	Size	Measurements	Color
	7559930	Right	7559937	Left	XS	12" – 14" / 30.5 – 35.5 cm	Black
	7559931	Right	7559938	Left	S	14" – 16" / 35.5 – 40.5 cm	Black
	7559932	Right	7559939	Left	M	16" – 18" / 40.5 – 46 cm	Black
	7559933	Right	7559940	Left	L	18" – 20" / 46 – 51 cm	Black
	7559934	Right	7559941	Left	XL	20" – 22" / 51 – 56 cm	Black
	7559935	Right	7559942	Left	XXL	22" – 24" / 56 – 61 cm	Black
	7559943	Right	7559936	Left	XXXL	24" – 26" / 61 – 66 cm	Black



# Actimove® Ankle Support

## Slim fit sleeve for shoes

- Enjoy sports, as pain relief and healing are supported through balanced warmth and medical compression.
- Great comfort with exceptional breathability through quick drying perforated performance material.
- Neoprene-free, made with COOLMAX® AIR technology.
- Fast & easy slip on style.
- Comfortable slim-fit with flat seams for shoes.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Available in four sizes.
- Can be worn left or right.
- Available in black or navy.
- For stock keeping, long shelf life of five years.

This product is intended to provide support to the ankle.


**Fields of application include:**

- Painful ankle injuries & prevention
- Mild sprains
- Weak ankles
- Swelling and tenderness



Perforated Performance Material



	Code	Color	Code	Color	Size	Measurements
	7560612	Navy	7560616	Black	S	8" - 9" / 20 - 23 cm
	7560613	Navy	7560617	Black	M	9" - 10 <sup>1/4</sup> " / 23 - 26 cm
	7560614	Navy	7560618	Black	L	10 <sup>1/4</sup> " - 11 <sup>1/2</sup> " / 26 - 29 cm
	7560615	Navy	7560619	Black	XL	11 <sup>1/2</sup> " - 12 <sup>1/2</sup> " / 29 - 32 cm
				7560626	Black	XXL

# Actimove® Ankle Support

## Adjustable

### Stays in place through conforming heel loop, adjustable for stability

- Enjoy sports, as pain relief and healing are supported through balanced warmth and medical compression.
- Comfort and good breathability through neoprene-free performance material.
- Stays in place with conforming elastic heel loop.
- Adjust the level of support, for comfort no matter the activity.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Universal sizing.
- Can be worn left or right.
- For stock keeping, long shelf life of five years.

This product is intended to provide support to the ankle.

**Fields of application include:**

- Painful ankle injuries & prevention
- Mild sprains
- Weak ankles
- Swelling and tenderness



Elastic Heel Loop



	Code	Size	Measurements	Color
	7562010	Universal	8" - 10 <sup>1/2</sup> " / 20 - 27 cm	Black



# Actimove® Ankle Support\*

## Elastic Wrap Around

### Ultra slim fit, extra thin, yet durable, adjustable for stability

- Enjoy sports, as pain relief and healing are supported through balanced warmth and medical compression.
- Increased comfort and good breathability with ventilated elastic strap design.
- Ultra slim fit due to extra thin but durable material.
- Adjust the level of support for comfort, no matter the activity.
- Available in three sizes.
- Can be worn left or right.
- Contains latex.
- Available in black.
- For stock keeping, long shelf life of five years.


This product is intended to provide support to the ankle.

**Fields of application include:**

- Painful ankle injuries & prevention
- Mild sprains
- Weak ankles
- Swelling and tenderness

\*Caution: This product contains natural rubber latex which may cause allergic reactions.



	Code	Size	Measurements	Color
	7341480	S	7" - 8" / 18 - 20.5 cm	Black
	7341481	M	8 <sup>1</sup> / <sub>4</sub> " - 9 <sup>1</sup> / <sub>2</sub> " / 21 - 24 cm	Black
	7341482	L	9 <sup>1</sup> / <sub>2</sub> " - 10 <sup>3</sup> / <sub>4</sub> " / 24 - 27 cm	Black
	7341492	XL	10 <sup>3</sup> / <sub>4</sub> " - 12" / 27 - 30.5 cm	Black





# Actimove® Ankle Stabilizer

## Criss-Cross Straps

### Simulating sports taping, increased stability for ankle and arch of foot

- Enjoy sports, as pain relief and healing are supported through balanced warmth and medical compression.
- Great comfort with exceptional breathability through quick drying performance material.
- Neoprene-free, made with COOLMAX® AIR technology.
- Increased stability around the ankle and arch of foot with criss-cross straps simulating sports taping.
- Stays in place during athletic activity with extra support from circular strap.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Universal sizing.
- Can be worn left or right.
- For stock keeping, long shelf life of five years.

This product is intended to provide support/stabilization to the ankle.


**Fields of application include:**

- Mild and moderate ankle sprains
- Sports injuries
- Helps to avoid re-injuries and chronic ankle instabilities



Criss-Cross Straps



	Code	Size	Measurements	Color
	7561130	Universal	8" - 10 <sup>1/2</sup> " / 20 - 27 cm	Black
	7561132	Universal Plus	10 <sup>1/2</sup> " - 14" / 27 - 30 cm	Black

# Actimove® Back Support, 4 Stays

## Adjustable Double Layer Compression

### Advanced support and pain relief with comfort

- Pain relief and healing are supported through balanced warmth and medical compression.
- Great comfort with exceptional breathability through quick drying perforated performance material.
- Neoprene-free, made with COOLMAX® AIR technology.
- Advanced support and pain relief from adjustable, circumferential double-layer compression.
- Comfortable to wear, even when sitting, as flexible stays maintain positioning.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Available in five sizes.
- For stock keeping, long shelf life of five years.

This product is intended to provide support and stabilization to the lumbar joints of the back.

**Fields of application include:**

- Acute and chronic back pain
- Lumbar strains
- Muscular insufficiency
- Provides support during sports



Adjustable Elastic Band

3D Diamond Lining with COOLMAX® AIR



	Code	Size	Measurements	Color
	7554120	S	28" – 33" / 71 – 84 cm	Black
	7554121	M	33" – 37" / 84 – 94 cm	Black
	7554122	L	37" – 41" / 94 – 104 cm	Black
	7554123	XL	41" – 45" / 104 – 114 cm	Black
	7554124	XXL	45" – 49" / 114 – 124 cm	Black
	7554125	XXXL	49" – 57" / 124 – 144 cm	Black

# Actimove® Back Stabilizer

Rigid Panel with Pressure Pads  
 SMART Easy-Closing-Pulley-System

Easy to apply and adjust, relieving massage for lower back pain

- Pain relief and healing are supported through balanced warmth and medical compression.
- Relieving massage effect through prominent, removable pressure pads.
- Comfort and good breathability through neoprene-free perforated performance material.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Easy and simple to apply and adjust the compression with just one thumb using the SMART Easy-Closing-Pulley-System.
- Available in four sizes.
- Additional strength and stability for the lumbar region through rigid back panel.
- For stock keeping, long shelf life of five years.

This product is intended to provide support and stabilization to the lumbar joints of the back.

Fields of application include:

- Acute and chronic lower back pain
- Degenerative back problems (e.g., spondylosis)
- Lumbar instabilities (e.g., muscular insufficiencies)
- Post-operative care and rehabilitation



Removable Pressure Pads



	Code	Size	Measurements	Color
	7553510	XS/S	29" - 35" / 74 - 89 cm	Black
	7553511	S/M	35" - 41" / 89 - 104 cm	Black
	7553512	L/XL	41" - 45" / 104 - 114 cm	Black
	7553513	2XL/3XL	45" - 49" / 114 - 124 cm	Black

# Actimove® Elbow Strap

## Hot/Cold Pack

### Extra pain relief and secure fit for tennis and golfers elbows

- Enjoy sports, as pain relief and healing are supported through balanced warmth and medical compression.
- Comfort and good breathability through neoprene-free performance material.
- Extra pain relief supported by a removable therapeutic hot/cold pack.
- Silicone dot lining for secure fit.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Universal sizing.
- Can be worn left or right.
- For stock keeping, long shelf life of five years.

This product is intended to support the elbow.

Fields of application include:

- Tennis, golfer elbow (epicondylitis)



Hot/Cold Pack



Silicone Dot Lining



Code	Size	Measurements	Color
7574230	Universal	11" - 15 <sup>3/4</sup> " / 28 - 40 cm	Black





# Actimove® Elbow Support

Adjustable

**Slip-on design with easy application and adjustment**

- Enjoy sports, as pain relief and healing are supported through balanced warmth and medical compression.
- Comfort and good breathability through neoprene-free performance material.
- Quick and easy to fit with slip-on design and individually adjustable straps.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Universal sizing.
- Can be worn left or right.
- For stock keeping, long shelf life of five years.

This product is intended to support the elbow.

Fields of application include:

- Painful elbow injuries & prevention



Slip-On Design



	Code	Size	Measurements	Color
	7561730	Universal	7 <sup>1</sup> / <sub>2</sub> " - 14 <sup>1</sup> / <sub>2</sub> " / 19 - 37 cm	Black
	7561734	Universal Plus	12 <sup>1</sup> / <sub>4</sub> " - 17" / 31 - 43 cm	Black



# Actimove® Shoulder Support

## Extra Pocket for Optional Hot/Cold Pack

### Adjustable support and protection for shoulder injuries and painful conditions

- Pain relief and healing are supported through balanced warmth and medical compression.
- Comfort and good breathability through neoprene-free performance material.
- Extra pain relief when optional hot/cold therapy pack is placed in build-in pocket.
- Supports and protects the shoulder while allowing flexibility.
- Adjustable chest and arm closure for a custom fit.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Available in three sizes.
- Can be worn left or right.
- For stock keeping, long shelf life of five years.

This product is intended to provide support to shoulder.


**Fields of application include:**

- Slight shoulder instabilities
- Degenerative shoulder conditions
- Tendomyopathies
- Shoulder sprains / strains
- After arthroscopy



Pocket for Optional Hot/Cold Pack



	Code	Size	Measurements	Color
	7557310	S	31 <sup>1/2</sup> " - 41 <sup>1/2</sup> " / 80 - 105 cm	Black
	7557311	M	41 <sup>1/2</sup> " - 51" / 105 - 130 cm	Black
	7557312	L	51" - 61" / 130 - 155 cm	Black

# Actimove® Thigh Support

Adjustable

**Adjustable, accommodating fit that stays in place for strains and contusions**

- Enjoy sports, as pain relief and healing are supported through balanced warmth and medical compression.
- Great comfort with exceptional breathability through quick drying performance material.
- Neoprene-free, made with COOLMAX® AIR technology.
- Feel secure as thigh is protected by an accommodating fit that stays in place.
- Adjust the level of support for comfort, no matter the activity.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Universal sizing.
- Can be worn left or right.
- For stock keeping, long shelf life of five years.

**This product is intended for thigh support.**

**Fields of application include:**

- For slight to moderate strains & contusions of the quadriceps or hamstring muscles



3D Diamond Lining with COOLMAX® AIR



Code	Size	Measurements	Color
7560202	Universal	17 <sup>3</sup> / <sub>4</sub> " - 31 <sup>1</sup> / <sub>2</sub> " / 45 - 80 cm	Black



# Actimove® EVERYDAY SUPPORTS

With age, patients frequently suffer from chronic degenerative joints diseases, impacting how they are able move. They are looking for stability to support their mobility.

Actimove® Everyday Supports provides medical compression that helps to reduce edema and aids pain relief. Its enduring comfort is combined in a range of reliable and durable supports.

- **Validated compression** properties provide support and stability which helps to reduce swelling and manage pain relief.
- The **unique helix design** is inspired by physiological taping and the supports are anatomically shaped for stability.
- Made from soft yarns, Actimove® Everyday Supports have special **comfort zones** for sensitive areas. Further, the products are suitable for those sensitive to latex, as this product is not made with natural rubber latex.



# Actimove® Wrist Stabilizer Carpal

## Pre-shaped Metal Stay

### Anatomical shape for correct positioning, full finger movement and hand dexterity

- Pain relief and healing are supported through balanced warmth and medical compression.
- Helps to manage pain by stabilizing wrist with anatomically shaped stay for correct positioning.
- Provides relief as hand can rest on palmar pad supported by rigid, angled stay.
- Full finger movement and hand dexterity from unique open design.
- Great comfort with exceptional breathability through quick drying perforated performance material.
- Neoprene-free, made with COOLMAX® AIR technology.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Universal sizing.
- Can be worn left or right.
- For stock keeping, long shelf life of five years.

This product is intended for wrist support.

Fields of application include:

- Symptomatic relief of carpal tunnel syndrome, e.g., while using a computer



Open Design for Dexterity



Code	Size	Measurements	Color
7571420	Universal	5 1/2" - 8" / 14 - 20 cm	Black

# Actimove® Mild Knee Support

Lightweight design, easy to apply


- Mild support and easy application due to lightweight design.
- Knitted design supports good fit while maintaining flexibility.
- Soft yarns for comfort.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Available in four sizes.
- Can be worn left or right.
- For stock keeping, long shelf life of five years.



This product is intended to provide support to the knee.

Fields of application include:

- Slightly swollen and aching knees

	Code	Size	Measurements	Color
	7558020	S	11 <sup>3</sup> / <sub>4</sub> " - 15" / 30 - 38 cm	White
	7558021	M	15" - 18" / 38 - 46 cm	White
	7558022	L	18" - 20 <sup>3</sup> / <sub>4</sub> " / 46 - 53 cm	White
	7558023	XL	20 <sup>3</sup> / <sub>4</sub> " - 24" / 53 - 61 cm	White



# Actimove® Knee Support

## Closed Patella

**Firm compression, reliable, durable support for everyday mobility, functional helix design**

- Firm compression to help manage pain from knee conditions.
- Reliable, durable support for everyday mobility.
- Seamless anatomical shape supported by a patterned helix design for a good fit.
- Soft, smooth Motion-Comfort-Zone in the popliteal area.
- Reinforced edges for long lasting wear.
- Soft and cotton-free yarns for enduring comfort.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Available in four sizes.
- Can be worn left or right.
- For stock keeping, long shelf life of five years.

**This product is intended to provide support to the knee.**

**Fields of application include:**

- Degenerative knee conditions
- Overuse or repetitive use syndromes
- Chronic knee pain



Helix Design

Motion-Comfort-Zone



	Code	Size	Measurements	Color
	7557536	S	12 <sup>1</sup> / <sub>4</sub> " – 14 <sup>1</sup> / <sub>4</sub> " / 31 – 36 cm	Beige
	7557537	M	14 <sup>1</sup> / <sub>4</sub> " – 16 <sup>1</sup> / <sub>4</sub> " / 36 – 41 cm	Beige
	7557538	L	16 <sup>1</sup> / <sub>4</sub> " – 18" / 41 – 46 cm	Beige
	7557539	XL	18" – 20 <sup>1</sup> / <sub>2</sub> " / 46 – 52 cm	Beige

# Actimove® Knee Support

Closed Patella, 2 Stays

**Firm compression, added stability, durable support, functional helix design**

- Firm compression to help manage pain from knee conditions.
- Added stability through two spiral stays.
- Reliable, durable support for everyday mobility.
- Anatomical shape supported by a patterned helix design for a good fit.
- Soft, smooth Motion-Comfort-Zone in the popliteal area.
- Reinforced edges for long lasting wear.
- Soft and cotton-free yarns for enduring comfort.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Available in four sizes.
- Can be worn left or right.
- For stock keeping, long shelf life of five years.

**This product is intended to provide support to the knee.**

**Fields of application include:**


- Degenerative knee conditions
- Overuse or repetitive use syndromes
- Chronic knee pain



Helix Design

Motion-Comfort-Zone



	Code	Size	Measurements	Color
	7557528	S	12 <sup>1</sup> / <sub>4</sub> " - 14 <sup>1</sup> / <sub>4</sub> " / 31 - 36 cm	Beige
	7557529	M	14 <sup>1</sup> / <sub>4</sub> " - 16 <sup>1</sup> / <sub>4</sub> " / 36 - 41 cm	Beige
	7557530	L	16 <sup>1</sup> / <sub>4</sub> " - 18" / 41 - 46 cm	Beige
	7557531	XL	18" - 20 <sup>1</sup> / <sub>2</sub> " / 46 - 52 cm	Beige

# Actimove® Knee Support

## Open Patella, 4 Stays

**Firm compression, extra stability, reduced patella pressure, durable support, functional helix design**

- Firm compression to help manage pain from knee conditions.
- Helps to relieve patella pain. Knee cap pressure is reduced by open design.
- Added stability through four spiral stays.
- Reliable, durable support for everyday mobility.
- Anatomical shape supported by a patterned helix design for a good fit.
- Soft, smooth Motion-Comfort-Zone in the popliteal area.
- Reinforced edges for long lasting wear.
- Soft and cotton-free yarns for enduring comfort.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Available in four sizes.
- Can be worn left or right.
- For stock keeping, long shelf life of five years.

This product is intended to provide support to the knee.

Fields of application include:


- Degenerative knee conditions
- Overuse or repetitive use syndromes
- Chronic knee pain



Motion-Comfort-Zone

Pattern Design



	Code	Size	Measurements	Color
	7557520	S	12 <sup>1</sup> / <sub>4</sub> " – 14 <sup>1</sup> / <sub>4</sub> " / 31 – 36 cm	Beige
	7557521	M	14 <sup>1</sup> / <sub>4</sub> " – 16 <sup>1</sup> / <sub>4</sub> " / 36 – 41 cm	Beige
	7557522	L	16 <sup>1</sup> / <sub>4</sub> " – 18" / 41 – 46 cm	Beige
	7557523	XL	18" – 20 <sup>1</sup> / <sub>2</sub> " / 46 – 52 cm	Beige

# Actimove® Mild Ankle Support

Lightweight design, easy to apply

- Mild support and easy application due to lightweight open heel design.
- Knitted design supports good fit while maintaining flexibility.
- Soft yarns for comfort.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Available in four sizes.
- Can be worn left or right.
- For stock keeping, long shelf life of five years.

This product is intended to provide support to the ankle.


Fields of application include:

- Slightly swollen and aching ankles.



Open Heel Design



	Code	Size	Measurements	Color
	7560320	S	7 <sup>3</sup> / <sub>4</sub> " - 9 <sup>1</sup> / <sub>2</sub> " / 20 - 24 cm	White
	7560321	M	9 <sup>1</sup> / <sub>2</sub> " - 11" / 24 - 28 cm	White
	7560322	L	11" - 12 <sup>1</sup> / <sub>2</sub> " / 28 - 32 cm	White
	7560323	XL	12 <sup>1</sup> / <sub>2</sub> " - 14 <sup>1</sup> / <sub>2</sub> " / 32 - 37 cm	White

# Actimove® Ankle Support

Firm compression, reliable, durable support for everyday mobility

- Firm compression to help manage pain from ankle conditions.
- Reliable, durable support for everyday mobility.
- Seamless anatomical shape supported by a patterned helix design for a good fit.
- Reinforced edges for long lasting wear.
- Soft, smooth Motion-Comfort-Zone in the instep.
- Soft and cotton-free yarns for enduring comfort.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Available in four sizes.
- Can be worn left or right.
- For stock keeping, long shelf life of five years.

This product is intended to provide support to the ankle.


Fields of application include:

- Degenerative ankle conditions
- Overuse or repetitive use syndromes
- Chronic ankle pain



Pattern Design



	Code	Size	Measurements	Color
	7560820	S	7" - 8 <sup>1</sup> / <sub>4</sub> " / 18 - 21 cm	Beige
	7560821	M	8 <sup>1</sup> / <sub>4</sub> " - 9 <sup>1</sup> / <sub>2</sub> " / 21 - 24 cm	Beige
	7560822	L	9 <sup>1</sup> / <sub>2</sub> " - 10 <sup>3</sup> / <sub>4</sub> " / 24 - 27 cm	Beige
	7560823	XL	10 <sup>3</sup> / <sub>4</sub> " - 11 <sup>3</sup> / <sub>4</sub> " / 27 - 30 cm	Beige



# Actimove® Lumbar Sacral Support

Low-profile support and compression to the lumbar region with asymmetrical wings for a contoured fit.

- Helps to reduce pain through support and compression to the lower back area.
- Contoured fit through its asymmetrical wings.
- Comfortable to wear even when sitting as flexible stays help to maintain positioning of the belt.
- Easy application and removal through wrap-around design.
- Enhances proprioception by supporting the back to keep an upright posture.
- Customized fit due to wide range of adjustability in a wrap-around design.
- Smooth, durable surface helps to prevent pilling.
- Suitable for daily routines due to its low-profile, anatomical, and breathable design which can be worn under most clothing.
- Enhanced wearing comfort thanks to its breathable and soft materials.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Easy to clean for better hygiene, machine washable at a gentle 40 °C / 105 °F wash.
- Suitable for both shorter and taller people as available in two heights.
- Unisex, fits both men and women.
- Available in six sizes.
- For stock keeping, long shelf life of five years.

**This product is intended to provide support and stabilization to the lumbar area.**

**Fields of application include:**

- Acute and chronic lower back pain
- Degenerative changes in the lumbar area (e.g., osteoporosis, osteochondrosis, spondylosis, spondylarthrosis, facet syndrome)
- Lumbar instabilities (insufficiencies of muscles and ligaments, sprains and strains)
- Lumbar disc protrusion
- Post-traumatic and post-operative use (e.g., after intervertebral disc operations)
- Prevention of repetitive strain injuries





	Code	Size	Measurements		Unit
			Circumference	Height	
	7552714	S	27 <sup>1</sup> / <sub>2</sub> " – 31 <sup>1</sup> / <sub>2</sub> " 70 – 80 cm	8" 20 cm	1 Lumbar Support
	7552715	M	31 <sup>1</sup> / <sub>2</sub> " – 35 <sup>1</sup> / <sub>2</sub> " 80 – 90 cm	8" 20 cm	1 Lumbar Support
	7552716	L	35 <sup>1</sup> / <sub>2</sub> " – 39 <sup>1</sup> / <sub>4</sub> " 90 – 100 cm	8" 20 cm	1 Lumbar Support
	7552717	XL	39 <sup>1</sup> / <sub>4</sub> " – 45 <sup>1</sup> / <sub>4</sub> " 100 – 115 cm	8" 20 cm	1 Lumbar Support
	7552718	XXL	45 <sup>1</sup> / <sub>4</sub> " – 51 <sup>1</sup> / <sub>4</sub> " 118 – 130 cm	8" 20 cm	1 Lumbar Support
	7552719	3XL	51 <sup>1</sup> / <sub>4</sub> " – 57" 130 – 145 cm	8" 20 cm	1 Lumbar Support
	7552720	S	27 <sup>1</sup> / <sub>2</sub> " – 31 <sup>1</sup> / <sub>2</sub> " 70 – 80 cm	10" 25 cm	1 Lumbar Support
	7552721	M	31 <sup>1</sup> / <sub>2</sub> " – 35 <sup>1</sup> / <sub>2</sub> " 80 – 90 cm	10" 25 cm	1 Lumbar Support
	7552722	L	35 <sup>1</sup> / <sub>2</sub> " – 39 <sup>1</sup> / <sub>4</sub> " 90 – 100 cm	10" 25 cm	1 Lumbar Support
	7552723	XL	39 <sup>1</sup> / <sub>4</sub> " – 45 <sup>1</sup> / <sub>4</sub> " 100 – 115 cm	10" 25 cm	1 Lumbar Support
	7552724	XXL	51 <sup>1</sup> / <sub>4</sub> " – 51 <sup>1</sup> / <sub>4</sub> " 118 – 130 cm	10" 25 cm	1 Lumbar Support
	7552725	3XL	51 <sup>1</sup> / <sub>4</sub> " – 57" 130 – 145 cm	10" 25 cm	1 Lumbar Support

# Actimove® Back Support

High-Density Foam Panel

Adjustable Double Layer Compression

## Relieving massage effect, stability for the lumbar region

- Pain relief and healing are supported through balanced warmth and medical compression.
- Comfort and good breathability through neoprene-free perforated performance material.
- Provides strength and stability for the lumbar region incorporating a tapered high-density foam panel.
- Relieving massage effect through prominent pressure pads.
- Easy to apply and adjustable compression, thanks to the finger pockets and loops.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Available in two sizes.
- For stock keeping, long shelf life of five years.

This product is intended to provide support and stabilization to the lumbar joints of the back.

**Fields of application include:**

- Acute and chronic lower back pain
- Degenerative back problems (e.g. spondylosis)
- Lumbar instabilities (e.g. muscular insufficiencies)
- Post-operative care and rehabilitation



High-Density Foam Panel

Finger Loops and Pockets



	Code	Size	Measurements	Color
	7345030	S/M	29" - 37" / 74 - 94 cm	Silver
	7345031	L/XL	37" - 45" / 94 - 114 cm	Silver
	7345036	XXL/XXXL	45" - 53" / 114 - 134 cm	Silver

# Actimove® Mild Elbow Support

Lightweight design, easy to apply


- Mild support and easy application due to lightweight design.
- Knitted design supports good fit while maintaining flexibility.
- Soft yarns for comfort.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Available in four sizes.
- Can be worn left or right.
- For stock keeping, long shelf life of five years.

This product is intended to support the elbow.

Fields of application include:

- Slightly swollen and aching elbows



	Code	Size	Measurements	Color
	7561620	S	8 <sup>1</sup> / <sub>4</sub> " - 9 <sup>3</sup> / <sub>4</sub> " / 21 - 25 cm	White
	7561621	M	9 <sup>3</sup> / <sub>4</sub> " - 11 <sup>1</sup> / <sub>2</sub> " / 25 - 29 cm	White
	7561622	L	11 <sup>1</sup> / <sub>2</sub> " - 13" / 29 - 33 cm	White
	7561623	XL	13" - 14 <sup>1</sup> / <sub>2</sub> " / 33 - 37 cm	White

# Actimove® Elbow Support

Firm compression, reliable, durable support for everyday mobility


- Firm compression to help manage pain from elbow conditions.
- Reliable, durable support for everyday mobility.
- Seamless anatomical shape supported by a patterned helix design for a good fit.
- Reinforced edges for long lasting wear.
- Soft and cotton-free yarns for enduring comfort.
- Soft, smooth Motion-Comfort-Zone in the elbow crease.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Available in four sizes.
- Can be worn left or right.
- For stock keeping, long shelf life of five years.

This product is intended to support the elbow.

Fields of application include:

- Degenerative elbow conditions
- Overuse or repetitive use syndromes
- Chronic elbow pain



	Code	Size	Measurements	Color
	7561430	S	7 <sup>3</sup> / <sub>4</sub> " - 9" / 20 - 23 cm	Beige
	7561431	M	9" - 10 <sup>1</sup> / <sub>4</sub> " / 23 - 26 cm	Beige
	7561432	L	10 <sup>1</sup> / <sub>4</sub> " - 11 <sup>1</sup> / <sub>2</sub> " / 26 - 29 cm	Beige
	7561433	XL	11 <sup>1</sup> / <sub>2</sub> " - 12 <sup>1</sup> / <sub>2</sub> " / 29 - 32 cm	Beige



# Actimove® Elbow Support

## Pressure Pads and Strap

**Firm compression, adjustable pressure, durable support for added pain relief**

- Firm compression to help manage pain from elbow conditions.
- Two shaped pressure pads create a massage effect for added pain relief.
- Targeted pain relief by adjusting pressure on forearm tendon by extra elastic strap.
- Reliable, durable support for everyday mobility.
- Seamless anatomical shape supported by a patterned helix design for a good fit.
- Soft, smooth Motion-Comfort-Zone in the elbow crease.
- Reinforced edges for long lasting wear.
- Soft and cotton-free yarns for enduring comfort.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Available in four sizes.
- Can be worn left or right.
- For stock keeping, long shelf life of five years.

**This product is intended to support the elbow.**

**Fields of application include:**

- Degenerative elbow conditions
- Overuse or repetitive use syndromes
- Chronic elbow pain



Pressure Pads for Massage Effect



	Code	Size	Measurements	Color
	7561438	S	7 <sup>3</sup> / <sub>4</sub> " - 9" / 20 - 23 cm	Beige
	7561439	M	9" - 10 <sup>1</sup> / <sub>4</sub> " / 23 - 26 cm	Beige
	7561440	L	10 <sup>1</sup> / <sub>4</sub> " - 11 <sup>1</sup> / <sub>2</sub> " / 26 - 29 cm	Beige
	7561441	XL	11 <sup>1</sup> / <sub>2</sub> " - 12 <sup>1</sup> / <sub>2</sub> " / 29 - 32 cm	Beige


# Actimove® ARTHRITIS CARE

Individuals suffering from chronic, painful osteoarthritis or rheumatoid arthritis are on a quest to find lasting pain-relief<sup>1</sup> to help them live their daily lives.

Actimove® Arthritis Care supports aid pain relief and increase joint mobility through soothing warmth and compression.

- Actimove® Arthritis Care are made with **special ceramic yarns** which retain and reflect body heat.
- **Validated compression** properties provide support and stability, helping to reduce swelling, increase joint mobility and help to manage pain.
- The **skin friendly materials** are soft and comfortable and suitable for those sensitive to latex, as this product is not made with natural rubber latex.

## CERAMIC YARNS RETAIN AND REFLECT BODY HEAT



<sup>1</sup> Essity Internal Shopper Study



# Actimove® Wrist Support

Soothing warmth for wrists from heat-retaining ceramic fibers


- Unique ceramic fibers retain and reflect body heat for soothing warmth to help relieve pain.
- Four-way stretch provides light, comfortable compression to reduce swelling for improved mobility.
- Therapeutic warming promotes circulation around the wrist for faster return to daily activities.
- Specially designed Flex-Zones offer higher elasticity for more comfort during joint movement.
- Flex-Zones aid proper placement of the support on the joint.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Available in four sizes.
- Can be worn left or right.
- For stock keeping, long shelf life of five years.

This product is intended to provide support and therapeutic warmth to the wrist.

Fields of application include:

- Painful arthritis of the wrist joint



	Code	Size	Measurements	Color
	7577920	S	4 <sup>3</sup> / <sub>4</sub> " - 6" / 12 - 15 cm	Beige
	7577921	M	6" - 6 <sup>3</sup> / <sub>4</sub> " / 15 - 17 cm	Beige
	7577922	L	6 <sup>3</sup> / <sub>4</sub> " - 7 <sup>1</sup> / <sub>2</sub> " / 17 - 19 cm	Beige
	7577923	XL	7 <sup>1</sup> / <sub>2</sub> " - 8 <sup>1</sup> / <sub>4</sub> " / 19 - 21 cm	Beige



# Actimove® Arthritis Gloves

Balanced warmth and medical compression, comfortable from outer seams

- Pain relief and healing are supported through balanced warmth and medical compression.
- Comfort and good breathability through neoprene-free performance material.
- Provides light compression to help reduce swelling for improved finger and hand mobility.
- Comfortable to wear because of premium soft material and outside seams.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Available in four sizes.
- Can be worn left or right.
- For stock keeping, long shelf life of five years.

This product is intended for arthritis of the hand.

Fields of application include:

- Arthritis
- Aching fingers and hands
- Overuse or repetitive use syndromes



Outside Seams



Inside of Glove



	Code	Size	Measurements	Color
	7578320	S	7" - 8" / 18 - 20 cm	Beige
	7578321	M	8" - 9 <sup>1</sup> / <sub>4</sub> " / 20 - 23.5 cm	Beige
	7578322	L	9 <sup>1</sup> / <sub>4</sub> " - 10 <sup>1</sup> / <sub>2</sub> " / 23.5 - 27 cm	Beige
	7578323	XL	10 <sup>1</sup> / <sub>2</sub> " - 11 <sup>3</sup> / <sub>4</sub> " / 27 - 30 cm	Beige



# Actimove® Knee Support

## Soothing warmth for knees from heat-retaining ceramic fibers


- Unique ceramic fibers retain and reflect body heat for soothing warmth to help relieve pain.
- Four-way stretch provides light, comfortable compression to reduce swelling for improved mobility.
- Therapeutic warming promotes circulation around the knee for faster return to daily activities.
- Specially designed Flex-Zones offer higher elasticity for more comfort during joint movement.
- Flex-Zones aid proper placement of the support on the joint.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Available in five sizes.
- Can be worn left or right.
- For stock keeping, long shelf life of five years.

This product is intended to provide support and therapeutic warmth to the knee.

Fields of application include:

- Painful arthritis of the knee joint



	Code	Size	Measurements	Color
	7578120	S	11" - 13" / 28 - 33 cm	Beige
	7578121	M	13" - 15" / 33 - 38 cm	Beige
	7578122	L	15" - 17" / 38 - 43 cm	Beige
	7578123	XL	17" - 19 <sup>1/4</sup> " / 43 - 49 cm	Beige
	7578124	XXL	19 <sup>1/4</sup> " - 22" / 49 - 56 cm	Beige

# Actimove® Ankle Support

Soothing warmth for ankles from heat-retaining ceramic fibers


- Unique ceramic fibers retain and reflect body heat for soothing warmth to help relieve pain.
- Four-way stretch provides light, comfortable compression to reduce swelling for improved mobility.
- Therapeutic warming promotes circulation around the ankle for faster return to daily activities.
- Specially designed Flex-Zones offer higher elasticity for more comfort during joint movement.
- Flex-Zones aid proper placement of the support on the joint.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Available in five sizes.
- Can be worn left or right.
- For stock keeping, long shelf life of five years.

This product is intended to provide support and therapeutic warmth to the ankle.

Fields of application include:

- Painful arthritis of the ankle joint



	Code	Size	Measurements	Color
	7578020	S	6 <sup>1</sup> / <sub>4</sub> " - 8 <sup>1</sup> / <sub>4</sub> " / 16 - 21 cm	Beige
	7578021	M	8 <sup>1</sup> / <sub>4</sub> " - 10 <sup>1</sup> / <sub>4</sub> " / 21 - 26 cm	Beige
	7578022	L	10 <sup>1</sup> / <sub>4</sub> " - 12 <sup>1</sup> / <sub>4</sub> " / 26 - 31 cm	Beige
	7578023	XL	12 <sup>1</sup> / <sub>4</sub> " - 15 <sup>1</sup> / <sub>4</sub> " / 31 - 39 cm	Beige
	7578024	XXL	15 <sup>1</sup> / <sub>4</sub> " - 19 <sup>1</sup> / <sub>4</sub> " / 39 - 43 cm	Beige

# Actimove® Elbow Support

Soothing warmth for elbows from heat-retaining ceramic fibers

- Unique ceramic fibers retain and reflect body heat for soothing warmth to help relieve pain.
- Four-way stretch provides light, comfortable compression to reduce swelling for improved mobility.
- Therapeutic warming promotes circulation around the elbow for faster return to daily activities.
- Specially designed Flex-Zones offer higher elasticity for more comfort during joint movement.
- Flex-Zones aid proper placement of the support on the joint.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Available in four sizes.
- Can be worn left or right.
- For stock keeping, long shelf life of five years.

This product is intended to provide support and therapeutic warmth to the elbow.

Fields of application include:

- Painful arthritis of the elbow joint



	Code	Size	Measurements	Color
	7578221	S	7 <sup>1</sup> / <sub>2</sub> " - 9" / 19 - 23 cm	Beige
	7578222	M	9" - 10 <sup>1</sup> / <sub>4</sub> " / 23 - 26 cm	Beige
	7578223	L	10 <sup>1</sup> / <sub>4</sub> " - 11 <sup>1</sup> / <sub>2</sub> " / 26 - 29 cm	Beige
	7578224	XL	11 <sup>1</sup> / <sub>2</sub> " - 12 <sup>3</sup> / <sub>4</sub> " / 29 - 32.5 cm	Beige



## Actimove® KIDS

Many children suffering from play or sport injuries may not like to wear supports, because most available neoprene made options may be too hot.

Actimove® Kids offers pain relief and healing through balanced warmth and medical compression. Kids will enjoy great comfort from skin-friendly materials offering excellent moisture management.

- Made with **COOLMAX® AIR<sup>1</sup>** technology, Actimove® Kids are engineered for quick drying and exceptional breathability to keep children cool and dry.
- The perforated **performance material** of Actimove® Kids is neoprene-free, friendly to the skin, and suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Furthermore, **validated compression** properties for support and stability help reduce swelling and manage pain relief.

<sup>1</sup> COOLMAX® AIR is a trademark of The LYCRA company

# Actimove® Wrist Stabilizer

Removable Metal Stay, Left or Right

**Anatomically shaped stay for secure stabilization, managing sprains and strains**

- Pain relief and healing are supported through balanced warmth and medical compression.
- Comfort and good breathability through neoprene-free performance material.
- Stabilizes wrist through an anatomically shaped removable stay for left and right hand application.
- Quick and easy to fit with slip-on design and individually adjustable straps.
- Special comfort in a premium soft material and seamless design.
- Can be easily worn under clothing due to low profile.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Pediatric sizes, ages 2–7.
- For stock keeping, long shelf life of five years.

This product is intended to provide support/stabilization to the wrist.

Fields of application include:

- Tendonitis
- Sprains, strains, soft tissue injuries
- Post-operative care
- After cast removal



Slip-On Design



	Code	Size	Measurements	Color
	7572930	Pediatric	4 <sup>1</sup> / <sub>2</sub> " – 5 <sup>1</sup> / <sub>2</sub> " / 11.5 – 14 cm	Black





# Actimove® Knee Support

## Open Patella

**Fast, easy and secure protection, relieving patella pain in minor knee injuries**

- Enjoy sports, as pain relief and healing are supported through balanced warmth and medical compression.
- Great comfort with exceptional breathability through quick drying perforated performance material.
- Neoprene-free, made with COOLMAX® AIR technology.
- Helps to relieve patella pain. Knee cap pressure is reduced by open design.
- Feel secure as knee is protected by a snug fit that stays in place.
- Fast & easy slip-on style.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Two sizes, suitable for ages 6-7 and 8-12.
- For stock keeping, long shelf life of five years.

**This product is intended to provide support to the knee.**

**Fields of application include:**

- Painful knee injuries and prevention
- Slight strains, sprains and contusions
- Joint swelling



3D Diamond Lining with COOLMAX® AIR



	Code	Size	Measurements	Color
	7558510	Pediatric	10" - 12" / 25 - 30.5 cm	Navy
	7558511	Youth	12" - 14" / 30.5 - 35.5 cm	Navy



# Actimove® Ankle Support

## Slim fit sleeve for shoes

- Enjoy sports, as pain relief and healing are supported through balanced warmth and medical compression.
- Great comfort with exceptional breathability through quick drying perforated performance material.
- Neoprene-free, made with COOLMAX® AIR technology.
- Fast & easy slip on style.
- Comfortable slim-fit with flat seams for shoes.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Two sizes, suitable for ages 6–10 and 11–15.
- For stock keeping, long shelf life of five years.

This product is intended to provide support to the ankle.

Fields of application include:

- Painful ankle injuries & prevention
- Mild sprains
- Weak ankles
- Swelling and tenderness



3D Diamond Lining with COOLMAX® AIR



	Code	Size	Measurements	Color
	7560610	Pediatric	6 <sup>1</sup> / <sub>2</sub> " – 7 <sup>1</sup> / <sub>4</sub> " / 16.5 – 19 cm	Navy
	7560611	Youth	7 <sup>1</sup> / <sub>4</sub> " – 8 <sup>1</sup> / <sub>2</sub> " / 18.5 – 22 cm	Navy



# Actimove® Mitella

## Arm Sling

### Colorful support for children

- Child-friendly designed shoulder, elbow and arm support limiting their range of motion with a rested, flexed elbow position.
- Helps to relieve pain by stabilizing the shoulder joint and keeping the flexed elbow in a pocket which limits the range of motion.
- Comfortable fit around the neck.
- Patient comfort thanks to the breathable pocket material.
- Appealing and child-friendly, unique, playful pocket design.
- Convenient for follow up exams, as product does not need to be removed prior to taking x-rays.
- Protects casts or bandaging from moisture and dirt through the pocket.
- Easy to clean for better hygiene, machine washable at 40 °C / 105 °F.
- Suitable for those sensitive to latex, this product is not made with natural rubber latex
- Available in a universal, pediatric size
- For stock keeping, long shelf life of five years.

The product is intended to provide support to the arm, elbow and shoulder.

#### Fields of application include:

- Postoperative use (e.g., after rotator cuff rupture)
- Fracture management (e.g., clavicle, scapular, humeral, and greater tuberosity fractures)
- Joint distortions (sprains) and dislocations (e.g., acromioclavicular & sternoclavicular dislocations, shoulder luxation)
- Inflammatory diseases (e.g., calcific tendinitis of the rotator cuff)
- Minor instability of the shoulder (e.g., dorsal instability)
- Support for hemiplegic patients
- Support while wearing a cast or splint
- Temporary arm stabilization in emergency/paramedic care



	Code	Size	Measurements	Color
	7456113	Pediatric	10 <sup>1/2</sup> " - 13 <sup>1/2</sup> " / 27 - 34 cm	Navy



# Actimove® Cervical

## Cervical Collar

### Child-friendly cervical support



- Holds the head in a stable position thanks to the pre-formed anatomical shape.
- Support of the cervical spine after circular application.
- Provides a precise fit and support thanks to the anatomical design.
- Reliable support due to the durable materials.
- For more convenience and child-friendliness, free of metal buckles and harsh closures.
- Comfortable wear thanks to the low-contour and child-friendly, soft, breathable foam.
- For odor control, made with anti-microbially treated soft foam.
- Easy to apply and adjust hook and loop closures.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Available in three sizes.
- For stock keeping, long shelf life of five years.

The product is intended to provide support of the cervical spine after circular application, to help hold the head in a stable position and to restrict the neck range of motion.

**Fields of application include:**

- Support/stabilization to the cervical area



	Code	Size	Measurements		Unit
			Circumference	Height	
	7555530	Infant	9" - 10 <sup>1/2</sup> " 23 - 27 cm	2" 5 cm	1 Lumbar Support
	7555531	(0 - 24 months)	10" - 12" 25 - 30 cm	2 <sup>1/2</sup> " 6 cm	1 Lumbar Support
	7555532	Pediatric	12" - 14" 30 - 36 cm	2 <sup>1/2</sup> " 6 cm	1 Lumbar Support

# Actimove® Elbow Support

## Adjustable

### Slip-on design with easy application and adjustment

- Enjoy sports, as pain relief and healing are supported through balanced warmth and medical compression.
- Comfort and good breathability through neoprene-free performance material.
- Quick and easy to fit with slip-on design and individually adjustable straps.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Pediatric sizes, ages 5–7.
- For stock keeping, long shelf life of five years.

This product is intended to support the elbow.

Fields of application include:

- Painful elbow injuries & prevention



Slip-On Design



Code	Size	Measurements	Color
7561731	Pediatric	6 <sup>1/2</sup> " - 7 <sup>1/2</sup> " / 16.5 - 19 cm	Black





# Actimove® PROFESSIONAL LINE

Solutions to help improve the life of your patients supporting you in your everyday work life, the Actimove Professional Line is designed to support the treatment of acute injuries, post-op use and rehabilitation of chronic, orthopedic issues.

A shopper study<sup>1</sup> has shown what the patients' concerns are. The most important one is pain relief – followed by quality and comfort.

## Pain Relief

Actimove Professional Line products provide the right functionality based on biomechanics to help support patients' pain relief.

## Quality

Actimove Professional Line products are tested thoroughly and constantly updated for a high standard of quality and longevity.

## Comfort

Well-fitting products made from carefully chosen materials, paired with the correct sizing ensure patient comfort when it is most needed.

The line comprises a great range of braces and supports – for most needs, ages, and body parts. Actimove pledges to improve its offering further and further through research and innovation. All to give you assurance and help manage your patients' pain and support their movement.

<sup>1</sup> Essity shopper study, data on file

# WRIST & THUMB



# Actimove® Rhizo Forte

## Thumb Orthosis

**The patented brace restricting thumb movement while maintaining the functionality of the hand during daily activities.**

- Helps to alleviate pain arising from osteoarthritis of the thumb.<sup>1</sup>
- Supports pain relief by placing the thumb in a functional position and immobilizing the carpometacarpal (CMC) and the first metacarpophalangeal (MCP1) joints.
- Improved design for optimized fit and comfort around the thumb.
- For everyday wear, as wrist and fingers are free to move while the thumb is effectively restricted.
- Provides excellent fit and high comfort by customizing anatomical shape.
- Individually adjustable to the palm by molding the aluminum core.
- Maximum wearing comfort thanks to the open, low-profile design.
- Measure the palm width for the correct size.
- Easy to close with one hand using hook & loop straps.
- Easy cleansing of the brace by using hand wash.
- For easy daily hygiene: Water-resistant brace dries quickly after hand washing.
- This orthosis can even be worn when in contact with food.<sup>2</sup>
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Protected by patent, see [www.essity.com/patents](http://www.essity.com/patents).
- Available in three sizes.
- Available in gray and black.
- For stock keeping, long shelf life of five years.


**This product is intended to provide immobilization of the thumb.**

**Fields of application include:**

- Thumb osteoarthritis
- Rheumatoid arthritis of the thumb
- Ulnar collateral ligament injuries of the MCP 1 joint (e.g., skier's thumb)
- Instability of the CMC joint
- Soft tissue injuries (sprains)
- Traumatic thumb dislocation
- Ligament instabilities
- Postsurgical and nonsurgical treatment





	Code		Size	Measurements	Color	Unit
	Left	Right				
	7623801	7623800	S	2 <sup>3</sup> / <sub>4</sub> " - 3" 6.8 - 8 cm	Gray	1 Thumb Orthosis
	7623803	7623802	M	3" - 3 <sup>1</sup> / <sub>2</sub> " 8 - 9 cm	Gray	1 Thumb Orthosis
	7623805	7623804	L	3 <sup>1</sup> / <sub>2</sub> " - 4" 9 - 10 cm	Gray	1 Thumb Orthosis
	7623819	7623818	S	2 <sup>3</sup> / <sub>4</sub> " - 3" 6.8 - 8 cm	Black	1 Thumb Orthosis
	7623821	7623820	M	3" - 3 <sup>1</sup> / <sub>2</sub> " 8 - 9 cm	Black	1 Thumb Orthosis
	7623823	7623822	L	3 <sup>1</sup> / <sub>2</sub> " - 4" 9 - 10 cm	Black	1 Thumb Orthosis

<sup>1</sup>Essity/BSN medical GmbH, 2021, data on file, submitted for publication 2 SGS group, 2021, data on file



# Actimove® ManuMotion

## Wrist Support

**Anatomically shaped wrist support featuring 3D knitting technology and award-winning helix design for medical compression providing stability with integrated pressure pad to help pain relief.**

- Supports pain relief and healing by stabilizing the wrist in a physiological position through the anatomically shaped, thermoplastic stay and pressure pad.
- Helps to support pain relief by stabilizing the wrist through the 3D flat-knitted structure providing compression.
- Added stability thanks to the circumferential, inelastic wrist band that is adjustable and fully removable.
- Contributes to healing by reducing edema, effusions, and hematoma, thanks to the anatomically shaped dorsal pressure pad for a targeted intermittent massage effect.
- For more wearing comfort especially designed soft, smooth Motion-Comfort-Zone helps to prevent wrinkling in the thumb opening area and allows it to move.
- Prevents further pressure peaks on bony prominences thanks to the cut-outs in the dorsal pressure pad.
- Anti-slip for pressure-reduced distal and proximal cuffs without an elastic band preventing constriction.
- Seamless look thanks to the fusion technology used for the integration of the pressure pad.
- Easily removable stay placed in a pocket using seamless fusion technology for fixation on the outer side of the wrist.
- Skin-friendly and breathable 3D knit using OEKO-TEX® certified combination of yarn materials.<sup>1</sup>
- Easier application thanks to the included donning aid.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Easy to clean for better hygiene, machine washable at 30 °C / 85 °F gentle wash cycle.
- Available in six sizes
- Available in white and charcoal.
- For stock keeping, long shelf life of five years.

**This product is intended for wrist support.**

**Fields of application include:**

- Painful arthrosis
- Rheumatoid arthritis of the wrist
- Distortions (sprains)
- Tendovaginitis
- Tendon injuries
- Ligament injuries
- Soft tissue injuries
- Posttraumatic irritations
- Postoperative use








Donning Aid included.

	Code		Size	Measurements	Color	Unit
	Left	Right				
	7349742	7349712	XS	5 <sup>1</sup> / <sub>2</sub> " – 6" 14 – 15 cm	White	1 Wrist Support
	7349743	7349713	S	6" – 6 <sup>1</sup> / <sub>4</sub> " 15 – 16 cm	White	1 Wrist Support
	7349744	7349714	M	6 <sup>1</sup> / <sub>4</sub> " – 6 <sup>3</sup> / <sub>4</sub> " 16 – 17 cm	White	1 Wrist Support
	7349745	7349715	L	6 <sup>3</sup> / <sub>4</sub> " – 7" 17 – 18 cm	White	1 Wrist Support
	7349746	7349716	XL	7" – 7 <sup>1</sup> / <sub>2</sub> " 18 – 19 cm	White	1 Wrist Support
	7349747	7349717	XXL	7 <sup>1</sup> / <sub>2</sub> " – 7 <sup>3</sup> / <sub>4</sub> " 19 – 20 cm	White	1 Wrist Support
	7349739	7349738	XS	5 <sup>1</sup> / <sub>2</sub> " – 6" 14 – 15 cm	Charcoal	1 Wrist Support
	7349761	7349760	S	6" – 6 <sup>1</sup> / <sub>4</sub> " 15 – 16 cm	Charcoal	1 Wrist Support
	7349763	7349762	M	6 <sup>1</sup> / <sub>4</sub> " – 6 <sup>3</sup> / <sub>4</sub> " 16 – 17 cm	Charcoal	1 Wrist Support
	7349765	7349764	L	6 <sup>3</sup> / <sub>4</sub> " – 7" 17 – 18 cm	Charcoal	1 Wrist Support
	7349767	7349766	XL	7" – 7 <sup>1</sup> / <sub>2</sub> " 18 – 19 cm	Charcoal	1 Wrist Support
	7349769	7349768	XXL	7 <sup>1</sup> / <sub>2</sub> " – 7 <sup>3</sup> / <sub>4</sub> " 19 – 20 cm	Charcoal	1 Wrist Support

<sup>1</sup> Tested for harmful substances according to OEKO-TEX® Standard 100. OEKO-TEX® certificate: 13.HTR.25112





	Code		Size	Measurements	Color	Unit
	Left	Right				
	7348202	7348201	XS	$\leq 5\frac{1}{2}''$ $\leq 14$ cm	Black	1 Wrist Brace
	7348204	7348203	S/M	$5\frac{1}{2}'' - 7\frac{1}{2}''$ 14 - 19 cm	Black	1 Wrist Brace
	7348206	7348205	L/XL	$7\frac{1}{2}'' - 9''$ 19 - 23 cm	Black	1 Wrist Brace

# Actimove® Wrist Splint 8"

Full-length wrist splint for secure stabilization, made of breathable open cell foam.

- Helps to reduce pain as wrist is supported in a functional position, preventing flexion and extension through the dorsal stay and the wide, pre-formed removable palmar stay.
- Aids pain relief as functional positioning of the wrist is provided by the wide, pre-formed and removable palmar stay.
- Helps to support pain relief through comfortable stabilization.
- Feel securely supported in a full 8" length splint.
- Additional support to the wrist from the lateral and medial stays.
- Comfortable to wear due to the open cell foam as a breathable alternative to neoprene.
- Individual fit and extra stabilization thanks to the three adjustable, sturdy straps.
- Increased comfort thanks to the large thumb opening.
- Full finger movement and dexterity thanks to the special open design.
- Easy one-hand application thanks to the tubular slip-on design.
- Durable thanks to the binding of the edges.
- Easy to close straps thanks to the curved D-rings.
- For hygiene, hand-washable at 30°C / 85 °F.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Available in five sizes.
- Available for left and right.

The product is intended to provide stabilization/support to the wrist.

Fields of application include:

- Tendovaginitis (US: tendinitis)
- Sprains and strains
- Soft tissue injuries
- Osteoarthritis (OA)
- Rheumatoid arthritis
- Post-operative care
- Care after cast removal
- Carpal tunnel syndrome (CTS)



	Code		Size	Measurements	Color	Unit
	Left	Right				
	7571755	7571750	XS	4 <sup>1</sup> / <sub>2</sub> " – 5 <sup>1</sup> / <sub>2</sub> " 11.5 – 14 cm	Black	1 Wrist Support
	7571756	7571751	S	5 <sup>1</sup> / <sub>2</sub> " – 6 <sup>1</sup> / <sub>2</sub> " 14 – 16.5 cm	Black	1 Wrist Support
	7571757	7571752	M	6 <sup>1</sup> / <sub>2</sub> " – 7 <sup>1</sup> / <sub>2</sub> " 16.5 – 19 cm	Black	1 Wrist Support
	7571758	7571753	L	7 <sup>1</sup> / <sub>2</sub> " – 8 <sup>1</sup> / <sub>2</sub> " 19 – 21.5 cm	Black	1 Wrist Support
	7571759	7571754	XL	> 8 <sup>1</sup> / <sub>2</sub> " > 21.5 cm	Black	1 Wrist Support

# Actimove® Wrist Splint 6"

Short-length wrist splint for comfortable stabilization, made of breathable open cell foam.


- Helps to reduce pain as wrist is supported in a functional position, preventing flexion and extension through the dorsal stay and the wide, pre-formed removable palmar stay.
- Aids pain relief as functional positioning of the wrist is provided by the wide, pre-formed and removable palmar stay.
- Helps to support pain relief through comfortable stabilization.
- Individual fit and extra stabilization thanks to the two adjustable, sturdy straps.
- Additional support to the wrist from the lateral and medial stays.
- More comfortable to wear as the short 6" length is less bulky.
- Comfortable to wear due to the open cell foam as a breathable alternative to neoprene.
- Full finger movement and dexterity thanks to the special open design.
- Increased comfort thanks to the large thumb opening.
- Easy one-hand application thanks to the tubular slip-on design.
- Easy to close straps thanks to the curved D-rings.
- Durable thanks to the binding of the edges.
- Neoprene-free.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- For hygiene, hand-washable at 30 °C / 85 °F.
- Available in four sizes.
- Available for left and right.

The product is intended to provide stabilization/support to the wrist.

Fields of application include:

- Tendovaginitis (US: tendinitis)
- Sprains and strains
- Soft tissue injuries
- Osteoarthritis (OA)
- Rheumatoid arthritis
- Post-operative care
- Care after cast removal
- Carpal tunnel syndrome (CTS)



	Code		Size	Measurements	Color	Unit
	Left	Right				
	7571534	7571530	S	5 <sup>1</sup> / <sub>2</sub> " – 6 <sup>1</sup> / <sub>2</sub> " 14 – 16.5 cm	Black	1 Wrist Brace
	7571535	7571531	M	6 <sup>1</sup> / <sub>2</sub> " – 7 <sup>1</sup> / <sub>2</sub> " 16.5 – 19 cm	Black	1 Wrist Brace
	7571536	7571532	L	7 <sup>1</sup> / <sub>2</sub> " – 8 <sup>1</sup> / <sub>2</sub> " 19 – 21.5 cm	Black	1 Wrist Brace
	7571537	7571533	XL	> 8 <sup>1</sup> / <sub>2</sub> " > 21.5 cm	Black	1 Wrist Brace



# Actimove® Wrist Brace

Suede Finish, Lace-Closure

Soft suede-finish wrist brace with easy single pull lace-closure

- Helps to reduce pain as wrist is supported in a functional position, preventing flexion and extension through the dorsal stay and the wide, pre-formed removable palmar stay.
- Aids pain relief as functional positioning of the wrist is provided by the wide, preformed, and removable palmar stay.
- Prevents flexion and extension through the anatomically designed palmar and dorsal stays.
- Extra comfortable to wear thanks to the perforated suede-finished material.
- Easy one-hand application thanks to the tubular slip-on design combined with the lace-pully closure system.
- Soft and even fit with conforming lace closure.
- Full finger movement and dexterity thanks to the special open design.
- Neoprene-free.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Available in five sizes.
- Available for left and right.
- Available in black and beige.

The product is intended to provide stabilization/support to the wrist.

Fields of application include:

- Tendovaginitis (US: tendinitis)
- Sprains and strains
- Soft tissue injuries
- Osteoarthritis (OA)
- Rheumatoid arthritis
- Post-operative care
- Care after cast removal
- Carpal tunnel syndrome (CTS)



	Code		Size	Measurements	Color	Unit
	Left	Right				
	7571685	7571680	XS	< 4 <sup>1</sup> / <sub>2</sub> " / < 11.5 cm	Black	1 Wrist Brace
	7571686	7571681	S	4 <sup>1</sup> / <sub>2</sub> " - 5 <sup>1</sup> / <sub>2</sub> " / 11.5 - 14 cm	Black	1 Wrist Brace
	7571687	7571682	M	5 <sup>1</sup> / <sub>2</sub> " - 6 <sup>1</sup> / <sub>2</sub> " / 14 - 16.5 cm	Black	1 Wrist Brace
	7571688	7571683	L	6 <sup>1</sup> / <sub>2</sub> " - 7 <sup>1</sup> / <sub>2</sub> " / 16.5 - 19 cm	Black	1 Wrist Brace
	7571689	7571684	XL	> 7 <sup>1</sup> / <sub>2</sub> " / > 19 cm	Black	1 Wrist Brace
	7571675	7571670	XS	< 4 <sup>1</sup> / <sub>2</sub> " / < 11.5 cm	Beige	1 Wrist Brace
	7571676	7571671	S	4 <sup>1</sup> / <sub>2</sub> " - 5 <sup>1</sup> / <sub>2</sub> " / 11.5 - 14 cm	Beige	1 Wrist Brace
	7571677	7571672	M	5 <sup>1</sup> / <sub>2</sub> " - 6 <sup>1</sup> / <sub>2</sub> " / 14 - 16.5 cm	Beige	1 Wrist Brace
	7571678	7571673	L	6 <sup>1</sup> / <sub>2</sub> " - 7 <sup>1</sup> / <sub>2</sub> " / 16.5 - 19 cm	Beige	1 Wrist Brace
	7571679	7571674	XL	> 7 <sup>1</sup> / <sub>2</sub> " / > 19 cm	Beige	1 Wrist Brace

# Actimove® Manus

## Wrist Stabilizer

Secure stabilization, removable stay for left or right hand.

- Pain relief and healing are supported by stabilization through anatomically shaped stay.
- Great comfort from soft material and contoured strap ends.
- Quick and easy to fit with slip on design and individually adjustable straps.
- Can be applied on left or right hand as stay is removable.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- For hygiene, hand washable at max. 30°C / 85°F.
- Available in four sizes.
- For stock keeping, long shelf life of five years.

The product is intended to provide support to the wrist.

Fields of application include:

- Tendonitis
- Sprains, strains
- Osteoarthritis
- Rheumatoid arthritis
- Post-operative care
- Post-traumatic care



	Code	Size	Measurements	Color	Unit
	7234850	S	5" - 6" 13 - 15 cm	Black	1 Wrist Stabilizer
	7234851	M	6" - 7" 15 - 17.5 cm	Black	1 Wrist Stabilizer
	7234852	L	7" - 8" 17.5 - 20 cm	Black	1 Wrist Stabilizer
	7234853	XL	8" - 9" 20 - 23 cm	Black	1 Wrist Stabilizer

# Actimove® Manus Forte Plus

## Wrist & Thumb Brace

**Highly comfortable, reliable brace for the immobilization of the wrist and thumb.**

- Supports pain relief by strongly and securely stabilizing the wrist in a functional position, preventing flexion and extension through the combination of the rigid aluminum stay, and the wide, circular and inelastic middle strap.
- Helps to alleviate pain as the thumb is immobilized by the rigid aluminum spica stay.
- Suitable for everyday wear, the secure-fit-design allows the fingers to move freely while fixating the wrist.
- Easy to apply, handle and remove with one hand due to the tubular slip-on design and the simple to adjust thumb closure.
- Additional overall stabilization of the wrist through two stays on the back of the hand.
- Individual fit by molding the anatomically shaped aluminum stay as required.
- Also usable in the late rehabilitation phase, when less stability is needed, by removing the palmar stay.
- For ease of use special security hook & loop closure to prevent straps from slipping out of the buckles.
- High wearing comfort and natural skin climate thanks to the soft, breathable spacer fabric.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Easy to clean for better hygiene, machine washable at 30 °C / 85 °F after palmar aluminum stay is removed.
- Available in three sizes.
- Available for left and right.
- For stock keeping, long shelf life of five years.


**The product is intended to provide stabilization/support to the wrist.**

**Fields of application include:**

- Osteoarthritis of the wrist and thumb (e.g., rhizarthrosis)
- Soft tissue injuries of the wrist and thumb (e.g. from repetitive motion activities, like texting or gaming)
- Skier's thumb
- Sprains and strains
- Rheumatoid arthritis
- Post-operative care
- Post-traumatic care
- Inflammatory diseases (e.g., De Quervain's tendosynovitis, tendovaginitis)
- Care after cast removal
- Carpal tunnel syndrome (CTS)





	Code		Size	Measurements	Color	Unit
	Left	Right				
	7349633	7349603	XS	$\leq 5\frac{1}{2}''$ $\leq 14$ cm	Black	1 Wrist Brace
	7349634	7349604	S/M	$5\frac{1}{2}'' - 7\frac{1}{2}''$ 14 - 19 cm	Black	1 Wrist Brace
	7349636	7349606	L/XL	$7\frac{1}{2}'' - 9''$ 19 - 23 cm	Black	1 Wrist Brace

# Actimove® Wrist Splint

## Abducted Thumb, Compression-Formed Composite

**Wrist splint with abducted thumb with extra protection through the composite closed cell foam thumb spica.**

- Helps to reduce pain through stabilization of wrist and holding the thumb in an abducted, favorable position.
- Helps to reduce pain as wrist is supported in a functional position, preventing flexion and extension through the dorsal stay and the wide, pre-formed removable palmar stay.
- Helps to reduce pain through comfortable compression thanks to soft foam material combined with three adjustable straps.
- Aids pain relief as functional positioning of the wrist is provided by the wide, pre-formed and removable palmar stay.
- Provides extra protection from impact through the durable, semi-rigid composite closed cell foam thumb spica.
- Extra stabilization through the three sturdy straps, combined with an elastic thumb strap.
- Full finger movement and dexterity thanks to the special open design.
- Comfortable to wear without pressure points around the MCP joint thanks to the seamless design.
- Comfortable wear thanks to the light, breathable inner and outer materials.
- Comfortable to wear due to the open cell foam as a breathable alternative to neoprene.
- Durable thanks to the binding of the edges.
- Easy one-hand application thanks to the tubular slip-on design.
- Easy to close straps thanks to the curved D-rings.
- For hygiene, hand-washable at 30 °C / 85 °F.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Available in five sizes.
- Available for left and right.

**The product is intended to provide stabilization/support to the wrist.**

**Fields of application include:**

- Osteoarthritis of the thumb (e.g., rhizarthrosis)
- Soft tissue injuries of the wrist and thumb (e.g., from repetitive motion activities, like texting or gaming)
- Skier's thumb
- Strains and sprains
- Rheumatoid arthritis
- Post-operative care
- Post-traumatic care
- Inflammatory diseases (e.g., De Quervain's tendosynovitis, tendovaginitis)
- Care after cast removal
- Carpal tunnel syndrome (CTS)





	Code		Size	Measurements	Color	Unit
	Left	Right				
	7571845	7571830	XS	4 <sup>1</sup> / <sub>2</sub> " – 5 <sup>1</sup> / <sub>2</sub> " 11.5 – 14 cm	Black	1 Wrist Support
	7571846	7571831	S	5 <sup>1</sup> / <sub>2</sub> " – 6 <sup>1</sup> / <sub>2</sub> " 14 – 16.5 cm	Black	1 Wrist Support
	7571847	7571832	M	6 <sup>1</sup> / <sub>2</sub> " – 7 <sup>1</sup> / <sub>2</sub> " 16.5 – 19 cm	Black	1 Wrist Support
	7571848	7571833	L	7 <sup>1</sup> / <sub>2</sub> " – 8 <sup>1</sup> / <sub>2</sub> " 19 – 21.5 cm	Black	1 Wrist Support
	7571849	7571834	XL	> 8 <sup>1</sup> / <sub>2</sub> " > 21.5 cm	Black	1 Wrist Support

# Actimove® Gauntlet

## Wrist & Thumb Stabilizer

Secure MCP joint stabilization, removable stay for left or right hand.

- Pain relief and healing are supported by stabilization through anatomically shaped stays.
- MCP joint stabilization through extra thumb stay while allowing full finger movement.
- Great comfort from soft material and contoured strap ends.
- Quick and easy to fit with slip on design and individually adjustable straps.
- Can be applied on left or right hand as palmar stay is removable.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- For hygiene, hand washable at max. 30 °C / 85 °F
- Available in four sizes.
- For stock keeping, long shelf life of five years.

The product is intended to provide support to the wrist and thumb.


Fields of application include:

- Osteoarthritis of the thumb
- Irritation of the thumb
- Skier's thumb
- Sprains, strains
- Rheumatoid arthritis
- Post-operative care
- Post-traumatic care
- DeQuervain's tendosynovitis



	Code	Size	Measurements	Color	Unit
	7285990	S	5" - 6" 13 - 15 cm	Black	1 Wrist & Thumb Stabilizer
	7285991	M	6" - 7" 15 - 17.5 cm	Black	1 Wrist & Thumb Stabilizer
	7285992	L	7" - 8" 17.5 - 20 cm	Black	1 Wrist & Thumb Stabilizer
	7285993	XL	8" - 9" 20 - 23 cm	Black	1 Wrist & Thumb Stabilizer

# KNEE





# Actimove® Genu Tri-Panel

## Knee Immobilizer

The knee immobilizer allows splint adjustments for a broad spectrum of patient leg circumferences.

- Three panels for ease in adjustment
- Maleable aluminum stays

The product is intended for knee support.

Fields of application include:

- Excellent for post-trauma, post surgical for knee immobilization





Code	Size	Measurements	Color	Unit
7281914	S / M 20" / 51 cm	< 25" / 63.5 cm	Blue	1 Knee Support
7281915	M / L 20" / 51 cm	> 25" / 63.5 cm	Blue	1 Knee Support
7281916	S / M 24" / 61 cm	< 25" / 63.5 cm	Blue	1 Knee Support
7281917	M / L 24" / 61 cm	> 25" / 63.5 cm	Blue	1 Knee Support



# Actimove® Genu Clips

Knee Immobilizer

# Actimove® Genu Eco

Knee Immobilizer

The knee immobilizer supports the knee and helps to reduce pain after trauma or surgery.

- Actimove® Genu Clips – clip closure for easy removal
- Actimove® Genu Eco – lightweight, maleable stays, velcro closure

The product is intended for knee support.

Fields of application include:

- Excellent for post-trauma, post-surgical and temporary knee immobilization





	Code	Size	Measurements	Color	Unit
Actimove® Genu Clips	7285945	Pediatric 12" / 30.5 cm	< 3'10" / < 12½" < 117 cm / < 32 cm	Blue	1 Knee Support
	7285946	Pediatric 16" / 40.5 cm	< 4'7" / < 16½" < 140 cm / < 42 cm	Blue	1 Knee Support
	7285947	S / M 20" / 51 cm	< 5'3" / < 20½" < 160 cm / < 52 cm	Blue	1 Knee Support
	7285948	M / L 20" / 51 cm	< 5'3" / 18 - 25½" < 160 cm / 46 - 65 cm	Blue	1 Knee Support
	7285949	S / M 24" / 61 cm	< 5'10" / < 20½" < 178 cm / < 52 cm	Blue	1 Knee Support
	7285950	M / L 24" / 61 cm	< 5'10" / 18 - 25½" < 178 cm / 46 - 65 cm	Blue	1 Knee Support
	7285951	S / M 28" / 71 cm	> 5'10" / < 20½" > 178 cm / < 52 cm	Blue	1 Knee Support
	7285952	M / L 28" / 71 cm	> 5'10" / 18 - 25½" > 178 cm / 46 - 65 cm	Blue	1 Knee Support
Actimove® Genu Eco	7281908	S / M 20" / 51 cm	< 5'3" / < 20½" < 160 cm / < 52 cm	Blue	1 Knee Support
	7281909	M / L 20" / 51 cm	< 5'3" / 18 - 25½" < 160 cm / 46 - 65 cm	Blue	1 Knee Support
	7281910	S / M 24" / 61 cm	< 5'10" / < 20½" < 178 cm / < 52 cm	Blue	1 Knee Support
	7281911	M / L 24" / 61 cm	< 5'10" / 18 - 25½" < 178 cm / 46 - 65 cm	Blue	1 Knee Support
	7281912	S / M 28" / 71 cm	> 5'10" / < 20½" > 178 cm / < 52 cm	Blue	1 Knee Support
	7281913	M / L 28" / 71 cm	> 5'10" / 18 - 25½" > 178 cm / 46 - 65 cm	Blue	1 Knee Support

# Actimove® GenuMotion

## Knee Support

Anatomically shaped knee support featuring 3D knitting technology and award-winning Helix design for gradient medical compression, providing stability through stays, and integrated pressure pad to help pain relief.

- Helps pain relief, as the support offers stability through the anatomical shape and the 3D flat-knitted structure, providing gradient medical compression and a secure fit.
- Contributes to pain relief and healing by reducing edema, through targeted pressure and an intermittent massage effect thanks to the anatomically shaped circular pressure pad.
- For more wearing comfort especially designed soft, smooth Motion-Comfort-Zone that helps to prevent wrinkling in the popliteal area.
- Anti-slip for pressure-reduced distal and proximal cuffs without an elastic band preventing constriction.
- Added stabilization with great mobility thanks to the lateral spiral stays.
- Skin-friendly and breathable 3D knit using OEKO-TEX® certified combination of yarn materials.<sup>1</sup>
- Seamless look thanks to the fusion technology used for the integration of the pressure pad and stays.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Easy to clean for better hygiene, machine washable at 30 °C / 85 °F gentle wash cycle.
- Available in 15 sizes, including six plus sizes.
- Available in white and charcoal.
- For stock keeping, long shelf life of five years.

The product is intended for knee support.

Fields of application include:

- Chondromalacia with or without recurrent joint swellings
- Gonarthrosis
- Arthritis
- Posttraumatic irritations
- Postoperative use
- Minor instabilities of the knee joint
- Minor meniscal damages
- Fibromyalgia (tendomyopathy)



	Code	Size	Measurements		Color	Unit
			Above Knee	Below Knee		
	7346810	XXS	13 <sup>3/4</sup> " - 15" 35 - 38 cm	9 <sup>7/8</sup> " - 11" 25 - 28 cm	White	1 Knee Support
	7346811	XS	15" - 16 <sup>1/8</sup> " 38 - 41 cm	11" - 12 <sup>1/4</sup> " 28 - 31 cm	White	1 Knee Support
	7346812	S	16 <sup>1/8</sup> " - 17 <sup>3/8</sup> " 41 - 44 cm	12 <sup>1/4</sup> " - 13 <sup>3/8</sup> " 31 - 34 cm	White	1 Knee Support
	7346813	M	17 <sup>3/8</sup> " - 18 <sup>1/2</sup> " 44 - 47 cm	13 <sup>3/8</sup> " - 14 <sup>5/8</sup> " 34 - 37 cm	White	1 Knee Support
	7346814	L	18 <sup>1/2</sup> " - 19 <sup>3/4</sup> " 47 - 50 cm	14 <sup>5/8</sup> " - 15 <sup>3/4</sup> " 37 - 40 cm	White	1 Knee Support
	7346815	XL	19 <sup>3/4</sup> " - 20 <sup>7/8</sup> " 50 - 53 cm	15 <sup>3/4</sup> " - 17" 40 - 43 cm	White	1 Knee Support
	7346816	XXL	20 <sup>7/8</sup> " - 22" 53 - 56 cm	17" - 18 <sup>1/8</sup> " 43 - 46 cm	White	1 Knee Support
	7346850	3XL	22" - 23 <sup>1/4</sup> " 56 - 59 cm	18 <sup>1/8</sup> " - 19 <sup>1/4</sup> " 46 - 49 cm	White	1 Knee Support
	7346851	4XL	23 <sup>1/4</sup> " - 24 <sup>3/8</sup> " 59 - 62 cm	19 <sup>1/4</sup> " - 20 <sup>1/2</sup> " 49 - 52 cm	White	1 Knee Support
	7346817	M+	18 <sup>1/2</sup> " - 19 <sup>3/4</sup> " 47 - 50 cm	13 <sup>3/8</sup> " - 14 <sup>5/8</sup> " 34 - 37 cm	White	1 Knee Support
	7346809	L+	19 <sup>3/4</sup> " - 20 <sup>7/8</sup> " 50 - 53 cm	14 <sup>5/8</sup> " - 15 <sup>3/4</sup> " 37 - 40 cm	White	1 Knee Support
	7346808	XL+	20 <sup>7/8</sup> " - 22" 53 - 56 cm	15 <sup>3/4</sup> " - 17" 40 - 43 cm	White	1 Knee Support
	7346807	XXL+	22" - 23 <sup>1/4</sup> " 56 - 59 cm	17" - 18 <sup>1/8</sup> " 43 - 46 cm	White	1 Knee Support
	7346852	3XL+	23 <sup>1/4</sup> " - 24 <sup>3/8</sup> " 59 - 62 cm	18 <sup>1/8</sup> " - 19 <sup>1/4</sup> " 46 - 49 cm	White	1 Knee Support
	7346853	4XL+	24 <sup>3/8</sup> " - 25 <sup>5/8</sup> " 62 - 65 cm	19 <sup>1/4</sup> " - 20 <sup>1/2</sup> " 49 - 52 cm	White	1 Knee Support
	7346843	XXS	13 <sup>3/4</sup> " - 15" 35 - 38 cm	9 <sup>7/8</sup> " - 11" 25 - 28 cm	Charcoal	1 Knee Support
	7346844	XS	15" - 16 <sup>1/8</sup> " 38 - 41 cm	11" - 12 <sup>1/4</sup> " 28 - 31 cm	Charcoal	1 Knee Support
	7346845	S	16 <sup>1/8</sup> " - 17 <sup>3/8</sup> " 41 - 44 cm	12 <sup>1/4</sup> " - 13 <sup>3/8</sup> " 31 - 34 cm	Charcoal	1 Knee Support
	7346846	M	17 <sup>3/8</sup> " - 18 <sup>1/2</sup> " 44 - 47 cm	13 <sup>3/8</sup> " - 14 <sup>5/8</sup> " 34 - 37 cm	Charcoal	1 Knee Support
	7346847	L	18 <sup>1/2</sup> " - 19 <sup>3/4</sup> " 47 - 50 cm	14 <sup>5/8</sup> " - 15 <sup>3/4</sup> " 37 - 40 cm	Charcoal	1 Knee Support
	7346848	XL	19 <sup>3/4</sup> " - 20 <sup>7/8</sup> " 50 - 53 cm	15 <sup>3/4</sup> " - 17" 40 - 43 cm	Charcoal	1 Knee Support
	7346849	XXL	20 <sup>7/8</sup> " - 22" 53 - 56 cm	17" - 18 <sup>1/8</sup> " 43 - 46 cm	Charcoal	1 Knee Support

<sup>1</sup> Tested for harmful substances according to OEKO-TEX® Standard 100. OEKO-TEX® certificate: 13.HTR.25112

# Actimove® Knee Brace

## Composite Polycentric Hinges

**Guided knee stabilization through polycentric hinges with built-in stops and additional lateral protection.**

- Helps to reduce pain through warmth and compression.
- Provides additional protection against lateral impacts during activity thanks to the contoured condyle shells.
- Guided stabilization of the knee with polycentric hinges closely tracking the natural range of motion.
- Protects the knee from hyperextension through the built-in hinge stop.
- Supports patella tracking thanks to the contoured double layer neoprene material.
- Long-lasting wearing comfort due to lightweight, durable composite hinge material and design.
- Comfortable and secure fit thanks to the anatomical, tapered neoprene sleeve design.
- For additional comfort, friction and irritation is reduced thanks to the popliteal opening.
- Stays in place with improved fit thanks to the four durable, adjustable support straps.
- Easy to apply and remove slip-on design.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Available in five sizes.
- Can be used both on left and right knee.

The product is intended to stabilize the knee.

Fields of application include:

- Intermediary support for mild to moderate ligament strains/sprains (ACL/PCL)
- Meniscal injuries
- Medial/lateral ligament instabilities and tears
- Hyperextension or chronic knee instability problems



	Code	Size	Measurements		Color	Unit
			Above Knee	Below Knee		
	7569130	S	15" - 18" 38 - 46 cm	11" - 13" 28 - 33 cm	Black	1 Knee Brace
	7569131	M	18" - 21" 46 - 53 cm	13" - 16" 33 - 41 cm	Black	1 Knee Brace
	7569132	L	21" - 25" 53 - 64 cm	16" - 19" 41 - 48 cm	Black	1 Knee Brace
	7569133	XL	25" - 29" 64 - 74 cm	19" - 22" 48 - 56 cm	Black	1 Knee Brace
	7569134	XXL	29" - 33" 74 - 84 cm	22" - 24" 56 - 61 cm	Black	1 Knee Brace



# ANKLE



# Actimove® TaloMotion

## Ankle Support

**Anatomically shaped ankle support featuring 3D knitting technology and award-winning Helix design for gradient medical compression, providing stability, and integrated pressure pads to help pain relief.**

- Helps pain relief, as the support offers stability through the anatomical shape and the 3D flat-knitted structure, providing gradient medical compression.
- Contributes to healing by reducing edema, helping to relieve pain through targeted pressure and an intermittent massage effect thanks to the two anatomically shaped medial and lateral pressure pads.
- For more wearing comfort especially designed soft Motion-Comfort-Zone for a smooth fit in the instep.
- Anti-slip for pressure-reduced distal and proximal cuffs without an elastic band preventing constriction.
- Easier application thanks to the included donning aid.
- Seamless look thanks to the fusion technology used for the integration of the pressure pads.
- For long lasting wear, the heel and sole area is reinforced and strengthened by a special knitting structure to help prevent material abrasion.
- Skin-friendly and breathable 3D knit using OEKO-TEX® certified combination of yarn materials.<sup>1</sup>
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Easy to clean for better hygiene, machine washable at 30 °C / 85 °F gentle wash cycle.
- Available in six sizes.
- Available in white and charcoal.
- For stock keeping, long shelf life of five years.

**This product is intended for ankle support.**


**Fields of application include:**

- Acute lateral ligament injury (ruptured ligaments in ankle joint)
- Acute ankle sprain
- Acute chronic and inflammatory edema (e.g., in arthrosis or arthritis)
- Chronic ankle instability
- Postoperative use
- Posttraumatic irritation (e.g., after distortions)





Donning Aid included.

	Code		Size	Measurements	Color	Unit
	Left	Right				
	7348742	7348712	XS	6 <sup>1</sup> / <sub>2</sub> " – 7 <sup>1</sup> / <sub>2</sub> " 17 – 19 cm	White	1 Ankle Support
	7348743	7348713	S	7 <sup>1</sup> / <sub>2</sub> " – 8 <sup>1</sup> / <sub>4</sub> " 19 – 21 cm	White	1 Ankle Support
	7348744	7348714	M	8 <sup>1</sup> / <sub>4</sub> " – 9" 21 – 23 cm	White	1 Ankle Support
	7348745	7348715	L	9" – 9 <sup>3</sup> / <sub>4</sub> " 23 – 25 cm	White	1 Ankle Support
	7348746	7348716	XL	9 <sup>3</sup> / <sub>4</sub> " – 10 <sup>1</sup> / <sub>2</sub> " 25 – 27 cm	White	1 Ankle Support
	7348737	7348707	XXL	10 <sup>1</sup> / <sub>2</sub> " – 11 <sup>1</sup> / <sub>2</sub> " 27 – 29 cm	White	1 Ankle Support
	7348761	7348760	XS	6 <sup>1</sup> / <sub>2</sub> " – 7 <sup>1</sup> / <sub>2</sub> " 17 – 19 cm	Charcoal	1 Ankle Support
	7348763	7348762	S	7 <sup>1</sup> / <sub>2</sub> " – 8 <sup>1</sup> / <sub>4</sub> " 19 – 21 cm	Charcoal	1 Ankle Support
	7348765	7348764	M	8 <sup>1</sup> / <sub>4</sub> " – 9" 21 – 23 cm	Charcoal	1 Ankle Support
	7348767	7348766	L	9" – 9 <sup>3</sup> / <sub>4</sub> " 23 – 25 cm	Charcoal	1 Ankle Support
	7348769	7348768	XL	9 <sup>3</sup> / <sub>4</sub> " – 10 <sup>1</sup> / <sub>2</sub> " 25 – 27 cm	Charcoal	1 Ankle Support
	7348771	7348770	XXL	10 <sup>1</sup> / <sub>2</sub> " – 11 <sup>1</sup> / <sub>2</sub> " 27 – 29 cm	Charcoal	1 Ankle Support

<sup>1</sup> Tested for harmful substances according to OEKO-TEX® Standard 100. OEKO-TEX® certificate: 13.HTR.25112

# Actimove® AchilloMotion

## Achilles Tendon Support

**Anatomically shaped Achilles tendon support featuring 3D knitting technology and award-winning Helix design for gradient medical compression, providing stability, and integrated pressure pads to help pain relief.**

- Helps pain relief, as the support offers stability through the anatomical shape and the 3D flat-knitted structure, providing gradient medical compression.
- Contributes to healing by reducing edema, helping to relieve pain through targeted pressure and an intermittent massage effect thanks to the two anatomically shaped pressure pads parallel to the Achilles tendon.
- For more wearing comfort especially designed soft Motion-Comfort-Zone smooth fit in the instep.
- Anti-slip for pressure-reduced distal and proximal cuffs without an elastic band preventing constriction.
- Easier application featuring a donning aid in combination with a slanted upper cuff.
- Seamless look thanks to the fusion technology used for the integration of the pressure pads.
- For long lasting wear, the heel and sole area is reinforced and strengthened by a special knitting structure to help prevent material abrasion.
- Additional relief to the Achilles tendon due to the firm, viscoelastic heel wedges that are re-usable after cleaning with water.
- Easy fixation in flat, closed shoes, as the heel wedges are self-adhesive.
- Skin-friendly and breathable 3D knit using OEKO-TEX® certified combination of yarn materials.<sup>1</sup>
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Easy to clean for better hygiene, machine washable at 30 °C / 85 °F gentle wash cycle.
- Available in five sizes.
- Available in white and charcoal.
- For stock keeping, long shelf life of five years.

**This product is intended for Achilles' tendon support.**

**Fields of application include:**

- Acute and chronic inflammation of the Achilles tendon
- Achillodynia
- Insertion tendinopathy
- Achillobursitis
- Haglund exostosis
- Degenerative changes of the Achilles tendon
- Postoperative use and posttraumatic irritations (e.g., after Achilles tendon repair)





Donning Aid included. Heel wedges included.

	Code	Size	Measurements	Color	Unit
	7623113	S	7 <sup>1</sup> / <sub>2</sub> " – 8 <sup>1</sup> / <sub>4</sub> " 19 – 21 cm	White	1 Achilles Tendon Support
	7623114	M	8 <sup>1</sup> / <sub>4</sub> " – 9"	White	1 Achilles Tendon Support
	7623115	L	9" – 9 <sup>3</sup> / <sub>4</sub> " 23 – 25 cm	White	1 Achilles Tendon Support
	7623116	XL	9 <sup>3</sup> / <sub>4</sub> " – 10 <sup>1</sup> / <sub>2</sub> " 25 – 27 cm	White	1 Achilles Tendon Support
	7623117	XXL	10 <sup>1</sup> / <sub>2</sub> " – 11 <sup>1</sup> / <sub>2</sub> " 27 – 29 cm	White	1 Achilles Tendon Support
	7623107	S	7 <sup>1</sup> / <sub>2</sub> " – 8 <sup>1</sup> / <sub>4</sub> " 19 – 21 cm	Charcoal	1 Achilles Tendon Support
	7623108	M	8 <sup>1</sup> / <sub>4</sub> " – 9"	Charcoal	1 Achilles Tendon Support
	7623109	L	9" – 9 <sup>3</sup> / <sub>4</sub> " 23 – 25 cm	Charcoal	1 Achilles Tendon Support
	7623110	XL	9 <sup>3</sup> / <sub>4</sub> " – 10 <sup>1</sup> / <sub>2</sub> " 25 – 27 cm	Charcoal	1 Achilles Tendon Support
	7623111	XXL	10 <sup>1</sup> / <sub>2</sub> " – 11 <sup>1</sup> / <sub>2</sub> " 27 – 29 cm	Charcoal	1 Achilles Tendon Support

<sup>1</sup> Tested for harmful substances according to OEKO-TEX® Standard 100. OEKO-TEX® certificate: 13.HTR.25112



# Actimove® TaloBrace

## Comfort Ankle Brace

**Comfortable, anatomically shaped stabilizer that secures the upper ankle joint. Special soft, flexible cushioned inner foam padding and elastic upper strap.**

- Supports pain relief through the stabilization of the upper ankle joint preventing inversion and eversion.
- High comfort from soft interior foam padding that can be worn directly on the skin.
- Better fit as the elastic upper hook and loop strap fits comfortably over the tibia.
- Great wearing comfort due to specially designed soft, flexible edges.
- Secure feeling due to the anatomically shaped shells offering a high rigidity for increased stability.
- Supports healing by protecting medial and lateral ankle ligaments.
- Allows for normal gait thanks to unrestricted flexion and plantarflexion.
- Easy for patients to apply correctly due to printed foot logo on the bottom of the heel strap.
- Long-lasting, extra durable, easily adjustable heel pad.
- Can be worn in most shoes thanks to the low-profile, light-weight shells.
- Good hygiene as removable foam liner can be washed at 40°C / 104°F. Shells can be wiped with a damp cloth.
- Available for left and right feet.
- Available in two sizes.
- For stock keeping, long shelf life of five years.

**This product is intended to stabilize the ankle joint.**

**Fields of application include:**

- Ligament and capsule injuries at the upper ankle joint (e.g., ligament ruptures, strains, distortions)
- Slight soft tissue injuries
- Chronic ankle instability
- Rehabilitation and prevention of ankle injuries
- Acute, chronic, and inflammatory edema



	Code	Side	Size	Measurements	Unit
	7311913	R	S/M Trainer	< 5'4" / < 165 cm	1 Ankle Brace
	7311915		L/XL Standard	> 5'4" / > 165 cm	1 Ankle Brace
	7311912	L	S/M Trainer	< 5'4" / < 165 cm	1 Ankle Brace
	7311914		L/XL Standard	> 5'4" / > 165 cm	1 Ankle Brace

# Actimove® TaloBrace Air

## Comfort Ankle Brace

**Comfortable, anatomically shaped stabilizer that secures the upper ankle joint, with adjustable air cells for compression and a customized fit. Special soft flexible edges, cushioned inner foam padding and elastic upper strap.**

- Supports pain relief through the stabilization of the upper ankle joint preventing inversion and eversion.
- Pre-inflated, adjustable air cells offer compression for a customized fit.
- Great wearing comfort due to specially designed soft, flexible edges.
- Better fit as the elastic upper hook and loop strap fits comfortably over the tibia.
- Secure feeling due to the anatomically shaped shells offering a high rigidity for increased stability.
- Supports healing by protecting medial and lateral ankle ligaments.
- Allows for normal gait thanks to unrestricted flexion and plantarflexion
- Easy for patients to apply correctly due to printed foot logo on the bottom of the heel strap.
- Long-lasting, extra durable, easily adjustable heel pad.
- Can be worn in most shoes thanks to the low-profile, light-weight shells.
- Good hygiene as liner and shells can be wiped with a damp cloth.
- Available for left and right feet.
- Available in two sizes.
- For stock keeping, long shelf life of five years.

**This product is intended to stabilize the ankle joint.**

**Fields of application include:**

- Ligament and capsule injuries at the upper ankle joint (e.g., ligament ruptures, strains, distortions)
- Slight soft tissue injuries
- Chronic ankle instability
- Rehabilitation and prevention of ankle injuries
- Acute, chronic, and inflammatory edema



	Code	Side	Size	Measurements	Unit
	7311917	R	S/M Trainer	< 5'4" / < 165 cm	1 Ankle Brace
	7311919		L/XL Standard	> 5'4" / > 165 cm	1 Ankle Brace
	7311916	L	S/M Trainer	< 5'4" / < 165 cm	1 Ankle Brace
	7311918		L/XL Standard	> 5'4" / > 165 cm	1 Ankle Brace

# Actimove® TaloCast Air

## Ankle Brace

**Anatomically shaped stabilizer that secures the upper ankle joint. Adjustable air cells for compression and a reliable and easy to adjust heel strap.**

- Supports pain relief through the stabilization of the upper ankle joint preventing inversion and eversion.
- Pre-inflated, adjustable air cells offer compression for a customized fit.
- Long-lasting, extra durable easily adjustable heel strap on the outside of the shell.
- Secure feeling due to the anatomically shaped shells offering a high rigidity for increased stability.
- Supports healing by protecting medial and lateral ankle ligaments.
- Allows for normal gait thanks to unrestricted flexion and plantarflexion.
- Can be worn in most shoes thanks to the low-profile, light-weight shells.
- Good hygiene as liner and shells can be wiped with a damp cloth.
- Fits left and right foot and simplifies stock keeping.
- Available in two sizes.
- For stock keeping, long shelf life of five years.

**This product is intended to stabilize the ankle joint.**

**Fields of application include:**

- Ligament and capsule injuries at the upper ankle joint (e.g., ligament ruptures, strains, distortions)
- Slight soft tissue injuries
- Chronic ankle instability
- Rehabilitation and prevention of ankle injuries
- Acute, chronic, and inflammatory edema



	Code	Size	Measurements	Unit
	7311904	S/M Trainer	< 5'4" / < 165 cm	1 Ankle Brace
	7311905	L/XL Standard	> 5'4" / > 165 cm	1 Ankle Brace



# Actimove® TaloCast AirGel

## Ankle Brace

**Anatomically shaped stabilizer that secures the upper ankle joint. Adjustable air cells for compression and optional cold therapy for extra pain relief.**

- Supports pain relief through the stabilization of the upper ankle joint preventing inversion and eversion.
- Pre-inflated, adjustable air cells offer compression for a customized fit.
- Extra pain management through integrated gel pad for cold therapy.
- Long-lasting, extra durable easily adjustable heel strap on the outside of the shell.
- Secure feeling due to the anatomically shaped shells offering a high rigidity for increased stability.
- Supports healing by protecting medial and lateral ankle ligaments.
- Allows for normal gait thanks to unrestricted flexion and plantarflexion.
- Can be worn in most shoes thanks to the low-profile, light-weight shells.
- Good hygiene as liner and shells can be wiped with a damp cloth.
- Fits left and right foot and simplifies stock keeping.
- Available in two sizes.
- For stock keeping, long shelf life of five years.

**This product is intended to stabilize the ankle joint.**

**Fields of application include:**

- Ligament and capsule injuries at the upper ankle joint (e.g., ligament ruptures, strains, distortions)
- Slight soft tissue injuries
- Chronic ankle instability
- Rehabilitation and prevention of ankle injuries
- Acute, chronic, and inflammatory edema



	Code	Size	Measurements	Unit
	7311900	S/M Trainer	< 5'4" / < 165 cm	1 Ankle Brace
	7311901	L/XL Standard	> 5'4" / > 165 cm	1 Ankle Brace

# Actimove® TaloCast

## Ankle Brace

**Anatomically shaped stabilizer that secures the upper ankle joint. Reliable and easy to adjust heel strap and cushioned inner foam padding.**


- Supports pain relief through the stabilization of the upper ankle joint preventing inversion and eversion.
- Good wearing comfort due to soft foam liner that can be worn directly on the skin.
- Long-lasting, extra durable easily adjustable heel strap on the outside of the shell.
- Secure feeling due to the anatomically shaped shells offering a high rigidity for increased stability.
- Supports healing by protecting medial and lateral ankle ligaments.
- Allows for normal gait thanks to unrestricted flexion and plantarflexion.
- Can be worn in most shoes thanks to the low-profile, light-weight shells.
- Good hygiene as removable foam liner can be washed at 40°C / 104°F. Shells can be wiped with a damp cloth.
- Fits left and right foot and simplifies stock keeping.
- Available in two sizes.
- For stock keeping, long shelf life of five years.

**This product is intended to stabilize the ankle joint.**


**Fields of application include:**

- Ligament and capsule injuries at the upper ankle joint (e.g., ligament ruptures, strains, distortions)
- Slight soft tissue injuries
- Chronic ankle instability
- Rehabilitation and prevention of ankle injuries
- Acute, chronic, and inflammatory edema



	Code	Size	Measurements	Unit
	7311908	S/M Trainer	< 5'4" / < 165 cm	1 Ankle Brace
	7311909	L/XL Standard	> 5'4" / > 165 cm	1 Ankle Brace





# Actimove® Standard Walker

## Ankle Walker

### Walker with circumferential strapping for increased stability.

- Designed to stabilize and protect the lower leg, ankle, and foot.
- Plastic uprights give additional stability by restricting motion of leg.
- Circumferential strapping for increased stability.
- Dense foam lining within the sole to provide comfort and support.
- Shock absorbing sole to minimize impact and improve comfort.
- Rocker bottom to simulate natural gait.
- Lightweight walker for high patient comfort.
- Washable inner sole.
- Available in High (knee-height) or Low (ankle-height).
- Available in five sizes.
- For stock keeping, long shelf life of five years.

The product is designed to stabilize and protect the lower leg, ankle and foot.

**Fields of application include:**

- Stabilization of fractures of the foot and ankle
- Soft tissue injuries of the lower extremity
- Severe ankle sprains
- Stress fractures of the lower leg
- Postoperative rehabilitation and trauma



Code	Size	Shoe Sizes		Color	Unit
		Male	Female		
7627030	XS - Low	< 4.5	< 6	Black	1 Ankle Walker
7627031	S - Low	4.5 - 7	6 - 8	Black	1 Ankle Walker
7627032	M - Low	7.5 - 10	8.5 - 11	Black	1 Ankle Walker
7627033	L - Low	10.5 - 12.5	11.5 - 13.5	Black	1 Ankle Walker
7627034	XL - Low	> 12.5	> 13.5	Black	1 Ankle Walker
7627035	XS - High	< 4.5	< 6	Black	1 Ankle Walker
7627036	S - High	4.5 - 7	6 - 8	Black	1 Ankle Walker
7627037	M - High	7.5 - 10	8.5 - 11	Black	1 Ankle Walker
7627038	L - High	10.5 - 12.5	11.5 - 13.5	Black	1 Ankle Walker
7627039	XL - High	> 12.5	> 13.5	Black	1 Ankle Walker

# Actimove® Standard Air Walker

## Ankle Walker

Walker with easily adjustable air bladders to help control swelling, with circumferential strapping for increased stability.

- Medial and Lateral Air bladders create a larger surface for contact pressure to help control swelling and ensure a snug fit.
- Designed to stabilize and protect the lower leg, ankle, and foot.
- Plastic uprights give additional stability by restricting motion of leg.
- Circumferential strapping for increased stability
- Dense foam lining within the sole to provide comfort and support.
- Shock absorbing sole to minimize impact and improve comfort
- Rocker bottom to simulate natural gait.
- Integrated pump and valve system designed to easily adjust desired compression level.
- Lightweight walker for high patient comfort.
- Washable inner sole.
- Available in High (knee-height) or Low (ankle-height).
- Available in five sizes.
- For stock keeping, long shelf life of five years.

The product is designed to stabilize and protect the lower leg, ankle and foot.

**Fields of application include:**

- Stabilization of fractures of the foot and ankle
- Soft tissue injuries of the lower extremity
- Severe ankle sprains
- Stress fractures of the lower leg
- Postoperative rehabilitation and trauma



Code	Size	Shoe Sizes		Color	Unit
		Male	Female		
7627230	XS - Low	< 4.5	< 6	Black	1 Ankle Walker
7627231	S - Low	4.5 - 7	6 - 8	Black	1 Ankle Walker
7627232	M - Low	7.5 - 10	8.5 - 11	Black	1 Ankle Walker
7627233	L - Low	10.5 - 12.5	11.5 - 13.5	Black	1 Ankle Walker
7627234	XL - Low	> 12.5	> 13.5	Black	1 Ankle Walker
7627235	XS - High	< 4.5	< 6	Black	1 Ankle Walker
7627236	S - High	4.5 - 7	6 - 8	Black	1 Ankle Walker
7627237	M - High	7.5 - 10	8.5 - 11	Black	1 Ankle Walker
7627238	L - High	10.5 - 12.5	11.5 - 13.5	Black	1 Ankle Walker
7627239	XL - High	> 12.5	> 13.5	Black	1 Ankle Walker

# Actimove® Closed Shell Air Walker

## Ankle Walker

**Walker with rigid outer shell and circumferential cuff for maximum stability and comfort.**

- Adjustable toe guard provides protection in the toe area.
- Designed to stabilize and protect the lower leg, ankle, and foot.
- Circumferential plastic cuff gives maximum stability and comfort.
- High quality rigid outer skeleton for protection.
- Circumferential strapping for increased stability.
- Dense foam lining within the sole to provide comfort and support.
- Shock absorbing sole to minimize impact and improve comfort.
- Rocker bottom to simulate natural gait.
- Medial and Lateral Air bladders create a larger surface for contact pressure to help control swelling and ensure a snug fit.
- Integrated pump and valve system designed to easily adjust desired compression level.
- Lightweight walker for high patient comfort.
- Washable inner sole.
- Optional, adjustable heel wedge for Achilles' tendon injuries available.
- Available in High (knee-height) or Low (ankle-height).
- Available in five sizes.
- For stock keeping, long shelf life of five years.

**The product is designed to stabilize and protect the lower leg, ankle and foot.**

**Fields of application include:**

- Stabilization of fractures of the foot and ankle
- Soft tissue injuries of the lower extremity
- Severe ankle sprains
- Stress fractures of the lower leg
- Achilles' tendon repairs (only in combination with separately available Actimove® Heel Wedge for Walkers)
- Postoperative rehabilitation and trauma



Optional Heel Wedge.



Code	Size	Shoe Sizes		Color	Unit
		Male	Female		
7627610	XS - Low	< 4.5	< 6	Black	1 Ankle Walker
7627611	S - Low	4.5 - 7	6 - 8	Black	1 Ankle Walker
7627612	M - Low	7.5 - 10	8.5 - 11	Black	1 Ankle Walker
7627613	L - Low	10.5 - 12.5	11.5 - 13.5	Black	1 Ankle Walker
7627614	XL - Low	> 12.5	> 12	Black	1 Ankle Walker
7627615	XS - High	< 4.5	< 6	Black	1 Ankle Walker
7627616	S - High	4.5 - 7	6 - 8	Black	1 Ankle Walker
7627617	M - High	7.5 - 10	8.5 - 11	Black	1 Ankle Walker
7627618	L - High	10.5 - 12.5	11.5 - 13.5	Black	1 Ankle Walker
7627619	XL - High	> 12.5	> 12	Black	1 Ankle Walker
7627794	One Size	Actimove Heel Wedges For Closed Shell Walker		Black	1 Heel Wedge



# BACK & LUMBAL





# Actimove® Lumbar Sacral Support Comfort

## Additional Support Belt

Comfortable lumbar sacral support in a double-layer design for additional lift and support to the lower back area, available in two heights.


- Helps to reduce pain through stabilization and compression to the lower back area.
- Pain relief through additional compression thanks to the double-layer design.
- Provides additional lift and support to the lower back area thanks to its double-layer design.
- Enhanced stability due to its firm material.
- Enhances proprioception by stabilizing the back to keep an upright posture.
- Contoured fit through its asymmetrical wings.
- Easily adjustable compression thanks to the double-layer design.
- Comfortable to wear even when sitting as four flexible stays help to maintain positioning of the belt.
- Easy application and removal through wrap-around design.
- Customized fit through quick and easy adjustability in a wrap-around design.
- Smooth, durable surface helps to prevent pilling.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Easy to clean for better hygiene, machine washable at a gentle 40°C/105°F wash.
- Unisex, fits both men and women.
- Available in six sizes, plus a universal size.
- For stock keeping, long shelf life of five years.

**This product is intended to provide support and stabilization to the lumbar area.**

**Fields of application include:**

- Acute and chronic lower back pain
- Degenerative changes in the lumbar area (e.g., osteoporosis, osteochondrosis, spondylosis, spondylarthrosis, facet syndrome)
- Lumbar instabilities (insufficiencies of muscles and ligaments, sprains and strains)
- Lumbar disc protrusion
- Post-traumatic and post-operative use (e.g., after intervertebral disc operations)
- Prevention of repetitive strain injuries





	Code	Size	Height	Measurements	Color	Unit
	7552949	S	10" 25 cm	27 <sup>1</sup> / <sub>2</sub> " - 31 <sup>1</sup> / <sub>2</sub> " 70 - 80 cm	White	1 Lumbar Support
	7552950	M	10" 25 cm	31 <sup>1</sup> / <sub>2</sub> " - 35 <sup>1</sup> / <sub>2</sub> " 80 - 90 cm	White	1 Lumbar Support
	7552951	L	10" 25 cm	35 <sup>1</sup> / <sub>2</sub> " - 39 <sup>1</sup> / <sub>4</sub> " 90 - 100 cm	White	1 Lumbar Support
	7552952	XL	10" 25 cm	39 <sup>1</sup> / <sub>4</sub> " - 45 <sup>1</sup> / <sub>4</sub> " 100 - 115 cm	White	1 Lumbar Support
	7552953	XXL	10" 25 cm	45 <sup>1</sup> / <sub>4</sub> " - 51 <sup>1</sup> / <sub>4</sub> " 115 - 130 cm	White	1 Lumbar Support
	7552954	3XL	10" 25 cm	51 <sup>1</sup> / <sub>4</sub> " - 57" 130 - 145 cm	White	1 Lumbar Support
	7552955	Universal	10" 25 cm	31 <sup>1</sup> / <sub>2</sub> " - 49 <sup>1</sup> / <sub>4</sub> " 80 - 125 cm	White	1 Lumbar Support

# Actimove® LumbaMotion

## Back Support with Pressure Pad

Back support featuring 3D knitting technology and award-winning Helix design that stabilizes the lower back area and supports pain relief through circumferential compression and removeable pressure pad.

- Supports pain relief by stabilizing the lower back through the specially designed 3D flat-knitted fabric, providing circumferential compression for increased intraabdominal pressure supporting the dorsal muscles.
- Added relief by providing a massage effect for better blood circulation and relief of muscular tension through the adjustable and removable lumbar silicone pressure pad.
- Enhances proprioception, supports correct posture, and provides warmth to the lower back area.
- Easy and secure to open, close and adjust thanks to the specially designed two-pocket closure system.
- Enhances patient compliance by providing a good anatomical fit.
- Good patient comfort through the skin-friendly and breathable 3D knit using OEKO-TEX® certified combination of yarn materials.<sup>1</sup>
- Ensures correct application as the care label indicates the top of the support.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Easy to clean for better hygiene, machine washable at 30°C/85°F gentle wash cycle after removal of the pressure pad.
- Available in five sizes.
- Available in white and charcoal.
- For stock keeping, long shelf life of five years.

This product is intended to provide stabilization to the lumbar area.

Fields of application include:

- Acute and chronic lower back pain
- Degenerative changes in the lumbar area (e.g. osteoporosis, osteochondrosis, spondylosis, spondylarthrosis, facet syndrome)
- Lumbar instabilities (insufficiencies of muscles and ligaments)
- Post-traumatic use
- Prevention of repetitive strain injuries







	Code	Size	Measurements	Color	Unit
	7342812	XS	25 <sup>1/2</sup> " – 29 <sup>1/2</sup> " 65 – 75 cm	White	1 Back Support
	7342813	S	29 <sup>1/2</sup> " – 34 <sup>1/4</sup> " 75 – 87 cm	White	1 Back Support
	7342814	M	34 <sup>1/4</sup> " – 39 <sup>1/4</sup> " 87 – 100 cm	White	1 Back Support
	7342815	L	39 <sup>1/4</sup> " – 44 <sup>3/4</sup> " 100 – 114 cm	White	1 Back Support
	7342816	XL	44 <sup>3/4</sup> " – 51 <sup>1/4</sup> " 114 – 130 cm	White	1 Back Support
	7342817	XS	25 <sup>1/2</sup> " – 29 <sup>1/2</sup> " 65 – 75 cm	Charcoal	1 Back Support
	7342818	S	29 <sup>1/2</sup> " – 34 <sup>1/4</sup> " 75 – 87 cm	Charcoal	1 Back Support
	7342819	M	34 <sup>1/4</sup> " – 39 <sup>1/4</sup> " 87 – 100 cm	Charcoal	1 Back Support
	7342820	L	39 <sup>1/4</sup> " – 44 <sup>3/4</sup> " 100 – 114 cm	Charcoal	1 Back Support
	7342821	XL	44 <sup>3/4</sup> " – 51 <sup>1/4</sup> " 114 – 130 cm	Charcoal	1 Back Support

<sup>1</sup> Tested for harmful substances according to OEKO-TEX® Standard 100. OEKO-TEX® certificate: 13.HTR.25112

# Actimove® LumbaMotion Lady

## Back Support with Pressure Pad

**Back support, especially designed for women, featuring 3D knitting technology and award-winning Helix design that stabilizes the lower back area and supports pain relief through circumferential compression and removeable pressure pad.**

- Supports pain relief by stabilizing the lower back through the specially designed 3D flat-knitted fabric, providing circumferential compression for increased intraabdominal pressure supporting the dorsal muscles.
- Special fit for women thanks to design that reflects the female waist-hip ratio.
- Added relief by providing a massage effect for better blood circulation and relief of muscular tension through the adjustable and removable lumbar silicone pressure pad.
- Enhances proprioception, supports correct posture, and provides warmth to the lower back area.
- Easy and secure to open, close and adjust thanks to the specially designed two-pocket closure system.
- Enhances patient compliance by providing a good anatomical fit.
- Good patient comfort through the skin-friendly and breathable 3D knit using OEKO-TEX® certified combination of yarn materials.<sup>1</sup>
- Ensures correct application as the care label indicates the top of the support.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Easy to clean for better hygiene, machine washable at 30 °C/85 °F gentle wash cycle after removal of the pressure pad.
- Available in five sizes.
- For stock keeping, long shelf life of five years.


**This product is intended to provide stabilization to the lumbar area.**

**Fields of application include:**

- Acute and chronic lower back pain
- Degenerative changes in the lumbar area (e.g., osteoporosis, osteochondrosis, spondylosis, spondylarthrosis, facet syndrome)
- Lumbar instabilities (insufficiencies of muscles and ligaments)
- Post-traumatic use
- Prevention of repetitive strain injuries





	Code	Size	Measurements	Color	Unit
	7342912	XS	29 <sup>1/2</sup> " - 33 <sup>1/2</sup> " 75 - 85 cm	White	1 Back Support
	7342913	S	33 <sup>1/2</sup> " - 38" 85 - 97 cm	White	1 Back Support
	7342914	M	38" - 43 <sup>1/4</sup> " 97 - 110 cm	White	1 Back Support
	7342915	L	43 <sup>1/4</sup> " - 48 <sup>3/4</sup> " 110 - 124 cm	White	1 Back Support
	7342916	XL	48 <sup>3/4</sup> " - 55" 124 - 140 cm	White	1 Back Support

<sup>1</sup> Tested for harmful substances according to OEKO-TEX® Standard 100. OEKO-TEX® certificate: 13.HTR.25112

# Actimove® Deluxe Lumbar Sacral Support

## Rigid Panel, For Women

**Stabilizes the lumbar region through circumferential compression, paired with a rigid anatomically-shaped back panel and stays, designed for women by women.**

- Helps to relieve lower back pain and supports healing through warmth and circumferential compression.
- Stabilizes the lumbar region thanks to circumferential compression and the rigid anatomically-shaped back panel combined with six rigid stays.
- Helps alleviate back pain by supporting abdominal muscles.
- For adjustable support, back panel can be removed.
- Designed for women by women.
- Premium comfort and good breathability through the extra soft fabric.
- Especially comfortable for women due to the female anatomical shape.
- Discreet wearing comfort under clothing due to its low-profile.
- Comfortable to wear, even when sitting, as flexible stays maintain positioning.
- Specially designed for women, with unique fashionable colors.
- Silky binding for a feminine touch.
- Easy to use thanks to the thumb loops and micro-hook and loop closure.
- Neoprene-free.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- For hygiene, hand-washable at 30 °C / 85 °F.
- Available in lavender and rose.
- Available in eight sizes.

**This product is intended to stabilize the lumbar sacral region to help relieve lower back pain and discomfort.**

**Fields of application include:**

- Muscle strain
- Lower back pain and discomfort
- Lumbar spine sprain/strain
- Mild lumbar instability
- Lumbar disc dysfunction





	Code	Size	Measurements	Color	Unit
	7279230	S	28" - 32" 71 - 81 cm	Lavender	1 Lumbar Support
	7279231	M	33" - 36" 84 - 91 cm	Lavender	1 Lumbar Support
	7279232	L	37" - 40" 94 - 102 cm	Lavender	1 Lumbar Support
	7279233	XL	41" - 44" 104 - 112 cm	Lavender	1 Lumbar Support
	7279234	XXL	45" - 48" 114 - 122 cm	Lavender	1 Lumbar Support
	7279243	3XL	49" - 52" 125 - 132 cm	Lavender	1 Lumbar Support
	7279244	4XL	53" - 56" 135 - 142 cm	Lavender	1 Lumbar Support
	7279245	5XL	57" - 59" 145 - 150 cm	Lavender	1 Lumbar Support
	7279130	S	28" - 32" 71 - 81 cm	Rose	1 Lumbar Support
	7279131	M	33" - 36" 84 - 91 cm	Rose	1 Lumbar Support
7279132	L	37" - 40" 94 - 102 cm	Rose	1 Lumbar Support	
7279133	XL	41" - 44" 104 - 112 cm	Rose	1 Lumbar Support	
7279134	XXL	45" - 48" 114 - 122 cm	Rose	1 Lumbar Support	
7279143	3XL	49" - 52" 125 - 132 cm	Rose	1 Lumbar Support	
7279144	4XL	53" - 56" 135 - 142 cm	Rose	1 Lumbar Support	
7279145	5XL	57" - 59" 145 - 150 cm	Rose	1 Lumbar Support	



# Actimove® LO•PLUS™ Lumbar Support

## Stabilizing Posterior Panel

**Lumbar Support with stabilizing, rigid posterior panel, with tapered elastic panels for compression and support.**

- Helps to relieve lower back pain and discomfort by supporting and stabilizing the lumbar region with the help of the rigid posterior panel.
- For compression and support, tapered elastic panels contour the body.
- For great comfort, neoprene-free, and made with soft padded lining.
- Easy, quick application, fully adjustable with simple hook-and-loop closure system.
- Neoprene-free.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- For hygiene, hand-washable at 30°C/85°F.
- Available in six sizes.

This product is intended to provide stability to the lumbar region to relieve lower back pain and discomfort.

Fields of application include:

- Occasional/chronic back pain
- Lumbar Sprains/Strains
- Lumbar Instability
- Muscle Weakness





	Code	Size	Measurements	Color	Unit
	7456330	S	24" - 28" 61 - 71 cm	Black	1 Lumbar Support
	7456331	M	29" - 33" 74 - 84 cm	Black	1 Lumbar Support
	7456332	L	34" - 39" 86 - 99 cm	Black	1 Lumbar Support
	7456333	XL	40" - 44" 102 - 112 cm	Black	1 Lumbar Support
	7456334	XXL	45" - 49" 114 - 124 cm	Black	1 Lumbar Support
	7456335	3XL	50" - 55" 127 - 140 cm	Black	1 Lumbar Support

# SHOULDER





# Actimove® Sling Comfort

## Arm Sling – Single Pack

Comfort arm sling for a stable and supported forearm, fully assembled and pre-cut for a fast application.

- Secure, one-size-fits-all forearm support in a customizable strap length.
- Good patient wearing comfort due to the complete, open construction and reduced neck pressure through the 5.5 cm wide and stable sling replacing a separate neck pad.
- Quick and easy to apply or remove with one hand thanks to the hook and loop Y-tabs, with one already attached.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Easy to clean for better hygiene, machine washable at 40 °C / 105 °F.
- For stock keeping, long shelf life of five years.

This product is intended to provide support to the arm, elbow and shoulder.

### Fields of application include:

- Postoperative use (e.g. after rotator cuff rupture)
- Fracture management (e.g. clavicle, scapular, humeral, and greater tuberosity fractures)
- Joint distortions (sprains) and dislocations (e.g. acromioclavicular & sternoclavicular dislocations, shoulder luxation)
- Inflammatory diseases (e.g. calcific tendinitis of the rotator cuff)
- Minor instability of the shoulder (e.g. dorsal instability)
- Support for hemiplegic patients
- Support while wearing a cast or splint
- Temporary arm stabilization in emergency/paramedic care



	Code	Measurements	Color	Unit
Actimove® Sling Comfort	7285918	2¼" x 13 yds 5.5 cm x 12 m	Blue	2 Rolls/Box
Actimove® Sling Comfort – single use	7285924	2¼" x 2 yds 5.5 cm x 1.9 m	Blue	1 Sling
Actimove® Y-Tab Hook Closure	7346600	-	Blue	1 Set

# Actimove® Mitella Comfort

## Comfort Arm Sling

**Comfort shoulder, elbow and arm support limiting the range of motion in the shoulder and elbow joint and providing rest in a flexed elbow position.**

- Easy to apply, adjust or remove with one hand thanks to the hook and loop Y-Tabs and open-and-close pocket design.
- For good patient wearing comfort, reduced neck pressure though the 5.5 cm wide and stable sling replacing a separate neck pad.
- Helps to relieve pain by stabilizing the shoulder joint and keeping the flexed elbow in a pocket which limits the range of motion.
- Secure to wear, as the hand positioning through the thumb loop limits the rotation of the forearm, keeping it safe in the pocket while still allowing full finger movement.
- Patient comfort thanks to the breathable pocket material.
- Convenient for follow up exams, as product does not need to be removed prior to taking x-rays.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Easy to clean for better hygiene, machine washable at 40 °C / 105 °F.
- Available in four sizes.
- For stock keeping, long shelf life of five years.

**This product is intended to provide support to the arm, elbow and shoulder.**

### Fields of application include:

- Postoperative use (e.g. after rotator cuff rupture)
- Fracture management (e.g. clavicle, scapular, humeral, and greater tuberosity fractures)
- Joint distortions (sprains) and dislocations (e.g. acromioclavicular & sternoclavicular dislocations, shoulder luxation)
- Inflammatory diseases (e.g. calcific tendinitis of the rotator cuff)
- Minor instability of the shoulder (e.g. dorsal instability)
- Support for hemiplegic patients
- Support while wearing a cast or splint
- Temporary arm stabilization in emergency/paramedic care







	Code	Size	Measurements	Color	Unit
	7281961	S	< 14 <sup>1/2</sup> " < 37 cm	Blue	1 Arm Sling
	7281962	M	14 <sup>1/2</sup> " - 15 <sup>3/4</sup> " 37 - 40 cm	Blue	1 Arm Sling
	7281963	L	15 <sup>3/4</sup> " - 17" 40 - 43 cm	Blue	1 Arm Sling
	7281964	XL	> 17" > 43 cm	Blue	1 Arm Sling

# Actimove® Umerus Comfort

## Comfort Shoulder Immobilizer

### Comfort shoulder and upper arm immobilizer for limited movements in the shoulder joint.

- Helps to relieve pain by immobilizing the shoulder joint and limiting the abduction, extension, internal and external rotation of the humerus.
- Improved shoulder stabilization through the waist strap.
- Secure, as the hand positioning through the thumb loop limits the rotation of the forearm, keeping it safe in the pocket while still allowing full finger movement.
- For an individual fit, suitable for most patients thanks to the adjustable waist and shoulder straps.
- Easy to apply, adjust or remove with one hand thanks to the hook and loop Y-Tabs and open-and-close pocket design.
- For good patient wearing comfort, reduced neck pressure through the 5.5cm wide and stable sling replacing a separate neck pad.
- Patient comfort thanks to the breathable pocket material.
- Convenient for follow up exams, as product does not need to be removed prior to taking x-rays.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Easy to clean for better hygiene, machine washable at 40 °C / 105 °F.
- Available in five sizes.
- For stock keeping, long shelf life of five years.

This product is intended to provide support to the arm, elbow and shoulder.

#### Fields of application include:

- Postoperative use (e.g. after rotator cuff rupture)
- Fracture management (e.g. clavicle, scapular, humeral, and greater tuberosity fractures)
- Joint distortions (sprains) and dislocations (e.g. acromioclavicular & sternoclavicular dislocations, shoulder luxation)
- Inflammatory diseases (e.g. calcific tendinitis of the rotator cuff)
- Minor instability of the shoulder (e.g. dorsal instability)
- Support for hemiplegic patients
- Support while wearing a cast or splint
- Temporary arm stabilization in emergency/paramedic care





	Code	Size	Measurements	Color	Unit
	7281994	XS	< 13" < 33 cm	Blue	1 Shoulder Immobilizer
	7281939	S	13" - 14 <sup>1/2</sup> " 33 - 37 cm	Blue	1 Shoulder Immobilizer
	7281940	M	14 <sup>1/2</sup> " - 15 <sup>3/4</sup> " 37 - 40 cm	Blue	1 Shoulder Immobilizer
	7281941	L	15 <sup>3/4</sup> " - 17" 40 - 43 cm	Blue	1 Shoulder Immobilizer
	7281942	XL	> 17" > 43 cm	Blue	1 Shoulder Immobilizer

# Actimove® Shoulder Abduction Sling

## Shoulder Immobilizer

### Easy to handle shoulder immobilization.

- Easy to apply and to remove, even with one hand, thanks to the stable snap-lock buckle connecting the strap with the pocket, and the quick release buckle attached to the waist strap.
- Designed to aid pain relief through shoulder immobilization at the upper-arm abduction at around 15°.
- Secure and tailored fit for most patients through the extra-long waist strap of which can be adjusted in length.
- Increased comfort through the adjustable padded neck strap reducing the pressure on the neck.
- For increased patient comfort during prolonged wear as the perforated material allows for breathability.
- Brings comfort during wear due to the soft outer and inner materials, edges and seams.
- Helps to maintain muscle strength through hand and arm muscle training exercises using the attached squeeze ball.
- Clothes are protected through the soft, supple and color fast material.
- Dirt and signs of wear are less visible due to the black inner lining of the arm pocket.
- For easy care and a hygienic wearing experience, can be machine washed at 40 °C.
- Convenient during medical check-ups through radiolucency.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Universal, right or left fit, eliminates the need for excess stock.
- Available in four sizes.
- For stock keeping, long shelf life of five years.

**This product is intended to provide support to the arm, elbow and shoulder.**

**Fields of application include:**

- Minor instability of the shoulder (e.g., dorsal instability)
- Postoperative use (e.g., after rotator cuff rupture)
- Joint dislocations (e.g., shoulder luxation)



	Code	Size	Measurements	Color	Unit
	7344604	S	< 14 <sup>1/2</sup> " < 37 cm	Blue	1 Immobilizer
	7344605	M	14 <sup>1/2</sup> " - 15 <sup>3/4</sup> " 37 - 40 cm	Blue	1 Immobilizer
	7344606	L	15 <sup>3/4</sup> " - 17" 40 - 43 cm	Blue	1 Immobilizer
	7344607	XL	> 17" > 43 cm	Blue	1 Immobilizer

# CERVICAL





## Actimove® Cervical

Cervical Collar Soft Density, Narrow Width

## Actimove® Cervical

Cervical Collar Soft Density, Regular Width

### Soft-density, comfortable cervical support

- Provides gentle head and neck support, restricting neck range of motion.
- Holds the head in a stable position thanks to the pre-formed anatomical shape.
- Support of the cervical spine after circular application.
- Reliable support due to the durable materials.
- Comfortable wear through the high-quality, soft-density foam.
- For added comfort covered with a soft lining.
- Easy to apply and adjust hook and loop closures.
- Suits most patients due to the variety in neck sizes covered.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Available in five sizes, including a universal size.


**The product is intended to provide support of the cervical spine after circular application, to help hold the head in a stable position and to restrict the neck range of motion.**

**Fields of application may include:**

- Degenerative disorders in the cervical area
- Minor cervical trauma
- Post-Operative Use
- Painful cervical syndrome
- Cervical migraine
- Minor ligament injury
- Cervical Strain
- Short-term use after whiplash.





	Code	Size	Measurements		Color	Unit
			Circumference	Height		
	<b>Narrow Width</b>					
	7554730	UNIVERSAL	13" - 19" 33 - 48 cm	2 <sup>1</sup> / <sub>2</sub> " 6.5 cm	White	1 Cervical Collar
	7554731	S	9" - 13" 23 - 33 cm	2 <sup>1</sup> / <sub>2</sub> " 6.5 cm	White	1 Cervical Collar
	7554732	M	11" - 15" 28 - 38 cm	2 <sup>1</sup> / <sub>2</sub> " 6.5 cm	White	1 Cervical Collar
	7554733	L	13" - 17" 33 - 43 cm	2 <sup>1</sup> / <sub>2</sub> " 6.5 cm	White	1 Cervical Collar
	7554734	XL	15" - 19" 38 - 48 cm	2 <sup>1</sup> / <sub>2</sub> " 6.5 cm	White	1 Cervical Collar
<b>Regular Width</b>						
7555131	UNIVERSAL	13" - 19" 33 - 48 cm	3 <sup>1</sup> / <sub>4</sub> " 8 cm	White	1 Cervical Collar	
7555134	S	9" - 13" 23 - 33 cm	3 <sup>1</sup> / <sub>4</sub> " 8 cm	White	1 Cervical Collar	
7555135	M	11" - 15" 28 - 38 cm	3 <sup>1</sup> / <sub>4</sub> " 8 cm	White	1 Cervical Collar	
7555136	L	13" - 17" 33 - 43 cm	3 <sup>1</sup> / <sub>4</sub> " 8 cm	White	1 Cervical Collar	
7555137	XL	15" - 19" 38 - 48 cm	3 <sup>1</sup> / <sub>4</sub> " 8 cm	White	1 Cervical Collar	

# Actimove® Cervical Comfort\*

## Comfort Cervical Collar

**Comfortable cervical collar in a pre-formed anatomical shape with a soft, navy foam cover.**

- Provides comfortable head and neck support, restricting range of motion.
- Easy to apply and adjust hook and loop closures.
- Holds the head in a stable position thanks to the pre-formed anatomical shape.
- Good patient comfort due to the soft foam cover.
- Maintains a cleaner look during wear thanks to its attractive navy-blue design.
- Suits most patients due to the variety in neck sizes covered.
- Caution: This product contains natural rubber latex which may cause allergic reactions.
- For extra hygiene and to keep the collar clean, the removable sleeve can be machine washed at 30°C/85°F.
- Available in six sizes, including a short height version in XL.
- For stock keeping, long shelf life of five years.

**The product is intended to provide support/stabilization to the cervical area.**

**Fields of application include:**

- Short-term use after whiplash
- Degenerative changes in the cervical area (e.g., osteochondrosis, spondylarthroses, minor cervical arthrosis)
- Cervical migraine
- Sprains
- Strains
- Acute and chronic painful cervical syndrome
- Muscular pain syndrome
- Minor cervical trauma
- Immediate treatment of acute cervical fractures
- Postoperative use





\*Caution: This product contains natural rubber latex which may cause allergic reactions.

	Code	Size	Measurements		Color	Unit
			Circumference	Height		
	7285936	XS	10" - 13" 26 - 33 cm	2 <sup>1</sup> / <sub>2</sub> " 6.5 cm	Blue	1 Collar
	7285937	S	11 <sup>3</sup> / <sub>4</sub> " - 14" 30 - 36 cm	3" 7.5 cm	Blue	1 Collar
	7285938	M	13" - 15 <sup>1</sup> / <sub>4</sub> " 33 - 39 cm	3 <sup>1</sup> / <sub>2</sub> " 9 cm	Blue	1 Collar
	7285939	L	14" - 16 <sup>1</sup> / <sub>2</sub> " 36 - 42 cm	4" 10 cm	Blue	1 Collar
	7285940	XL	15" - 19" 38 - 48 cm	4 <sup>1</sup> / <sub>4</sub> " 11 cm	Blue	1 Collar
	7285941	XL-SHORT	15 <sup>1</sup> / <sub>4</sub> " - 17 <sup>3</sup> / <sub>4</sub> " 39 - 45 cm	2 <sup>3</sup> / <sub>4</sub> " 7 cm	Blue	1 Collar

# ABDOMEN & RIBS







# Actimove® Abdominal Binder

## Three-Panel

# Actimove® Abdominal Binder

## Four-Panel

**Abdominal binder that helps to reduce pain and supports wound healing through constant individual compression around the abdomen after surgery or trauma, breathable, available in two heights.**


- Helps to reduce pain through support to the injured area post-operatively.
- Helps to reduce pain through individually adjustable compression.
- Provides support and relief to the abdomen after surgery and trauma through compression.
- Helps to support wound healing through constant individual compression around the abdominal area.
- Helps to perform pulmonary exercises through the elastic material.
- For better comfort, made from breathable materials.
- More flexible and conforming fit through special, soft connectors between the panels.
- Comfortable against the skin thanks to its inner, soft, plush material.
- Customized fit due to wide range of adjustability in a wrap-around design.
- Easy application and removal through wrap-around design.
- Anatomical fit through V-shape.
- Helps to avoid skin maceration through non-absorbing, cotton-free materials.
- Supports a wide range of medical applications (e.g., drainage tubes, ostomy bags) as materials will not unravel when cut.
- Easy wound inspection due to fast opening and closure of belt.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Easy to clean for better hygiene, machine washable at a gentle 40 °C/105 °F wash.
- Unisex, fits both men and women.
- Actimove Abdominal Binder, Four-Panel: Suitable for taller patients, or longline surgeries, thanks to high four-panel design.
- Available in four sizes, the Actimove Abdominal Binder Three-Panel at a height of 23 cm / 9" and the Actimove Abdominal Binder Four-Panel at a height of 30 cm / 12".
- For stock keeping, long shelf life of five years.

**This product is intended for abdominal support.**

**Fields of application include:**

- Postoperative and posttraumatic use
- Relief of the abdomen and thorax



	Product	Code	Size	Measurements	Unit
	Actimove Abdominal Binder Three-Panel	7551306	XS	21 <sup>3</sup> / <sub>4</sub> " – 33 <sup>1</sup> / <sub>2</sub> " 55 – 85 cm	1 Abdominal Support
		7551307	S	33 <sup>1</sup> / <sub>2</sub> " – 45 <sup>1</sup> / <sub>4</sub> " 85 – 115 cm	1 Abdominal Support
		7551308	M	45 <sup>1</sup> / <sub>4</sub> " – 61" 115 – 155 cm	1 Abdominal Support
		7551309	L	61" – 76 <sup>3</sup> / <sub>4</sub> " 155 – 195 cm	1 Abdominal Support
	Actimove Abdominal Binder Four-Panel	7551310	XS	21 <sup>3</sup> / <sub>4</sub> " – 33 <sup>1</sup> / <sub>2</sub> " 55 – 85 cm	1 Abdominal Support
		7551311	S	33 <sup>1</sup> / <sub>2</sub> " – 45 <sup>1</sup> / <sub>4</sub> " 85 – 115 cm	1 Abdominal Support
		7551312	M	45 <sup>1</sup> / <sub>4</sub> " – 61" 115 – 155 cm	1 Abdominal Support
		7551313	L	61" – 76 <sup>3</sup> / <sub>4</sub> " 155 – 195 cm	1 Abdominal Support

# Actimove® Abdominal Binder Comfort

## Soft Pad

**Comfortable abdominal binder with soft foam pad for extra cushioning to the injured area. Helps to reduce pain and supports wound healing through constant individual compression around the abdomen after surgery or trauma, available in two heights.**

- Helps to reduce pain through support to the injured area post-operatively.
- Helps to reduce pain through individually adjustable compression.
- Provides support and relief to the abdomen after surgery and trauma through compression.
- Helps to support wound healing through constant individual compression around the abdominal area.
- Helps to perform pulmonary exercises through the elastic material.
- Enhanced wearing comfort due to extra cushioning from the soft foam pad to help protect the injured area.
- Enhanced wearing comfort due to extra cushioning of the injured area through terry-cloth lined, soft foam pad.
- Enhanced wearing comfort thanks to its breathable and soft materials.
- Easy application and removal through wrap-around design.
- Easily applied and removed thanks to the velvety, hook-compatible outer surface of pad.
- Customized fit due to wide range of adjustability in a wrap-around design.
- Helps to avoid skin maceration through non-absorbing, cotton-free materials.
- Supports a wide range of medical applications (e.g., drainage tubes, ostomy bags) as materials will not unravel when cut.
- Easy wound inspection due to fast opening and closure of belt.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Easy to clean for better hygiene, machine washable at a gentle 40 °C/105 °F wash.
- Unisex, fits both men and women.
- Suitable for both shorter and taller patients, or longline surgeries, as available in two heights.
- Available in five sizes.
- For stock keeping, long shelf life of five years.

**This product is intended for abdominal support.**

**Fields of application include:**

- Postoperative and posttraumatic use
- Relief of the abdomen and thorax







	Code	Size	Height	Measurements	Unit
	7551413	XS	9" 23 cm	21 <sup>3</sup> / <sub>4</sub> " – 33 <sup>1</sup> / <sub>2</sub> " 55 – 85 cm	1 Abdominal Support
	7551414	S	9" 23 cm	33 <sup>1</sup> / <sub>2</sub> " – 45 <sup>1</sup> / <sub>4</sub> " 85 – 115 cm	1 Abdominal Support
	7551415	M	9" 23 cm	45 <sup>1</sup> / <sub>4</sub> " – 61" 115 – 155 cm	1 Abdominal Support
	7551416	L	9" 23 cm	61" – 76 <sup>3</sup> / <sub>4</sub> " 155 – 195 cm	1 Abdominal Support
	7551417	XS	12" 30 cm	21 <sup>3</sup> / <sub>4</sub> " – 33 <sup>1</sup> / <sub>2</sub> " 55 – 85 cm	1 Abdominal Support
	7551418	S	12" 30 cm	33 <sup>1</sup> / <sub>2</sub> " – 45 <sup>1</sup> / <sub>4</sub> " 85 – 115 cm	1 Abdominal Support
	7551419	M	12" 30 cm	45 <sup>1</sup> / <sub>4</sub> " – 61" 115 – 155 cm	1 Abdominal Support
	7551420	L	12" 30 cm	61" – 76 <sup>3</sup> / <sub>4</sub> " 155 – 195 cm	1 Abdominal Support



# Actimove® Rib Belt

## Two-Panel, For Men

**Low-profile, breathable rib belt providing external support and individually adjustable compression with an anatomical fit.**


- Helps to reduce pain through individually adjustable compression.
- Helps to reduce pain giving external support to the injured rib area.
- For pain management, elastic material provides compression through limiting expansion of the rib cage.
- For better comfort, made from breathable materials.
- More flexible and conforming fit through special, soft connectors between the panels.
- Comfortable against the skin thanks to its inner, soft, plush material.
- Easy application and removal through wrap-around design.
- Customized fit due to wide range of adjustability in a wrap-around design.
- Inconspicuous wear due to its low-profile design.
- Helps to avoid skin maceration through non-absorbing, cotton-free materials.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Easy to clean for better hygiene, machine washable at a gentle 40 °C/105 °F wash.
- Available in two sizes, including a universal size.
- For stock keeping, long shelf life of five years.

This product is intended to support and stabilize the rib area.

Fields of application include:

- Rib fracture
- Post-traumatic conditions in the rib area (e.g., bruising)



	Code	Size	Measurements	Unit
	7550804	S	23 <sup>1</sup> / <sub>2</sub> " – 35 <sup>1</sup> / <sub>2</sub> " 60 – 90 cm	1 Rib Belt
	7550806	UN/L	< 51 <sup>1</sup> / <sub>4</sub> " < 130 cm	1 Rib Belt

# Actimove® Rib Belt

## Two-Panel, For Women

**Low-profile, breathable rib belt providing external support and individually adjustable compression with an anatomical, female fit.**


- Helps to reduce pain through individually adjustable compression.
- Helps to reduce pain giving external support to the injured rib area.
- For pain management, elastic material provides compression through limiting expansion of the rib cage.
- For better comfort, made from breathable materials.
- More flexible and conforming fit through special, soft connectors between the panels.
- Comfortable against the skin thanks to its inner, soft, plush material.
- Easy application and removal through wrap-around design.
- Customized fit due to wide range of adjustability in a wrap-around design.
- Enhanced comfort and improved fit through the bust line contour.
- Discreet wear thanks to its low-profile design.
- Helps to avoid skin maceration through non-absorbing, cotton-free materials.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Easy to clean for better hygiene, machine washable at a gentle 40 °C/105 °F wash.
- Available in two sizes, including a universal size.
- For stock keeping, long shelf life of five years.

**This product is intended to support and stabilize the rib area.**

**Fields of application include:**

- Rib fracture
- Post-traumatic conditions in the rib area (e.g., bruising)



	Code	Size	Measurements	Unit
	7550805	S	23 <sup>1</sup> / <sub>2</sub> " – 35 <sup>1</sup> / <sub>2</sub> " 60 – 90 cm	1 Rib Belt
	7550807	UN/L	< 51 <sup>1</sup> / <sub>4</sub> " <130 cm	1 Rib Belt

ELBOW







# Actimove® EpiMotion

## Elbow Support


**Anatomically shaped elbow support featuring 3D knitting technology and award-winning Helix design for gradient medical compression, providing stability, and integrated pressure pads to help pain relief.**

- Helps pain relief, as the support offers stability through the anatomical shape and the 3D flat-knitted structure, providing gradient medical compression.
- Contributes to pain relief and healing by reducing edema, through targeted pressure and an intermittent massage effect thanks to the two anatomically shaped ulnar and radial pressure pads.
- Prevents further pressure peaks on bony prominences thanks to the epicondyle cut-outs in the pressure pads.
- For more wearing comfort especially designed soft Motion-Comfort-Zone smooth fit in the elbow crease.
- Anti-slip for pressure-reduced distal and proximal cuffs without an elastic band preventing constriction.
- Seamless look thanks to the fusion technology used for the integration of the pressure pads.
- Skin-friendly and breathable 3D knit using OEKO-TEX® certified combination of yarn materials.<sup>1</sup>
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Easy to clean for better hygiene, machine washable at 30°C/85°F gentle wash cycle.
- Available in seven sizes.
- For stock keeping, long shelf life of five years.

**This product is intended for elbow support.**

**Fields of application include:**

- Tennis elbow (lateral epicondylitis)
- Golfers Elbow (medial epicondylitis)
- Arthritis (joint inflammation)
- Arthrosis
- Rheumatoid joint effusion
- Bursitis olecrani (student's elbow)
- Insertion tendopathy of the triceps brachii
- Posttraumatic irritations
- Postoperative use







	Code	Size	Measurements	Color	Unit
<p>12 cm = 4<sup>3</sup>/<sub>4</sub>"</p>	7347711	XXS	6 <sup>1</sup> / <sub>2</sub> " – 7 <sup>1</sup> / <sub>2</sub> " 17 – 19 cm	White	1 Elbow Support
	7347712	XS	7 <sup>1</sup> / <sub>2</sub> " – 8 <sup>1</sup> / <sub>4</sub> " 19 – 21 cm	White	1 Elbow Support
	7347713	S	8 <sup>1</sup> / <sub>4</sub> " – 9" 21 – 23 cm	White	1 Elbow Support
	7347714	M	9" – 9 <sup>3</sup> / <sub>4</sub> " 23 – 25 cm	White	1 Elbow Support
	7347715	L	9 <sup>3</sup> / <sub>4</sub> " – 10 <sup>1</sup> / <sub>2</sub> " 25 – 27 cm	White	1 Elbow Support
	7347716	XL	10 <sup>1</sup> / <sub>2</sub> " – 11 <sup>1</sup> / <sub>2</sub> " 27 – 29 cm	White	1 Elbow Support
	7347717	XXL	11 <sup>1</sup> / <sub>2</sub> " – 12 <sup>1</sup> / <sub>4</sub> " 29 – 31 cm	White	1 Elbow Support

<sup>1</sup> Tested for harmful substances according to OEKO-TEX® Standard 100. OEKO-TEX® certificate: 13.HTR.25112

# Actimove® EpiSport

## Elbow Clasp, Adjustable

**Low-profile anatomically shaped elbow clasp for adjustable, targeted muscle and tendon pressure with a secure and comfortable fit.**

- Relieves pain through individual and targeted dual compression on the radial and ulnar elbow tendons achieved by the anatomically shaped clasp and the adjustable VELCRO® brand<sup>1</sup> hook and loop strap.
- Optimal fit as the anatomically shaped clasp aligns to the natural form and allows physiological movements of the elbow.
- Inconspicuous to wear, as the low-profile design allows a smooth silhouette under shirt sleeves.
- Stays in place and comfortable to wear during sports, like tennis and golf, as well as everyday activities, thanks to the soft, non-slip inner liner.
- Easy one-handed application thanks to the elastic VELCRO® brand straps and buckle closure.
- Easy to keep clean for better hygiene, as the VELCRO® brand strap can be removed for cleaning and the clasp wiped down.
- Long-lasting, made from elastic, durable and flexible materials.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Available in five sizes.
- For stock keeping, long shelf life of five years.

**This product is intended to support the elbow.**

**Fields of application include:**

- Tennis elbow (lateral epicondylitis)
- Golfers' elbow (medical epicondylitis)



	Code	Size	Measurements	Color	Unit
	7347005	XS	7 <sup>3</sup> / <sub>4</sub> " – 9" 20 – 23 cm	Black	1 Elbow Clasp
	7347006	S	9" – 9 <sup>3</sup> / <sub>4</sub> " 23 – 25 cm	Black	1 Elbow Clasp
	7347007	M	9 <sup>3</sup> / <sub>4</sub> " – 11" 25 – 28 cm	Black	1 Elbow Clasp
	7347008	L	11" – 11 <sup>3</sup> / <sub>4</sub> " 28 – 30 cm	Black	1 Elbow Clasp
	7347009	XL	11 <sup>3</sup> / <sub>4</sub> " – 13 <sup>1</sup> / <sub>2</sub> " 30 – 34 cm	Black	1 Elbow Clasp

<sup>1</sup> VELCRO® is a registered trademark of Velcro BVBA. Used with permission.

# THIGH & CALF



# Actimove® PowerMotion

## Thigh Support

**Patented thigh support providing circular, gradient compression with a unique Y-shape strap design simulating professional taping techniques.**


- Aids pain relief and enhances healing by stimulating blood circulation through the unique Y-shaped strap design providing additional adjustable compression.
- Helps to reduce edema and enhances pain reduction by simulating taping techniques normally used for strains and fiber ruptures of the calf muscle.
- Promotes injured thigh muscle regeneration by controlling its elongation through targeted compression and Y-shaped straps.
- Comfortable to wear over the highly sensitive injured area thanks to the specially designed 3D flat-knitted Pressure-Relief-Zone.
- Secure fit due its anatomical shape and adjustable straps. Fast and easy to apply for the patient when compared to professional taping techniques, thanks to the individually adjustable closure strapping system.
- Provides patient comfort thanks to the breathable materials.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Easy to clean for better hygiene, machine washable at 40 °C/105 °F gentle wash cycle.
- Available in five sizes.
- Protected by patent, for more info visit: [www.essity.com/patents](http://www.essity.com/patents)
- For stock keeping, long shelf life of five years.

**This product is intended for thigh support.**

**Fields of application include:**

- Muscle strain and impact bruising
- Minor muscle fiber tears
- Myogelosis of the hamstring or quadriceps muscle



	Code	Size	Measurements	Unit
	7286706	S	16 <sup>1</sup> / <sub>8</sub> " – 17 <sup>3</sup> / <sub>4</sub> " 41 – 45 cm	1 Thigh Muscle Support
	7286707	M	17 <sup>3</sup> / <sub>4</sub> " – 19 <sup>3</sup> / <sub>4</sub> " 45 – 50 cm	1 Thigh Muscle Support
	7286708	L	19 <sup>3</sup> / <sub>4</sub> " – 22 <sup>1</sup> / <sub>8</sub> " 50 – 56 cm	1 Thigh Muscle Support
	7286709	XL	22 <sup>1</sup> / <sub>8</sub> " – 24 <sup>7</sup> / <sub>8</sub> " 56 – 63 cm	1 Thigh Muscle Support
	7286710	XXL	24 <sup>7</sup> / <sub>8</sub> " – 27 <sup>1</sup> / <sub>2</sub> " 63 – 70 cm	1 Thigh Muscle Support

# Actimove® PowerMotion

## Calf Support

**Patented calf support providing circular, gradient compression with a unique Y-shape strap design simulating professional taping techniques.**

- Aids pain relief and enhances healing by stimulating blood circulation through the unique Y-shaped strap design providing additional adjustable compression.
- Helps to reduce edema and enhances pain reduction by simulating taping techniques normally used for strains and fiber ruptures of the calf muscle.
- Promotes injured calf muscle regeneration by controlling its elongation through targeted compression and Y-shaped straps.
- Comfortable to wear over the highly sensitive injured area thanks to the specially designed 3D flat-knitted Pressure-Relief-Zone.
- Secure fit due its anatomical shape and adjustable straps.
- Fast and easy to apply for the patient when compared to professional taping techniques, thanks to the individually adjustable closure strapping system.
- Provides patient comfort thanks to the breathable materials.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Easy to clean for better hygiene, machine washable at 40 °C/105 °F gentle wash cycle.
- Available in three sizes.
- Protected by patent, for more info visit: [www.essity.com/patents](http://www.essity.com/patents)
- For stock keeping, long shelf life of five years.

**This product is intended for calf support.**

**Fields of application include:**

- Muscle strain and impact bruising
- Minor muscle fiber tears
- Myogelosis of the hamstring or quadriceps muscle



	Code	Size	Measurements	Unit
	7286604	S	11 <sup>1</sup> / <sub>2</sub> " – 13 <sup>3</sup> / <sub>8</sub> " 29 – 34 cm	1 Calf Muscle Support
	7286605	M	13 <sup>3</sup> / <sub>8</sub> " – 15 <sup>3</sup> / <sub>4</sub> " 34 – 40 cm	1 Calf Muscle Support
	7286606	L	15 <sup>3</sup> / <sub>4</sub> " – 18 <sup>1</sup> / <sub>2</sub> " 40 – 47 cm	1 Calf Muscle Support



# SPECIALITY SUPPORTS

# Actimove® Athletic Supporter

Comfortable protection for men during athletic activities.

- Comfortable support and protection during athletic, strenuous activities thanks to the stretchy pouch material.
- For more comfort, made with a wide, elastic waistband and soft, non-binding leg straps.
- For hygiene, hand-washable at 30°C/85°F.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Available in four sizes.

**Purpose:**

- All day protection and support. Suitable during sports and strenuous activities.



	Code	Size	Measurements	Color	Unit
	7551930	S	28" – 32" 71 – 81 cm	White	1 Athletic Supporter
	7551931	M	32" – 36" 81 – 91 cm	White	1 Athletic Supporter
	7551932	L	36" – 40" 91 – 102 cm	White	1 Athletic Supporter
	7551933	XL	40" – 46" 102 – 117 cm	White	1 Athletic Supporter

# Actimove® Hernia Support

## Support Belt

Single or double hernia support, wearable over undergarments.

- Helps to relieve pain as specially designed pressure pads provide targeted support to reduce inguinal hernias.
- Suitable for single or double hernias thanks to the removable right and left pads.
- Flexible, removable foam cushions designed to provide gentle, comfortable compression on the hernias.
- Easily adjustable for desired compression and fit thanks to its elastic waist band.
- For patient comfort, can be worn on top of underwear.
- For more comfort, adjustable design without metal buckles, nor snaps.
- For skin comfort, made with a soft, plush lining.
- For individual fit, made with adjustable leg straps.
- Inconspicuous wear under most clothing thanks to its neutral beige color.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Easy hygiene, hand, or machine washable in cold water with mild soap.
- Available in three sizes.

The product is intended to provide adjustable compression to inguinal hernias.

Fields of application include:

- Reducible single or double inguinal hernias



	Code	Size	Measurements	Color	Unit
	7455910	S	30" – 35" 76 – 89 cm	Beige	1 Hernia Support Belt
	7455911	M	35" – 41" 89 – 104 cm	Beige	1 Hernia Support Belt
	7455912	L	41" – 46" 104 – 117 cm	Beige	1 Hernia Support Belt

# Actimove® Hernia Support

## Support Brief

Single or double hernia support, wearable like underwear.


- Helps to relieve pain as specially designed pressure pads provide targeted support to reduce inguinal hernias.
- Flexible, removable foam cushions designed to provide gentle, comfortable compression on the hernias.
- Ideal amount of focused compression can easily be increased or decreased thanks to the adjustable straps.
- Suitable for single or double hernias thanks to the removable right and left pads.
- Individual support due to easily adjustable compression.
- For patient comfort and ease of wear combines underwear and hernia support in one.
- Extra wearing comfort from the soft, plush padding of the waistband.
- Comfortable, customized fit by adjusting the hook and loop straps.
- For comfort, made with soft, smooth materials.
- Inconspicuous wear under most clothing thanks to its low-profile design.
- For hygiene, hand-washable at 30°C / 85°F.
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Available in five sizes.

The product is intended to provide adjustable compression to inguinal hernias.

Fields of application include:

- Reducible single or double inguinal hernias



	Code	Size	Measurements	Color	Unit
	7551830	S	28" - 30" / 71 - 76 cm	White	1 Hernia Support Brief
	7551831	M	32" - 34" / 82 - 86 cm	White	1 Hernia Support Brief
	7551832	L	36" - 38" / 91 - 97 cm	White	1 Hernia Support Brief
	7551833	XL	40" - 42" / 102 - 107 cm	White	1 Hernia Support Brief
	7551834	XXL	44" - 46" / 112 - 117 cm	White	1 Hernia Support Brief

Code	Description	Size	Color	Page
7234850	Actimove® Manus – Wrist Stabilizer	S	Black	75
7234851	Actimove® Manus – Wrist Stabilizer	M	Black	75
7234852	Actimove® Manus – Wrist Stabilizer	L	Black	75
7234853	Actimove® Manus – Wrist Stabilizer	XL	Black	75
7245300	Actimove® Knee Stabilizer – Adjustable Horseshoe and Stays	XS	Black	20
7245301	Actimove® Knee Stabilizer – Adjustable Horseshoe and Stays	S	Black	20
7245302	Actimove® Knee Stabilizer – Adjustable Horseshoe and Stays	M	Black	20
7245303	Actimove® Knee Stabilizer – Adjustable Horseshoe and Stays	L	Black	20
7245304	Actimove® Knee Stabilizer – Adjustable Horseshoe and Stays	XL	Black	20
7245305	Actimove® Knee Stabilizer – Adjustable Horseshoe and Stays	XXL	Black	20
7245306	Actimove® Knee Stabilizer – Adjustable Horseshoe and Stays	XXXL	Black	20
7245307	Actimove® Knee Stabilizer – Adjustable Horseshoe and Stays	XXXXL	Black	20
7279130	Actimove® Deluxe Lumbar Sacral Support – Rigid Panel, For Women	S	Rose	113
7279131	Actimove® Deluxe Lumbar Sacral Support – Rigid Panel, For Women	M	Rose	113
7279132	Actimove® Deluxe Lumbar Sacral Support – Rigid Panel, For Women	L	Rose	113
7279133	Actimove® Deluxe Lumbar Sacral Support – Rigid Panel, For Women	XL	Rose	113
7279134	Actimove® Deluxe Lumbar Sacral Support – Rigid Panel, For Women	XXL	Rose	113
7279143	Actimove® Deluxe Lumbar Sacral Support – Rigid Panel, For Women	XXXL	Rose	113
7279144	Actimove® Deluxe Lumbar Sacral Support – Rigid Panel, For Women	XXXXL	Rose	113
7279145	Actimove® Deluxe Lumbar Sacral Support – Rigid Panel, For Women	XXXXXL	Rose	113
7279230	Actimove® Deluxe Lumbar Sacral Support – Rigid Panel, For Women	S	Lavender	113
7279231	Actimove® Deluxe Lumbar Sacral Support – Rigid Panel, For Women	M	Lavender	113
7279232	Actimove® Deluxe Lumbar Sacral Support – Rigid Panel, For Women	L	Lavender	113
7279233	Actimove® Deluxe Lumbar Sacral Support – Rigid Panel, For Women	XL	Lavender	113
7279234	Actimove® Deluxe Lumbar Sacral Support – Rigid Panel, For Women	XXL	Lavender	113
7279243	Actimove® Deluxe Lumbar Sacral Support – Rigid Panel, For Women	XXXL	Lavender	113
7279244	Actimove® Deluxe Lumbar Sacral Support – Rigid Panel, For Women	XXXXL	Lavender	113
7279245	Actimove® Deluxe Lumbar Sacral Support – Rigid Panel, For Women	XXXXXL	Lavender	113
7281908	Actimove® Genu Eco – Knee Immobilizer	S/M – 20"	Blue	85
7281909	Actimove® Genu Eco – Knee Immobilizer	M/L – 20"	Blue	85
7281910	Actimove® Genu Eco – Knee Immobilizer	S/M – 24"	Blue	85
7281911	Actimove® Genu Eco – Knee Immobilizer	M/L – 24"	Blue	85
7281912	Actimove® Genu Eco – Knee Immobilizer	S/M – 28"	Blue	85
7281913	Actimove® Genu Eco – Knee Immobilizer	M/L – 28"	Blue	85
7281914	Actimove® Genu Tri-Panel – Knee Immobilizer	S/M – 20"	Blue	83
7281915	Actimove® Genu Tri-Panel – Knee Immobilizer	M/L – 20"	Blue	83
7281916	Actimove® Genu Tri-Panel – Knee Immobilizer	S/M – 24"	Blue	83
7281917	Actimove® Genu Tri-Panel – Knee Immobilizer	M/L – 24"	Blue	83
7281939	Actimove® Umerus Comfort – Comfort Shoulder Immobilizer	S	Blue	121
7281940	Actimove® Umerus Comfort – Comfort Shoulder Immobilizer	M	Blue	121
7281941	Actimove® Umerus Comfort – Comfort Shoulder Immobilizer	L	Blue	121
7281942	Actimove® Umerus Comfort – Comfort Shoulder Immobilizer	XL	Blue	121
7281961	Actimove® Mitella Comfort – Comfort Arm Sling	S	Blue	119
7281962	Actimove® Mitella Comfort – Comfort Arm Sling	M	Blue	119
7281963	Actimove® Mitella Comfort – Comfort Arm Sling	L	Blue	119
7281964	Actimove® Mitella Comfort – Comfort Arm Sling	XL	Blue	119
7281994	Actimove® Umerus Comfort – Comfort Shoulder Immobilizer	XS	Blue	121
7285918	Actimove® Sling Comfort – Arm Sling – Single Pack	Universal	Blue	117
7285924	Actimove® Sling Comfort – Arm Sling – Single Pack – Single use	Universal	Blue	117
7285936	Actimove® Cervical Comfort – Comfort Cervical Collar	XS	Blue	127
7285937	Actimove® Cervical Comfort – Comfort Cervical Collar	S	Blue	127
7285938	Actimove® Cervical Comfort – Comfort Cervical Collar	M	Blue	127



Code	Description	Size	Color	Page
7285939	Actimove® Cervical Comfort – Comfort Cervical Collar	L	Blue	127
7285940	Actimove® Cervical Comfort – Comfort Cervical Collar	XL	Blue	127
7285941	Actimove® Cervical Comfort – Comfort Cervical Collar	XL – Short	Blue	127
7285945	Actimove® Genu Clips – Knee Immobilizer	Pediatric – 12”	Blue	85
7285946	Actimove® Genu Clips – Knee Immobilizer	Pediatric – 16”	Blue	85
7285947	Actimove® Genu Clips – Knee Immobilizer	S/M – 20”	Blue	85
7285948	Actimove® Genu Clips – Knee Immobilizer	M/L – 20”	Blue	85
7285949	Actimove® Genu Clips – Knee Immobilizer	S/M – 24”	Blue	85
7285950	Actimove® Genu Clips – Knee Immobilizer	M/L – 24”	Blue	85
7285951	Actimove® Genu Clips – Knee Immobilizer	S/M – 28”	Blue	85
7285952	Actimove® Genu Clips – Knee Immobilizer	M/L – 28”	Blue	85
7285990	Actimove® Gauntlet – Wrist & Thumb Stabilizer	S	Black	80
7285991	Actimove® Gauntlet – Wrist & Thumb Stabilizer	M	Black	80
7285992	Actimove® Gauntlet – Wrist & Thumb Stabilizer	L	Black	80
7285993	Actimove® Gauntlet – Wrist & Thumb Stabilizer	XL	Black	80
7286604	Actimove® PowerMotion – Calf Support	S	Black	143
7286605	Actimove® PowerMotion – Calf Support	M	Black	143
7286606	Actimove® PowerMotion – Calf Support	L	Black	143
7286706	Actimove® PowerMotion – Thigh Support	S	Black	142
7286707	Actimove® PowerMotion – Thigh Support	M	Black	142
7286708	Actimove® PowerMotion – Thigh Support	L	Black	142
7286709	Actimove® PowerMotion – Thigh Support	XL	Black	142
7286710	Actimove® PowerMotion – Thigh Support	XXL	Black	142
7311800	Actimove® Knee Brace – Wrap Around, Polycentric Hinges, Condyle Pads	XS	Black	22
7311801	Actimove® Knee Brace – Wrap Around, Polycentric Hinges, Condyle Pads	S	Black	22
7311802	Actimove® Knee Brace – Wrap Around, Polycentric Hinges, Condyle Pads	M	Black	22
7311803	Actimove® Knee Brace – Wrap Around, Polycentric Hinges, Condyle Pads	L	Black	22
7311804	Actimove® Knee Brace – Wrap Around, Polycentric Hinges, Condyle Pads	XL	Black	22
7311805	Actimove® Knee Brace – Wrap Around, Polycentric Hinges, Condyle Pads	XXL	Black	22
7311806	Actimove® Knee Brace – Wrap Around, Polycentric Hinges, Condyle Pads	XXXL	Black	22
7311828	Actimove® Knee Brace – Wrap Around, Polycentric Hinges, Condyle Pads	XXXXL	Black	22
7311900	Actimove® TaloCast AirGel – Ankle Brace	S/M Trainer	Black	97
7311901	Actimove® TaloCast AirGel – Ankle Brace	L/XL Standard	Black	97
7311904	Actimove® TaloCast Air – Ankle Brace	S/M Trainer	Black	96
7311905	Actimove® TaloCast Air – Ankle Brace	L/XL Standard	Black	96
7311908	Actimove® TaloCast – Ankle Brace	S/M Trainer	Black	98
7311909	Actimove® TaloCast – Ankle Brace	L/XL Standard	Black	98
7311912	Actimove® TaloBrace – Comfort Ankle Brace	Left Side S/M Trainer	Black	94
7311913	Actimove® TaloBrace – Comfort Ankle Brace	Right Side S/M Trainer	Black	94
7311914	Actimove® TaloBrace – Comfort Ankle Brace	Left Side L/XL Standard	Black	94
7311915	Actimove® TaloBrace – Comfort Ankle Brace	Right Side L/XL Standard	Black	94
7311916	Actimove® TaloBrace Air – Comfort Ankle Brace	Left Side S/M Trainer	Black	95
7311917	Actimove® TaloBrace Air – Comfort Ankle Brace	Right Side S/M Trainer	Black	95
7311918	Actimove® TaloBrace Air – Comfort Ankle Brace	Left Side L/XL Standard	Black	95
7311919	Actimove® TaloBrace Air – Comfort Ankle Brace	Right Side L/XL Standard	Black	95
7341480	Actimove® Ankle Support – Elastic Wrap Around	S	Black	27

Code	Description	Size	Color	Page
7341481	Actimove® Ankle Support – Elastic Wrap Around	M	Black	27
7341482	Actimove® Ankle Support – Elastic Wrap Around	L	Black	27
7341492	Actimove® Ankle Support – Elastic Wrap Around	XL	Black	27
7341680	Actimove® Wrist Support – Elastic Wrap Around	Universal	Black	12
7342812	Actimove® LumbaMotion – Back Support with Pressure Pad	XS	White	109
7342813	Actimove® LumbaMotion – Back Support with Pressure Pad	S	White	109
7342814	Actimove® LumbaMotion – Back Support with Pressure Pad	M	White	109
7342815	Actimove® LumbaMotion – Back Support with Pressure Pad	L	White	109
7342816	Actimove® LumbaMotion – Back Support with Pressure Pad	XL	White	109
7342817	Actimove® LumbaMotion – Back Support with Pressure Pad	XS	Charcoal	109
7342882	Actimove® LumbaMotion – Back Support with Pressure Pad	S	Charcoal	109
7342819	Actimove® LumbaMotion – Back Support with Pressure Pad	M	Charcoal	109
7342820	Actimove® LumbaMotion – Back Support with Pressure Pad	L	Charcoal	109
7342821	Actimove® LumbaMotion – Back Support with Pressure Pad	XL	Charcoal	109
7342912	Actimove® LumbaMotion Lady – Back Support with Pressure Pad	XS	White	111
7342913	Actimove® LumbaMotion Lady – Back Support with Pressure Pad	S	White	111
7342914	Actimove® LumbaMotion Lady – Back Support with Pressure Pad	M	White	111
7342915	Actimove® LumbaMotion Lady – Back Support with Pressure Pad	L	White	111
7342916	Actimove® LumbaMotion Lady – Back Support with Pressure Pad	XL	White	111
7344604	Actimove® Shoulder Abduction Sling – Shoulder Immobilizer	S	Blue	122
7344605	Actimove® Shoulder Abduction Sling – Shoulder Immobilizer	M	Blue	122
7344606	Actimove® Shoulder Abduction Sling – Shoulder Immobilizer	L	Blue	122
7344607	Actimove® Shoulder Abduction Sling – Shoulder Immobilizer	XL	Blue	122
7345030	Actimove® Back Support – High-Density Foam Panel Adjustable Double Layer Compression	S/M	Silver	46
7345031	Actimove® Back Support – High-Density Foam Panel Adjustable Double Layer Compression	L/XL	Silver	46
7345036	Actimove® Back Support – High-Density Foam Panel Adjustable Double Layer Compression	XXL/XXXL	Silver	46
7346600	Actimove® Y-Tab Hook Closure	-	Blue	117
7346807	Actimove® GenuMotion – Knee Support	XXL+	White	87
7346808	Actimove® GenuMotion – Knee Support	XL+	White	87
7346809	Actimove® GenuMotion – Knee Support	L+	White	87
7346810	Actimove® GenuMotion – Knee Support	XXS	White	87
7346811	Actimove® GenuMotion – Knee Support	XS	White	87
7346812	Actimove® GenuMotion – Knee Support	S	White	87
7346813	Actimove® GenuMotion – Knee Support	M	White	87
7346814	Actimove® GenuMotion – Knee Support	L	White	87
7346815	Actimove® GenuMotion – Knee Support	XL	White	87
7346816	Actimove® GenuMotion – Knee Support	XXL	White	87
7346817	Actimove® GenuMotion – Knee Support	M+	White	87
7346843	Actimove® GenuMotion – Knee Support	XXS	Charcoal	87
7346844	Actimove® GenuMotion – Knee Support	XS	Charcoal	87
7346845	Actimove® GenuMotion – Knee Support	S	Charcoal	87
7346846	Actimove® GenuMotion – Knee Support	M	Charcoal	87
7346847	Actimove® GenuMotion – Knee Support	L	Charcoal	87
7346848	Actimove® GenuMotion – Knee Support	XL	Charcoal	87
7346849	Actimove® GenuMotion – Knee Support	XXL	Charcoal	87
7346850	Actimove® GenuMotion – Knee Support	XXXL	White	87
7346851	Actimove® GenuMotion – Knee Support	XXXXL	White	87
7346852	Actimove® GenuMotion – Knee Support	XXXL+	White	87
7346853	Actimove® GenuMotion – Knee Support	XXXXL+	White	87

Code	Description	Size	Color	Page
7347005	Actimove® EpiSport – Elbow Clasp, Adjustable	XS	Black	140
7347006	Actimove® EpiSport – Elbow Clasp, Adjustable	S	Black	140
7347007	Actimove® EpiSport – Elbow Clasp, Adjustable	M	Black	140
7347008	Actimove® EpiSport – Elbow Clasp, Adjustable	L	Black	140
7347009	Actimove® EpiSport – Elbow Clasp, Adjustable	XL	Black	140
7347711	Actimove® EpiMotion – Elbow Support	XXS	White	139
7347712	Actimove® EpiMotion – Elbow Support	XS	White	139
7347713	Actimove® EpiMotion – Elbow Support	S	White	139
7347714	Actimove® EpiMotion – Elbow Support	M	White	139
7347715	Actimove® EpiMotion – Elbow Support	L	White	139
7347716	Actimove® EpiMotion – Elbow Support	XL	White	139
7347717	Actimove® EpiMotion – Elbow Support	XXL	White	139
7348201	Actimove® Manus Forte – Wrist Brace	Right XS	Black	71
7348202	Actimove® Manus Forte – Wrist Brace	Left XS	Black	71
7348203	Actimove® Manus Forte – Wrist Brace	Right S/M	Black	71
7348204	Actimove® Manus Forte – Wrist Brace	Left S/M	Black	71
7348205	Actimove® Manus Forte – Wrist Brace	Right L/XL	Black	71
7348206	Actimove® Manus Forte – Wrist Brace	Left L/XL	Black	71
7348707	Actimove® TaloMotion – Ankle Support	Right XXL	White	91
7348712	Actimove® TaloMotion – Ankle Support	Right XS	White	91
7348713	Actimove® TaloMotion – Ankle Support	Right S	White	91
7348714	Actimove® TaloMotion – Ankle Support	Right M	White	91
7348715	Actimove® TaloMotion – Ankle Support	Right L	White	91
7348716	Actimove® TaloMotion – Ankle Support	Right XL	White	91
7348737	Actimove® TaloMotion – Ankle Support	Left XXL	White	91
7348742	Actimove® TaloMotion – Ankle Support	Left XS	White	91
7348743	Actimove® TaloMotion – Ankle Support	Left S	White	91
7348744	Actimove® TaloMotion – Ankle Support	Left M	White	91
7348745	Actimove® TaloMotion – Ankle Support	Left L	White	91
7348746	Actimove® TaloMotion – Ankle Support	Left XL	White	91
7348760	Actimove® TaloMotion – Ankle Support	Right XS	Charcoal	91
7348761	Actimove® TaloMotion – Ankle Support	Left XS	Charcoal	91
7348762	Actimove® TaloMotion – Ankle Support	Right S	Charcoal	91
7348763	Actimove® TaloMotion – Ankle Support	Left S	Charcoal	91
7348764	Actimove® TaloMotion – Ankle Support	Right M	Charcoal	91
7348765	Actimove® TaloMotion – Ankle Support	Left M	Charcoal	91
7348766	Actimove® TaloMotion – Ankle Support	Right L	Charcoal	91
7348767	Actimove® TaloMotion – Ankle Support	Left L	Charcoal	91
7348768	Actimove® TaloMotion – Ankle Support	Right XL	Charcoal	91
7348769	Actimove® TaloMotion – Ankle Support	Left XL	Charcoal	91
7348770	Actimove® TaloMotion – Ankle Support	Right XXL	Charcoal	91
7348771	Actimove® TaloMotion – Ankle Support	Left XXL	Charcoal	91
7349603	Actimove® Manus Forte Plus – Wrist & Thumb Brace	XS	Black	77
7349604	Actimove® Manus Forte Plus – Wrist & Thumb Brace	S/M	Black	77
7349606	Actimove® Manus Forte Plus – Wrist & Thumb Brace	L/XL	Black	77
7349633	Actimove® Manus Forte Plus – Wrist & Thumb Brace	XS	Black	77
7349634	Actimove® Manus Forte Plus – Wrist & Thumb Brace	S/M	Black	77
7349636	Actimove® Manus Forte Plus – Wrist & Thumb Brace	L/XL	Black	77
7349712	Actimove® ManuMotion – Wrist Support	Right XS	White	69
7349713	Actimove® ManuMotion – Wrist Support	Right S	White	69
7349714	Actimove® ManuMotion – Wrist Support	Right M	White	69
7349715	Actimove® ManuMotion – Wrist Support	Right L	White	69

Code	Description	Size	Color	Page
7349716	Actimove® ManuMotion – Wrist Support	Right XL	White	69
7349717	Actimove® ManuMotion – Wrist Support	Right XXL	White	69
7349738	Actimove® ManuMotion – Wrist Support	Right XS	Charcoal	69
7349739	Actimove® ManuMotion – Wrist Support	Left XS	Charcoal	69
7349742	Actimove® ManuMotion – Wrist Support	Left XS	White	69
7349743	Actimove® ManuMotion – Wrist Support	Left S	White	69
7349744	Actimove® ManuMotion – Wrist Support	Left M	White	69
7349745	Actimove® ManuMotion – Wrist Support	Left L	White	69
7349746	Actimove® ManuMotion – Wrist Support	Left XL	White	69
7349747	Actimove® ManuMotion – Wrist Support	Left XXL	White	69
7349760	Actimove® ManuMotion – Wrist Support	Right S	Charcoal	69
7349761	Actimove® ManuMotion – Wrist Support	Left S	Charcoal	69
7349762	Actimove® ManuMotion – Wrist Support	Right M	Charcoal	69
7349763	Actimove® ManuMotion – Wrist Support	Left M	Charcoal	69
7349764	Actimove® ManuMotion – Wrist Support	Right L	Charcoal	69
7349765	Actimove® ManuMotion – Wrist Support	Left L	Charcoal	69
7349766	Actimove® ManuMotion – Wrist Support	Right XL	Charcoal	69
7349767	Actimove® ManuMotion – Wrist Support	Left XL	Charcoal	69
7349768	Actimove® ManuMotion – Wrist Support	Right XXL	Charcoal	69
7349769	Actimove® ManuMotion – Wrist Support	Left XXL	Charcoal	69
7455910	Actimove® Hernia Support – Support Belt	S	Beige	146
7455911	Actimove® Hernia Support – Support Belt	M	Beige	146
7455912	Actimove® Hernia Support – Support Belt	L	Beige	146
7456113	Actimove® Mitella – Arm Sling	Pediatric	Navy	61
7456330	Actimove® LO•PLUS™ Lumbar Support – Stabilizing Posterior Panel	S	Black	115
7456331	Actimove® LO•PLUS™ Lumbar Support – Stabilizing Posterior Panel	M	Black	115
7456332	Actimove® LO•PLUS™ Lumbar Support – Stabilizing Posterior Panel	L	Black	115
7456333	Actimove® LO•PLUS™ Lumbar Support – Stabilizing Posterior Panel	XL	Black	115
7456334	Actimove® LO•PLUS™ Lumbar Support – Stabilizing Posterior Panel	XXL	Black	115
7456335	Actimove® LO•PLUS™ Lumbar Support – Stabilizing Posterior Panel	XXXL	Black	115
7550010	Actimove® Knee Brace – Wrap Around, Simple Hinges, Condyle Pads	XS	Black	21
7550011	Actimove® Knee Brace – Wrap Around, Simple Hinges, Condyle Pads	S	Black	21
7550012	Actimove® Knee Brace – Wrap Around, Simple Hinges, Condyle Pads	M	Black	21
7550013	Actimove® Knee Brace – Wrap Around, Simple Hinges, Condyle Pads	L	Black	21
7550014	Actimove® Knee Brace – Wrap Around, Simple Hinges, Condyle Pads	XL	Black	21
7550015	Actimove® Knee Brace – Wrap Around, Simple Hinges, Condyle Pads	XXL	Black	21
7550016	Actimove® Knee Brace – Wrap Around, Simple Hinges, Condyle Pads	XXXL	Black	21
7550804	Actimove® Rib Belt – Two-Panel, For Men	S	White	134
7550805	Actimove® Rib Belt – Two-Panel, For Women	S	White	135
7550806	Actimove® Rib Belt – Two-Panel, For Men	UN/L	White	134
7550807	Actimove® Rib Belt – Two-Panel, For Women	UN/L	White	135
7551306	Actimove® Abdominal Binder – Three-Panel	XS	White	132
7551307	Actimove® Abdominal Binder – Three-Panel	S	White	132
7551308	Actimove® Abdominal Binder – Three-Panel	M	White	132
7551309	Actimove® Abdominal Binder – Three-Panel	L	White	132
7551310	Actimove® Abdominal Binder – Four-Panel	XS	White	132
7551311	Actimove® Abdominal Binder – Four-Panel	S	White	132
7551312	Actimove® Abdominal Binder – Four-Panel	M	White	132
7551313	Actimove® Abdominal Binder – Four-Panel	L	White	132
7551413	Actimove® Abdominal Binder Comfort – Soft Pad	XS	White	133
7551414	Actimove® Abdominal Binder Comfort – Soft Pad	S	White	133
7551415	Actimove® Abdominal Binder Comfort – Soft Pad	M	White	133

Code	Description	Size	Color	Page
7551416	Actimove® Abdominal Binder Comfort – Soft Pad	L	White	133
7551417	Actimove® Abdominal Binder Comfort – Soft Pad	XS	White	133
7551418	Actimove® Abdominal Binder Comfort – Soft Pad	S	White	133
7551419	Actimove® Abdominal Binder Comfort – Soft Pad	M	White	133
7551420	Actimove® Abdominal Binder Comfort – Soft Pad	L	White	133
7551830	Actimove® Hernia Support – Support Brief	S	White	147
7551831	Actimove® Hernia Support – Support Brief	M	White	147
7551832	Actimove® Hernia Support – Support Brief	L	White	147
7551833	Actimove® Hernia Support – Support Brief	XL	White	147
7551834	Actimove® Hernia Support – Support Brief	XXL	White	147
7551930	Actimove® Athletic Supporter	S	White	145
7551931	Actimove® Athletic Supporter	M	White	145
7551932	Actimove® Athletic Supporter	L	White	145
7551933	Actimove® Athletic Supporter	XL	White	145
7552714	Actimove® Lumbar Sacral Support	S	White	45
7552715	Actimove® Lumbar Sacral Support	M	White	45
7552716	Actimove® Lumbar Sacral Support	L	White	45
7552717	Actimove® Lumbar Sacral Support	XL	White	45
7552718	Actimove® Lumbar Sacral Support	XXL	White	45
7552719	Actimove® Lumbar Sacral Support	XXXL	White	45
7552720	Actimove® Lumbar Sacral Support	S	White	45
7552721	Actimove® Lumbar Sacral Support	M	White	45
7552722	Actimove® Lumbar Sacral Support	L	White	45
7552723	Actimove® Lumbar Sacral Support	XL	White	45
7552724	Actimove® Lumbar Sacral Support	XXL	White	45
7552725	Actimove® Lumbar Sacral Support	XXXL	White	45
7552943	Actimove® Lumbar Sacral Support Comfort – Additional Support Belt	S	White	107
7552949	Actimove® Lumbar Sacral Support Comfort – Additional Support Belt	S	White	107
7552950	Actimove® Lumbar Sacral Support Comfort – Additional Support Belt	M	White	107
7552951	Actimove® Lumbar Sacral Support Comfort – Additional Support Belt	L	White	107
7552952	Actimove® Lumbar Sacral Support Comfort – Additional Support Belt	XL	White	107
7552953	Actimove® Lumbar Sacral Support Comfort – Additional Support Belt	XXL	White	107
7552954	Actimove® Lumbar Sacral Support Comfort – Additional Support Belt	XXXL	White	107
7552955	Actimove® Lumbar Sacral Support Comfort – Additional Support Belt	Universal	White	107
7553510	Actimove® Back Stabilizer – Rigid Panel with Pressure Pads SMART Easy-Closing-Pulley-System	XS/S	Black	31
7553511	Actimove® Back Stabilizer – Rigid Panel with Pressure Pads SMART Easy-Closing-Pulley-System	S/M	Black	31
7553512	Actimove® Back Stabilizer – Rigid Panel with Pressure Pads SMART Easy-Closing-Pulley-System	L/XL	Black	31
7553513	Actimove® Back Stabilizer – Rigid Panel with Pressure Pads SMART Easy-Closing-Pulley-System	XXL/XXXL	Black	31
7554120	Actimove® Back Support, 4 Stays – Adjustable Double Layer Compression	S	Black	30
7554121	Actimove® Back Support, 4 Stays – Adjustable Double Layer Compression	M	Black	30
7554122	Actimove® Back Support, 4 Stays – Adjustable Double Layer Compression	L	Black	30
7554123	Actimove® Back Support, 4 Stays – Adjustable Double Layer Compression	XL	Black	30
7554124	Actimove® Back Support, 4 Stays – Adjustable Double Layer Compression	XXL	Black	30
7554125	Actimove® Back Support, 4 Stays – Adjustable Double Layer Compression	XXXL	Black	30
7554730	Actimove® Cervical – Cervical Collar Soft Density, Narrow Width	Universal	White	125
7554731	Actimove® Cervical – Cervical Collar Soft Density, Narrow Width	S	White	125
7554732	Actimove® Cervical – Cervical Collar Soft Density, Narrow Width	M	White	125
7554733	Actimove® Cervical – Cervical Collar Soft Density, Narrow Width	L	White	125
7554734	Actimove® Cervical – Cervical Collar Soft Density, Narrow Width	XL	White	125



Code	Description	Size	Color	Page
7555131	Actimove® Cervical – Cervical Collar Soft Density, Regular Width	Universal	White	125
7555134	Actimove® Cervical – Cervical Collar Soft Density, Regular Width	S	White	125
7555135	Actimove® Cervical – Cervical Collar Soft Density, Regular Width	M	White	125
7555136	Actimove® Cervical – Cervical Collar Soft Density, Regular Width	L	White	125
7555137	Actimove® Cervical – Cervical Collar Soft Density, Regular Width	XL	White	125
7555530	Actimove® Cervical – Cervical Collar	Infant	White	62
7555531	Actimove® Cervical – Cervical Collar	(0 – 24 months)	White	62
7555532	Actimove® Cervical – Cervical Collar	Pediatric	White	62
7557310	Actimove® Shoulder Support – Extra Pocket for Optional Hot/Cold Pack	S	Black	34
7557311	Actimove® Shoulder Support – Extra Pocket for Optional Hot/Cold Pack	M	Black	34
7557312	Actimove® Shoulder Support – Extra Pocket for Optional Hot/Cold Pack	L	Black	34
7557520	Actimove® Knee Support – Open Patella, 4 Stays	S	Beige	41
7557521	Actimove® Knee Support – Open Patella, 4 Stays	M	Beige	41
7557522	Actimove® Knee Support – Open Patella, 4 Stays	L	Beige	41
7557523	Actimove® Knee Support – Open Patella, 4 Stays	XL	Beige	41
7557528	Actimove® Knee Support – Closed Patella, 2 Stays	S	Beige	40
7557529	Actimove® Knee Support – Closed Patella, 2 Stays	M	Beige	40
7557530	Actimove® Knee Support – Closed Patella, 2 Stays	L	Beige	40
7557531	Actimove® Knee Support – Closed Patella, 2 Stays	XL	Beige	40
7557536	Actimove® Knee Support – Closed Patella	S	Beige	39
7557537	Actimove® Knee Support – Closed Patella	M	Beige	39
7557538	Actimove® Knee Support – Closed Patella	L	Beige	39
7557539	Actimove® Knee Support – Closed Patella	XL	Beige	39
7558020	Actimove® Mild Knee Support	S	White	38
7558021	Actimove® Mild Knee Support	M	White	38
7558022	Actimove® Mild Knee Support	L	White	38
7558023	Actimove® Mild Knee Support	XL	White	38
7558510	Actimove® Knee Support – Open Patella	Pediatric	Navy	59
7558511	Actimove® Knee Support – Open Patella	Youth	Navy	59
7558512	Actimove® Knee Support – Open Patella	XS	Navy	18
7558513	Actimove® Knee Support – Open Patella	S	Navy	18
7558514	Actimove® Knee Support – Open Patella	M	Navy	18
7558515	Actimove® Knee Support – Open Patella	L	Navy	18
7558516	Actimove® Knee Support – Open Patella	XL	Navy	18
7558517	Actimove® Knee Support – Open Patella	XXL	Navy	18
7558518	Actimove® Knee Support – Open Patella	XXXL	Navy	18
7558519	Actimove® Knee Support – Open Patella	XS	Black	18
7558520	Actimove® Knee Support – Open Patella	S	Black	18
7558521	Actimove® Knee Support – Open Patella	M	Black	18
7558522	Actimove® Knee Support – Open Patella	L	Black	18
7558523	Actimove® Knee Support – Open Patella	XL	Black	18
7558524	Actimove® Knee Support – Open Patella	XXL	Black	18
7558525	Actimove® Knee Support – Open Patella	XXXL	Black	18
7558548	Actimove® Knee Support – Open Patella	XXXXL	Black	18
7558620	Actimove® Knee Support – Closed Patella	XS	Navy	17
7558621	Actimove® Knee Support – Closed Patella	S	Navy	17
7558622	Actimove® Knee Support – Closed Patella	M	Navy	17
7558623	Actimove® Knee Support – Closed Patella	L	Navy	17
7558624	Actimove® Knee Support – Closed Patella	XL	Navy	17
7558625	Actimove® Knee Support – Closed Patella	XXL	Navy	17
7558626	Actimove® Knee Support – Closed Patella	XXXL	Navy	17
7558627	Actimove® Knee Support – Closed Patella	XS	Black	17
7558628	Actimove® Knee Support – Closed Patella	S	Black	17

Code	Description	Size	Color	Page
7558629	Actimove® Knee Support – Closed Patella	M	Black	17
7558630	Actimove® Knee Support – Closed Patella	L	Black	17
7558631	Actimove® Knee Support – Closed Patella	XL	Black	17
7558632	Actimove® Knee Support – Closed Patella	XXL	Black	17
7558633	Actimove® Knee Support – Closed Patella	XXXL	Black	17
7558920	Actimove® Patella Strap – Adjustable	Universal	Black	15
7559110	Actimove® Dual Knee Strap – Adjustable	S	Black	16
7559111	Actimove® Dual Knee Strap – Adjustable	M	Black	16
7559112	Actimove® Dual Knee Strap – Adjustable	L	Black	16
7559113	Actimove® Dual Knee Strap – Adjustable	XL	Black	16
7559310	Actimove® Knee Support – Open Patella Adjustable	Universal	Black	19
7559313	Actimove® Knee Support – Open Patella Adjustable	Universal Plus	Black	19
7559930	Actimove® PF Knee Brace – Lateral Support Simple Hinges, Condyle Pads, J-shaped Buttress	Right XS	Black	24
7559931	Actimove® PF Knee Brace – Lateral Support Simple Hinges, Condyle Pads, J-shaped Buttress	Right S	Black	24
7559932	Actimove® PF Knee Brace – Lateral Support Simple Hinges, Condyle Pads, J-shaped Buttress	Right M	Black	24
7559933	Actimove® PF Knee Brace – Lateral Support Simple Hinges, Condyle Pads, J-shaped Buttress	Right L	Black	24
7559934	Actimove® PF Knee Brace – Lateral Support Simple Hinges, Condyle Pads, J-shaped Buttress	Right XL	Black	24
7559935	Actimove® PF Knee Brace – Lateral Support Simple Hinges, Condyle Pads, J-shaped Buttress	Right XXL	Black	24
7559936	Actimove® PF Knee Brace – Lateral Support Simple Hinges, Condyle Pads, J-shaped Buttress	Left XXXL	Black	24
7559937	Actimove® PF Knee Brace – Lateral Support Simple Hinges, Condyle Pads, J-shaped Buttress	Left XS	Black	24
7559938	Actimove® PF Knee Brace – Lateral Support Simple Hinges, Condyle Pads, J-shaped Buttress	Left S	Black	24
7559939	Actimove® PF Knee Brace – Lateral Support Simple Hinges, Condyle Pads, J-shaped Buttress	Left M	Black	24
7559940	Actimove® PF Knee Brace – Lateral Support Simple Hinges, Condyle Pads, J-shaped Buttress	Left L	Black	24
7559941	Actimove® PF Knee Brace – Lateral Support Simple Hinges, Condyle Pads, J-shaped Buttress	Left XL	Black	24
7559942	Actimove® PF Knee Brace – Lateral Support Simple Hinges, Condyle Pads, J-shaped Buttress	Left XXL	Black	24
7559943	Actimove® PF Knee Brace – Lateral Support Simple Hinges, Condyle Pads, J-shaped Buttress	Right XXXL	Black	24
7560202	Actimove® Thigh Support – Adjustable	Universal	Black	35
7560320	Actimove® Mild Ankle Support	S	White	42
7560321	Actimove® Mild Ankle Support	M	White	42
7560322	Actimove® Mild Ankle Support	L	White	42
7560323	Actimove® Mild Ankle Support	XL	White	42
7560610	Actimove® Ankle Support	Pediatric	Navy	60
7560611	Actimove® Ankle Support	Youth	Navy	60
7560612	Actimove® Ankle Support	S	Navy	25
7560613	Actimove® Ankle Support	M	Navy	25
7560614	Actimove® Ankle Support	L	Navy	25
7560615	Actimove® Ankle Support	XL	Navy	25
7560616	Actimove® Ankle Support	S	Black	25
7560617	Actimove® Ankle Support	M	Black	25
7560618	Actimove® Ankle Support	L	Black	25
7560619	Actimove® Ankle Support	XL	Black	25

Code	Description	Size	Color	Page
7560626	Actimove® Ankle Support	XXL	Black	25
7560820	Actimove® Ankle Support	S	Beige	43
7560821	Actimove® Ankle Support	M	Beige	43
7560822	Actimove® Ankle Support	L	Beige	43
7560823	Actimove® Ankle Support	XL	Beige	43
7561130	Actimove® Ankle Stabilizer – Criss-Cross Straps	Universal	Black	29
7561132	Actimove® Ankle Stabilizer – Criss-Cross Straps	Universal Plus	Black	29
7561430	Actimove® Elbow Support	S	Beige	48
7561431	Actimove® Elbow Support	M	Beige	48
7561432	Actimove® Elbow Support	L	Beige	48
7561433	Actimove® Elbow Support	XL	Beige	48
7561438	Actimove® Elbow Support – Pressure Pads and Strap	S	Beige	49
7561439	Actimove® Elbow Support – Pressure Pads and Strap	M	Beige	49
7561440	Actimove® Elbow Support – Pressure Pads and Strap	L	Beige	49
7561441	Actimove® Elbow Support – Pressure Pads and Strap	XL	Beige	49
7561620	Actimove® Mild Elbow Support	S	White	47
7561621	Actimove® Mild Elbow Support	M	White	47
7561622	Actimove® Mild Elbow Support	L	White	47
7561623	Actimove® Mild Elbow Support	XL	White	47
7561730	Actimove® Elbow Support – Adjustable	Universal	Black	33
7561731	Actimove® Elbow Support – Adjustable	Pediatric	Black	63
7561734	Actimove® Elbow Support – Adjustable	Universal Plus	Black	33
7562010	Actimove® Ankle Support – Adjustable	Universal	Black	26
7562610	Actimove® Wrist Support – Adjustable	Universal	Black	11
7563630	Actimove® Thumb Stabilizer – Extra Stays	S/M	Black	14
7563631	Actimove® Thumb Stabilizer – Extra Stays	L/XL	Black	14
7563632	Actimove® Thumb Stabilizer – Extra Stays	S/M	Beige	14
7563633	Actimove® Thumb Stabilizer – Extra Stays	L/XL	Beige	14
7568750	Actimove® Knee Brace – Adjustable Horseshoe, Simple Hinges, Condyle Pads	XS	Black	23
7568751	Actimove® Knee Brace – Adjustable Horseshoe, Simple Hinges, Condyle Pads	S	Black	23
7568752	Actimove® Knee Brace – Adjustable Horseshoe, Simple Hinges, Condyle Pads	M	Black	23
7568753	Actimove® Knee Brace – Adjustable Horseshoe, Simple Hinges, Condyle Pads	L	Black	23
7568754	Actimove® Knee Brace – Adjustable Horseshoe, Simple Hinges, Condyle Pads	XL	Black	23
7568755	Actimove® Knee Brace – Adjustable Horseshoe, Simple Hinges, Condyle Pads	XXL	Black	23
7568756	Actimove® Knee Brace – Adjustable Horseshoe, Simple Hinges, Condyle Pads	XXXL	Black	23
7568757	Actimove® Knee Brace – Adjustable Horseshoe, Simple Hinges, Condyle Pads	XXXXL	Black	23
7569130	Actimove® Knee Brace – Composite Polycentric Hinges	S	Black	88
7569131	Actimove® Knee Brace – Composite Polycentric Hinges	M	Black	88
7569132	Actimove® Knee Brace – Composite Polycentric Hinges	L	Black	88
7569133	Actimove® Knee Brace – Composite Polycentric Hinges	XL	Black	88
7569134	Actimove® Knee Brace – Composite Polycentric Hinges	XXL	Black	88
7571420	Actimove® Wrist Stabilizer Carpal – Pre-shaped Metal Stay	Universal	Black	37
7571530	Actimove® Wrist Splint 6"	Right S	Black	73
7571531	Actimove® Wrist Splint 6"	Right M	Black	73
7571532	Actimove® Wrist Splint 6"	Right L	Black	73
7571533	Actimove® Wrist Splint 6"	Right XL	Black	73
7571534	Actimove® Wrist Splint 6"	Left S	Black	73
7571535	Actimove® Wrist Splint 6"	Left M	Black	73
7571536	Actimove® Wrist Splint 6"	Left L	Black	73
7571537	Actimove® Wrist Splint 6"	Left XL	Black	73
7571670	Actimove® Wrist Brace – Suede Finish, Lace-Closure	Right XS	Beige	74
7571671	Actimove® Wrist Brace – Suede Finish, Lace-Closure	Right S	Beige	74

Code	Description	Size	Color	Page
7571672	Actimove® Wrist Brace – Suede Finish, Lace-Closure	Right M	Beige	74
7571673	Actimove® Wrist Brace – Suede Finish, Lace-Closure	Right L	Beige	74
7571674	Actimove® Wrist Brace – Suede Finish, Lace-Closure	Right XL	Beige	74
7571675	Actimove® Wrist Brace – Suede Finish, Lace-Closure	Left XS	Beige	74
7571676	Actimove® Wrist Brace – Suede Finish, Lace-Closure	Left S	Beige	74
7571677	Actimove® Wrist Brace – Suede Finish, Lace-Closure	Left M	Beige	74
7571678	Actimove® Wrist Brace – Suede Finish, Lace-Closure	Left L	Beige	74
7571679	Actimove® Wrist Brace – Suede Finish, Lace-Closure	Left XL	Beige	74
7571680	Actimove® Wrist Brace – Suede Finish, Lace-Closure	Right XS	Black	74
7571681	Actimove® Wrist Brace – Suede Finish, Lace-Closure	Right S	Black	74
7571682	Actimove® Wrist Brace – Suede Finish, Lace-Closure	Right M	Black	74
7571683	Actimove® Wrist Brace – Suede Finish, Lace-Closure	Right L	Black	74
7571684	Actimove® Wrist Brace – Suede Finish, Lace-Closure	Right XL	Black	74
7571685	Actimove® Wrist Brace – Suede Finish, Lace-Closure	Left XS	Black	74
7571686	Actimove® Wrist Brace – Suede Finish, Lace-Closure	Left S	Black	74
7571687	Actimove® Wrist Brace – Suede Finish, Lace-Closure	Left M	Black	74
7571688	Actimove® Wrist Brace – Suede Finish, Lace-Closure	Left L	Black	74
7571689	Actimove® Wrist Brace – Suede Finish, Lace-Closure	Left XL	Black	74
7571750	Actimove® Wrist Splint 8"	Right XS	Black	72
7571751	Actimove® Wrist Splint 8"	Right S	Black	72
7571752	Actimove® Wrist Splint 8"	Right M	Black	72
7571753	Actimove® Wrist Splint 8"	Right L	Black	72
7571754	Actimove® Wrist Splint 8"	Right XL	Black	72
7571755	Actimove® Wrist Splint 8"	Left XS	Black	72
7571756	Actimove® Wrist Splint 8"	Left S	Black	72
7571757	Actimove® Wrist Splint 8"	Left M	Black	72
7571758	Actimove® Wrist Splint 8"	Left L	Black	72
7571759	Actimove® Wrist Splint 8"	Left XL	Black	72
7571830	Actimove® Wrist Splint – Abducted Thumb, Compression-Formed Composite	Right XS	Black	79
7571831	Actimove® Wrist Splint – Abducted Thumb, Compression-Formed Composite	Right S	Black	79
7571832	Actimove® Wrist Splint – Abducted Thumb, Compression-Formed Composite	Right M	Black	79
7571833	Actimove® Wrist Splint – Abducted Thumb, Compression-Formed Composite	Right L	Black	79
7571834	Actimove® Wrist Splint – Abducted Thumb, Compression-Formed Composite	Right XL	Black	79
7571845	Actimove® Wrist Splint – Abducted Thumb, Compression-Formed Composite	Left XS	Black	79
7571846	Actimove® Wrist Splint – Abducted Thumb, Compression-Formed Composite	Left S	Black	79
7571847	Actimove® Wrist Splint – Abducted Thumb, Compression-Formed Composite	Left M	Black	79
7571848	Actimove® Wrist Splint – Abducted Thumb, Compression-Formed Composite	Left L	Black	79
7571849	Actimove® Wrist Splint – Abducted Thumb, Compression-Formed Composite	Left XL	Black	79
7572930	Actimove® Wrist Stabilizer – Removable Metal Stay, Left or Right	Pediatric	Black	58
7572931	Actimove® Wrist Stabilizer – Removable Metal Stay, Left or Right	S	Black	13
7572932	Actimove® Wrist Stabilizer – Removable Metal Stay, Left or Right	M	Black	13
7572933	Actimove® Wrist Stabilizer – Removable Metal Stay, Left or Right	L	Black	13
7574230	Actimove® Elbow Strap – Hot/Cold Pack	Universal	Black	32
7577920	Actimove® Wrist Support	S	Beige	52
7577921	Actimove® Wrist Support	M	Beige	52
7577922	Actimove® Wrist Support	L	Beige	52
7577923	Actimove® Wrist Support	XL	Beige	52
7578020	Actimove® Ankle Support	S	Beige	55
7578021	Actimove® Ankle Support	M	Beige	55
7578022	Actimove® Ankle Support	L	Beige	55
7578023	Actimove® Ankle Support	XL	Beige	55
7578024	Actimove® Ankle Support	XXL	Beige	55

Code	Description	Size	Color	Page
7578120	Actimove® Knee Support	S	Beige	54
7578121	Actimove® Knee Support	M	Beige	54
7578122	Actimove® Knee Support	L	Beige	54
7578123	Actimove® Knee Support	XL	Beige	54
7578124	Actimove® Knee Support	XXL	Beige	54
7578221	Actimove® Elbow Support	S	Beige	56
7578222	Actimove® Elbow Support	M	Beige	56
7578223	Actimove® Elbow Support	L	Beige	56
7578224	Actimove® Elbow Support	XL	Beige	56
7578320	Actimove® Arthritis Gloves	S	Beige	53
7578321	Actimove® Arthritis Gloves	M	Beige	53
7578322	Actimove® Arthritis Gloves	L	Beige	53
7578323	Actimove® Arthritis Gloves	XL	Beige	53
7623107	Actimove® AchilloMotion – Achilles Tendon Support	S	Charcoal	93
7623108	Actimove® AchilloMotion – Achilles Tendon Support	M	Charcoal	93
7623109	Actimove® AchilloMotion – Achilles Tendon Support	L	Charcoal	93
7623110	Actimove® AchilloMotion – Achilles Tendon Support	XL	Charcoal	93
7623111	Actimove® AchilloMotion – Achilles Tendon Support	XXL	Charcoal	93
7623113	Actimove® AchilloMotion – Achilles Tendon Support	S	White	93
7623114	Actimove® AchilloMotion – Achilles Tendon Support	M	White	93
7623115	Actimove® AchilloMotion – Achilles Tendon Support	L	White	93
7623116	Actimove® AchilloMotion – Achilles Tendon Support	XL	White	93
7623117	Actimove® AchilloMotion – Achilles Tendon Support	XXL	White	93
7623800	Actimove® Rhizo Forte – Thumb Orthosis	Right S	Gray	67
7623801	Actimove® Rhizo Forte – Thumb Orthosis	Left S	Gray	67
7623802	Actimove® Rhizo Forte – Thumb Orthosis	Right M	Gray	67
7623803	Actimove® Rhizo Forte – Thumb Orthosis	Left M	Gray	67
7623804	Actimove® Rhizo Forte – Thumb Orthosis	Right L	Gray	67
7623805	Actimove® Rhizo Forte – Thumb Orthosis	Left L	Gray	67
7623818	Actimove® Rhizo Forte – Thumb Orthosis	Right S	Black	67
7623819	Actimove® Rhizo Forte – Thumb Orthosis	Left S	Black	67
7623820	Actimove® Rhizo Forte – Thumb Orthosis	Right M	Black	67
7623821	Actimove® Rhizo Forte – Thumb Orthosis	Left M	Black	67
7623822	Actimove® Rhizo Forte – Thumb Orthosis	Right L	Black	67
7623823	Actimove® Rhizo Forte – Thumb Orthosis	Left L	Black	67
7627030	Actimove® Standard Walker – Ankle Walker	XS Low	Black	100
7627031	Actimove® Standard Walker – Ankle Walker	S Low	Black	100
7627032	Actimove® Standard Walker – Ankle Walker	M Low	Black	100
7627033	Actimove® Standard Walker – Ankle Walker	L Low	Black	100
7627034	Actimove® Standard Walker – Ankle Walker	XL Low	Black	100
7627035	Actimove® Standard Walker – Ankle Walker	XS High	Black	100
7627036	Actimove® Standard Walker – Ankle Walker	S High	Black	100
7627037	Actimove® Standard Walker – Ankle Walker	M High	Black	100
7627038	Actimove® Standard Walker – Ankle Walker	L High	Black	100
7627039	Actimove® Standard Walker – Ankle Walker	XL High	Black	100
7627230	Actimove® Standard Air Walker – Ankle Walker	XS Low	Black	101
7627231	Actimove® Standard Air Walker – Ankle Walker	S Low	Black	101
7627232	Actimove® Standard Air Walker – Ankle Walker	M Low	Black	101
7627233	Actimove® Standard Air Walker – Ankle Walker	L Low	Black	101
7627234	Actimove® Standard Air Walker – Ankle Walker	XL Low	Black	101
7627235	Actimove® Standard Air Walker – Ankle Walker	XS High	Black	101
7627236	Actimove® Standard Air Walker – Ankle Walker	S High	Black	101



Code	Description	Size	Color	Page
7627237	Actimove® Standard Air Walker – Ankle Walker	M High	Black	101
7627238	Actimove® Standard Air Walker – Ankle Walker	L High	Black	101
7627239	Actimove® Standard Air Walker – Ankle Walker	XL High	Black	101
7627610	Actimove® Closed Shell Air Walker – Ankle Walker	XS Low	Black	103
7627611	Actimove® Closed Shell Air Walker – Ankle Walker	S Low	Black	103
7627612	Actimove® Closed Shell Air Walker – Ankle Walker	M Low	Black	103
7627613	Actimove® Closed Shell Air Walker – Ankle Walker	L Low	Black	103
7627614	Actimove® Closed Shell Air Walker – Ankle Walker	XL Low	Black	103
7627615	Actimove® Closed Shell Air Walker – Ankle Walker	XS High	Black	103
7627616	Actimove® Closed Shell Air Walker – Ankle Walker	S High	Black	103
7627617	Actimove® Closed Shell Air Walker – Ankle Walker	M High	Black	103
7627618	Actimove® Closed Shell Air Walker – Ankle Walker	L High	Black	103
7627619	Actimove® Closed Shell Air Walker – Ankle Walker	XL High	Black	103
7627794	Actimove® Closed Shell Air Walker – Heel Wedge	One Size	Black	103



Actimove<sup>®</sup>,  
an Essity Company



[essity.com](http://essity.com) | [actimove.com](http://actimove.com)

BSN Medical Inc., an Essity Company  
5825 Carnegie Blvd., Charlotte NC 28209-4633  
Tel. (+1) 704 554 9933 Fax (+1) 800 835 4325  
To order toll-free Tel (+1) 800 552 1157