

# Hospital product catalogue





# Forward

It's just not a direction.  
**It's an attitude.**

It's about not resting on our laurels,  
**but leading the way.**



**Forward-thinking is what makes us a leader in health, hygiene, and medical solutions.**

Our extensive portfolio of leading therapeutic brands reflects the spirit of innovation we bring to everything we do.


We have the courage to lead, to break barriers, and to always push for better.

**Better Products.**

**Better Service.**

**Better Value.**

**Better Outcomes.**



Our wide range of therapeutic solutions, clinical expertise and world-class brands **mean one partner to support your needs.**

Wound Care: **Leukoplast®** **Cutimed®**

Orthopedics: **Actimove®**

Compression: **JOBST®**

## Your advantage

**Wide range of therapeutic solutions means one partner to support your needs.**

### The value of higher quality products

Do more with less  
Lower overall costs  
Improve patient outcomes

### Clinical Expertise

Strong clinical in-servicing through nurses, field staff, advisory and web-guides.  
Keep current with clinical training and webinars

### Innovative Solutions

Proven product solutions you can count on to help you better manage medical conditions.

## Our goal

**Being your ally in better outcomes for you and those in your care.  
Now and for the future.**



# Contents

## Acute Wound Care 6

Introduction	7
Portfolio overview	7
Leukomed® Control	8
Leukomed® Sorbact®	8
Leukomed® T Plus	9
Leukomed®	9
Leukomed® Skin Sensitive	10
Leukomed® T Plus Skin Sensitive	10
Leukomed® T Skin Sensitive	11
Leukomed® T	11
Leukomed® LA	12
Leukosan® Adhesive	12
Leukoplast® Wound Closure Strips	13

## Advanced Wound Care 14

Introduction	15
Portfolio overview	15
Cutimed® Debriclean	16
Cutimed® Gel	16
Cutimed® Sorbact® Compress	17
Cutimed® Sorbact® Ribbon Gauze	17
Cutimed® Sorbact® Gel Dressing	18
Cutimed® Sorbact® Absorption Dressing	18
Cutimed® Sorbact® Foam Dressing	19
Cutimed® Siltec Sorbact®	19
Cutimed® Sorbion Sorbact®	20
Cutimed® Sorbion Sachet S	20
Cutimed® Sorbion Sachet Multi-Star	21
Cutimed® Sorbion Sachet XL	21
Cutimed® Sorbion Sachet S Drainage	22
Cutimed® Sorbion Sana Gentle	22
Cutimed® Siltec	23
Cutimed® Siltec L	23
Cutimed® Siltec B	24
Cutimed® Siltec Heel 2D	24
Cutimed® Siltec Heel 3D	25
Cutimed® Alginate	25
Cutimed® Hydrocontrol	26
Cuticell® Contact	26
Cuticell®	27
Cuticell® Classic	27

## Surgical Tapes & Bandages 28

Introduction	
Portfolio overview	29
Fixomull® Stretch	30
Hypafix®	30
Fixomull® skin sensitive	31
Hypafix® skin sensitive	31
Fixomull® Transparent	32
Hypafix® Transparent	32
Leukoplast® Skin sensitive	33
Leukoplast® Standard	33
Leukoplast® elastic	34
Leukoplast® Waterproof	34
Leukofix®	35
Leukoflex®	35
Leukosilk®	36
Leukopor®	36
Leukoplast® Sleek	37
Leukoplast® Sleek LF	37
Easifix®	38
Easifix® K	38
Handyband® Conforming	39
Easifix® Cohesive	39
Handy®Gauze Cohesive	40
Surgifix®	40
Tricofix®	41
Elastolite®	41
Handy®Crepe Hospital	42
Handy®Crepe Medium Weight	42
Tensocrepe® Medium Weight	43
Tensogrip®	43
Co-Plus® LF	44
Co-Plus®	44
Tensocrepe® Heavy Weight BP	45
Tensocrepe® Hospital Heavy Weight	45
Handy®Crepe Heavy	46
Tensoplast®	46
Acrylastic®	47
JOBST® Comprilan®	47
Tensopress®	48
Eloflex®	48
JOBST® Compri2	49

## First Aid Dressings 50

Introduction	51
Portfolio overview	51
Leukoplast® Visual	52
Leukoplast® Kids	52
Leukoplast® Universal	53
Leukoplast® Elastic	53
Leukoplast® Strong	54
Leukoplast® Soft White	55
Leukoplast® Detectable	55
Leukoplast® Barrier	56
Leukoplast® Aqua Pro	56

*Always follow the directions for use. These products are not available for purchase by the general public*
















## Acute Wound Care



Medical experience since 1901 and professional, high-quality products make Leukoplast a strong and trusted partner in wound care - whether used by healthcare experts or by consumers at home. Our comprehensive portfolio is designed to serve the needs of your patients in everyday wound management to a professional standard. From first aid dressings for everyday use to specialized products for more complex wounds.

**Leukoplast. Wound care in best hands.**

## Portfolio overview

Dressings	Post Op Dressings				
	Leukomed® Control	Leukomed® Sorbact	Leukomed® T Plus	Leukomed®	
	Skin Sensitive Dressings				
	Leukomed® skin sensitive	Leukomed® T skin sensitive	Leukomed® T Plus skin sensitive		
Film Dressings					
Leukomed® T					
Low Adherent Dressings					
Leukomed® LA					
Wound Closure	Skin Glue				
	Leukosan® Adhesive				
Closure Strips					
Leukosan® Strips					

## Leukomed® Control

### Fully Transparent Post-Op Dressing

#### Product benefits

- Ultra thin and conformable film for a second skin feeling
- Absorbent hydropolymer gel wound pad
- Completely transparent composition which allows for early detection of wound infection
- Waterproof with up to 14 day wear time



Leukomed® Control			
Code	Size	Pad Size	Pack Qty
73230-00	5 cm x 7 cm	2cm x 4cm	10/box
73230-01	7 cm x 10 cm	4cm x 7cm	10/box
73230-02	8 cm x 15 cm	4.5cm x 11.5cm	10/box
73230-03	10 cm x 24 cm	6cm x 20cm	5/box
73230-05	10 cm x 35 cm	6cm x 31cm	5/box

#### Indications

Dry and low-exuding acute wounds including surgical wounds, superficial and partial thickness burns, skin donor sites, lacerations and abrasions. Ideal for postoperative care.

## Leukomed® Sorbact®

### Antimicrobial Post-Op Dressing

#### Product benefits

- Innovative post-operative dressing to help to reduce wound bacteria colonisation by its physical mode of action
- Unique bacteria binding Sorbact®
- Highly breathable and conformable
- Sterile and waterproof barrier
- Absorbent, cushioning, low adherent wound pad



Leukomed® Sorbact®			
Code	Size	Pad Size	Pack Qty
76199-26	5cm x 7.2cm	3cm x 4cm	100/box
76199-01	8cm x 10cm	4cm x 6cm	20/box
76199-02	8cm x 15cm	4cm x 11cm	20/box
76199-03	10cm x 20cm	5cm x 16cm	20/box
76199-04	10cm x 25cm	5cm x 20.5cm	20/box
76199-05	10cm x 30cm	5cm x 25cm	20/box
76199-06	10cm x 35cm	5cm x 30 cm	20/box

#### Indications

Post-operative and traumatic wounds with dry to low exudate levels

- Surgical incisions
- Lacerations, cuts, abrasions

## Leukomed® T Plus

### Film Island Dressing

#### Product benefits

- Ultra thin, conformable film
- Bacterial and waterproof barrier
- Unique application technique with the distinctive red strip
- Protects skin integrity
- Absorbent wound pad



Leukomed® T Plus			
Code	Size	Pad Size	Pack Qty
72382-00	5 cm x 7.2 cm	2.8cm x 4cm	50/box
72382-01	8 cm x 10 cm	4cm x 6cm	50/box
72382-02	8 cm x 15 cm	4cm x 11cm	50/box
72382-03	10 cm x 25 cm	5cm x 20cm	50/box
72382-05	10 cm x 35 cm	5cm x 30cm	50/box

#### Indications

To cover minor to medium size post-operative and traumatic wounds

## Leukomed®

### Non-Woven Post-Op Dressing

#### Product benefits

- Reliable and gently adherent
- Highly conformable and breathable
- Sterile and protective
- Secure absorption



Leukomed®			
Code	Size	Pad Size	Pack Qty
72380-00	5 cm x 7.2 cm	2.8cm x 4cm	50/box
72380-01	8 cm x 10 cm	4cm x 6cm	50/box
72380-02	8 cm x 15 cm	4cm x 11cm	50/box
72380-04	10 cm x 25 cm	5cm x 20cm	50/box
72380-06	10cm x 35cm	5cm x 30cm	50/box

#### Indications

Surgical or traumatic wounds

## Leukomed® Skin Sensitive

Absorbent dressing for wound protection on fragile or sensitive skin

### Product benefits

- Skin friendly silicone adhesive
- Reduced pain at dressing removal
- Covers and protects wounds
- Absorbs blood and exudate



Leukomed® Skin Sensitive			
Code	Size	Pad Size	Pack Qty
71239-00	5 cm x 7.2 cm	2.8cm x 4cm	20/box
71239-01	8 cm x 10 cm	4cm x 6cm	20/box
71239-02	8 cm x 15 cm	4cm x 11cm	20/box
71239-03	10 cm x 25 cm	5cm x 20cm	20/box
71239-04	5 cm x 7.2 cm	2.8cm x 4cm	5/box
71239-05	8 cm x 10 cm	4cm x 6cm	5/box

### Indications

Covers and protects wounds, such as closed surgical wounds, cuts and abrasions

## Leukomed® T Plus Skin Sensitive

Film island dressing for fragile or sensitive skin

### Product benefits

- Secure wound protection on fragile skin
- Skin friendly silicone adhesive
- Bacteria barrier and waterproof protection
- Atraumatic removal



Leukomed® T Plus Skin Sensitive			
Code	Size	Pad Size	Pack Qty
76178-00	5 cm x 7.2 cm	2.8cm x 4cm	5/box
76178-01	8 cm x 10 cm	4cm x 6cm	5/box
76178-02	8 cm x 15 cm	4cm x 11cm	5/box
76178-03	10 cm x 25 cm	5cm x 20cm	5/box

### Indications

- Minor traumatic wounds, surgical wounds and low to moderate exuding skin lesions
- Fragile and compromised skin

## Leukomed® T Skin Sensitive

Transparent film dressing for fragile or sensitive skin

### Product benefits

- Secure wound protection and fixation on fragile skin
- Skin friendly silicone adhesive
- Bacteria barrier and waterproof protection



Leukomed® T Plus		
Code	Size	Pack Qty
76464-00	5 cm x 7.2 cm	5/box
76464-01	8 cm x 10 cm	5/box
76464-02	8 cm x 15 cm	5/box

### Indications

- Sterile fixation of primary dressings and compresses
- Additional fixation of devices
- Protection of non-exuding wounds from friction and shear forces

## Leukomed® T

Film Dressing

### Product benefits

- Ultra thin, conformable film
- Gently and almost invisibly protects freshly epithelised skin against shear forces, contaminants, bacteria and viruses
- Unique application technique with the distinctive red strip



Leukomed® T		
Code	Size	Pack Qty
72381-00	5 cm x 7.2 cm	50/box
72381-01	8 cm x 10 cm	50/box
72381-08	10 cm x 12.5 cm	50/box

### Indications

- For cover and protection of post-operative and superficial traumatic wounds
- For sterile fixation of primary dressings and compresses

## Leukomed® LA

*Low adherent absorbent dressing*

### Product benefits

- Reliable, non-stick absorbent cotton dressing
- Quick and easy to apply
- High wearing comfort
- Reduces the risk of skin maceration, blistering and ensures atraumatic dressing changes

Leukomed® LA		
Code	Size	Pack Qty
73289-00	5cm x 5cm	50/box
73289-01	10cm x 10cm	50/box
73289-02	10cm x 20cm	50/box
73289-03	53cm x 60cm	40/box
73289-04	53cm x 120cm	20/box
73289-05	53cm x 180cm	20/box



### Indications

- Low to moderate exuding wounds including post-operative wounds
- Clean and sutured wounds
- Abrasions, lacerations and minor burns

## Leukosan® Adhesive

*Topical skin adhesive for precision bonding*

### Product benefits

- **Express polymerisation** Quick wound closure, lasting adaption
- **Ultra-high viscosity** Quick, easy and safe to apply
- **Action flexibility** Flexes with skin, reliable adaption and protection
- **No heat generation** Pain-free application for high patient comfort, child-friendly
- **2 in 1 applicator tip** Optional precise and/or protective application

Leukosan® Adhesive		
Code	Size	Pack Qty
72541-00	0.7ml tube	10/box



### Indications

To be used topically by trained healthcare professionals primarily to close surgical incisions and minor traumatic wounds with easily opposable skin edges in areas of low skin tension.

## Leukoplast® Wound Closure Strips

*Wound closure strips*

### Product benefits

- **Intelligent threads** Optimal elasticity for adaptability to successive healing phases
- **Reliable, skin-friendly adhesion** Keeps the strips securely in place
- **No fraying or rolling up of edges** High wearing comfort
- **Breathable and permeable material** Prevents maceration and exudate strikethrough
- **Optimal durability** Minimises aftercare, time and cost savings

### Indications

- Ideal for closing facial lacerations or incisions without additional trauma and best cosmetic results
- Additional tension relief for wounds closed with intracutaneous sutures
- Early replacement of staples or sutures to improve cosmetic results.



Leukoplast® Wound Closure Strips					
Code	Size	Colour	Strips/Pouch	Pouches/Box	Total Strips
72740-00	6mm x 38mm	White	6/pouch	50/box	300
72740-01	3mm x 75mm	White	5/pouch	50/box	250
72740-02	6mm x 75mm	White	3/pouch	50/box	150
72740-03	6mm x 100mm	White	10/pouch	50/box	500
72740-04	12mm x 100mm	White	6/pouch	50/box	300
72740-05	25mm x 100mm	White	4/pouch	25/box	100
72740-06	6mm x 38mm	Tan	6/pouch	50/box	300
72740-07	6mm x 75mm	Tan	3/pouch	50/box	150
72740-08	12mm x 100mm	Tan	6/pouch	50/box	300



# Cutimed®

Innovative wound management for complex or chronic wounds



## Advanced Wound Care

Wound healing complications are common in diabetes or venous disorder, but efficient and effective products help patients get back to their lives sooner. Especially developed for the treatment of chronic and complex wounds, the Cutimed Advanced Wound Care range covers all stages of the healing process.

With innovative wound management solutions and the perfect balance of advanced technology and user-friendly concepts, Cutimed empowers professionals to help their patients heal quickly and effectively for the best possible patient outcome.

### Wound depth

- Superficial
- Deep
- Superficial + deep

### Wound phase

- Necrotic
- Infected
- Sloughy
- Granulating
- Epithelialized

### Exudate level

- Dry to low
- Low to moderate
- Moderate to high
- High to excessive

## Portfolio overview

Clean & Debride	 Cutimed® DebriClean	 Cutimed® Gel				
Manage Infection	 Cutimed® Sorbact® Compress	 Cutimed® Sorbact® Ribbon Gauze	 Cutimed® Sorbact® Gel Dressing	 Cutimed® Sorbact® Absorption Dressing		
Manage Infection & Exudate	 Cutimed® Siltec Foam Dressing	 Cutimed® Siltec Sorbact®	 Cutimed® Sorbion Sorbact®			
Manage Exudate	 Cutimed® Sorbion Sachet S	 Cutimed® Sorbion Sachet Multi-Star	 Cutimed® Sorbion Sachet XL	 Cutimed® Sorbion Sachet S Drainage	 Cutimed® Sorbion Sana Gentle	
Manage Exudate	 Cutimed® Siltec	 Cutimed® Siltec L	 Cutimed® Siltec B	 Cutimed® Siltec Heel 2D	 Cutimed® Siltec Heel 3D	 Cutimed® Alginate
Manage Tissue Boost	 Cutimed® Hydrocontrol					
Complementary	 Cuticell® Contact	 Cuticell®	 Cuticell® Classic			



## Cutimed® DebriClean

### Monofilament fibre mechanical debridement pad

#### Product benefits

- Fast and effective wound cleansing and debridement
- A dual action debridement pad with soft and abrasive looped fibres
- Effectively cleans wounds of slough, biofilm and microorganisms
- Easy and comfortable to use



#### Indications

Aids in the debridement of superficial wounds and surrounding skin in, for example:

- diabetic ulcers
- venous and arterial ulcers
- pressure sores
- postoperative wounds and wounds healing by secondary intention
- superficial burns
- the absorption of exudate, slough and dry skin flakes



Code	Size	Pack Qty
76170-00	9cm x 13cm	5/box

## Cutimed® Gel

### Amorphous hydrogel

#### Product Benefits

- A clear, amorphous hydrogel which can be used to debride necrotic and sloughy tissue
- Easy to apply with the sterile nozzle provided, a spatula or directly from the tube
- Compatible with any type of secondary dressing, dependent upon exudate levels
- Preservative-free for less risk of allergic reactions



Code	Size	Pack Qty
72610-00	8g	10/box
72610-01	15g	10/box
72610-02	25g	10/box

#### Indications

- Management of dry to low exuding wounds and effective autolytic debridement of slough and eschar
- Provides moisture to dry wounds, even in later stages of wound healing

## Cutimed® Sorbact® Compress

### Bacteria and fungi binding wound contact layer

#### Product Benefits

- Non-absorbent bacteria and fungi binding wound contact layer, based on Sorbact Technology, that allows passage of wound exudate into a secondary dressing
- Prevents and treats wound infections and facilitates the wound healing process
- Can be used in compression therapy and in conjunction with NPWT



Code	Size	Pack Qty
72165-00	7cm x 9cm	40/box
72164-00	4cm x 6cm	40/box
72693-08	NPWT 10cm x 100cm	10/box

#### Indications

For use in management of clean, contaminated, colonised or infected exuding wounds, such as surgical wounds, traumatic wounds, pressure ulcers, diabetic foot ulcers and leg ulcers. Can be used on both superficial and deep wounds.

## Cutimed® Sorbact® Ribbon Gauze

### Bacteria and fungi binding ribbon gauze

#### Product Benefits

- Infection treatment and prevention for cavity wounds
- Effectively binds and removes both bacteria and fungi
- Allows the passage of wound exudate into a secondary dressing
- Soft cotton gauze prevents shearing and skin to skin contact in skin folds



Code	Size	Pack Qty
72166-09	1cm x 50cm	20/box
72166-10	10cm x 200cm	10/box
72167-00	5cm x 200cm	10/box
72166-00	2cm x 50cm	20/box

#### Indications

For use in management of clean, contaminated, colonised or infected exuding wounds, such as cavity wounds and fistulas. Also intended for use in fungal infections in skin folds.

## Cutimed® Sorbact® Gel Dressing

### Bacteria and fungi binding dressing covered with gel

#### Product Benefits

- Combines Sorbact wound contact layer with a high moisture donation water-based gel
- Donates moisture, enables a moist wound environment and facilitates autolytic debridement
- Two in one: easy to use
- Manages and prevents wound infection



Code	Size	Pack Qty
72611-00	7.5cm x 7.5cm	10/box
72611-01	7.5cm x 15cm	10/box
72764-00	3cm x 15cm	10/box
72611-16	15cm x 20cm	5/box
72764-01	20cm x 25cm	5/box

#### Indications

For use in management of clean, contaminated, colonised or infected wounds with dry to low exudate levels, such as traumatic wounds, burns, cavity wounds, fistulas, pressure ulcers, diabetic ulcers and foot and leg ulcers.

## Cutimed® Sorbact® Absorption Dressing

### Bacteria and fungi binding wound dressing

#### Product Benefits

- An absorbent dressing pad covered with a Sorbact technology coated acetate swab
- Infection treatment and prevention for low to moderately exuding wounds



Code	Size	Pack Qty
72162-16	10cm x 10cm	20/box
72163-09	10cm x 20cm	10/box
72161-14	7cm x 9cm	20/box

#### Indications

For use in management of clean, contaminated, colonised or infected wounds with moderate exudate levels, such as surgical wounds, traumatic wounds, pressure ulcers, diabetic ulcers and foot and leg ulcers.

#### Manage Exudate & Infection

## Cutimed® Sorbact® Foam Dressing

### Bacteria and fungi binding foam dressing

#### Product Benefits

- Foam dressing with Sorbact technology for managing and preventing infection in low to moderately exuding wounds
- Absorbs and retains exudate, thereby reducing the risk of maceration and enabling a moist wound environment
- Semi-permeable backing film allows excess fluid to evaporate



Code	Size	Pack Qty
72763-00	10cm x 10cm	10/box
72763-01	15cm x 15cm	10/box
72763-02	10cm x 20cm	10/box

#### Indications

For use in the management of clean, contaminated, colonised or infected wounds with low to moderate exudate levels, such as surgical wounds, traumatic wounds, pressure ulcers, diabetic ulcers and foot and leg ulcers.

## Cutimed® Siltec Sorbact®

### Bacteria and fungi binding silicone foam dressing

#### Product Benefits

- Foam dressing with silicone adhesive border, integrated with Sorbact technology for managing and preventing infection
- Super-absorbent foam for vertical absorption of exudate reducing the risk of maceration
- The moisture vapor transmission rate (MVTR) allows for dressing permeability while also maintaining a moisture balance



Code	Size	Pad Size	Pack Qty
73251-30	7.5cm x 7.5cm	4cm x 4cm	10/box
73251-31	10cm x 10cm	6cm x 6cm	10/box
73251-32	15cm x 15cm	10cm x 10cm	10/box

#### Indications

For use in management of clean, contaminated, colonised or infected wounds with moderate exudate levels, such as surgical wounds, traumatic wounds, pressure ulcers, diabetic ulcers and foot and leg ulcers.

## Cutimed® Sorbion Sorbact®

### Bacteria and fungi binding super absorbent dressing

#### Product Benefits

- Binds and removes bacteria with a physical mode of action
- Absorbs large amount of exudate
- Retains exudate under compression, reducing the risk of maceration
- Helps maintain a moist wound healing environment



Code	Size	Pad Size	Pack Qty
72698-00	10cm x 10cm	8cm x 9cm	10/box
72698-01	10cm x 10cm	7.7cm x 17.8cm	10/box
72698-02	20cm x 20cm	17.8cm x 17.8cm	10/box
72698-03	20cm x 20cm	17.5cm x 27.5cm	10/box

#### Indications

For use in the management of clean, contaminated, colonised or infected wounds with high to excessive exudate levels, such as surgical wounds, traumatic wounds, pressure ulcers, diabetic ulcers and foot and leg ulcers.

#### Manage Exudate

## Cutimed® Sorbion Sachet S

### Non-adhesive superabsorbent dressing

#### Product Benefits

- Versatile superabsorbent dressing with high wearing comfort for primary wound contact
- Absorbs wound exudate and debris and locks it in
- Binds debris and inflammatory factors like matrix metalloproteinases MMPs safely inside the absorbent core
- Maintains a moist wound environment
- Can be used under compression



Code	Size	Pack Qty
73232-00	7.5cm x 7.5cm	10/box
73232-06	10cm x 10cm	10/box
73232-09	20cm x 10cm	10/box
73232-15	20cm x 20cm	10/box
73232-18	30cm x 20cm	10/box

#### Indications

Suitable for moderately to highly exuding wounds, including:

- acute and chronic wounds
- leg, foot and pressure ulcers
- wounds healing by secondary intention

## Cutimed® Sorbion Sachet Multi-Star

### Shaped non-adhesive superabsorbent dressing

#### Product Benefits

- Versatile superabsorbent dressing with high wearing comfort for primary wound contact
- Absorbs wound exudate and debris and locks it in
- Allows flexible applications for use on sites such as the toes, heel, foot, sacrum, chest, axilla, elbows or wound cavities



Code	Size	Pack Qty
73237-01	ø 8cm	10/box
73237-03	ø 14cm	10/box

#### Indications

Suitable for moderately to highly exuding wounds, including:

- acute and chronic wounds
- leg, foot and pressure ulcers
- wounds healing by secondary intention

## Cutimed® Sorbion Sachet XL

### Shaped non-adhesive superabsorbent dressing

#### Product Benefits

- Versatile superabsorbent dressing for primary wound contact
- Absorbs wound exudate and debris and locks it in
- Specially designed for larger wounds covering the entire lower leg
- High wearing comfort on areas where application is difficult, such as the lower leg, sacral region, thorax, thigh and lower part of the abdomen



Code	Size	Pack Qty
73240-00	45cm x 25cm	10/box

#### Indications

Suitable for moderately to highly exuding wounds, including:

- acute and chronic wounds,
- leg, foot and pressure ulcers,
- wounds healing by secondary intention.



## Cutimed® Sorbion Sachet S Drainage

### Non-adhesive superabsorbent dressing for catheters and drainage

#### Product Benefits

- Versatile superabsorbent dressing for primary wound contact
- Absorbs wound exudate and debris and locks it in
- The unique split design ensures an especially good fit at the exit points of tubes e.g. in PEG (percutaneous endoscopic gastrostomy), probes, tracheostomy



Code	Size	Pack Qty
73241-00	10cm x 10cm	10/box

#### Indications

Suitable for moderately to highly exuding wounds, including:

- acute and chronic wounds
- leg, foot and pressure ulcers
- wounds healing by secondary intention

## Cutimed® Sorbion Sana Gentle

### Gentle non-adhesive superabsorbent dressing

#### Product Benefits

- Versatile superabsorbent dressing for primary wound contact
- Absorbs wound exudate and debris and locks it in
- Integrated wound contact layer made of polyethylene to protect granulation tissue
- Atraumatic dressing change for sensitive wounds and skin
- Also available in Multi-Star shape



Code	Size	Pack Qty
73233-19	8.5cm x 8.5cm	10/box
73233-22	12cm x 12cm	10/box
73233-25	22cm x 12cm	10/box
73233-14	22cm x 22cm	10/box
73233-16	32cm x 22cm	10/box
73238-01	Multi Star Ø 11cm	10/box
73238-03	Multi Star Ø 17cm	10/box

#### Indications

For use on chronic or traumatic wounds with moderate to high exudation such as pressure ulcers, leg ulcers, diabetic foot syndrome and after burns: donor sites, fragile skin, and epidermolysis bullosa (EB).

## Cutimed® Siltec

### Non-bordered silicone foam dressing

#### Product Benefits

- Superabsorbent hotmelt stripes for additional absorption capacity
- Highly breathable
- Maintains critical moisture balance
- Skin protection with low adherence while reducing risk for maceration
- Excellent exudate retention capacity, even under compression.



Code	Size	Pack Qty
73285-00	5cm x 6cm	10/box
73285-01	10cm x 10cm	10/box
73285-02	10cm x 20cm	10/box
73285-03	15cm x 15cm	10/box

#### Indications

For the treatment of exuding wounds with low to high exudate levels such as venous and arterial ulcers, pressure ulcers, diabetic foot ulcers, surgical incisions, skin grafts and donor sites, lacerations or abrasions.

## Cutimed® Siltec L

### Non-bordered light silicone foam dressing

#### Product Benefits

- Thinner version of Cutimed Siltec
- Soft and conformable wound dressing
- Protective dressing for compromised and fragile skin
- Excellent exudate retention capacity, even under compression.



Code	Size	Pack Qty
73283-00	5cm x 6cm	10/box
73283-01	10cm x 10cm	10/box
73283-02	15cm x 15cm	10/box

#### Indications

For exuding wounds such as: venous and arterial ulcers, pressure ulcers, diabetic gangrenes, surgical incisions, skin grafts and donor sites, lacerations or abrasions. Recommended for wounds with low to moderate exudate levels.

## Cutimed® Siltec B

### Silicone foam dressing with border

#### Product Benefits

- Fast application without secondary fixation
- Cushioning effect that may help to prevent pressure ulcers
- Superabsorbent hotmelt stripes for additional absorption capacity
- Highly breathable
- Maintains critical moisture balance



Code	Size	Pack Qty
73284-00	7.5cm x 7.5cm	10/box
73284-05	10cm x 10cm	10/box
73284-01	12.5cm x 12.5cm	10/box
73284-06	10cm x 22.5cm	10/box
73284-14	7cm x 10cm Oval	10/box
73287-00	17.5 cm x 17.5 cm Sacrum	5/box
73287-01	23 cm x 23 cm Sacrum	5/box

#### Indications

For the treatment of exuding wounds with low to high exudate levels such as venous and arterial ulcers, pressure ulcers, diabetic foot ulcers, surgical incisions, skin grafts and donor sites, lacerations or abrasions.

## Cutimed® Siltec Heel 2D

### Silicone foam dressing

#### Product Benefits

- The dressing is designed to fit the anatomical shape of the heel and is folded around the heel at application
- Superabsorbent hotmelt stripes for additional absorption capacity
- Maintains critical moisture balance
- Reduces pressure peaks, provides cushioning



Code	Size	Pack Qty
73288-05	16cm x 24cm	5/box

#### Indications

For exuding wounds such as: venous and arterial ulcers, pressure ulcers, diabetic gangrenes, surgical incisions, skin grafts and donor sites, lacerations or abrasions. Recommended for wounds with low to moderate exudate levels.

## Cutimed® Siltec Heel 3D

### Silicone foam dressing

#### Product Benefits

- Especially designed to provide an optimal fit to every heel shape during ulcer management
- Anatomical shape provides cushioning and pressure relief from the heel to the Achilles tendon
- The extra thick foam of the dressing allows for reliable absorption of exudate
- Ready-to-use, easy application



Code	Size	Pack Qty
73286-01	Heel 3D	5/box

#### Indications

For exuding wounds such as: venous and arterial ulcers, pressure ulcers, diabetic gangrenes, surgical incisions, skin grafts and donor sites, lacerations or abrasions. Recommended for wounds with low to moderate exudate levels.

## Cutimed® Alginate

### Calcium alginate dressing

#### Product Benefits

- Haemostatic effect to control minor bleeding in superficial wounds
- Helps maintain a moist environment
- Fast gelling
- Highly absorbent
- Retains the exudate in a gel
- Excellent wet integrity



Code	Size	Pack Qty
72634-08	5cm x 5cm	10/box
72634-09	10cm x 10cm	10/box
72634-10	10cm x 20cm	10/box
72634-11	2.5cm x 30cm rope	5/box

#### Indications

For moderate to heavily exuding, superficial and deep wounds in the exudation and granulation phase:

- Ulcers of various origins (venous, arterial)
- Leg ulcers
- Pressure ulcers
- Diabetic ulcers
- Cavity wounds
- Acute wounds such as skin donor sites, post-operative surgical wounds, superficial and partial thickness burns, lacerations and abrasions

# Cutimed® Hydrocontrol

## Moisture-balancing, absorbent hydropolymer dressing

### Product Benefits

- Removal of necrotic tissue and moist wound healing
- Promotes wound stimulation potential due to MMP binding and growth factor protection
- Excellent exudate management
- Effective moisture balance
- Support of pain relief
- Cooling effect



Code	Size	Pack Qty
73231-00	4.5cm x 4.5cm	10/box
73231-01	7.5cm x 7.5cm	10/box
73231-02	10cm x 10cm	10/box
73231-03	15cm x 15cm	5/box

### Indications

For dry to moderately exuding chronic and secondary healing wounds, such as leg ulcers (venous, arterial and mixed aetiology), pressure ulcers and diabetic foot ulcers. Suitable for use under compression therapy.

### Complementary

# Cuticell® Contact

## Soft silicone wound contact layer

### Product benefits

- Single-sided adherence without sticking to wound bed
- Especially flexible and conformable
- Tissue-friendly, atraumatic dressing changes
- Minimised risk of maceration
- Easy wound inspection due to transparency
- Allows quick transition of highly viscous exudate



Code	Size	Pack Qty
72680-00	5cm x 7.5cm	5/box
72680-01	7.5cm x 10cm	5/box
72680-02	10cm x 18cm	5/box
72680-03	15cm x 25cm	5/box

### Indications

For acute and chronic wounds such as:

- Superficial and partial thickness burns
- Donor and recipient skin graft sites
- Lacerations and abrasions
- Leg and pressure ulcers
- Surgical incisions
- Blistering
- Skin tears

# Cuticell®

## Ointment dressing

### Product Benefits

- Small mesh sized, smooth acetate fabric
- Does not charge the wound with excess ointment
- Tissue-friendly, atraumatic dressing changes
- Helps to maintain a moist environment
- Does not fray when cut to size
- Prevents occlusive effects



Code	Size	Pack Qty
72539-02	7.5cm x 20cm	10/box
72539-04	20cm x 40cm	25/box

### Indications

For the treatment of superficial exuding wounds such as:

- Superficial and partial thickness burns
- Donor and recipient skin graft sites
- Lacerations and abrasions
- Radiation injuries

# Cuticell® Classic

## Paraffin gauze dressing

### Product Benefits

- Open weave 100% cotton fabric with interlocking threads
- Low-adherent and soothing
- Tissue-friendly, atraumatic dressing changes
- Helps to maintain a moist environment
- Does not fray when cut to size
- Permeable to air and exudates



Code	Size	Pack Qty
72538-00	5cm x 5cm	5/box
72538-01	5cm x 5cm	50/box
72538-02	10cm x 10cm	10/box
72538-03	10cm x 10cm	100/box
72538-06	10cm x 7m	1/box

### Indications

For the treatment of superficial exuding wounds such as:


- Superficial and partial thickness burns
- Donor and recipient skin graft sites
- Lacerations and abrasions
- Radiation injuries





# Portfolio overview

**Leukoplast®** **Leuko® - Tapes Guide**



Dressing type	Leukoplast	Leukoplast elastic	Leukoplast Waterproof	Leukopor	Leukosilk	Leukofix	Leukoflex	Leukoplast Fixomull	Leukoplast Fixomull Transparent	Leukoplast Skin Sensitive	Leukoplast Fixomull Skin Sensitive
Normal skin	○	○	○		○			○	○		
Sensitive skin				○		○	○			○	○
Damaged skin				○						○	○
Very old / young				○						○	○
Active, mobile people	○	○	○		○	○	○	○	○		
Injuries on joints		○					○	○			
Wet, dirty conditions			○				○		○		
Cosmetic areas						○	○				
Frequent dressing changes				○	○					○	○

## Surgical Tapes & Bandages A solution for all situations

Surgical tapes and bandages are a key modality for the management of wide range of medical conditions including soft tissue injuries, wound healing, and the fixation of medical devices. Essity has an extensive range of options to meet the varying demands of all these challenges. Leukoplast® and Handy® brands have seen a number of new developments such as Leukoplast® antimicrobial spools to minimize the risk of infection.

The Leukoplast® skin sensitive range of silicone based adhesive tapes are ideal for the customers with fragile & sensitive skin and those who tape frequently. With the addition of the skin sensitive range, Leukoplast® now have a tape for all situations: secure latex option, skin-friendly polyacrylate tapes and the gentle silicone adhesive solution.

Our Handy® bandaging range allows the choice of medium and heavy weight crepe bandages, light weight conforming and cohesive retention bandages and firm compression bandaging solutions.

### Bandage Guide



Dressing type	Retention				Light Support	Strong Support	Compression	
	Handyband Conforming	Handy Gauze Cohesive	Fixomull Stretch	Fixomull Transparent	Handy Crepe Medium	Handy Crepe Heavy	Comprilan	Eloflex
Keeps dressing in place	📦	📦	📦	📦	📦	📦		
Ongoing support for strains and sprains					📦	📦		
Restricts movement						📦		
Used for light compression					📦	📦		
Used for RICER / heavy compression							📦	📦
Sticks to itself / non-slip		📦						
Washable	📦				📦	📦	📦	📦
Used for under taping			📦	📦				
Skin friendly adhesive			📦	📦				
Waterproof				📦				

Wide area fixation

## Fixomull® Stretch

### Low-allergy retention tape

**Product Benefits**

- Porous – allowing for the evaporation of perspiration
- Polyacrylate adhesive mass is very well tolerated by the skin and alleviates skin trauma when taping



Code	Size	Pack Qty
02033-01	15cm X 2m	1/box
02036-00	5cm X 10m	1/box
02037-00	10cm X 10m	1/box
02038-00	15cm X 10m	1/box
02039-00	20cm x 10m	1/box
02040-00	30cm X 10m	1/box

**Indications**

- Device and dressing retention
- Regular shoulder strapping and taping
- For patellofemoral pain (taping of the knee)

## Hypafix®

### Conformable retention tape

**Product Benefits**

- Combines security with comfort
- Highly conformable
- Allows joints to be dressed without constraining movement
- Gentle secure adhesive ensures excellent dressing retention



Code	Size	Pack Qty
71443-00	2.5cm X 10m	1/box
71443-01	5cm X 10m	1/box
71443-02	10cm X 10m	1/box
71443-03	15cm X 10m	1/box
71443-04	20cm X 10m	1/box
71443-05	30cm X 10m	1/box

**Indications**

- Device and dressing retention

Wide area fixation: skin sensitive

## Fixomull® skin sensitive

### Retention tape with silicone adhesive

**Product Benefits**

- Breathable/no skin macerations
- Skin-friendly/atraumatic removal
- High wearing comfort
- Residue free removal



Code	Size	Pack Qty
79965-01	5cm X 5m	1/box
79965-02	10cm X 5m	1/box
79965-03	15cm X 5m	1/box
79965-04	10cm X 2m	1/box

**Indications**

- Fixation of primary wound dressings and additional securement of devices, such as catheters, tubes, probes or drainage tubes

## Hypafix® skin sensitive

### Retention tape with silicone adhesive

**Product Benefits**

- Reliable adhesion with low impact on fragile and compromised skin
- Atraumatic removal, protecting skin integrity
- Prevents skin stripping and tearing
- Breathable material prevents maceration
- Can be repositioned without stressing the skin or losing adhesive strength



Code	Size	Pack Qty
79966-01	5cm X 5m	1/box
79966-02	10cm X 5m	1/box

**Indications**

- Fixation of primary wound dressings and additional securement of devices, such as catheters, tubes, probes or drainage tubes

## Fixomull® Transparent

### Transparent wide area fixation tape

#### Product Benefits

- Transparent for easy inspection
- Reduced risk of tension blisters
- Breathable
- Residue-free removal
- Protects against friction
- Can be worn in the shower



Code	Size	Pack Qty
72217-00	10cm x 2m	1/box
72216-00	5cm x 10m	1/box
72216-01	10cm x 10m	1/box

#### Indications

- securing small and large wound dressings
- can be used in conjunction with waterproof bandages, for excellent wound cover

## Hypafix® Transparent

### Transparent wide area fixation tape

#### Product Benefits

- Transparent for easy inspection
- Breathable
- Residue-free removal
- Protects against friction
- Waterproof with reduced risk of contamination



Code	Size	Pack Qty
72378-00	5cm x 10m	1/box
72378-01	10cm x 10m	1/box

#### Indications

- securing small and large wound dressings, devices or instruments

## Leukoplast® Skin sensitive

### Skin sensitive silicone tape

#### Product Benefits

- Re-positionable
- Antimicrobial
- Skin-friendly/atraumatic removal
- Residue free removal
- Easy to tear by hand
- Residue free removal



Code	Size	Pack Qty
76173-00	1.25cm x 2.6m snap spool	24/box
76173-01	2.5cm x 2.6mtr snap spool	12/box
76175-00	2.5cm x 2.6mtr	12/box

#### Indications

- Fixation of wound dressings and additional securement of devices such as catheters, probes or drainage tubes
- Fragile & compromised skin
- Secure wound dressings & devices

## Leukoplast® Standard

### Rigid, tan, zinc oxide tape

#### Product Benefits

- Tan rigid tape
- Strong adhesion for long-lasting wear
- Porous, allowing skin to breathe
- Water-repellent
- Well tolerated by normal skin



Code	Size	Pack Qty
01691-01	1.25cm x 9.2m	24/drum
01692-01	2.5cm x 9.2m	12/drum
01690-00	3.8cm x 9.2m	8/drum
01693-01	5cm x 9.2m	6/drum
01521-00*	1.25cm x 5m	24/box
01522-00*	2.5cm x 5m	12/box
01524-00*	5cm x 5m	6/box
47725-00*	7.5cm x 5m	4/box

#### Indications

- Securing dressings
- Heavy duty and rough work
- Universal applications

\*Antimicrobial spools



## Leukoplast® elastic

### Elastic, tan, zinc oxide tape

#### Product Benefits

- Tan coloured
- Elastic for flexible area dressing retention
- Elastic, enabling easier application
- Flexible to conform to body movement
- Strong zinc oxide adhesive for long lasting wear
- Well tolerated by normal skin



Code	Size	Pack Qty
01061-00*	2.5cm x 1m	12/box
01062-00*	5cm x 1m	12/box
01071-00*	2.5cm x 2.5m	12/box
01072-00*	5cm x 2.5m	12/box
01073-00*	7.5cm x 2.5m	6/box

#### Indications

- Securing dressings
- Non-sensitive skin types

## Leukoplast® Waterproof

### Waterproof tape

#### Product Benefits

- Waterproof and dirt-repellent
- Strong zinc oxide adhesive for long lasting wear
- Well tolerated by normal skin
- White



Code	Size	Pack Qty
02321-00*	1.25cm x 5m	24/box
02322-00*	2.5cm x 5m	12/box
02324-00*	5cm x 5m	6/box

#### Indications

- Securing dressings
- Wet or dirty environments
- Heavy duty or rough work conditions
- Non-sensitive types

\*Antimicrobial spools

## Leukofix®

### Hypoallergenic transparent tape

#### Product Benefits

- Clear, hypoallergenic tape
- Hand-tearable in any direction
- Low-irritant and gentle on sensitive areas
- Water-repellent



Code	Size	Pack Qty
02136-00	1.25cm x 9.2m	24/box
02137-00	2.5cm x 9.2m	12/box
02138-00	5cm x 9.2m	6/box
02121-00*	1.25cm x 5m	24/box
02122-00*	2.5cm x 5m	12/box
02124-00*	5cm x 5m	6/box

#### Indications

- Securing dressings
- Sensitive skin
- Face, neck and cosmetic areas
- Active patients requiring strong adhesion

## Leukoflex®

### Hypoallergenic waterproof tape

#### Product Benefits

- Sheer, flexible, stretchable waterproof tape
- Hypoallergenic and gentle on sensitive skin
- Conforms to body movement
- Cosmetically appealing



Code	Size	Pack Qty
01121-00*	1.25cm x 5m	24/box
01122-00*	2.5cm x 5m	12/box
01124-00*	5cm x 5m	6/box
01131-00*	1.25cm x 5m	24/box
01132-00*	2.5cm x 5m	12/box
01134-00*	5cm x 5m	6/box

#### Indications

- Securing dressings
- Sensitive skin
- Hands which may be exposed to water or oil
- Injuries on or around joints
- Cosmetically sensitive areas
- Sealing in topical medications

## Leukosilk®

### High strength hypoallergenic tape

#### Product Benefits

- White, hypoallergenic cloth tape
- Strong adhesion, yet gentle on the skin
- Hand-tearable
- Water-repellent
- Porous, allowing skin to breathe



Code	Size	Pack Qty
09566-00*	1.25cm x 9.2m	24/box
09567-00*	2.5cm x 9.2m	12/box
09569-00*	5cm x 9.2m	6/box
01021-00*	1.25cm x 5m	24/box
01022-00*	2.5cm x 5m	12/box
01024-00*	5cm x 5m	6/box

#### Indications

- Securing dressings
- Sensitive skin
- Injuries on or around joints
- Cosmetically sensitive areas

## Leukopor®

### Hypoallergenic paper tape

#### Product Benefits

- Off-white, light weight, hypoallergenic paper tape
- Porous – allowing skin to breathe
- Water-repellent
- Easily hand-tearable



Code	Size	Pack Qty
02453-00	1.25cm x 9.2m	24/box
02454-00	2.5cm x 9.2m	12/box
02455-00	5cm x 9.2m	6/box
02456-00	7.5cm x 9.2m	4/box
02471-00*	1.25cm x 5m	24/box
02472-00*	2.5cm x 5m	12/box
02474-00*	5cm x 5m	6/box

#### Indications

- Elderly or young children
- Patients with sensitive skin
- Securing dressings
- Frequent dressing changes
- Injuries on delicate areas or damaged skin

\*Antimicrobial spools

## Leukoplast® Sleek

### Waterproof surgical tape

#### Product Benefits

- PVC substrate is waterproof and resistant to water, oil and grease
- Low profile tan colour
- The zinc oxide adhesive ensures secure fixation in difficult conditions
- Highly conformable ensuring easy application to varying limb contours



Code	Size	Pack Qty
72379-02	2.5cm X 5m	12/box
72379-03	5cm X 5m	12/box
72379-04	7.5cm x 5m	12/box

#### Indications

- Fixation of dressings and appliances where they are likely to be in contact with dirt or water
- Suitable for use in the workplace

## Leukoplast® Sleek LF

### Latex-free waterproof surgical tape

#### Product Benefits

- Latex free adhesive
- PVC substrate is waterproof and resistant to water, oil and grease
- Low profile tan colour
- Highly conformable ensuring easy application to varying limb contours



Code	Size	Pack Qty
72359-02*	2.5cm x 5m	12/box
72359-03*	5cm x 5m	12/box
72359-04*	7.5cm x 5m	12/box

#### Indications

- Fixation of dressings and appliances where they are likely to be in contact with dirt or water
- Suitable for use in the workplace

**Easifix®****Conforming retention bandage****Product Benefits**

- Highly conformable for bandaging difficult areas
- Safe high stretch low regain avoids limb constriction
- Cool light construction ensures excellent breathability



Code	Size	Pack Qty
71961-01	2.5cm x 1.75m	12/box
71429-10	5cm x 1.75m	12/box
71430-03	7.5cm x 1.75m	12/box
71428-04	10cm x 1.75m	12/box
71429-11	15cm x 1.75m	12/box

**Indications**

- Dressing retention on high-body movement areas such as joints

**Easifix® K****Bi conformable retention bandage****Product Benefits**

- Two-way stretch knit construction is highly conformable
- Open design improves breathability
- Finished edges prevent fraying



Code	Size	Pack Qty
72617-00	2.5cm x 2.4m	20/box
72617-01	5cm x 2.4m	20/box
72617-02	7.5cm x 2.4m	20/box
72617-03	10cm x 2.4m	20/box
72617-04	15cm x 2.4m	20/box

**Indications**

- Dressing retention on high-body movement areas such as joints

**Handyband® Conforming****Conforming gauze bandage****Product Benefits**

- Highly conformable, lightweight retention bandage
- Woven from viscose and polyamide yarns, this bandage is soft and porous, permitting air circulation and avoiding the risk of limb constriction
- High stretch and limited regain properties
- The carefully finished edges also minimise fraying



Code	Size	Pack Qty
92562-00	2.5cm x 1.5m	12/box
92562-01	5cm x 1.5m	12/box
92562-02	7.5cm x 1.5m	12/box
92562-03	10cm x 1.5m	12/box
92562-04	15cm x 1.5m	12/box

**Indications**

- Bandage retention on high-body movement areas such as joints
- Use where very light, comfortable dressing is required
- Active patients and prolonged use

**Easifix® Cohesive****Self adhesive retention bandage****Product Benefits**

- Economical 10m rolls can be cut so you only use what you need
- Cohesive construction sticks to itself - not to skin or clothing
- Non slip, highly conformable construction allows bandaging of difficult areas



Code	Size	Pack Qty
71439-21	4cm x 2m	10/box
71439-19	6cm x 2m	10/box
71439-16	8cm x 2m	10/box
71439-17	10cm x 2m	10/box
71439-20	6cm x 10m	10/box
71439-18	8cm x 10m	10/box

**Indications**

- Dressing retention on high-body movement areas such as joints
- Ideal for dressing bony prominences without impairing movement



# Handy® Gauze Cohesive

## Self adhesive retention bandage

### Product Benefits

- White, low bulk bandage
- Sticks to itself, not skin or clothing
- Air-permeable
- Economical: can be cut to any length

Code	Size	Pack Qty
92557-00	2.5cm x 2m	2/box
92557-01	6cm x 2m	1/box
92557-02	8cm x 2m	1/box
92557-03	10cm x 2m	1/box
92557-04	6cm x 10m	1/box
92557-05	8cm x 10m	1/box
92557-06	10cm x 10m	1/box



### Indications

- Non-slip dressing retention
- Difficult to dress and high body movement areas
- Self application

### Retention Bandages: Tubular

# Surgifix®

## Tubular elastic net bandage

### Product Benefits

- The elasticated knitted design allows easy application to the most difficult areas
- Roll format can be cut ensuring you only use what you need.
- The light, cool design is breathable and comfortable

Code	Size	Unstretched Length	Pack Qty
73172-03	Size 0	9m	1/box
73172-04	Size 1	9m	1/box
73172-05	Size 2	9m	1/box
73172-06	Size 3	9m	1/box
73172-07	Size 4	9m	1/box
73172-08	Size 5	9m	1/box
73172-12	Size 5.5	9m	1/box
73172-09	Size 6	9m	1/box
73172-10	Size 7	9m	1/box



### Indications

- Ideal for dressing retention on heads, fingers and other awkward sites
- Dressing retention in chronic wound treatment.
- Versatile bandage to fit a broad range of indications

# Tricofix®

## Tubular gauze bandage

### Product Benefits

- Light weight cotton design is comfortable on the skin
- Low regain characteristics make it ideal for use in sub bandage applications

Code	Size	Pack Qty
02191-00	1.5cm x 20m (Size A/1)	1/box
02192-00	2.5cm x 20m (Size B/2)	1/box
02193-00	4cm x 20m (Size C/3)	1/box
02195-00	6cm x 20m (Size D/5)	1/box
02197-00	8cm x 20m (Size E/6)	1/box
02198-00	10cm x 20m (Size F/7)	1/box
02199-00	12cm x 20m (Size G/9)	1/box
02201-00	16cm x 10m (Size K)	1/box
02202-00	21cm x 10m (Size L)	1/box



### Indications

- Protecting the skin surface under compression, zinc gel or adhesive bandages
- Can be used as a lining for casting or compression bandaging

### Light Support Bandages: Non-Adhesive

# Elastolite®

## Light support crepe bandage

### Product Benefits

- High twist cotton yarns ensure excellent stretch and regain characteristics thus permitting free movement while maintaining effective therapy
- Selvedged edges prevent fraying
- Non bleached natural colour
- Washable

Code	Size	Pack Qty
NON-STERILE		
72092-25	5cm x 1.5m	12/box
72092-28	7.5cm x 1.5m	12/box
72092-26	10cm x 1.5m	12/box
72092-27	15cm x 1.5m	12/box
STERILE		
72092-58	5cm x 1.5m	160/box
72092-56	7.5cm x 1.5m	136/box
72092-59	10cm x 1.5m	90/box
72092-57	15cm x 1.5m	18/box



### Indications

- Ideal for use on awkward contours of the body

## Handy® Crepe Hospital

### Light support crepe bandage

#### Product Benefits

- 100% high twist cotton yarns provide good stretch and regain
- Clinical white finish
- Selvedged edges prevent fraying

Code	Size	Pack Qty
73072-17	5cm x 1.5m	12/box
73072-18	7.5cm x 1.5m	12/box
73072-21	10cm x 1.5m	12/box
73072-24	15cm x 1.5m	12/box



#### Indications

- Dressing retention where light compression is required
- Post-operative use

## Handy® Crepe Medium Weight

### Light support crepe bandage

#### Product Benefits

- White
- Cool and comfortable bandaging
- Non-fraying

Code	Size	Pack Qty
73057-03	5cm x 1.6m	12/box
73057-00	7.5cm x 1.6m	12/box
73057-01	10cm x 1.6m	12/box
73057-05	15cm x 1.6m	12/box



#### Indications

- Dressing retention where light compression is required
- Post-operative use

## Tensocrepe® Medium Weight

### Light support crepe bandage

#### Product Benefits

- Premium quality white bandage
- Economical – can be washed and reused
- Cool and absorbent construction makes it comfortable for the patient
- Provides good soft tissue support

Code	Size	Pack Qty
72082-25	5cm x 1.5m	12/box
72082-26	7.5cm x 1.5m	12/box
72082-22	10cm x 1.5m	12/box
72082-23	15cm x 1.5m	12/box



#### Indications

- Dressing retention where light compression is required
- Post-operative use

#### Light Support Bandages: Tubular

## Tensogrip®

### Tubular support bandage

#### Product Benefits

- Easily applied, circular knitted construction
- Design ensures consistent even pressure is applied
- Available in tan or white
- Economical – cut off only what you need
- Washable, re-usable and durable

Code	Size	Colour	Pack Qty
71513-00	Size A 4.5cm x 10m	White	1/box
71514-00	Size B 6.25cm x 10m	White	1/box
71519-00	Size C 6.75cm x 10m	White	1/box
71515-00	Size D 7.5cm x 10m	White	1/box
71520-00	Size E 8.75cm x 10m	White	1/box
71516-00	Size F 10cm X 10m	White	1/box
71521-00	Size G 12cm x 10m	White	1/box
71517-00	Size J 17.5cm x 10m	White	1/box
71523-00	Size K 21.5cm x 10m	White	1/box
71522-00	Size L 32.5cm x 10m	White	1/box
71514-01	Size B 6.25cm x 10m	Tan	1/box
71519-01	Size C 6.75cm x 10m	Tan	1/box
71515-01	Size D 7.5cm x 10m	Tan	1/box
71520-01	Size E 8.75cm x 10m	Tan	1/box
71516-01	Size F 10cm X 10m	Tan	1/box
71521-01	Size G 12cm x 10m	Tan	1/box



#### Indications

- Mild strains and sprains
- Mild oedema
- Firm device retention
- For easy easy/fast application and reapplication

## Co-Plus® LF

### Latex free, hand tearable self adhesive bandage

#### Product Benefits

- Latex free
- Hand tearable in both directions
- Flexible, non slip substrate adheres to itself, not to skin, hair or clothing
- Provides good tissue support and some compression if required
- Does not shrink sideways – avoids limb constriction
- Wide range of colours to suit individual needs



Code	Size	Colour	Pack Qty
72100-14	2.5cm x 2m	Tan	60/box
72100-16	5cm X 2m	Tan	36/box
72100-17	7.5cm X 2m	Tan	24/box
72100-18	10cm X 2m	Tan	18/box
72100-19	15cm X 2m	Tan	12/box
72100-20	7.5cm X 2m	Colour*	24/box
72100-21	10cm x 2m	Colour*	18/box

#### Indications

- For bandaging awkward anatomical areas
- Active patients with a low tolerance of adhesives
- Mild strains and sprains
- Mild oedema
- Firm device retention

\* mixed colour pack containing: light blue, dark blue, red, pink, green, purple

#### Light Support Bandages: Cohesive

## Co-Plus®

### Self adhesive light support bandage

#### Product Benefits

- Adheres to itself, not to skin, hair or clothing
- Flexible non slip substrate allows the bandage to conform to awkward contours
- Provides good tissue support and some compression if required
- Does not shrink sideways – avoids limb constriction
- Wide range of colours to suit individual needs



Code	Size	Colour	Pack Qty
72100-06	2.5cm x 3m	Tan	60/box
72100-08	5cm x 3m	Tan	36/box
72100-09	7.5cm x 3m	Tan	24/box
72100-04	10cm x 3m	Tan	18/box
72100-05	15cm x 3m	Tan	12/box
72102-03	7.5cm x 3m	Colour*	24/box
72102-02	10cm x 3m	Colour*	18/box

#### Indications

- For bandaging awkward anatomical areas
- Active patients with a low tolerance of adhesives
- Mild strains and sprains
- Mild oedema
- Firm device retention

\* mixed colour pack containing: light blue, dark blue, red, pink, green, purple

#### Strong Support Bandages: Non-Adhesive

## Tensocrepe® Heavy Weight BP

### Strong support crepe bandage

#### Product Benefits

- BP specification – ensures reliable performance
- 100% cotton high twist yarns deliver firm support
- Central yellow line facilitates application
- Washable and reusable



Code	Size	Colour	Pack Qty
72082-31	5cm x 2.3m	Tan	12/box
72082-32	7.5cm x 2.3m	Tan	12/box
72082-29	10cm x 2.3m	Tan	12/box
72082-30	15cm x 2.3m	Tan	12/box
72082-20	5cm x 2.3m	White	12/box
72082-21	7.5cm x 2.3m	White	12/box
72082-18	10cm x 2.3m	White	12/box
72082-19	15cm x 2.3m	White	12/box

#### Indications

- Dressing retention where light, peripheral compression is required
- On-going support for minor sprains and strains
- Medium to long-term bandaging
- Restricting movement

## Tensocrepe® Hospital Heavy Weight

### Strong support crepe bandage

#### Product Benefits

- Economical cotton and polyamide construction
- Cool breathable fabric
- Soft selvedged edges are comfortable and prevent fraying
- Washable and reusable



Code	Size	Pack Qty
72082-15	5cm x 2.3m	12/box
72082-17	7.5cm x 2.3m	12/box
72082-16	10cm x 2.3m	12/box
72082-14	15cm x 2.3m	12/box

#### Indications

- Dressing retention where light, peripheral compression is required
- On-going support for minor sprains and strains
- Medium to long-term bandaging
- Restricting movement



## Handy® Crepe Heavy

### Strong support crepe bandage

#### Product Benefits

- Tan coloured
- Centre guide for easy application
- Washable and reusable
- Non-fraying



Code	Size	Colour	Pack Qty
73050-00	5cm x 2.3m	100% cotton	12/pack
73050-02	7.5cm x 2.3m	100% cotton	12/pack
73050-03	10cm x 2.3m	100% cotton	12/pack
73050-04	15cm x 2.3m	100% cotton	12/pack
73050-10	5cm x 2.3m	Viscose/polyamide	12/pack
73050-12	7.5cm X 2.3m	Viscose/polyamide	12/pack
73050-14	10cm x 2.3m	Viscose/polyamide	12/pack
73050-16	15cm x 2.3m	Viscose/polyamide	12/pack

#### Indications

- Dressing retention where light, peripheral compression is required
- On-going support for minor sprains and strains
- Medium to long-term bandaging
- Restricting movement

## Acrylastic®

### Hypoallergenic elastic adhesive bandage

#### Product Benefits

- 100% cotton fabric is comfortable and cool
- The polyacrylate adhesive is latex free and suitable for use on sensitive skin
- Low profile tan colour



Code	Size	Pack Qty
02697-00	6cm x 2.5m	1/box
02698-00	8cm x 2.5m	1/box
02699-00	10cm x 2.5m	1/box

#### Indications

- Treatment of all conditions which require limitation of movement, e.g. sprains and strains, sport injuries
- Prophylactic use in active sport
- Strong support and consistent control for injured ligaments or joints

#### Strong Support Bandages: Adhesive

## Tensoplast®

### Elastic adhesive bandage

#### Product Benefits

- Very secure zinc oxide adhesive ensures very secure adhesion
- Central yellow line facilitates application
- Feathered selvedged edges improve comfort and stop fraying
- Air blown adhesive allows the skin to breathe



Code	Size	Pack Qty
72050-51	5cm x 2.5m	12/box
72050-44	7.5cm x 2.5m	12/box
72050-47	10cm x 2.5m	12/box

#### Indications

- Treatment of all conditions which require limitation of movement, e.g. sprains and strains, sport injuries
- Prophylactic use in active sport
- Strong support and consistent control for injured ligaments or joints

#### Compression Bandages: Non-Adhesive Single Layer

## JOBST® Comprilan®

### Short stretch compression bandage

#### Product Benefits

- The interlocking weave provides the short stretch performance to help ensure consistent graduated compression
- Low resting pressure, well tolerated by patients at night
- 100% cotton bandage is comfortable to the skin
- Easy care and durable can be washed and reused
- Low profile tan colour



Code	Size	Pack Qty
01026-00	6cm x 2.6m	1/box
01027-00	8cm x 2.6m	1/box
01028-00	10cm x 2.6m	1/box
01029-00	12cm x 2.6m	1/box

#### Indications

- Sustained compression for the management of vascular and lymphatic conditions
- Suitable for RICER (soft tissue, injury management)
- Oedema control



# Tensopress®

## Elastic compression bandage

### Product Benefits

- High stretch substrate provides therapeutic compression for both ambulatory and non ambulatory patients
- Central green line aids in application
- Washable and durable
- Low profile – can be worn with existing footwear



Code	Size	Pack Qty
71723-01	7.5cm x 3m	10/box
71723-00	10cm x 3m	10/box

### Indications

- Sustained compression for the management of vascular and lymphatic conditions
- Suitable for RICER (soft tissue, injury management)
- Oedema control

# Eloflex®

## High stretch compression bandage

### Product Benefits

- High stretch bandage offers strong compression
- Ideal for sporting related injuries where oedema is present
- Washable and reusable



Code	Size	Pack Qty
02476-00	6cm x 1.25m	1/box
02478-00	8cm x 1.25m	1/box
02488-00	8cm x 3.5m	1/box
02490-00	10cm x 3.5m	1/box
02492-00	12cm x 3.5m	1/box

### Indications

- Sustained compression for the management of vascular and lymphatic conditions
- Suitable for RICER (soft tissue, injury management)
- Oedema control

## Compression Bandages: 2-layer

# JOBST® Compriz

## 2-layer compression bandage system

### Product Benefits

- Provides therapeutic graduated compression
- Printed ellipses change to circles when correct tension has been applied for ease of application
- Can be left in place for up to 7 days



Code	Compression Level	Ankle Size	Pack Qty
76271-00	Regular 40mmHg	18cm - 25cm	1/box
76271-01	Regular 40mmHg	25cm - 32cm	1/box
76271-02	Lite 30mmHg	18cm - 25cm	1/box
76271-03	Lite 30mmHg	25cm - 32cm	1 box

### Indications

- First line compression therapy treatment of venous leg ulcers
- Also suitable for Lymphoedema, Lipedema, night compression, heavy & swollen legs and thrombosis prophylaxis



# Portfolio overview

Dressings	Everyday Wound Dressings				
		Leukoplast® Visual	Leukoplast® Kids	Leukoplast® Universal	Leukoplast® Elastic
	Everyday Wound Dressings				
		Leukoplast® Strong	Leukoplast® Soft White		
	Detectable Dressings				
		Leukoplast® Detectable			
	Waterproof Dressings				
		Leukoplast® Barrier	Leukoplast® Aqua Pro		

## First Aid Dressings



Caring for wounds on fingers, fingertips, or even moving body parts: Leukoplast provides a choice of first aid dressings that meets all needs! Whether it is the type and size of the wound, the skin condition, allergy problems, an infection-sensitive environment - there are Leukoplast first aid dressings for all conditions, which no first aid kit should be without!

## Leukoplast® Visual

Blue sterile bandage strips and adhesive dressings for minor wounds

### Product benefits

- Highly visible and easily found if lost, with its contrasting signal blue colour
- Waterproof
- Bacterial barrier
- Highly absorbent, low adherent wound pad



Leukoplast® Visual		
Code	Size	Pack Qty
76451-00	22 x 72mm Strips Sterile	100/box
76451-01	22 X 72mm Strips Sterile	50/box

### Indications

- Minor wounds such as cuts, lacerations and abrasions
- Food handling and catering environments

## Leukoplast® Kids

Water- and dirt-repellent dressings specially for children to treat minor wounds

### Product benefits

- Water- and dirt-repellent
- Skin-friendly
- Reliable adhesion
- Fun kids motifs make it easy to smile again after a small accident



Leukoplast® Kids		
Code	Size	Pack Qty
76458-05	Strips Assortment Non-Sterile 8 x 19x63 mm 4 x 38x63 mm	12/box

### Indications

- Minor cuts and grazes

## Leukoplast® Universal

Finger dressings for short-term use on minor wounds, lacerations or punctures

### Product benefits

- Well suited to sensitive skin
- Keeps wound clean
- Flexible and easy to use
- Secure, reliable adhesion



Leukoplast® Universal		
Code	Size	Pack Qty
76461-02	19 x 72 mm Strips Non-Sterile	20/box
76462-00	19 x 72 mm Strips Sterile	50/box
76462-01	19 x 72 mm Strips Sterile	100/box
76462-02	19 x 72 mm Strips Sterile	10/box
76462-03	19 x 72 mm Strips Sterile	50/box
76462-04	19 x 72 mm Strips Sterile	20/box

### Indications

Short-term use on minor wounds, lacerations, abrasions or puncture sites

## Leukoplast® Elastic

Flexible first aid dressings for minor wounds, lacerations, abrasions or punctures

### Product benefits

- Highly elastic
- Low-adherent wound pad
- Quick treatment of minor wounds
- Reliable adhesive
- Available in different sizes



Leukoplast® Elastic		
Code	Size	Pack Qty
76455-23	Ø 2.2 cm Non-Sterile	100/box

### Indications

- For minor wounds, lacerations, abrasions or puncture sites.
- Ideal for dressing fingers or joints and other frequently moved areas



## Leukoplast® Strong

Highly adhesive water-repellent bandage for short-term use on minor wounds

### Product benefits

- Secure fixation
- Highly absorbent wound pad
- Breathable
- Water repellent



Leukoplast® Strong		
Code	Size	Pack Qty
76459-07	4cm x 1m Dressing Non-Sterile	1 Roll/Box
76459-05	6cm x 1m Dressing Non-Sterile	1 Roll/Box
76459-01	6cm x 5m Dressing Non-Sterile	1 Roll/Box
76459-02	8cm x 5m Dressing Non-Sterile	1 Roll/Box
76459-10	22 x 63mm Strips Non-Sterile	20/box
76460-05	22 x 38mm Strips Sterile	100/box
76460-06	22 x 63mm Strips Sterile	100/box
76460-07	22 x 72mm Strips Sterile	100/box
76460-08	22 x 72mm Strips Sterile	50/box
76460-10	38 x 75mm Finger Strips Sterile	100/box
76460-11	50 x 63mm Finger Strips Sterile	100/box
76460-12	45 x 50mm Finger Strips Sterile	50/box
76459-09	Classic Strips Assortment 22 x 36 mm 50 x 63 mm 38 x 75 mm 50 x 72 mm	20/box

### Indications

Short-term use on minor wounds, lacerations or puncture sites

## Leukoplast® Soft White

Breathable, soft and highly absorbent wound dressing with reliable adhesion

### Product benefits

- Skin friendly
- Breathable
- Highly absorbent
- Conformable and soft



Leukoplast® Soft White		
Code	Size	Pack Qty
7645026	1.9 cm x 4 cm Injection Dressing Non-Sterile	500/box
7645027	6 cm x 10 cm Strips Non-Sterile	10/box

### Indications

For minor wounds, lacerations, abrasions or puncture sites.

## Leukoplast® Detectable

X-Ray, electromagnetically and visually detectable waterproof dressing

### Product benefits

- Detectable and visible (contains metal and x-ray detectable materials)
- Waterproof
- Breathable
- Bacteria barrier



Leukoplast® Detectable		
Code	Size	Pack Qty
76453-03	22 x 72mm Strips Sterile	100/box
76453-01	28 x 72mm Strips Sterile	50/box
76453-04	50 x 72mm Strips Sterile	100/box
76453-05	44 x 50mm Fingerstrips Sterile	50/box
76453-02	Finger Strips Assortment Sterile 44 x 50 mm 38 x 75 mm	20/box

### Indications

- Minor wounds such as cuts, lacerations and abrasions
- Food handling and catering environments



## Leukoplast® Barrier

*Waterproof yet breathable dressing for minor wounds helps to prevent infection*

### Product benefits

- Waterproof yet breathable
- Absorbent yet non-stick wound pad
- Flexible and conformable
- Protects against contamination



Leukoplast® Barrier		
Code	Size	Pack Qty
76452-00	22 x 38mm Strips Sterile	100/box
76452-01	22 x 63mm Strips Sterile	100/box
76452-08	22 x 72mm Strips Sterile	100/box
76452-02	38 x 38mm Strips Sterile	100/box
76452-04	38 x 63mm Strips Sterile	100/box
76452-05	50 x 72mm Strips Sterile	100/box
76452-07	44 x 50mm Fingerstrips Sterile	50/box

### Indications

Minor wounds including cuts and grazes

## Leukoplast® Aqua Pro

*Waterproof yet breathable wound dressing protecting against water and contaminants*

### Product benefits

- Waterproof
- Secure fixation
- Protective microbial barrier
- Breathable
- Available in different sizes
- Skin-friendly



Leukoplast® Aqua Pro		
Code	Size	Pack Qty
76457-08	Strips Assortment Non-Sterile ø 24 mm 19 x 72 mm 38 x 63 mm	20/box

### Indications

Short term use on minor wounds including lacerations and puncture sites



Essity Health & Medical  
PO Box 337, Mount Waverley,  
Victoria, Australia 3149  
[medical.essity.com.au](http://medical.essity.com.au)

Customer Service Australia  
T 1300 276 633 F 1300 998 830  
Customer Service New Zealand  
T 0508 276 111 F 0508 998 830