

JOBST® READY- TO-WEAR COMPRESSION



DEFY GRAVITY



Ready-to-wear compression garments for venous disorders and post venous surgery (including vascular ablation, sclerotherapy and vein stripping)

JOBST.CO.UK



JOBST® OPAQUE





Ready-to-Wear RAL Compression Hosiery

JOBST® Opaque is a range of unisex ready-to-wear, circular-knit hosiery for the treatment of mild to moderate oedema, venous conditions, and post venous surgery.

JOBST® Opaque garments come in pairs and can be worn for up to six-months*.

UNISEX AND ADAPTABLE

JOBST® Opaque is a unisex garment offering a wide range of product options for both male and female patients, helping to cater for a variety of patient needs and conditions.

SIZES AND OPTIONS / MATERIALS					
Six sizes	Two lengths (Petite and Regular)	Wide circ.	Made-to-measure**	Open Toe	Closed Toe
✓	✓	✓	✓	✓	✓
<div><div> Knitted top band</div><div> Lace silicone top band</div><div> Dotted silicone top band</div><div> Soft silicone top band</div></div>					

COMPRESSION CLASS

Class 1 - 18-21 mmHg

Class 2 - 23-32 mmHg

Class 3 - 34-46 mmHg**

COLOURS



Sand



Black



Caramel



Bronze



Navy

STYLES



Knee
High



Thigh
High



Tights



*Six month guarantee applicable for garment failure caused by manufacturing or material defect. Please refer to the instructions for use for 'Care Instructions'

**CCL3 is only available as a custom fit / made to measure garment

JOBST® ULTRASHEER




Ready-to-Wear RAL Compression Hosiery

JOBST® UltraSheer is a range of ready-to-wear, circular-knit hosiery for the treatment of mild to moderate oedema, venous conditions, and post venous surgery.

JOBST UltraSheer garments come in pairs and can be worn for up to six-months*.

ELEGANT AND STYLISH

JOBST® UltraSheer provides outstanding sheerness and a silky look for an elegant and stylish appearance inspired by womens' dress tights.

SIZES AND OPTIONS / MATERIALS				
Six sizes	Two lengths (Petite and Regular)	Wide circ.	Open Toe	Closed Toe
✓	✓	✓	✓	✓
<div><div> Lace silicone top band</div><div> Soft silicone top band</div><div> Decorative Dotted</div></div>				

COMPRESSION CLASS

Class 1 - 18-21 mmHg

Class 2 - 23-32 mmHg

COLOURS



STYLES



*Six month guarantee applicable for garment failure caused by manufacturing or material defect. Please refer to the instructions for use for 'Care Instructions'

JOBST® FORMEN AMBITION AND JOBST® FORMEN EXPLORE

Ready-to-Wear RAL Compression Hosiery

JOBST® ForMen is a range of ready-to-wear, circular-knit hosiery for the treatment of mild to moderate oedema, venous conditions, and post venous surgery.

JOBST ForMen garments come in pairs and can be worn for up to six-months*.

AMBITION - STYLISH DRESS SOCK

Inspired by a man's dress sock, the timeless ribbed design is ideal for work and formal occasions.

EXPLORE - DURABLE WITH ADDED COMFORT

Reinforced foot zones offer increased strength and durability for active wear, combined with a high cotton content for added comfort and reduced potential for skin irritations.

SIZES AND OPTIONS / MATERIALS

Six sizes	Two lengths (Regular and Long)	Closed Toe
✓	✓	✓

COMPRESSION CLASS

Class 1 - 18-21 mmHg

Class 2 - 23-32 mmHg

COLOURS - JOBST® FORMEN AMBITION



Black Dark Grey Brown Navy Khaki

COLOURS - JOBST® FORMEN EXPLORE



Black Dark Grey Navy Khaki

STYLES



Knee
High



*Six month guarantee applicable for garment failure caused by manufacturing or material defect. Please refer to the instructions for use for 'Care Instructions'





Find ordering and prescribing information fast with JOBST Code Finder.

No more searching through a catalogue, use our web app on your mobile device to find the information you need.

CODEFINDER.JOBST.COM



CHECK OUT OUR RANGE OF SUPPORT TOOLS TO HELP MAKE CHOOSING JOBST EVEN EASIER



Download printable copies of JOBST order forms

You can download copies of the order forms for the full JOBST range that can be printed on standard office printers at jobst.co.uk/professionals/resources/measuring-forms/ or scan the QR to access.



Measuring and donning help videos

Scan the QR code to access our UK playlist with help videos on measuring and donning JOBST ready-to-wear garments