

Cutimed® Siltec® Plus

Cutimed® Siltec® Plus is a silicone coated, sterile, absorbent polyurethane foam dressing for atraumatic dressing changes that contains super absorbent stripes which absorb and lock wound exudate.

Product indication

Cutimed® Siltec® Plus is intended for the treatment of exuding wounds with low to high exudate levels such as venous and arterial ulcers, pressure ulcers, diabetic foot ulcers, surgical incisions, skin grafts and donor sites, lacerations or abrasions.

Product benefits



Reliable exudate management

The vertical absorption makes sure exudate is taken up reliably and kept inside the dressing - even under compression. Super absorbent stripes additionally support the absorption capacity.^{1,4}



Gentle with high adherence

The perforated silicone wound contact layer - with higher adherence than Cutimed® Siltec® - prevents adherence to newly formed tissue, protects wound margins and ensures atraumatic, pain-reduced dressing changes.^{1,2,3,5}



Skin-friendly

The dressing has a low potential for skin irritations and allergies.¹



Water-repellent and breathable backing

The outer film is bacteria proof, water-repellent and ensures proper fluid evaporation for changing stages of saturation.




Can be cut to size

The dressing can also be cut to size, making it as comfortable for the patient as possible.



Applications



1 Seckam A., A multicentre, observational evaluation of the product characteristics of two absorbent foam dressings. *Br J Nurs.* Jun 27, 2019;28(12):S10-S17.


2 King B., A clinical evaluation of 20 patients when using a new absorbent silicone foam wound dressing: Cutimed Siltec B. *Wounds UK.* 2018; Vol 14. No 3.

3 Zimmermann S., *In vitro* exudate management performance of a new hydrophilic foam dressing range. *Poster EWMA.* 2017; Amsterdam.

4 Bateman SD., 150 patient experiences with a soft silicone foam dressing. *Br J Nurs.* Jun 25-Jul 8, 2015;24(12):S16, S18-23.

















5 Bateman, S.D., The challenge of skin tears and lacerations in an emergency environment. *Wounds UK.* 2014; 10. 94-99.

Cutimed® Siltec® Plus

Ordering Information	Size	Product Code	Pack size	Nappi Code
Cutimed® Siltec® Plus	5 x 6 cm	73288-00	10	278976-001
Highly absorbent	10 x 10 cm	73288-01	10	278977-001
Silicone foam dressing	10 x 20 cm	73288-02	10	278978-001
with stronger initial adhesion	15 x 15 cm	73288-03	10	278980-001
	20 x 20 cm	73288-04	5	278983-001

**LATEX
FREE
FORMULA**

Icon overview

Improved guidance with symbols	
Exudate	 Dry  Low  Moderate  High  Excessive  Low to Excessive
Wound depth	 superficial  deep  superficial and deep
Wound state	 Infected  Sloughy  Granulating  Epithelializing  Necrotic
Both side application	 Not on pack (in IFU only)
Secondary dressing needed	

Find icons here

