

Cutimed® Sorbion® Sachet Multi Star

Cutimed® Sorbion® Sachet Multi Star is a sterile, hydroactive, gel-forming, super absorbent wound dressing.

Product indication

Cutimed® Sorbion® Sachet Multi Star is intended for the treatment of exuding wounds with moderate to excessive exudate levels, such as decubitus, ulcer cruris, diabetic foot ulcers and other wounds healing through secondary intention, oozing laparotomy wounds, fistulas to the skin, exulcerating carcinomas, postoperative wound dehiscence and similar exuding wounds.

Product benefits



Outstanding absorption capacity

Outstanding absorption capacity¹ and high retention performance. It is capable of absorbing up to 44.7 % more than comparable super absorbent dressings.

Protection against maceration through safe, vertical absorption.



Time- and cost-saving

Time- and cost-saving thanks to reduced frequency of dressing changes and therapy efforts.²



Flexible and conformable

It is especially suitable for application on contoured body parts, such as the toes or heels.



Wound cleansing effect

Creates optimal wound healing conditions and provides soft debridement, minimizing the treatment effort.^{3,4}



Apply either side down

Can be used on both sides for easy and comfortable application.



Applications



1 BSN medical GmbH. In vitro study, Cutimed® Sorbion® is capable of absorbing and retaining more exudate than comparable super absorbent dressings. Testing was performed by the independent laboratory Surgical Material Testing Laboratory (S.M.T.L.). 2017.

2 Panca Met *et al.* Clinical and cost-effectiveness of absorbent dressings in the treatment of highly exuding VLU. *J Wound Care.* Mar 2013;22(3):109-10, 112-8.

3 Romanelli M, Vowden K, Weir D. Exudate management made easy. *Wounds Int.* 2010; 1(2):1-6.

4 Romanelli M *et al.* A pilot study evaluating the wound and skin care performances of the Hydration Response Technology dressing: a new concept of debridement. *Journal of Wound Technology.* 2009; Vol. 5.

