

EVERYDAY Ready-To-Wear

JOBST® Ultrasheer



Name and Style	S		M		L		XL		Colour
	Open Toe	Closed Toe	Open Toe	Closed Toe	Open Toe	Closed Toe	Open Toe	Closed Toe	
Knee High	75165-30	75165-25	75166-30	75166-25	75167-30	75167-25	75168-30	75168-25	● Natural Regular
		75408-02		75409-02		75410-02		75415-02	● Natural Petite
		75165-04		75166-04		75167-04		75168-04	● Black
		73633-00		73633-01		73633-02		73633-03	● Diamond Pattern Black
		75201-19		75202-19		75203-19		75204-19	● Natural Regular
		75435-02		75436-02		75437-02		75438-02	● Natural Petite
		75201-04		75202-04		75203-04		75204-04	● Black
		73634-00		73634-01		73634-02		73634-03	● Diamond Pat Black
		75241-21		75242-21		75243-21		75244-21	● Natural
		75241-04		75242-04		75243-04		75244-04	● Black
Thigh High		75169-22		75170-22		75171-22		75172-22	● Natural
		75169-04		75170-04		75171-04		75172-04	● Black
		73635-00		73635-01		73635-02		73635-03	● Diamond Pattern Black
		75205-23		75206-23		75207-23		75208-23	● Natural
		75205-04		75206-04		75207-04		75208-04	● Black
		73636-00		73636-01		73636-02		73636-03	● Diamond Pattern Black
		75245-20		75246-20		75247-20		75248-20	● Natural
		75245-04		75246-04		75247-04		75248-04	● Black
Sensitive Thigh High		77686-12		77686-13		77686-14		77686-15	● Natural
		77686-04		77686-05		77686-06		77686-07	● Black
		77687-12		77687-13		77687-14		77687-15	● Natural
		77687-04		77687-05		77687-06		77687-07	● Black
		77688-12		77688-13		77688-14		77688-15	● Natural
		77688-04		77688-05		77688-06		77688-07	● Black
Waist High		75173-12		75174-12		75175-12		75176-12	● Natural
		75173-04		75174-04		75175-04		75176-04	● Black
		73637-00		73637-01		73637-02		73637-03	● Diamond Pattern Black
		75209-13		75210-13		75211-13		75212-13	● Natural
		75209-04		75210-04		75211-04		75212-04	● Black
		73638-00		73638-01		73638-02		73638-03	● Diamond Pattern Black
		75249-13		75250-13		75251-13		75252-13	● Natural
		75249-04		75250-04		75251-04		75252-04	● Black
Maternity Waist High		75177-04		75178-04		75179-04		75180-04	● Natural
		75177-01		75178-01		75179-01		75180-01	● Black
		75213-06		75214-06		75215-06			● Natural
	75213-03		75214-03		75215-03			● Black	

SITUATIONAL Ready-To-Wear

JOBST® Travel Sock



Length	Size 1	Size 2	Size 3	Size 4	Size 5	Colour
	Closed Toe	Closed Toe	Closed Toe	Closed Toe	Closed Toe	
Knee High	78844-07	78845-07	78846-07	78847-07	79272-07	● Black
	78844-12	78845-12	78846-12	78847-12	79272-12	● Beige

JOBST® Travel Sock measurement table

Size	1	2	3	4	5
Ankle Circumference	18 - 21 cm	20 - 23 cm	22 - 25 cm	24 - 27 cm	26 - 29 cm
Calf Circumference	27 - 37 cm	30 - 40 cm	33 - 43 cm	36 - 46 cm	39 - 49 cm
AUS Women's Shoe	4.5 - 6.5	7 - 9	9.5 - 11	11.5 - 13.5	-
AUS Men's Shoe	3.5 - 5.5	6 - 8	8.5 - 10	10.5 - 12	12.5 - 14
EU Shoe	37 - 38	39 - 40	41 - 42	43 - 44	45 - 46
UK Shoe	3 - 5	5 - 7	7 - 9	9 - 11	11 - 13

JOBST® Activewear



Length	S	M	L	XL	Colour
	Closed Toe	Closed Toe	Closed Toe	Closed Toe	
Knee High	75147-03	75148-03	75149-03	75150-03	● Black
	75147-02	75148-02	75149-02	75150-02	○ White
	75155-03	75156-03	75157-03	75158-03	● Black
	75155-02	75156-02	75157-02	75158-02	○ White
	75411-03	75412-03	75413-03	75414-03	● Black
	75411-02	75412-02	75413-02	75414-02	○ White

JOBST® Sport



Length	S	M	L	XL	Colour
	Closed Toe	Closed Toe	Closed Toe	Closed Toe	
Knee High	75289-70	75289-71	75289-72	75289-73	● Black
	75289-80	75289-81	75289-82	75289-83	● Royal Blue
	75289-90	75289-91	75289-92	75289-93	● Pink

JOBST® forMen



Length	S	M	L	XL	Colour
	Closed Toe	Closed Toe	Closed Toe	Closed Toe	
Knee High	75253-01	75254-01	75255-01	75256-01	● Black
	75257-01	75258-01	75259-01	75260-01	● Black
	75265-01	75266-01	75267-01	75268-01	● Black
Thigh High	75285-00	75286-00	75287-00	75288-00	● Black
	75261-00	75262-00	75263-00	75264-00	● Black
	75269-00	75270-00	75271-00	75272-00	● Black

JOBST® forMen Casual



Length	S	M	L	XL	Colour
	Closed Toe	Closed Toe	Closed Toe	Closed Toe	
Knee High	75476-02	75477-02	75478-02	75479-02	● Black
	75476-03	75477-03	75478-03	75479-03	● Khaki Beige
	75484-02	75485-02	75486-02	75487-02	● Black
	75484-03	75485-03	75486-03	75487-03	● Khaki Beige
	75492-02	75493-02	75494-02	75495-02	● Black
	75492-03	75493-03	75494-03	75495-03	● Khaki Beige

JOBST® Relief



S	S	M	M	L	L	L	L	XL	XL	XL	XL	Colour
Open Toe	Closed Toe	Open Toe	Closed Toe	Open Toe	Closed Toe	Open Toe Full Calf	Closed Toe Full Calf	Open Toe	Closed Toe	Open Toe Full Calf	Closed Toe Full Calf	
Relief Garments Knee High												
78649-00	75339-12	78649-01	75340-12	78649-02	75341-12	78649-04	75343-12	78649-03	75342-12	78649-05	75344-12	● Black
75339-10	75339-11	75340-10	75340-11	75341-10	75341-11	75343-10	75343-11	75342-10	75342-11	75344-10	75344-11	● Beige
78052-20	78052-32	78052-21	78052-33	78052-22	78052-34	78052-24	78052-36	78052-23	78052-35	78052-25	78052-37	● Black
78052-26	78052-38	78052-27	78052-39	78052-28	78052-40	78052-30	78052-42	78052-29	78052-41	78052-31	78052-43	● Black Petite
78050-04	78050-15	78051-04	78051-15	78052-04	78052-15	77975-15	77975-16	78053-04	78053-15	77976-15	77976-16	● Beige
78054-04	78054-15	78055-04	78055-15	78056-04	78056-15	77977-15	77977-16	78057-04	78057-15	77978-15	77978-16	● Beige
Relief Garments Thigh High												
78649-20	75345-02	78649-21	75346-02	78649-22	75347-02			78649-23	75348-02			● Black
78650-20	78042-07	78650-21	78043-07	78650-22	78044-07			78650-23	78045-07			● Black
78042-08	78042-09	78043-08	78043-09	78044-08	78044-09			78045-08	78045-09			● Beige
78059-08	78059-09	78047-08	78047-09	78048-08	78048-09			78049-08	78049-09			● Beige
Relief Garments Waist High												
78476-10	78476-11	78477-10	78477-11	78478-10	78478-11			78479-10	78479-11			● Beige
78480-10		78481-10		78482-10				78483-10				● Beige
Relief Garments Chaps												
78031-05		78032-05		78033-05				78034-05				● Beige L
78031-06		78032-06		78033-06				78034-06				● Beige R
78035-05		78036-05		78037-05				78041-05				● Beige L
78035-06		78036-06		78037-06				78041-06				● Beige R

JOBST Compression Classes

- CCL 1: 15-20 mmHg
- CCL 2: 20-30 mmHg
- CCL 3: 30-40 mmHg

THE RIGHT FIT. LOOK NO FURTHER



Trusted JOBST® Partners

You can purchase JOBST Ready-To-Wear Compression garments from a range of online stockists and local pharmacies. Scan the QR code for a list of all the online stockists



JOBST® How to Measure Videos

Please scan the QR code to view videos explaining how to measure for the following garments- Knee High, Thigh High and Waist High.



JOBST®, an Essity brand

Always read the label and follow the directions for use. ESS0876 v2 12/2024

Essity Health & Medical
PO Box 337, Mount Waverley,
Victoria, Australia 3149
medical.essity.com.au
jobst.com.au

JOBST® Customer Service AU
T 1300 998 810 F 1300 998 820
JOBST® Customer Service NZ
T 0508 998 810 F 0508 998 820

JOBST® LOWER LIMB COMPRESSION GARMENTS



JOBST® DEFY GRAVITY

JOBST® Medical Grade, Ready-To-Wear garments with graduated compression help to improve blood flow and reduce oedema

What type of compression does my patient need?

1 Select the required compression

- CCL 1: 15-20 mmHg**
 - Tired aching legs
 - Minor ankle, leg and foot swelling
 - Minor varicoses during pregnancy
 - Long distance travel
- CCL 2: 20-30 mmHg**
 - Severe varicose veins
 - Lymphatic oedema
 - varicoses during pregnancy
 - Post vein surgery/post sclerotherapy
 - Moderate tendency of oedema
 - Manages manifestations of Post Thrombotic Syndrome
 - Helps prevent recurrence of Superficial thrombophlebitis (PTS), Chronic Venous Insufficiency CVI/CVD
 - Aids in prevention of thrombotic syndrome
- CCL 3: 30-40 mmHg**
 - Severe varicose veins
 - Lymphatic oedema
 - varicoses during pregnancy
 - Post vein surgery/post sclerotherapy
 - Manages manifestations of Post Thrombotic Syndrome
 - Helps prevent recurrence of Superficial thrombophlebitis (PTS), Chronic Venous Insufficiency CVI/CVD
 - Aids in prevention of thrombotic syndrome

2 Select garment style

Some considerations when selecting the most appropriate style for your patient:

- Indication being treated/managed
- Lifestyle, mobility and age
- Ability to don and tolerate level of compression

Open toe available

3 Measure for JOBST® Medical LegWear

Ultrasheer, Relief, for Men, for Women Casual, ActiveWear and Sport

Measuring positions

- CH - Hip** Circumference at largest dimension around buttocks
- CG - Thigh** Circumference at largest dimension
- CC - Calf** Circumference at largest dimension
- CB - Ankle** Circumference at smallest dimension

Provide the following measurements

- Knee High - CB & CC
- Thigh High - CB, CC & CG
- Waist High - CB, CG & CH

NOTE Each leg is measured independently and measure in the early morning when possible as swelling often increases later in the day.

Ultrasheer: Petite Knee High Length 34 - 39 cm, Regular Knee High Length 39 - 45 cm

Style	Measurement	S	M	L	L Full Calf	XL	XL Full Calf
MATERNITY	Hip	81 - 132 cm	86 - 142 cm	91 - 152 cm	112 - 180 cm		
	Thigh	40 - 62 cm	46 - 70 cm	54 - 78 cm	60 - 81 cm		
	Ankle	18 - 21 cm	21 - 25 cm	25 - 29 cm	29 - 33 cm		
UNISEX	Hip	71 - 117 cm	76 - 127 cm	81 - 137 cm	102 - 166 cm		
	Thigh	40 - 62 cm	46 - 70 cm	54 - 78 cm	60 - 81 cm		
	Calf	28 - 38 cm	30 - 42 cm	32 - 46 cm	34 - 50 cm		
	Ankle	18 - 21 cm	21 - 25 cm	25 - 29 cm	29 - 33 cm		