

JOBST® ELVAREX SOFT

Ready to Wear

Flat knit compression garments



- JOBST® Elvarex® Soft Knee High garments are the perfect alternative for patients where management with circular knitted garments are no longer suitable due to uncontrolled oedema
- Indicated for the management of mild Lymphoedema and vascular insufficiency
- Soft and Skin Friendly: The soft knitted yarn combines skin protecting softness with excellent breathability, making it an ideal solution for those with sensitive and fragile skin
- Therapeutically effective: Combining the right level of compression with the right fabric stiffness, JOBST® Elvarex Soft offers an extremely comfortable and effective treatment for low to moderate oedema



Sizing & Compression Classes

CCL 2: 20-30 mmHg

CCL 3: 30-40 mmHg



JOBST® Elvarex Soft Regular Calf

JOBST® Elvarex Soft 20-30mmHg Silicone Band - short 2.5 cm, long 5 cm						
Size	I	II	III	IV	V	VI
Regular (ID AD= 40-45)	79559-012	79559-013	79559-014	79559-015	79559-016	79559-017
Petite (ID AD= 34-39)	79558-012	79558-013	79558-014	79558-015	79558-016	79558-017
JOBST® Elvarex Soft 30-40mmHg Silicone Band - short 2.5 cm, long 5 cm						
Size	I	II	III	IV	V	VI
Regular (ID AD= 40-45)	73552-012	73552-013	73552-014	73552-015	73552-016	73552-017
Petite (ID AD= 34-39)	73579-012	73579-013	73579-014	73579-015	73579-016	73579-017
Size	I	II	III	IV	V	VI
cD	28 - 32 cm	30 - 34 cm	32 - 36 cm	34 - 38 cm	36 - 40 cm	38 - 42 cm
cC	29 - 33 cm	31 - 35 cm	33 - 37 cm	35 - 39 cm	37 - 41 cm	39 - 43 cm
cB1	23 - 26 cm	25 - 28 cm	27 - 30 cm	29 - 32 cm	31 - 34 cm	33 - 36 cm
cB	18 - 20 cm	20 - 22 cm	22 - 24 cm	24 - 26 cm	26 - 28 cm	28 - 30 cm
cY	26 - 31 cm	28 - 33 cm	39 - 35 cm	31 - 37 cm	32 - 38 cm	33 - 40 cm
cA	17 - 22 cm	19 - 24 cm	21 - 26 cm	23 - 29 cm	25 - 32 cm	27 - 34 cm

JOBST® Elvarex Soft PLUS Calf

JOBST® Elvarex Soft 20-30mmHg Silicone Band - short 2.5 cm, long 5 cm						
Size- PLUS	I Plus	II Plus	III Plus	IV Plus	V Plus	VI Plus
Regular (ID AD= 40-45)	79559-018	79559-019	79559-020	79559-021	79559-022	79559-023
Petite (ID AD= 34-39)	79558-018	79558-019	79558-020	79558-021	79558-022	79558-023
JOBST® Elvarex Soft 30-40mmHg Silicone Band - short 2.5 cm, long 5 cm						
Size- PLUS	I Plus	II Plus	III Plus	IV Plus	V Plus	VI Plus
Regular (ID AD= 40-45)	73552-018	73552-019	73552-020	73552-021	73552-022	73552-023
Petite (ID AD= 34-39)	73579-018	73579-019	73579-020	73579-021	73579-022	73579-023
Size chart	I Plus	II Plus	III Plus	IV Plus	V Plus	VI Plus
cD Plus	32 - 36 cm	34 - 38 cm	36 - 40 cm	38 - 42 cm	40 - 44 cm	42 - 46 cm
cC plus	33 - 37 cm	35 - 39 cm	37 - 41 cm	39 - 43 cm	41 - 45 cm	43 - 47 cm
cB1 plus	24 - 27 cm	26 - 29 cm	28 - 31 cm	30 - 33 cm	32 - 35 cm	34 - 37 cm
cB	18 - 20 cm	20 - 22 cm	22 - 24 cm	24 - 26 cm	26 - 28 cm	28 - 30 cm
cY	26 - 31 cm	28 - 33 cm	39 - 35 cm	31 - 37 cm	32 - 38 cm	33 - 40 cm
cA	17 - 22 cm	19 - 24 cm	21 - 26 cm	23 - 29 cm	25 - 32 cm	27 - 34 cm



JOBST®,
an Essity brand

Always read the label and
follow the directions for use.

ESS0971
10/2024

Essity Health & Medical

PO Box 337, Mount Waverley,
Victoria, Australia 3149
medical.essity.com.au
jobst.com.au

JOBST® Customer Service AU
T 1300 998 810 F 1300 998 820

JOBST® Customer Service NZ
T 0508 998 810 F 0508 998 820