

DEFY GRAVITY



THE EASY-TO-USE TWO LAYER BANDAGE SYSTEM



JOBST® COMPRI2



JOBST® COMPRI2

EFFECTIVE OEDEMA REDUCTION. EASY TO USE.

The JOBST® Compri2 - 2 Layer Bandage System is especially designed for **easy application**. It provides **graduated and reliable compression** for up to seven days to help achieve **effective oedema reduction**. The 1st layer - a nonwoven padding is **tearable and stretchable** in two dimensions which makes it easy to add and adjust padding to vulnerable pressure points. The 2nd layer - a short-stretch bandage with a **simple visual compression indicator**, supports an easy, safe and appropriate application of the compression level.

Scannable **QR codes** and a handy, easy to read IFU, offer precise instructions for **achieving correct application**, to deliver high working pressure and low resting pressure, 24 hours a day.

1ST LAYER - NONWOVEN PADDING



- Offers even pressure distribution
- Stretchability in both dimensions simplifies application
- Tearable padding layer bandage allows additional protection to sensitive areas and lowers risk of skin damage due to pressure marks
- PU layer lamination means no slip on the skin and easy application





2ND LAYER - SHORT-STRETCH BANDAGE

- Cohesive short-stretch material delivers high working pressure and low resting pressure
- Allows re-adjustment during application
- Visual indicator allowing simple and effective application of correct compression level
- Two compression levels; standard and lite, offers compression to suit individual needs

ORDERING INFORMATION:

QR code leads to simple application video



	Ankle size	REF No.	Kits per box	HIMF Code	
Regular compression (ABPI >0.8) (at least 40 mm Hg)					
JOBST® Compri2	18 - 25 cm 7 1/8" - 9 3/4"	76271-00	1		
		76271-04	6		
	25 - 32 cm 9 3/4" - 12 1/2"	76271-01	1	710909	
		76271-05	6		
Lite compression (ABPI 0.5 - 0.8) (reduced average pressure of 30 mm Hg)					
JOBST® Compri2 lite	18 - 25 cm 7 1/8" - 9 3/4"	76271-02	1		
		76271-06	6		
	25 - 32 cm 9 3/4" - 12 1/2"	76271-03	1	791235	
		76271-07	6		

Always read the label and follow the directions for use.



JOBST®,
an Essity brand

ESS0974
11/2024

Essity Health & Medical

PO Box 337, Mount Waverley,
Victoria, Australia 3149
medical.essity.com.au
jobst.com.au

JOBST® Customer Service AU

T 1300 276 633 E customerservice.AU@essity.com

JOBST® Customer Service NZ

T 1300 276 111 E customerservice.NZ@essity.com