

# SUPPORTED SELF-CARE: HOW MUCH COULD YOU SAVE?

JOBST® FARROWWRAP® &  
JOBST® ULCERCARE



COMPRESSION FOR EMPOWERING PATIENT SELF-CARE  
[JOBST.CO.UK](http://JOBST.CO.UK)

# REDUCE YOUR CLINICAL WASTE TO LANDFILL

Do you know how much more waste is generated using two-layer compression bandages instead of JOBST® FarrowWrap® or JOBST UlcerCare? Let us put it into perspective.

The below scenarios are based on **only one patient with one venous leg ulcer (VLU) on one leg**, treated over six-months<sup>1-9</sup>.

Based on three bandage changes per week, up to **72 bandage changes** are needed<sup>1-9</sup>. Double this amount (up to 144 bandage changes) if both legs are bandaged...

The combined height<sup>10</sup> (m) of two-layer compression bandages\* used is up to **914.4 METRES**



**THAT'S UP TO  
2.8 EIFFEL TOWERS  
IN HEIGHT OF  
CLINICAL WASTE  
TO LANDFILL  
PER LEG<sup>10</sup>**



Two-layer compression bandages are single use only and cannot be recycled. This clinical waste must go to a landfill.

**How does this impact your Trust's sustainability targets?**

# NURSING TIME IMPACT

Compression bandages require specialist training to apply, meaning **up to 72 nurse visits** would be needed to heal a VLU over a six month period<sup>3</sup>.

Each visit takes on average 45 minutes<sup>11</sup>, this would total **up to 54 hours of nursing time**.

- By switching to JOBST FarrowWrap or JOBST UlcerCare<sup>3</sup>:**
- Patients can apply and adjust the garments themselves as part of supported self-care
  - If the patient cannot apply the garments, carers or family members can apply (specialist nursing training is not needed)

**UP TO 54 HOURS  
OF NURSE VISITS**

**UP TO 16.5 TIMES  
MORE COST  
INCURRED**

## ADDITIONAL COST OF BANDAGES<sup>12,13</sup>

Compression bandages need to be changed up to three times per week<sup>1-10</sup>, meaning a significantly higher cost of overall treatment.



	JOBST Compr2	JOBST FarrowWrap 4000	JOBST UlcerCare
NHS Drug Tariff price per unit	£8.34	£111.32	£34.29
Total cost of treatment over six months	£600.48 (up to 72 bandage changes <sup>1-10</sup> )	£111.32	£34.29
<b>Potential savings (6 monthly) by switching to JOBST FarrowWrap or JOBST UlcerCare Up to £489.16 - £566.19</b>			

# MOVE PATIENTS TO SUPPORTED SELF-CARE

**One JOBST FarrowWrap or JOBST UlcerCare** provides the same six month's worth of treatment per patient currently using compression bandages.

JOBST FarrowWrap and JOBST UlcerCare have a six month guarantee against garment failure (caused by manufacturing or material defect\*\*).



\*Comparison is based on JOBST Compr2 two-layer compression bandaging stretched.  
\*\*Please refer to the instructions for use for 'Care Instructions'

**Disclaimer**  
The patient treatment plan utilising compression bandaging is indicated following an in-depth patient assessment and diagnosis. The frequency of bandage changes per week varies. The reported average is two to three times per week<sup>2-6</sup>. During the decongestion stage for lymphoedema management bandage changes may be required daily<sup>7-9</sup>.

**References**  
Available on request



## JOBST® FARROWWRAP® 4000

Suitable for patients with **moderate to severe** venous conditions **with or without** a **venous leg ulcer** present.


- Designed for easy application with just four bands
- Made with an inner sleeve to aid donning
- Helps promote healing of venous ulceration
- Comes with a JOBST FarrowHybrid compression sock (20-30mmHg wide)

Lower limb styles available on Drug Tariff



Legpiece

Colours



Tan

Black

Compression level

30-40mmHg

When you order a JOBST FarrowWrap 4000 you will receive 1 x JOBST FarrowHybrid (20-30mmHg)

JOBST FarrowHybrid	X-small	Small	Medium	Large
Size received	Medium-wide	Medium-wide	Large-wide	Large-wide*

Note: Open wounds must be covered with an appropriate dressing before the compression garment is applied.

\*Do not apply JOBST FarrowHybrid if calf circumference exceeds 60 cm, consider a JOBST FarrowWrap footpiece in this circumstance. When a Large JOBST FarrowWrap 4000 is ordered, an additional non-compressive liner will be provided to accommodate calf circumference greater than 60 cm.



## JOBST® FARROWHYBRID®

Treatment of **mild to moderate** oedema in the **foot only**.


- Non-compressive above the ankle
- Ability to wear normal shoes
- Available in a regular or wide foot
- Supplied in pairs

Lower limb styles available on Drug Tariff



JOBST FarrowHybrid

Colours



Tan

Black

Compression level

20-30mmHg

UK	Small	Medium	Large	X-Large
Shoe size	2-5	6-8	9-11	12-13

Regular - **A'** = 20-27 cm

Wide - **A'** = 28-35 cm



## JOBST® FARROWWRAP® LITE

Suitable for patients with **mild to moderate** oedema when **reduced compression** is required.

- Can support patients with sensitive skin
- Liner included

Lower limb styles available on Drug Tariff



Footpiece



Legpiece



Kneepiece



Thighpiece\*

Colours



Tan

Compression level

20-30mmHg

\*Thighpiece comes with knee piece



# JOBST® FARROWWRAP® STRONG

Suitable for patients with **moderate to severe oedema** with minimal shape distortion and minor skin folds

- Reliable hold for fluctuating oedema
- Liner included

## Lower limb styles available on Drug Tariff



Footpiece



Legpiece



Kneepiece



Thighpiece\*

### Colours



Tan

### Compression level

30-40mmHg

\*Thighpiece comes with knee piece



# JOBST® FARROWWRAP® CLASSIC

Made from more rigid material than JOBST FarrowWrap Strong. Suitable for patients with **moderate to severe oedema** with irregular shaped limbs and deep skin folds.

- Rigid fabric with reinforced spine helps with oedema reduction
- Designed to sit flat against skin folds and not dig into the skin
- Liner included

## Lower limb styles available on Drug Tariff



Footpiece



Legpiece



Kneepiece



Thighpiece\*

### Colours



Tan

### Compression level

30-40mmHg

\*Thighpiece comes with knee piece



# JOBST® FARROW® TOE CAP

Suitable for patients with **mild to moderate** oedema in the toes. Can be worn in conjunction with JOBST® FarrowWrap® footpiece.

- Can support patients with sensitive skin
- Fine soft seams for greater wearing comfort
- Durable and soft double laminated fabric
- Can be trimmed to patient requirements

## Lower limb styles available on Drug Tariff



Toe Cap

### Colours



Tan



Black

### Compression level

15-20mmHg

20-30mmHg





# JOBST® ULCERCARE

2-in-1 compression system for the **treatment of active venous leg ulcers (VLU)**.

- Kit provides 40mmHg of compression at the ankle, as recommended by the National Wound Care Strategy Programme (NWCSP) for the healing of VLUs<sup>1</sup>
- Liners provide 15-20mmHg of mild compression, ideal for Immediate and Necessary Care (NWCSP)<sup>1</sup>
- Liners contain silk for ease of donning
- Only two measurements required
- Available in seven sizes as a Ready-to-Wear (RtW) range. Also available as a Custom-Fit (MtM) garment
- Ability to wear normal shoes helping to reduce the risk of falls compared to compression bandages from an improved range of ankle motion (ROAM)<sup>14</sup>
- Available with or without a zipper

## Styles available on Drug Tariff



Legpiece

### Colours

JOBST UlcerCare Kit



Tan



Black

JOBST UlcerCare Liner



White

### Compression level

40mmHg

JOBST UlcerCare Kit

15-20mmHg

JOBST UlcerCare Liner



1. National Wound Care Strategy Programme: (2024) Recommendations for Leg Ulcers.

14. Atkin, L., Stephenson, J., Parfitt, G. et al. An investigation to assess ankle mobility in healthy individuals from the application of multi-component compression bandages and compression hosiery. J Foot Ankle Res 9, 18 (2016). <https://doi.org/10.1186/s13047-016-0151-8>



# CHECK OUT OUR RANGE OF SUPPORT TOOLS TO HELP MAKE CHOOSING JOBST EVEN EASIER



## Download printable copies of JOBST order forms

You can download copies of the order forms for the full JOBST range that can be printed on standard office printers at [jobst.co.uk/professionals/resources/measuring-forms/](https://jobst.co.uk/professionals/resources/measuring-forms/) or scan the QR to access.



## Measuring and donning help videos


Scan the QR code to access our UK playlist with help videos on measuring and donning JOBST FarrowWrap and JOBST UlcerCare garments



Find ordering and prescribing information fast with JOBST Code Finder.

No more searching through a catalogue, use our web app on your mobile device to find the information you need.

[CODEFINDER.JOBST.COM](https://codefinder.jobst.com)





# STEP

Sustainable  
Treatment  
Economical  
Practice



To learn more about the benefits of moving patients to supported self-care by choosing re-usable compression therapy solutions, scan the QR code or visit: <https://www.jobst.co.uk/professionals/resources/step-sustainable-treatment-economical-practice>

