

JOBST® Relief

15-20mmHg compression garment

Ready-to-wear order form and prescription request


Email order to Customer Services on: compression.uk@jobst.com

Date: _____ Purchase Order No.: _____ Patient Name: _____

Measured By: _____ Tel: _____ Email: _____

Delivery Address: _____

Post Code: _____



Measuring instructions

1. Measure narrowest ankle **B** and widest calf **C** circumferences
2. Measure leg length from floor to **D** (two finger widths below the crease of the knee)
3. Determine size (**S-XLWC**) and length (**Regular**) using chart below

	Small	Medium	Large	X Large	L Wide Calf	XL Wide Calf
C	28 - 38cm	30 - 42cm	32 - 46cm	34 - 50cm	46 - 61cm	46 - 61cm
B	18 - 21cm	21 - 25cm	25 - 29cm	29 - 33cm	25 - 30cm	30 - 36cm
ID Regular	> 38cm					

☐ **JOBST® Relief**

Size: ☐ Small ☐ Medium ☐ Large ☐ X Large
☐ Large Wide Calf ☐ X Large Wide Calf

Length **ID** ☒ Regular

Colour: ☐ Black ☐ Beige

Toe: ☐ Open Toe ☐ Closed Toe

PIP Code: -

Product Code: -

Quantity Required:

Repeat prescription required every months

JOBST Relief Order Codes

JOBST® Relief (15-20mmHg)

Product Code	Size	PIP Code
Beige - Closed Toe		
75339-11	S	436-2885
75340-11	M	436-2893
75341-11	L	436-2901
75342-11	XL	436-2919
75343-11	L Wide Calf	436-2927
75344-11	XL Wide Calf	436-2992
Black - Closed Toe		
75339-12	S	436-2844
75340-12	M	436-2851
75341-12	L	436-2869
75342-12	XL	436-2877
75343-12	L Wide Calf	436-2935
75344-12	XL Wide Calf	436-2984

Product Code	Size	PIP Code
Beige - Open Toe		
75339-10	S	436-2786
75340-10	M	436-2794
75341-10	L	436-2802
75342-10	XL	436-2810
75343-10	L Wide Calf	436-2943
75344-10	XL Wide Calf	436-2968
Black - Open Toe		
75339-20	S	436-2760
75340-20	M	436-2778
75341-20	L	436-2836
75342-20	XL	436-2828
78649-40	L Wide Calf	436-2950
78649-41	XL Wide Calf	436-2976

Find ordering and prescribing information fast with JOBST Code Finder

No more searching through a catalogue, use our web app on your mobile device to find the information you need.

Scan the QR code or search for codefinder.jobst.com on your mobile device.



Go digital and order using JOBST Online

Bringing you an easier way to order JOBST. Order direct with payment by account or on prescription by selecting JOBST Delivered* as the dispenser on JOBST Online.

Scan QR code to register or visit order.jobst.com/UK

*JOBST Delivered is operated by Daylong Direct, 10 Cossall Industrial Estate, Ilkeston, Derby, DE7 5UG

Privacy Policy: www.daylong.co.uk

