

# JOBST® Confidence

Knee High Custom-Fit Order Form

Fax order to customer services on: **0845 122 3450**

Email order to customer services on: **compression.uk@jobst.com**

Date: \_\_\_\_\_ Purchase Order No.: \_\_\_\_\_ Patient Name: \_\_\_\_\_ DoB: \_\_\_\_\_  
 Measured By: \_\_\_\_\_ Tel: \_\_\_\_\_ Email: \_\_\_\_\_  
 Delivery Address: \_\_\_\_\_ Invoice Address: \_\_\_\_\_  
 \_\_\_\_\_




<input type="checkbox"/> <b>JOBST® Confidence</b>	
<b>Compression Class (RAL)</b>	<b>Quantity</b>
<input type="checkbox"/> CCL 1 (18-21mmHg)	Left <input type="text"/>
<input type="checkbox"/> CCL 2 (23-32mmHg)	Right <input type="text"/>
<input type="checkbox"/> CCL 3 (34-46mmHg)	
<b>Style and Options</b>	
<input type="checkbox"/> Slipform	<input type="checkbox"/> Ankle functional zone
<input type="checkbox"/> Deco line	<input type="checkbox"/> Initials (max 2 letters) <input type="text"/>
<input type="checkbox"/> Hallux valgus (bunion relief zone)*	
<input type="checkbox"/> Closed toe	*Only available with slant foot
<b>Colour</b>	
<input type="checkbox"/> Beige	<input type="checkbox"/> Black
<input type="checkbox"/> Caramel	<input type="checkbox"/> Red heather
<input type="checkbox"/> Anthracite heather	<input type="checkbox"/> Jeans heather
Note: the ankle functional zone on red, anthracite and jeans heather coloured garments is more visible	
<b>Top band option</b>	
<input type="checkbox"/> 5cm SoftFit	
<b>Remarks</b>	

### Left Leg Length Measurements (cms)

### Circumference Measurements (cms)

### Right Leg Length Measurements (cms)




### Foot Length - Left

For open toe \_\_\_\_\_ cm  
 For closed toe \_\_\_\_\_ cm (longest toe)


### Foot Length - Right

For open toe \_\_\_\_\_ cm  
 For closed toe \_\_\_\_\_ cm (longest toe)

### Left Foot Slant



### Right Foot Slant



Note: medial = inside  
 lateral = outside

# JOBST® Confidence Ordering Information

Compression Class / Option    Code    Pieces

## JOBST® Confidence Knee High



### Basic style

CCL 1 (18-21mmHg)	L1-16-04	1
CCL 2 (23-32mmHg)	L2-17-04	1
CCL 3 (34-46mmHg)	L3-18-04	1

### Options


Closed toe	L-A001C	1
SoftFit	L-A016C	1
Slipform	L-A011C	1
Ankle functional zone	L-A010C	1

## Options



### Ankle Functional Zone

Circumferential functional zone allows greater movement and flexibility at the ankle



### Slipform

Slipform provides a higher shape to the outside of the garment, at the top, conforming to the limb shape



### Deco Line

A decorative seam line can be added to the front of the garment on request



### Hallux Valgus (bunion relief zone)\*

Creates a looser knitting pattern around the ball of the foot to provide relief to painful bunions



### SoftFit Top Band

Soft and secure top band that has been uniquely knitted into the garment top band for a more comfortable secure fit



### Initials

Patient initials can be knitted into the garment below the JOBST logo for a personal touch

\*Only available with slant foot