

JOBST® FarrowWrap® Lite


Armpiece & Gauntlet

Ready-to-wear order form and prescription request

Email order to Customer Services on: **compression.uk@jobst.com**
 Fax order to Customer Services on: **0345 122 3450**

Date: _____ Purchase Order No.: _____ Patient Name: _____ DoB: _____
 Measured By: _____ Tel: _____ Email: _____
 Delivery Address: _____ Invoice Address: _____

Armpiece



Measuring instructions

Armpiece

1. Measure narrowest part of the wrist **C** circumference
2. Measure elbow **E** circumference, ensuring arm is bent slightly
3. Measure axilla **G** circumference (top of the arm before it meets the armpit)
4. Measure arm **C-G** length from the outside of the wrist to the axilla, ensuring arm is bent slightly
5. Determine size (**S-L**), length (**Short to Long**) & colour (**Tan or Black**)

Gauntlet

1. Measure the circumference of the palm at the base of the fingers **A**
2. Measure the circumference of the palm at the thumb fold **B**
3. Determine size (**XS-L**)

Check out our measuring guide videos

Scan the QR code for help with measuring for JOBST® FarrowWrap® Lite



Armpiece

	Small	Medium	Large
G Axilla	22 - 31cm	29 - 39cm	32 - 45cm
E Elbow	20 - 27cm	25 - 34cm	30 - 40cm
C Wrist	14 - 18cm	16 - 21cm	19 - 25cm
C-G Long	48 - 51cm		
C-G Regular	44 - 47cm		
C-G Short	40 - 43cm		

Gauntlet

	X Small	Small	Medium	Large
B Palm at fold of thumb	< 19cm	19 - 21cm	21 - 23cm	> 24cm
A Palm at base of fingers	< 18cm	18 - 20cm	20 - 22cm	> 22cm

Armpiece

Size: Small Medium Large
Length: Short Regular Long
Arm: Left Right
Colour: Tan Black

PIP Code: -
 Product Code: -
 Quantity required:
 Repeat prescription required every months

Gauntlet

Size: X Small Small Medium Large
Colour: Tan Black

PIP Code: -
 Product Code: -
 Quantity required:
 Repeat prescription required every months

JOBST® FarrowWrap® Lite Order Codes

JOBST® FarrowWrap® Lite (20-30mmHg) Armpiece & Gauntlet

For patients with **mild to moderate** oedema caused by venous or lymphatic conditions, where reduced compression is required

Product Code	Option	Size	Length	Colour	PIP Code
Armpiece Left					
76650-11	Left	Small	Short	Tan	392-6953
76650-05	Left	Small	Regular	Tan	392-7019
76650-35	Left	Small	Long	Tan	392-7076
76650-13	Left	Medium	Short	Tan	392-6961
76650-07	Left	Medium	Regular	Tan	392-7027
76650-01	Left	Medium	Long	Tan	392-7084
76650-15	Left	Large	Short	Tan	392-6979
76650-09	Left	Large	Regular	Tan	392-7035
76650-03	Left	Large	Long	Tan	392-7092
76650-29	Left	Small	Short	Black	392-6987
76650-23	Left	Small	Regular	Black	392-7043
76650-13	Left	Small	Long	Black	392-7100
76650-31	Left	Medium	Short	Black	392-6995
76650-25	Left	Medium	Regular	Black	392-7050
76650-19	Left	Medium	Long	Black	392-7118
76650-33	Left	Large	Short	Black	392-7001
76650-27	Left	Large	Regular	Black	392-7068
76650-21	Left	Large	Long	Black	392-7126

Product Code	Option	Size	Length	Colour	PIP Code
Armpiece Right					
76650-12	Right	Small	Short	Tan	392-6771
76650-06	Right	Small	Regular	Tan	392-6839
76650-00	Right	Small	Long	Tan	392-6896
76650-14	Right	Medium	Short	Tan	392-6789
76650-08	Right	Medium	Regular	Tan	392-6847
76650-02	Right	Medium	Long	Tan	392-6904
76650-16	Right	Large	Short	Tan	392-6797
76650-10	Right	Large	Regular	Tan	392-6854
76650-04	Right	Large	Long	Tan	392-6912
76650-30	Right	Small	Short	Black	392-6805
76650-24	Right	Small	Regular	Black	392-6862
76650-18	Right	Small	Long	Black	392-6920
76650-32	Right	Medium	Short	Black	392-6813
76650-26	Right	Medium	Regular	Black	392-6870
76650-20	Right	Medium	Long	Black	392-6938
76650-34	Right	Large	Short	Black	392-6821
76650-28	Right	Large	Regular	Black	392-6888
76650-22	Right	Large	Long	Black	392-6946

Product Code	Option	Size	Length	Colour	PIP Code
Gauntlet					
76653-00	Ambidextrous	Extra Small	-	Tan	404-5647
76653-01	Ambidextrous	Small	-	Tan	404-5654
76653-02	Ambidextrous	Medium	-	Tan	404-5662
76653-03	Ambidextrous	Large	-	Tan	404-5670

JOBST Code Finder: find ordering and prescribing info easily

No more searching through a catalogue, use our web app to find the information you need.

Scan the QR code or visit
[jobst.co.uk/
professionals/
ordering/find-
ordering-codes](http://jobst.co.uk/professionals/ordering/find-ordering-codes)



Go digital and order using JOBST Online

Bringing you an easier way to order JOBST compression garments. Order direct with payment by account or on prescription by selecting JOBST Delivered* as the dispenser. Scan QR code to register or visit order.jobst.com/UK

*JOBST Delivered is operated by Daylong Direct, 10 Cossall Industrial Estate, Ilkeston, Derby, DE7 5UG



JOBST Delivered: fast prescription delivery

To use this dispensing service, please send this form and prescription to:
**JOBST Delivered, FREEPOST,
10 Cossall Industrial Estate,
Ilkeston, Derby, DE7 5UG.**

If using an Electronic Prescription Service (EPS), please select code **FJ708** to nominate Daylong Direct as the dispenser or order online at jobstonline.co.uk and select JOBST Delivered as the dispenser.

For envelopes, please contact:
jobstdelivered@daylongdirect.co.uk
or **0808 196 6250**