

# JOBST® FarrowWrap®


## A-D Below Knee

Ready-to-wear order form and prescription request

Email order to Customer Services on: **compression.uk@jobst.com**  
 Fax order to Customer Services on: **0345 122 3450**

Date: \_\_\_\_\_ Purchase Order No.: \_\_\_\_\_ Patient Name: \_\_\_\_\_ DoB: \_\_\_\_\_  
 Measured By: \_\_\_\_\_ Tel: \_\_\_\_\_ Email: \_\_\_\_\_  
 Delivery Address: \_\_\_\_\_ Invoice Address: \_\_\_\_\_

### Legpiece



### Measuring instructions

#### Legpiece

1. Measure narrowest ankle **B** and widest calf **C** circumferences
2. Measure length posteriorly (back of leg) following the contour from **A** (floor) to **D** (two finger widths below the crease of the knee)
3. Determine size (**XS-XL**) & length (**Regular or Tall**) using chart below

A pair of appropriately sized liners will be provided free of charge with each JOBST FarrowWrap legpiece

#### Footpiece

1. Measure **A1** circumference
2. Measure length **X** from heel to base of little toe
3. Determine size (**XS-XL**) & length (**Regular or Long**) using chart below

#### JOBST FarrowHybrid (Liners)


1. JOBST FarrowHybrid (Liners) can be used instead of a footpiece, providing calf circumference does not exceed 60cm
2. Measure **A1** circumference
3. Record **shoe size**
4. Determine size (**S-XL**), width (**Regular or Wide**) & colour (**Tan or Black**)

#### Check out our measuring guide videos

Scan the QR code for help with measuring for JOBST FarrowWrap



### Footpiece



	X Small	Small	Medium	Large	X Large
<b>Legpiece</b>					
<b>C</b> Mid-Calf	36-43cm	42-50cm	48-58cm	53-63cm	58-68cm
<b>B</b> Ankle	21-25cm	25-30cm	30-36cm	36-42cm	42-50cm
<b>A-D</b> Regular	33-37cm	35-39cm	37-41cm	39-43cm	39-43cm
<b>A-D</b> Tall	38-41cm	40-43cm	42-45cm	44-47cm	44-47cm
NOTE: for <b>A-D</b> length measure posteriorly (back of leg) and follow contours					
<b>Footpiece</b>					
<b>A1</b> Mid-foot	22-24cm	25-27cm	28-30cm	31-34cm	35-40cm
<b>X</b> Regular	16-17cm	17-18cm	19-20cm	20-21cm	22-23cm
<b>X</b> Long	18-19cm	19-20cm	21-22cm	22-23cm	24-25cm
NOTE: for <b>X</b> length measure from the heel to the base of the little toe					

Note: 3cm allowance on either side of measurement

**Legpiece**

**Size:**  X Small  Small  Medium  Large  X Large

**Length **A-D****  Regular  Tall

**Fabric:**  Lite  Strong  Classic

PIP Code:  -

Product Code:  -

Quantity Required:

Repeat prescription required every  months

**Footpiece**

**Size:**  X Small  Small  Medium  Large  X Large

**Length **X****  Regular  Long

**Fabric:**  Lite  Strong  Classic

PIP Code:  -

Product Code:  -

Quantity Required:

Repeat prescription required every  months

**JOBST® FarrowHybrid® (Liners)**

**Size:**  Small (2-5)  Medium (6-8)  Large (9-11)  X Large (12-13) (shoe size)

**Width **A1****  Regular (20-27cm)  Wide (28-35cm)

**Colour:**  Tan  Black (available in Wide variant only)

PIP Code:  -

Product Code:  -

Quantity Required:

Repeat prescription required every  months

# JOBST FarrowWrap Order Codes

## JOBST® FarrowWrap® Lite (20-30mmHg)

For patients with **mild to moderate** oedema caused by venous or lymphatic conditions, where reduced compression is required

Product Code	Size	Length	Colour	PIP Code
<b>Legpiece</b>				
76654-00	Extra Small	Regular	Tan	360-6233
76654-05	Extra Small	Tall	Tan	360-6282
76654-01	Small	Regular	Tan	360-6241
76654-06	Small	Tall	Tan	360-6290
76654-02	Medium	Regular	Tan	360-6258
76654-07	Medium	Tall	Tan	360-6308
76654-03	Large	Regular	Tan	360-6266
76654-08	Large	Tall	Tan	360-6316
76654-04	Extra Large	Regular	Tan	360-6274
76654-09	Extra Large	Tall	Tan	360-6324



Product Code	Size	Length	Colour	PIP Code
<b>Footpiece</b>				
76652-05	Extra Small	Regular	Tan	360-6332
76652-00	Extra Small	Long	Tan	360-6381
76652-06	Small	Regular	Tan	360-6340
76652-01	Small	Long	Tan	360-6399
76652-07	Medium	Regular	Tan	360-6357
76652-02	Medium	Long	Tan	360-6407
76652-08	Large	Regular	Tan	360-6365
76652-03	Large	Long	Tan	360-6415
76652-09	Extra Large	Regular	Tan	360-6373
76652-04	Extra Large	Long	Tan	360-6423



## JOBST® FarrowWrap® Strong (30-40mmHg)

For patients with **moderate to severe** oedema caused by venous or lymphatic conditions, suitable for skin folds and shape distortion

Product Code	Size	Length	Colour	PIP Code
<b>Legpiece</b>				
76648-00	Extra Small	Regular	Tan	382-0073
76648-05	Extra Small	Tall	Tan	382-0024
76648-01	Small	Regular	Tan	382-0107
76648-06	Small	Tall	Tan	382-0057
76648-02	Medium	Regular	Tan	382-0099
76648-07	Medium	Tall	Tan	382-0040
76648-03	Large	Regular	Tan	382-0081
76648-08	Large	Tall	Tan	382-0032
76648-04	Extra Large	Regular	Tan	382-0065
76648-09	Extra Large	Tall	Tan	382-0016



Product Code	Size	Length	Colour	PIP Code
<b>Footpiece</b>				
76646-04	Extra Small	Regular	Tan	381-9976
76646-09	Extra Small	Long	Tan	381-9927
76646-05	Small	Regular	Tan	382-0008
76646-00	Small	Long	Tan	381-9950
76646-06	Medium	Regular	Tan	381-9992
76646-01	Medium	Long	Tan	381-9943
76646-07	Large	Regular	Tan	381-9984
76646-02	Large	Long	Tan	381-9935
76646-08	Extra Large	Regular	Tan	381-9968
76646-03	Extra Large	Long	Tan	381-9919



## JOBST® FarrowWrap® Classic (30-40mmHg)

For patients with **moderate to severe** oedema caused by venous or lymphatic conditions, suitable for irregular shaped limbs and deep skin folds

Product Code	Size	Length	Colour	PIP Code
<b>Legpiece</b>				
76657-00	Extra Small	Regular	Tan	360-6589
76657-05	Extra Small	Tall	Tan	360-6639
76657-01	Small	Regular	Tan	360-6597
76657-06	Small	Tall	Tan	360-6647
76657-02	Medium	Regular	Tan	360-6605
76657-07	Medium	Tall	Tan	360-6654
76657-03	Large	Regular	Tan	360-6613
76657-08	Large	Tall	Tan	360-6662
76657-04	Extra Large	Regular	Tan	360-6621
76657-09	Extra Large	Tall	Tan	360-6670

Product Code	Size	Length	Colour	PIP Code
<b>Footpiece</b>				
76656-05	Extra Small	Regular	Tan	360-6688
76656-00	Extra Small	Long	Tan	360-6738
76656-06	Small	Regular	Tan	360-6696
76656-01	Small	Long	Tan	360-6746
76656-07	Medium	Regular	Tan	360-6704
76656-02	Medium	Long	Tan	360-6753
76656-08	Large	Regular	Tan	360-6712
76656-03	Large	Long	Tan	360-6761
76656-09	Extra Large	Regular	Tan	360-6720
76656-04	Extra Large	Long	Tan	360-6779

## JOBST® FarrowHybrid® (20-30mmHg)

For the treatment of **mild to moderate** oedema in the foot and ankle

Product Code	Size	Colour	PIP Code
76664-00	Small - Regular (shoe size 2-5)	Tan	408-1626
76664-01	Medium - Regular (shoe size 6-8)	Tan	408-1634
76664-02	Large - Regular (shoe size 9-11)	Tan	408-1642
76664-03	Extra Large - Regular (shoe size 12-13)	Tan	408-1659
76664-12	Small - Wide (shoe size 2-5)	Tan	408-1667
76664-13	Medium - Wide (shoe size 6-8)	Tan	408-1675
76664-14	Large - Wide (shoe size 9-11)	Tan	408-2384
76664-15	Extra Large - Wide (shoe size 12-13)	Tan	408-1683
76664-16	Small - Wide (shoe size 2-5)	Black	408-6765
76664-17	Medium - Wide (shoe size 6-8)	Black	408-6757
76664-18	Large - Wide (shoe size 9-11)	Black	408-6740
76664-19	Extra Large - Wide (shoe size 12-13)	Black	408-6732

Do not apply JOBST FarrowHybrid (Liners) if calf circumference exceeds 60cm, consider a JOBST FarrowWrap footpiece in this circumstance



### Find ordering and prescribing information fast with JOBST Code Finder

No more searching through a catalogue, use our web app on your mobile device to find the information you need.

Scan the QR code or search for **codefinder.jobst.com** on your mobile device.



### Fast and free home delivery of prescription garments with JOBST Delivered\*

To use this dispensing service, send this form and prescription to: **JOBST Delivered, FREEPOST, 10 Cossall Industrial Estate, Ilkeston, Derby, DE7 5UG** or order online at [order.jobst.com/UK](http://order.jobst.com/UK) and select JOBST Delivered as the dispenser (note: prescription to be posted to JOBST Delivered unless using EPS).

If the patient is using an Electronic Prescription Service (EPS), select code **FJ708** to nominate Daylong Direct as their dispenser of choice.

For envelopes, please contact: **jobstdelivered@daylongdirect.co.uk** or **0808 196 6250**

\*JOBST Delivered is operated by Daylong Direct, 10 Cossall Industrial Estate, Ilkeston, Derby, DE7 5UG  
Privacy Policy: [www.daylong.co.uk](http://www.daylong.co.uk)

### Go digital and order using JOBST Online

Bringing you an easier way to order JOBST. Order direct with payment by account or on prescription by selecting JOBST Delivered\* as the dispenser on JOBST Online.

Scan QR code to register or visit [order.jobst.com/UK](http://order.jobst.com/UK)



\*JOBST Delivered is operated by Daylong Direct, 10 Cossall Industrial Estate, Ilkeston, Derby, DE7 5UG

Privacy Policy: [www.daylong.co.uk](http://www.daylong.co.uk)



JOBST®, an Essity brand



Essity  
T/A BSN medical Limited  
PO Box 258 • Willerby • Hull • HU10 6WT

www.jobst.co.uk  
Tel: 01482 670100 • Fax: 01482 670111  
E-mail: [compression.uk@jobst.com](mailto:compression.uk@jobst.com)



By completing this order form, you are confirming that you are aware of your obligation to obtain consent from the patient on the processing of their data for the production of their JOBST® compression garment in accordance with the Data Protection Act 2018, the UK's implementation of the (EU) General Data Protection Regulation [GDPR]. For more information and our Privacy Policy, visit [www.jobst.co.uk](http://www.jobst.co.uk).