

## Custom-fit measuring form

Email order to Customer Services on: [compression.uk@jobst.com](mailto:compression.uk@jobst.com)  
 Fax order to Customer Services on: **0345 122 3450**

Date: \_\_\_\_\_ Purchase Order No.: \_\_\_\_\_ Patient Name: \_\_\_\_\_ DoB: \_\_\_\_\_  
 Measured By: \_\_\_\_\_ Tel: \_\_\_\_\_ Email: \_\_\_\_\_  
 Delivery Address: \_\_\_\_\_ Invoice Address: \_\_\_\_\_

### Length

**E-G**

Left: \_\_\_\_\_  
 Right: \_\_\_\_\_

**D-E**

Left: \_\_\_\_\_  
 Right: \_\_\_\_\_


**A-D**

Left: \_\_\_\_\_  
 Right: \_\_\_\_\_

**X** Straight distance

Left: \_\_\_\_\_  
 Right: \_\_\_\_\_

### Circumferences



### Measuring instructions

#### Thighpiece

1. Measure **E** (5cm above the patella), **F** (midpoint of the thigh) and **G** (2cm below gluteal fold) circumferences
2. Measure **E-G** length medially (inside leg) following the contour
3. Measure **D-E** length posteriorly (back of leg) following the contour from **D** (two finger widths below the crease of the knee) to **E**

#### Legpiece

1. Measure the following circumferences:
  - **B** = narrowest part of the ankle
  - **B<sup>1</sup>** = base of the calf/insertion of the Achilles tendon
  - **B<sup>2</sup>** = midpoint between **B<sup>1</sup>** and **C**
  - **C** = widest part of the calf
2. Measure **A-D** length posteriorly (back of leg) following the contour from **A** (floor) to **D** (two finger widths below the crease of the knee)

#### Footpiece

1. Measure **A<sup>1</sup>** circumference
2. Measure **X** length from heel to base of little toe (straight)

**Custom-fit**

**Thighpiece** (Kneepiece included)  Lite  Strong  Classic  Quantity Required:

**Legpiece**  Lite  Strong  Classic

**Footpiece**  Lite  Strong  Classic

**Remarks**

---

A pair of appropriately sized liners will be provided free of charge with each JOBST® FarrowWrap® lower limb piece.