

JOBST® forMen Ambition and JOBST® forMen Explore

Ready-to-Wear RAL Compression Hosiery

Available on
Drug Tariff











JOBST® forMen is a range of ready-to-wear, circular-knit hosiery for the treatment of mild to moderate lymphatic and venous disorders.

Ambition - Stylish and effective

Inspired by a man's dress sock the timeless ribbed design, available in a variety of colours, offers good patient compliance. The balanced mix of soft yarns and high quality cotton ensures high elasticity for ease of donning and increased comfort. For the effective management of lymphatic and venous disorders and helps prevent the recurrence of leg ulcers.

Explore - Durable and effective

Reinforced foot zones offer increased strength and durability for active wear. The high cotton content of JOBST® forMen Explore - the classic of all clothing fabrics - ensures a natural feel against the skin and reduces the potential for skin irritations. For the effective management of lymphatic and vascular disorders and helps prevent the recurrence of leg ulcers.

Compression Classes				
RAL CCL 1 = 18-21mmHg		RAL CCL 2 = 23-32mmHg		
✓		✓		
Sizes and Options / Materials				
Six sizes	Regular	Long	Closed Toe	
✓	✓	✓	✓	
JOBST® forMen Ambition - Colours				
 Black	 Dark Grey	 Brown	 Navy	 Khaki
JOBST® forMen Explore - Colours				
 Black	 Dark Grey	 Navy	 Khaki	
Style				
				
AD Knee High				

For ordering information, please refer to the Compression Therapy Catalogue, or contact our Customer Services department on 0845 122 3600



Luxury compression hosiery for men & women



His and Hers

Ready-to-wear compression garments for lymphoedema and venous disorders



JOBST®,
an Essity brand

4437/319

Essity
T/A BSN medical Limited
PO Box 258 • Willerby • Hull • HU10 6WT

www.bsnmedical.co.uk
Tel: 01482 670100 • Fax: 01482 670111
E-mail: compression.uk@essity.com



@JOBSTUK

JOBST® Opaque













Ready-to-Wear RAL Compression Hosiery

Available on
Drug Tariff

JOBST Opaque is a range of ready-to-wear, circular-knit hosiery for the treatment of mild to moderate lymphoedema and venous disorders.

Soft and comfortable

JOBST® Advanced Comfort combines high-quality yarns and a superior knit structure, to make JOBST® Opaque even softer than ever and creates unbelievable wearing comfort. The Nilit® Aquarius fibres effectively wick away moisture and the double-covered inlay yarns make JOBST® Opaque incredibly easy to don.

Compression Classes						
RAL CCL 1 = 18-21mmHg			RAL CCL 2 = 23-32mmHg			
✓			✓			
Sizes and Options / Materials						
Six sizes	Regular	Petite	Wider circ.	Made-to-measure	Open Toe	Closed Toe
✓	✓	✓	✓	✓	✓	✓
<div><div></div><div></div><div></div><div></div></div> <div>Knitted top bandLace silicone top bandDotted silicone top bandSoft silicone top band</div>						
Colours						
 Sand	 Black	 Caramel	 Bronze	 Navy		
Styles						
						
AD Knee High	AG Thigh High		AT Tights			

For ordering information, please refer to the Compression Therapy Catalogue, or contact our Customer Services department on 0845 122 3600

JOBST® UltraSheer













Ready-to-Wear RAL Compression Hosiery

Available on
Drug Tariff

JOBST® UltraSheer is a range of ready-to-wear, circular-knit hosiery designed to offer highly effective, therapeutic compression for the management of lymphatic and venous disorders.

Elegant and stylish

JOBST® UltraSheer provides outstanding sheerness and a silky look for an elegant and stylish appearance that encourages good patient compliance. The fine yarns deliver superior breathability and moisture management making the garments more comfortable to wear. With six fashionable colours for increased patient choice, JOBST® UltraSheer is suitable for any wardrobe.

Compression Classes					
RAL CCL 1 = 18-21mmHg			RAL CCL 2 = 23-32mmHg		
✓			✓		
Sizes and Options / Materials					
Six sizes	Regular	Petite	Wider circ.	Open Toe	Closed Toe
✓	✓	✓	✓	✓	✓
					
Lace silicone top band		Soft silicone top band		Decorative Dotted	
Colours					
 Natural	 Caramel	 Bronze	 Café	 Anthracite	 Black
Styles					
					
AD Knee High		AG Thigh High		Tights	

For ordering information, please refer to the Compression Therapy Catalogue, or contact our Customer Services department on 0845 122 3600

JOBST® Bellavar®








Ready-to-Wear RAL Compression Hosiery

Available on
Drug Tariff

JOBST® Bellavar® is a range of ready-to-wear, circular-knit hosiery for the treatment of moderate lymphatic and venous disorders.

Firm feel fabric

JOBST® Bellavar® is manufactured using a firmer feel fabric with stiff resistance, to apply a stronger level of compression. The reinforced toe area and robust elastic heel provides high durability with a comfortable fit.

Compression Classes				
RAL CCL 2 = 23-32mmHg			RAL CCL 3 = 34-46mmHg	
✓			✓	
Sizes and Options / Materials - Knitted or dotted silicone top band for thigh high				
Six sizes	Regular	Petite	Open Toe	Closed Toe
✓	✓	✓	✓	✓
<div><div></div><div></div></div> <div>Knitted top bandDotted silicone top band</div>				
Colours				
 Sahara			 Pearl	
Styles				
				
AD Knee High	AG Thigh High		AT Tights	

For ordering information, please refer to the Compression Therapy Catalogue, or contact our Customer Services department on 0845 122 3600