



DEFY GRAVITY






JOBST® COMPRESSION

Lower Limb Compression Garments

CONTENTS

Lower Limb Compression Garments

	EVERYDAY Ready-To-Wear	18
	JOBST® Ultrasheer	
	JOBST® forMen	
	JOBST® forMenCasual	
	JOBST® Relief	
	SITUATIONAL Ready-To-Wear	24
	JOBST® Travel	
	JOBST® Activewear	
	JOBST® Sport	
	JOBST® soSoft	
	MEDICAL Ready-To-Wear	32
	JOBST® Elvarex® Soft	
	JOBST® UlcerCARE	
	JOBST® Comprinet pro	
	Accessories & Care Guidelines	39

WHAT TYPE OF COMPRESSION DOES MY PATIENT NEED?

1 Select the required compression

CCL 1: 15-20 mmHg

- Tired aching legs
- Minor ankle, leg and foot swelling
- Minor varicose during pregnancy
- Long distance travel

CCL 2: 20-30 mmHg

- Tired aching legs
- Moderate to severe varicose during pregnancy
- Post vein surgery/post sclerotherapy
- Moderate tendency of oedema
- Superficial thrombophlebitis
- Helps prevent recurrence of ulcerations

CCL 3: 30-40 mmHg


- Severe varicose veins
- Lymphatic oedema
- Post vein surgery/post sclerotherapy
- Venous ulcers
- Manages manifestations of Post Thrombotic Syndrome (PTS), Chronic Venous Insufficiency CVI/CVD
- Aids in prevention of thrombotic syndrome

JOBST® Medical LegWear should not be worn and is contraindicated for the following conditions:

- Ischemia (eg. Advanced arterial disease)*
- Uncontrolled congestive heart failure
- Untreated septic phlebitis
- Phlegmasia cerulea dolens

*As part of a comprehensive assessment, a Doppler Ultrasound is used in the vascular screening process to clarify the vascular aetiology. A hand-held Doppler will assist in determining the degree, if any, of arterial disease. Using the Doppler the Ankle/Brachial Index (ABPI) can be calculated.

2 Select Garment style



* Open toe available

Some considerations when selecting the most appropriate style for your patient:

- Indication being treated/managed
- Lifestyle, mobility and age
- Ability to don and tolerate level of compression
- Skin Integrity
- Open toe or closed toe requirements
- Maternity Style only available in JOBST® UltraSheer

3

Measure for JOBST® Medical LegWear

Measuring positions

cH - Hip

Circumference at largest dimension around buttocks

cH

cG - Thigh

Circumference at largest dimension

cG

cC - Calf


Circumference at largest dimension

cC

cB - Ankle

Circumference at smallest dimension

cB










Provide the following measurements

Knee High - cB & cC

Thigh High - cB, cC & cG

Waist high - cB, cG & cH

NOTE Each leg is measured independently and measure in the early morning when possible as swelling often increases later in the day.

Style	Measurement	S	M	L	L Full Calf	XL	XL Full Calf
UNISEX	Ankle 	18 - 21 cm	21 - 25 cm	25 - 29 cm	25 - 30 cm	29 - 33 cm	30 - 36 cm
	Calf 	28 - 38 cm	30 - 42 cm	32 - 46 cm	46 - 61 cm	34 - 50 cm	46 - 61 cm
	Thigh 	40 - 62 cm	46 - 70 cm	54 - 78 cm		60 - 81 cm	
	Hip 	71 - 117 cm	76 - 127 cm	81 - 137 cm		102 - 166 cm	
MATERNITY	Ankle 	18 - 21 cm	21 - 25 cm	25 - 29 cm		29 - 33 cm	
	Thigh 	40 - 62 cm	46 - 70 cm	54 - 78 cm		60 - 81 cm	
	Hip 	81 - 132 cm	86 - 142 cm	91 - 152 cm		112 - 180 cm	

UltraSheer: Petite Knee High Length 34 - 39 cm, Regular Knee High Length 39 - 45 cm

JOBST® ULTRASHEER

Circular Knit


Ready-To-Wear (RTW)

- Graduated compression helps improve blood flow and reduce swelling
- Medical compression stockings to help manage specific venous and lymphatic conditions
- Fine high quality yarns with an open mesh like knitting structure allows air to pass for greater breathability
- Excellent moisture wicking for better patient compliance
- Shimmering appearance for a glossy, stylish look
- Silicone top bands that grip in style
- Maternity style: Elastic SoftFlex Zone adapts to the growing belly



Sizing & Compression Classes

- CCL 1: 15-20 mmHg**
- CCL 2: 20-30 mmHg**
- CCL 3: 30-40 mmHg**



- Sheerness for beautiful legs
- Allows air to pass
- For a glossy look

Name and Style	S		M		L		XL		Colour
	Open Toe	Closed Toe	Open Toe	Closed Toe	Open Toe	Closed Toe	Open Toe	Closed Toe	
UltraSheer Stockings Knee High	75165-30	75165-25	75166-30	75166-25	75167-30	75167-25	75168-30	75168-25	● Natural Regular
		75408-02		75409-02		75410-02		75415-02	● Natural Petite
		75165-04		75166-04		75167-04		75168-04	● Black
		73633-00		73633-01		73633-02		73633-03	● Diamond Pattern Black
		75201-19		75202-19		75203-19		75204-19	● Natural Regular
		75435-02		75436-02		75437-02		75438-02	● Natural Petite
		75201-04		75202-04		75203-04		75204-04	● Black
		73634-00		73634-01		73634-02		73634-03	● Diamond Pattern Black
		75241-21		75242-21		75243-21		75244-21	● Natural
		75241-04		75242-04		75243-04		75244-04	● Black
UltraSheer Stockings Thigh High		75169-22		75170-22		75171-22		75172-22	● Natural
		75169-04		75170-04		75171-04		75172-04	● Black
		73635-00		73635-01		73635-02		73635-03	● Diamond Pattern Black
		75205-23		75206-23		75207-23		75208-23	● Natural
		75205-04		75206-04		75207-04		75208-04	● Black
		73636-00		73636-01		73636-02		73636-03	● Diamond Pattern Black
		75245-20		75246-20		75247-20		75248-20	● Natural
	75245-04		75246-04		75247-04		75248-04	● Black	
UltraSheer SENSITIVE Stockings Thigh High		77686-12		77686-13		77686-14		77686-15	● Natural
		77686-04		77686-05		77686-06		77686-07	● Black
		77687-12		77687-13		77687-14		77687-15	● Natural
		77687-04		77687-05		77687-06		77687-07	● Black
		77688-12		77688-13		77688-14		77688-15	● Natural
		77688-04		77688-05		77688-06		77688-07	● Black
UltraSheer Stockings Waist High		75173-12		75174-12		75175-12		75176-12	● Natural
		75173-04		75174-04		75175-04		75176-04	● Black
		73637-00		73637-01		73637-02		73637-03	● Diamond Pattern Black
		75209-13		75210-13		75211-13		75212-13	● Natural
		75209-04		75210-04		75211-04		75212-04	● Black
		73638-00		73638-01		73638-02		73638-03	● Diamond Pattern Black
		75249-13		75250-13		75251-13		75252-13	● Natural
		75249-04		75250-04		75251-04		75252-04	● Black
UltraSheer Stockings Maternity		75177-04		75178-04		75179-04		75180-04	● Natural
		75177-01		75178-01		75179-01		75180-01	● Black
		75213-06		75214-06		75215-06			● Natural
		75213-03		75214-03		75215-03			● Black

JOBST® FORMEN

Circular knit

Ready-To-Wear (RTW)

- Graduated compression helps improve blood flow and reduce swelling
- A unique knit structure
- Ribbed design offers a fashionable men's look dress sock
- JOBST® forMen offers therapy and style with unparalleled comfort and durability




Sizing & Compression Classes

CCL 1: 15-20 mmHg

CCL 2: 20-30 mmHg

CCL 3: 30-40 mmHg



Durable

For long-term use



Odor Control

Helps to keep feet fresh



Moisture Management

Helps to keep feet dry

Length	S	M	L	XL	Colour
	Closed Toe	Closed Toe	Closed Toe	Closed Toe	
Knee High	75253-01	75254-01	75255-01	75256-01	● Black
Thigh High	75285-00	75286-00	75287-00	75288-00	● Black
Knee High	75257-01	75258-01	75259-01	75260-01	● Black
Thigh High	75261-00	75262-00	75263-00	75264-00	● Black
Knee High	75265-01	75266-01	75267-01	75268-01	● Black
Thigh High	75269-00	75270-00	75271-00	75272-00	● Black

JOBST® FORMEN CASUAL

Circular knit


Ready-To-Wear (RTW)

- Graduated compression helps improve blood flow and reduce swelling
- A soft, contemporary, relaxed design offering graduated compression
- Dri-release technology delivers superior moisture wicking technology to keep feet dry and cool
- With its reinforced heel and strong yarns, JOBST® forMen Casual Socks are able to resist daily wear and tear
- Contains a wool blend



Sizing & Compression Classes

- CCL 1: 15-20 mmHg
- CCL 2: 20-30 mmHg
- CCL 3: 30-40 mmHg



Helps to keep feet dry



Helps to keep feet fresh



Wool fibres promote softness


S	M	L	XL	Colour
Closed Toe	Closed Toe	Closed Toe	Closed Toe	
75476-02	75477-02	75478-02	75479-02	● Black
75476-03	75477-03	75478-03	75479-03	● Khaki Beige
75484-02	75485-02	75486-02	75487-02	● Black
75484-03	75485-03	75486-03	75487-03	● Khaki Beige
75492-02	75493-02	75494-02	75495-02	● Black
75492-03	75493-03	75494-03	75495-03	● Khaki Beige

JOBST® RELIEF

Circular knit


Ready-To-Wear (RTW)

- Graduated compression helps improve blood flow and reduce swelling
- Universal graduated compression in a classical stocking appearance
- Reinforced areas in the foot make Relief stockings especially durable
- Circular knit and seamless providing comfort and support
- Available in open and closed toe for choice
- A wide top band promotes all day comfort



Sizing & Compression Classes

- CCL 1: 15-20 mmHg
- CCL 2: 20-30 mmHg
- CCL 3: 30-40 mmHg



S	S	M	M	L	L	L	L	XL	XL	XL	XL	Colour
Open Toe	Closed Toe	Open Toe	Closed Toe	Open Toe	Closed Toe	Open Toe Full Calf	Closed Toe Full Calf	Open Toe	Closed Toe	Open Toe Full Calf	Closed Toe Full Calf	
Relief Garments Knee High												
78649-00	75339-12	78649-01	75340-12	78649-02	75341-12	78649-04	75343-12	78649-03	75342-12	78649-05	75344-12	● Black
75339-10	75339-11	75340-10	75340-11	75341-10	75341-11	75343-10	75343-11	75342-10	75342-11	75344-10	75344-11	● Beige
78052-20	78052-32	78052-21	78052-33	78052-22	78052-34	78052-24	78052-36	78052-23	78052-35	78052-25	78052-37	● Black
78052-26	78052-38	78052-27	78052-39	78052-28	78052-40	78052-30	78052-42	78052-29	78052-41	78052-31	78052-43	● Black Petite
78050-04	78050-15	78051-04	78051-15	78052-04	78052-15	77975-15	77975-16	78053-04	78053-15	77976-15	77976-16	● Beige
78054-04	78054-15	78055-04	78055-15	78056-04	78056-15	77977-15	77977-16	78057-04	78057-15	77978-15	77978-16	● Beige
Relief Garments Thigh High												
78649-20	75345-02	78649-21	75346-02	78649-22	75347-02			78649-23	75348-02			● Black
78650-20	78042-07	78650-21	78043-07	78650-22	78044-07			78650-23	78045-07			● Black
78042-08	78042-09	78043-08	78043-09	78044-08	78044-09			78045-08	78045-09			● Beige
78059-08	78059-09	78047-08	78047-09	78048-08	78048-09			78049-08	78049-09			● Beige
Relief Garments Waist High												
78476-10	78476-11	78477-10	78477-11	78478-10	78478-11			78479-10	78479-11			● Beige
78480-10		78481-10		78482-10				78483-10				● Beige
Relief Garments Chaps												
78031-05		78032-05		78033-05				78034-05				● Beige L
78031-06		78032-06		78033-06				78034-06				● Beige R
78035-05		78036-05		78037-05				78041-05				● Beige L
78035-06		78036-06		78037-06				78041-06				● Beige R

JOBST® TRAVEL

Circular knit


Ready-To-Wear (RTW)

- Graduated compression helps improve blood flow and reduce swelling
- Promotes leg health during long hours of inactivity
- Helps reduce swelling and aching legs, ankles and feet due to long periods of sitting or standing
- 15-20mmHg compression helps reduce the risk of deep vein thrombosis (Blood Clot)
- Improves circulation and energises legs



Sizing & Compression Classes

CCL 1: 15-20 mmHg



- Improves blood circulation
Helps Minimize Swelling
- Promotes leg health
Travel
- Allows air to pass
Breathable

Size 1	Size 2	Size 3	Size 4	Size 5	Colour
Closed Toe	Closed Toe	Closed Toe	Closed Toe	Closed Toe	
78844-07	78845-07	78846-07	78847-07	79272-07	● Black
78844-12	78845-12	78846-12	78847-12	79272-12	● Beige

JOBST® Travel Socks measurements table

Size	1	2	3	4	5
Ankle Circumference	18 - 21 cm	20 - 23 cm	22 - 25cm	24 - 27 cm	26 - 29 cm
Calf Circumference	27 - 37 cm	30 - 40 cm	33 - 43 cm	36 - 46 cm	39 - 49 cm
US Women's Shoe*	4.5 - 6.5	7 - 9	9.5 - 11	11.5 - 13.5	-
US Men's Shoe*	3.5 - 5.5	6 - 8	8.5 - 10	10.5 - 12	12.5 - 14
EU Shoe	37 - 38	39 - 40	41 - 42	43 - 44	45 - 46
UK Shoe	3 - 5	5 - 7	7 - 9	9 - 11	11 - 13

*Australian shoe sizing is the same as US women's and men's sizing



JOBST® ACTIVEWEAR

Circular knit


Ready-To-Wear (RTW)

- Graduated compression helps improve blood flow and reduce swelling
- Soft and latex free
- Unisex
- Easy to apply
- Cushioned foot and seamless toe help protect feet from shoe seams
- Dri-release technology: providing a dry and cooling effect



Sizing & Compression Classes

- CCL 1: 15-20 mmHg
- CCL 2: 20-30 mmHg
- CCL 3: 30-40 mmHg



-  Cushioning
Cushioned foot part provides comfort
-  Moisture Management
Helps to keep feet dry
-  Soft Touch
Blend of soft yarns

S	M	L	XL	Colour
Closed Toe	Closed Toe	Closed Toe	Closed Toe	
75147-03	75148-03	75149-03	75150-03	● Black
75147-02	75148-02	75149-02	75150-02	○ White
75155-03	75156-03	75157-03	75158-03	● Black
75155-02	75156-02	75157-02	75158-02	○ White
75411-03	75412-03	75413-03	75414-03	● Black
75411-02	75412-02	75413-02	75414-02	○ White



JOBST® SPORT

Circular knit


Ready-To-Wear (RTW)

- Graduated compression helps improve blood flow and reduce swelling
- Soft fabric with a breathable knit provides a high level of comfort for sporting activities
- JOBST® Moisture management system for moisture control of the feet
- Anatomically cushioned heel and sole support
- Activated carbon fibers for effective odour control
- Extra flat toe seam for extended wearing comfort



Sizing & Compression Classes

CCL 1: 15-20 mmHg



Support zone for tendons



Helps to keep feet fresh



Helps keep feet dry

S	M	L	XL	Colour
75289-70	75289-71	75289-72	75289-73	Black
75289-80	75289-81	75289-82	75289-83	Royal Blue
75289-90	75289-91	75289-92	75289-93	Pink



JOBST® soSoft

Circular knit


Ready-To-Wear (RTW)

- JOBST® soSoft compression socks provide graduated compression to help support circulation
- With gentle compression of 8-15mmHg, JOBST® soSoft compression socks can help tired and achy legs
- Suitable for longer periods of sitting and standing
- Provide relief from minor oedema



Sizing & Compression Classes

**8-15 mmHg
Supportive Compression**



S	M	L	XL	Colour
73679-11	73679-21	73679-31	73673-05	● Black
73680-12	73680-22	73680-32	73673-10	● Sand

Sizing	Ankle cB	Calf cC	
	S	18 - 21 cm	28 - 38 cm
	M	20 - 24 cm	30 - 42 cm
	L	22 - 26 cm	32 - 46 cm
	XL	24 - 28 cm	34 - 50 cm

US Sizing Women		
	S	4.5 - 6.5
	M	7 - 9
	L	9.5 - 11
	XL	11.5 - 13

Shoe size alone is not an indicator for correct sizing. Measurements of calf and ankle circumference are mandatory for the proper performance of this product



JOBST® ELVAREX SOFT

Ready to Wear

Flat knit compression garments



- JOBST® Elvarex® Soft Knee High garments are the perfect alternative for patients where management with circular knitted garments are no longer suitable due to uncontrolled oedema
- Indicated for the management of mild Lymphoedema and vascular insufficiency
- Soft and Skin Friendly: The soft knitted yarn combines skin protecting softness with excellent breathability, making it an ideal solution for those with sensitive and fragile skin
- Therapeutically effective: Combining the right level of compression with the right fabric stiffness, JOBST® Elvarex Soft offers an extremely comfortable and effective treatment for low to moderate oedema



Sizing & Compression Classes

CCL 2: 20-30 mmHg

CCL 3: 30-40 mmHg



JOBST® Elvarex Soft Regular Calf

JOBST® Elvarex Soft 20-30mmHg Silicone Band - short 2.5 cm, long 5 cm						
Size	I	II	III	IV	V	VI
Regular (ID AD= 40-45)	79559-012	79559-013	79559-014	79559-015	79559-016	79559-017
Petite (ID AD= 34-39)	79558-012	79558-013	79558-014	79558-015	79558-016	79558-017
JOBST® Elvarex Soft 30-40mmHg Silicone Band - short 2.5 cm, long 5 cm						
Size	I	II	III	IV	V	VI
Regular (ID AD= 40-45)	73552-012	73552-013	73552-014	73552-015	73552-016	73552-017
Petite (ID AD= 34-39)	73579-012	73579-013	73579-014	73579-015	73579-016	73579-017
Size	I	II	III	IV	V	VI
cD	28 - 32 cm	30 - 34 cm	32 - 36 cm	34 - 38 cm	36 - 40 cm	38 - 42 cm
cC	29 - 33 cm	31 - 35 cm	33 - 37 cm	35 - 39 cm	37 - 41 cm	39 - 43 cm
cB1	23 - 26 cm	25 - 28 cm	27 - 30 cm	29 - 32 cm	31 - 34 cm	33 - 36 cm
cB	18 - 20 cm	20 - 22 cm	22 - 24 cm	24 - 26 cm	26 - 28 cm	28 - 30 cm
cY	26 - 31 cm	28 - 33 cm	39 - 35 cm	31 - 37 cm	32 - 38 cm	33 - 40 cm
cA	17 - 22 cm	19 - 24 cm	21 - 26 cm	23 - 29 cm	25 - 32 cm	27 - 34 cm

JOBST® Elvarex Soft PLUS Calf

JOBST® Elvarex Soft 20-30mmHg Silicone Band - short 2.5 cm, long 5 cm						
Size- PLUS	I Plus	II Plus	III Plus	IV Plus	V Plus	VI Plus
Regular (ID AD= 40-45)	79559-018	79559-019	79559-020	79559-021	79559-022	79559-023
Petite (ID AD= 34-39)	79558-018	79558-019	79558-020	79558-021	79558-022	79558-023
JOBST® Elvarex Soft 30-40mmHg Silicone Band - short 2.5 cm, long 5 cm						
Size- PLUS	I Plus	II Plus	III Plus	IV Plus	V Plus	VI Plus
Regular (ID AD= 40-45)	73552-018	73552-019	73552-020	73552-021	73552-022	73552-023
Petite (ID AD= 34-39)	73579-018	73579-019	73579-020	73579-021	73579-022	73579-023
Size chart	I Plus	II Plus	III Plus	IV Plus	V Plus	VI Plus
cD Plus	32 - 36 cm	34 - 38 cm	36 - 40 cm	38 - 42 cm	40 - 44 cm	42 - 46 cm
cC plus	33 - 37 cm	35 - 39 cm	37 - 41 cm	39 - 43 cm	41 - 45 cm	43 - 47 cm
cB1 plus	24 - 27 cm	26 - 29 cm	28 - 31 cm	30 - 33 cm	32 - 35 cm	34 - 37 cm
cB	18 - 20 cm	20 - 22 cm	22 - 24 cm	24 - 26 cm	26 - 28 cm	28 - 30 cm
cY	26 - 31 cm	28 - 33 cm	39 - 35 cm	31 - 37 cm	32 - 38 cm	33 - 40 cm
cA	17 - 22 cm	19 - 24 cm	21 - 26 cm	23 - 29 cm	25 - 32 cm	27 - 34 cm

JOBST® ULCERCARE




JOBST® UlcerCARE™ enables effective management of hard to heal wounds on the lower leg while offering ease of use for patients and care givers - the perfect combination for providing patients with an improved quality of life.

JOBST® UlcerCARE™ Compression System:

- Provides graduated compression of 40 mmHg, the recommended pressure for the healing of these types of wounds
- Is a 2-in-1 stocking system designed for patient's self-management to facilitate easy changing of the wound dressing and to allow for improved personal hygiene (e.g. bathing and skin care)
- JOBST® UlcerCARE Liner provides mild compression (17mmHg) for the non-ambulatory patient to help manage lower leg oedema during bed rest
- Outer stocking comes with a zipper on the left or right side for comfortable closure on the opposite side of the wound. It is removed at night and applied first thing in the morning

Sizing & Compression Classes

Size			JOBST® UlcerCARE Zipper Stocking and Liners (1 stocking/2 liners)		
	Ankle	Ankle	Zipper on Left Side	Zipper on Right Side	JOBST® UlcerCARE Replacement Liners (3 liners)
Small	18 - 21 cm	29 - 35 cm	73631-21	73631-61	73632-21
Medium	21 - 24 cm	34 - 41 cm	73631-22	73631-62	73632-22
Large	24 - 27 cm	40 - 47 cm	73631-23	73631-63	73632-23
X-Large	27 - 30 cm	46 - 54 cm	73631-24	73631-64	73632-24
XX-Large	30 - 33 cm	48 - 56 cm	73631-25	73631-65	73632-25
XXX-Large	33 - 36 cm	50 - 58 cm	73631-26	73631-66	73632-26
XXXX-Large	36 - 39 cm	52 - 60 cm	73631-27	73631-67	73632-27



JOBST® COMPRINET® PRO

Graduated compression anti-embolism stockings

- Indicated for thrombosis prophylaxis in **surgical patients** and patients requiring prolonged periods of immobilisation
- Comprinet® pro anti-embolism stockings can help reduce the risk of thrombosis in immobile and surgical patients
- Comprinet® pro exerts pressure on the vessels from outside. This narrows the veins and, although you are lying in bed, the blood flows at a faster rate back to the heart




Sizing & Compression Classes

The colour-coded sizing system allows rapid allocation of the right size.

The colours displayed at the inspection window correspond to the width (inner colour) and the length (outer colour) of the stocking.

Select the correct size garment based upon measurement guide on this page.



Comprinet® pro Knee High

Size	Width	Colour Code	Product Code
Small	23 - 30 cm		77026-02
Medium	30 - 38 cm		77027-02
Large	38 - 46 cm		77028-02
X-Large	46 - 53 cm		77124-05
XX-Large	53 - 66 cm		77133-03

Comprinet® pro Thigh High

Size	Width	Length	Colour Code	Product Code
Small Regular	23 - 30 cm	64 - 80 cm		46335-12
Small Long	23 - 30 cm	80v97 cm		46336-02
Medium Regular	30 - 38 cm	64v80 cm		46337-02
Medium Long	30 - 38 cm	80 - 97 cm		46338-02
Large Regular	38 - 46 cm	64 - 80 cm		46388-02
Large Long	38 - 46 cm	80 - 97 cm		46389-02
X-Large Regular	46 - 53 cm	64 - 80 cm		77125-02
X-Large Long	46 - 53 cm	80 - 97 cm		77128-02
XX-Large Regular	53 - 66 cm	64 - 80 cm		77130-01
XX-Large Long	53 - 66 cm	80 - 97 cm		77131-01

JOBST® ACCESSORIES AND CARE GUIDELINES

JOBST® DONNING AND DOFFING SOLUTIONS



Application of knee high stockings



Application of thigh high stockings



Doffing of thigh high stockings instructions

JOBST® Application



1 Put the hand inside the stocking and grip the heel.



2 Turn the stocking inside-out, whilst still holding onto the heel.



3 Pull the inside-out of the stocking over the front of the foot as far as the heel.



4 Gradually pull the stocking onto the calf, like a normal stocking.



5 Push the overlapping part of the stocking over the heel.



6 Fold over the stocking and pull without using force.



7 Make sure that the heel of the stocking is free of creases.



8 Pull the stocking up to the knee, thigh or waist; let it slide over the thumbs. Knee high socks or stockings should finish 2 finger width's below the knee crease at the back of the leg.

JOBST® ACCESSORIES



JOBST® Stocking Donner

Product Code: 79668-00

Due to the nature of compression stockings, people may find putting them on to be initially a challenging task. The **JOBST® Stocking Donner**, which features a **lightweight frame** and foam grip handles, is designed to aid in the application of all JOBST® Medical LegWear.



Place your JOBST® compression stocking inside the semi-circle frame and pull the top of the stocking down over the semi-circle. The heel of the stocking must face the back of your JOBST® Stocking Donner.

Continue to slide the stocking down over the semi-circle until the heel protrudes from the frame, and is centred and even.



Preferably whilst seated, grasp the padded handles and gently pull the stocking donner upwards whilst moving your leg slightly outwards, until the stocking is above the calf.



Whilst using a support for balance, insert your foot into the stocking until your foot is completely on the floor.

Helpful Tip: Sit down and hold the heel pocket against the frame as you slide your toes into the stocking.



Push the JOBST® Stocking Donner backward and down to free the stocking.

Be sure the heel of the stocking is properly placed on the foot. Then adjust the length and smooth out any wrinkles by stroking with the palm of your hands.

The knee band in Knee High socks or stockings should stop two finger widths below the crease at the back of the knee.



ARION DONNING AND DOFFING AIDS

Slide solutions for medical compression garments

Product benefits

Easy to Use



- ✓ **Minimise friction** between compression garment and skin
- ✓ **Reduce physical strain** for both compression wearers and caregivers
- ✓ Promote **independence** of compression wearers

Safe to Apply
























- ✓ **Reduce risk** of skin damage; can be especially beneficial in case of high skin sensitivity or leg ulcers
- ✓ Contribute to a **longer garment lifetime** by preventing damage from fingers and nails
- ✓ Help **maintaining the original compression longer** and avoiding garments from overstretching and distortion

Simple Solutions



- ✓ Offer application aids for all types of compression garments: stockings, tights, arm sleeves, and gloves
- ✓ Designed for different shapes and sizes of limbs
- ✓ Support the concordance to compression therapy

Arion Application Aids for your medical compression garments

	Donning and Doffing Aids for medical compression stockings and tights		Donning aid for arm sleeves																																							
	 <p>MD CE</p>	 <p>MD CE</p>	 <p>MD CE</p>																																							
Open toe			—																																							
Closed toe			—																																							
Donning (On)/ Doffing (Off)	 <p>ON</p>																																									
Sizes	<table border="1"> <thead> <tr> <th>Code</th> <th>Size</th> <th>Shoe size EU</th> </tr> </thead> <tbody> <tr> <td>79658-03</td> <td>Size XS</td> <td>≤ 33</td> </tr> <tr> <td>79658-04</td> <td>Size S</td> <td>34 - 36</td> </tr> <tr> <td>79658-02</td> <td>Size M</td> <td>37 - 38,5</td> </tr> <tr> <td>79659-02</td> <td>Size L</td> <td>39 - 41</td> </tr> <tr> <td>79660-01</td> <td>Size XL</td> <td>≥ 42</td> </tr> </tbody> </table>	Code	Size	Shoe size EU	79658-03	Size XS	≤ 33	79658-04	Size S	34 - 36	79658-02	Size M	37 - 38,5	79659-02	Size L	39 - 41	79660-01	Size XL	≥ 42	<table border="1"> <thead> <tr> <th>Code</th> <th>Size</th> <th>Shoe size EU</th> </tr> </thead> <tbody> <tr> <td>79571-00</td> <td>Size M</td> <td>≤ 38</td> </tr> <tr> <td>79571-01</td> <td>Size L</td> <td>39 - 43</td> </tr> <tr> <td>79571-02</td> <td>Size XL</td> <td>≥ 44</td> </tr> </tbody> </table> <p>XL also suitable for extra wide feet (circumference up to 28 cm) </p>	Code	Size	Shoe size EU	79571-00	Size M	≤ 38	79571-01	Size L	39 - 43	79571-02	Size XL	≥ 44	<table border="1"> <thead> <tr> <th>Code</th> <th>Size</th> <th>Arm Length*</th> </tr> </thead> <tbody> <tr> <td>79661-02</td> <td>Size M</td> <td>59 - 67 cm</td> </tr> <tr> <td>75100-01</td> <td>Size L</td> <td>> 67 cm</td> </tr> </tbody> </table> <p>*Arm length measured from the tip of the middle finger to the armpit.</p>	Code	Size	Arm Length*	79661-02	Size M	59 - 67 cm	75100-01	Size L	> 67 cm
	Code	Size	Shoe size EU																																							
	79658-03	Size XS	≤ 33																																							
	79658-04	Size S	34 - 36																																							
	79658-02	Size M	37 - 38,5																																							
79659-02	Size L	39 - 41																																								
79660-01	Size XL	≥ 42																																								
Code	Size	Shoe size EU																																								
79571-00	Size M	≤ 38																																								
79571-01	Size L	39 - 43																																								
79571-02	Size XL	≥ 44																																								
Code	Size	Arm Length*																																								
79661-02	Size M	59 - 67 cm																																								
75100-01	Size L	> 67 cm																																								
Video instructions: scan the QR code to watch the video	 	 	 																																							
Care instructions	 <p>Wash by hand (30°C)</p>  <p>Do not Bleach</p>	 <p>Tumble dry</p>  <p>Dry with a towel</p>	 <p>Do not Iron</p>  <p>Do not dry clean</p>																																							



JOBST® DONNING GLOVES

JOBST® Donning Gloves help create friction and provide a better grip when you put your garments on. They support an even distribution of the fabric of your compression garments on the leg. Wearing donning gloves also reduces the risk of damage to the fabric.

Extra Tip

Besides donning aids and accessories, there are helpful tips and tricks for easier application. Make sure that your skin is dry before putting on your garments. Avoid using any greasy ointments or body lotions prior to applying compression stockings.



DO NOT:



Tumble dry



Bleach



Iron



Use softeners



Dry clean

XS	S	M	L	XL	Colour
77400-27	77400-28	77400-29	77400-30	77400-31	● Blue

JOBST® ROLL-ON-ADHESIVE

An adhesive lotion to help hold garments in place

Suitable for daily use.

A single application with the roll-on-ball on the exposed skin around the arm or leg is sufficient. Then the edge of the garment is folded up and pressed against the skin.

For garments with a silicone edge, the adhesive lotion should not be applied directly under it.

Cleaning instructions: Remnants of the adhesive lotion can be easily washed off the skin or garment with water after removing the garment.

Product code: 79274-11





JOBST®
GENERAL
CARE
INSTRUCTIONS

Washing Instructions

- Hand or gentle machine warm wash at 40 degrees celsius (Laundry bag recommended)
- Use mild liquid detergent
- Do not use powder detergent as it clogs elastic fibres
- Rinse well to remove all traces of detergent
- Do not bleach

For garments with silicone band:

The use of skin care products (such as skin creams, body lotion or powder) as well as various skin types (such as dry or peeling skin) and hair growth, may prevent the silicone band from keeping the thigh high in place. Therefore, moisturisers should be well absorbed before applying stockings.

We recommend the following steps to maximise the grip performance:

1. Turn the garment inside out before washing
2. Before the main wash, thoroughly hand wash the silicone band using a delicate fabric detergent. Concentrate on removing any lotions, hair or foreign material from the silicone.

Drying Instructions

- Dry JOBST® UltraSheer, JOBST® Travel Socks and JOBST® UlcerCARE™ liner either indoors or in the shade, laid flat
- JOBST® forMen Casual, ActiveWear, Sport, soSoft, forMen, Relief and UlcerCARE™ outer stocking can all be tumble dried at low temperature
- Do not wring



Do not bleach



Do not iron



Do not dry clean

Please refer to the product label and / or package insert for full instructions on the safe use of these products.

Always read the label and follow the instructions for use. If symptoms persist, see your doctor / healthcare professional.

