

FIBER CONTENT

Product ¹ Produit	Polyamide (Nylon) ² Polyamide (Nylon)	Elastane (Spandex) ³ Élasthanne (Spandex)	Silicone ⁴ Silicone	Polyester ⁵ Polyester	Other ⁶ Autre
Active Wear 15-20/20-30	20%	9%	0%	60%	11% Cotton ⁷
Active Wear 30-40	20%	12%	0%	58%	10% Cotton ⁷
Style SoftFit 15-20	53%	12%	2%	0%	16% Pima Cotton ⁸ + 17% MicroModal ⁹
Style SoftFit 20-30	53%	17%	2%	0%	14% Pima Cotton ⁸ + 14% MicroModal ⁹
Style SoftFit 30-40	45%	26%	2%	0%	14% Pima Cotton ⁸ + 13% MicroModal ⁹
forMen 15-20 KH	72%	26%	2%	0%	–
forMen 15-20 TH	75%	20%	5%	0%	–
forMen 20-30 KH	69%	29%	2%	0%	–
forMen 20-30 TH	70%	26%	4%	0%	–
forMen 30-40 KH	67%	31%	2%	0%	–
forMen 30-40 TH	66%	30%	4%	0%	–

COMPOSITION

Product ¹ Produit	Polyamide (Nylon) ² Polyamide (Nylon)	Elastane (Spandex) ³ Élasthanne (Spandex)	Silicone ⁴ Silicone	Polyester ⁵ Polyester	Other ⁶ Autre
Opaque 15-20/20-30 KH	76%	24%	0%	0%	–
Opaque 15-20 TH	70%	20%	10%	0%	–
Opaque 15-20 THW (Sensitive)	74%	23%	3%	0%	–
Opaque 15-20 WH	70%	30%	0%	0%	–
Opaque 20-30 TH	67%	19%	14%	0%	–
Opaque 20-30 THW (Sensitive)	73%	24%	3%	0%	–
Opaque 20-30 WH	71%	29%	0%	0%	–
Opaque 30-40 KH	64%	36%	0%	0%	–
Opaque 30-40 TH	56%	34%	10%	0%	–
Opaque 30-40 THW (Sensitive)	60%	37%	3%	0%	–
Opaque 30-40 WH	62%	38%	0%	0%	–
Opaque SoftFit 15-20/20-30	70%	26%	4%	0%	–
Opaque SoftFit 30-40	60%	36%	4%	0%	–
Relief 15-20 KH	80%	20%	0%	0%	–
Relief 15-20 TH / w / Band	73%	19%	8%	0%	–
Relief 20-30 KH CT / OT	73%	27%	0%	0%	–
Relief 20-30 KH w / Band	70%	24%	6%	0%	–
Relief 20-30 TH CT / OT	73%	27%	0%	0%	–
Relief 20-30 TH w / Band Relief 20-30 TH w / Band Petite	69%	24%	7%	0%	–
Relief 20-30 WH CT / OT, CHAP OT	77%	23%	0%	0%	–

FIBER CONTENT

Product ¹ Produit	Polyamide (Nylon) ² Polyamide (Nylon)	Elastane (Spandex) ³ Elasthanne (Spandex)	Silicone ⁴ Silicone	Polyester ⁵ Polyester	Other ⁶ Autre
Relief 30-40 WH CT/OT, Chap OT	70%	30%	0%	0%	-
Relief 30-40 KH CT/OT	70%	30%	0%	0%	-
Relief 30-40 KH w/Band OT	66%	29%	5%	0%	-
Relief 30-40 TH CT	65%	35%	0%	0%	-
Relief 30-40 TH OT	66%	34%	0%	0%	-
Relief 30-40 TH w/Band	64%	30%	6%	0%	-
Sport 15-20/20-30	54%	10%	0%	36	-
UltraSheer/Pattern 15-20 KH/MAT	67%	33%	0%	0%	-
UltraSheer/Pattern 15-20 TH (Dot)	60%	25%	15%	0%	-
UltraSheer 15-20 TH (Lace)	54%	31%	15%	0%	-
UltraSheer 15-20 THW (Sensitive)	69%	27%	4%	0%	-
UltraSheer/Pattern 15-20 WH	65%	35%	0%	0%	-
UltraSheer/Pattern 20-30 KH	71%	29%	0%	0%	-
UltraSheer 20-30 MAT	73%	27%	0%	0%	-
UltraSheer/Pattern 20-30 TH (Dot)	61%	27%	12%	0%	-
UltraSheer 20-30 TH (Lace)	56%	32%	12%	0%	-
UltraSheer 20-30 THW (Sensitive)	69%	28%	3%	0%	-
UltraSheer 20-30 WH/Pattern	70%	30%	0%	0%	-

COMPOSITION

Product ¹ Produit	Polyamide (Nylon) ² Polyamide (Nylon)	Elastane (Spandex) ³ Elasthanne (Spandex)	Silicone ⁴ Silicone	Polyester ⁵ Polyester	Other ⁶ Autre
UltraSheer 30-40 KH	52%	48%	0%	0%	-
UltraSheer 30-40 TH (Dot)	48%	42%	10%	0%	-
UltraSheer 30-40 TH (Lace)	47%	45%	8%	0%	-
UltraSheer 30-40 THW (Sensitive)	56%	41%	3%	0%	-
UltraSheer 30-40 WH	51%	49%	0%	0%	-
UltraSheer SoftFit 15-20/20-30	60%	34%	6%	0%	-
UltraSheer SoftFit 30-40	50%	46%	4%	0%	-
Maternity Opaque 15-20 KH	84%	16%	0%	0%	-
Maternity Opaque 15-20 TH	75%	15%	10%	0%	-
Maternity Opaque 15-20 WH	78%	22%	0%	0%	-
Maternity Opaque 20-30 KH	81%	19%	0%	0%	-
Maternity Opaque 20-30 TH	73%	18%	9%	0%	-
Maternity Opaque 20-30 WH	76%	24%	0%	0%	-