


# JOBST® Bella™ Lite

## Upper Limb

Ready-to-wear order form and prescription request

Email order to Customer Services on: [compression.uk@jobst.com](mailto:compression.uk@jobst.com)  
Please Note: PDF is the preferred file format.

Date: \_\_\_\_\_ Purchase Order No.: \_\_\_\_\_ Patient Name: \_\_\_\_\_  
 Measured By: \_\_\_\_\_ Tel: \_\_\_\_\_ Email: \_\_\_\_\_  
 Delivery Address: \_\_\_\_\_  
 Post Code: \_\_\_\_\_



**Length**

Left Right

**C-G**

--	--

**Note:**  
Length measurement taken on inside of the arm from wrist (C) to axilla (G)

**Circumferences**

Left Right

Mid-Upper

--	--

Mid-Lower

--	--

Wrist

--	--

Palm

--	--

**JOBST Bella Lite**  
Armsleeve

**Size**

Small  
 Medium  
 Large

**Compression Class**

CCL 1 (15-20mmHg)  
 CCL 2 (20-30mmHg)

**Length**

Regular  
 Long

**Top Band Option**

Knitted  
 Silicone

**Colour**

Beige  
 Black (CCL2 knitted band, long only)

PIP Code: [ ] [ ] [ ] - [ ] [ ] [ ] [ ]  
 Product Code: [ ] [ ] [ ] [ ] [ ] - [ ] [ ]  
 Quantity Required: [ ] [ ] [ ]  
 Repeat prescription required every [ ] months

**JOBST Bella Lite**  
Combined Armsleeve

**Size**

Small  
 Medium  
 Large

**Compression Class**

CCL 1 (15-20mmHg)  
 CCL 2 (20-30mmHg)

**Length**

Regular  
 Long

**Top Band Option**

Knitted  
 Silicone

**Colour**

Beige

PIP Code: [ ] [ ] [ ] - [ ] [ ] [ ] [ ]  
 Product Code: [ ] [ ] [ ] [ ] [ ] - [ ] [ ]  
 Quantity Required: [ ] [ ] [ ]  
 Repeat prescription required every [ ] months

**JOBST Bella Lite**  
Gauntlet

**Size**

Small  
 Medium  
 Large

**Compression Class**

CCL 1 (15-20mmHg)  
 CCL 2 (20-30mmHg)

**Colour**

Beige  
 Black (CCL2 only)

PIP Code: [ ] [ ] [ ] - [ ] [ ] [ ] [ ]  
 Product Code: [ ] [ ] [ ] [ ] [ ] - [ ] [ ]  
 Quantity Required: [ ] [ ] [ ]  
 Repeat prescription required every [ ] months

**Remarks**

		Armsleeve			Combined Armsleeve			Gauntlet	
		C Wrist	D Mid-Lower	F Mid-Upper	B Palm	D Mid-Lower	F Mid-Upper	B Palm	C Wrist
Circumference	Small	14cm - 20cm	18cm - 23cm	22cm - 40cm	15cm - 19cm	18cm - 23cm	22cm - 40cm	15cm - 19cm	14cm - 20cm
	Medium	16cm - 22cm	23cm - 29cm	27cm - 45cm	19cm - 22cm	23cm - 29cm	27cm - 45cm	19cm - 22cm	16cm - 22cm
	Large	18cm - 24cm	29cm - 34cm	33cm - 51cm	22cm - 26cm	29cm - 34cm	33cm - 51cm	22cm - 26cm	18cm - 24cm
Length	Regular	up to (≤) 43cm						N/A	
	Long	over (>) 43 cm							

# JOBST® Bella™ Lite



## Armsleeve Ordering Information

Product Code	CCL	Option	Size	Length	Colour	PIP Code
75273-02	CCL1	Knitted Band	Small	Regular	Beige	382-1196
75273-00	CCL1	Knitted Band	Small	Long	Beige	325-7011
75274-02	CCL1	Knitted Band	Medium	Regular	Beige	382-1188
75274-00	CCL1	Knitted Band	Medium	Long	Beige	325-7037
75275-02	CCL1	Knitted Band	Large	Regular	Beige	382-1170
75275-00	CCL1	Knitted Band	Large	Long	Beige	325-7029
75273-03	CCL1	Silicone Band	Small	Regular	Beige	382-1238
75273-01	CCL1	Silicone Band	Small	Long	Beige	325-6955
75274-03	CCL1	Silicone Band	Medium	Regular	Beige	382-1220
75274-01	CCL1	Silicone Band	Medium	Long	Beige	325-6971
75275-03	CCL1	Silicone Band	Large	Regular	Beige	382-1212
75275-01	CCL1	Silicone Band	Large	Long	Beige	325-6963
75276-05	CCL2	Knitted Band	Small	Regular	Beige	382-1261
75276-00	CCL2	Knitted Band	Small	Long	Beige	325-6989
75277-05	CCL2	Knitted Band	Medium	Regular	Beige	382-1204
75277-00	CCL2	Knitted Band	Medium	Long	Beige	325-7003
75278-05	CCL2	Knitted Band	Large	Regular	Beige	382-1246
75278-00	CCL2	Knitted Band	Large	Long	Beige	325-6997
75276-06	CCL2	Silicone Band	Small	Regular	Beige	382-1287
75276-01	CCL2	Silicone Band	Small	Long	Beige	325-6898
75277-06	CCL2	Silicone Band	Medium	Regular	Beige	382-1279
75277-01	CCL2	Silicone Band	Medium	Long	Beige	325-6914
75278-06	CCL2	Silicone Band	Large	Regular	Beige	382-1253
75278-01	CCL2	Silicone Band	Large	Long	Beige	325-6906
75276-04	CCL2	Knitted Band	S	Long	Black	363-7568
75277-04	CCL2	Knitted Band	M	Long	Black	363-7576
75278-04	CCL2	Knitted Band	L	Long	Black	363-7550

## Combined Armsleeve Ordering Information

Product Code	CCL	Option	Size	Length	Colour	PIP Code
77678-00	CCL1	Knitted Band	Small	Regular	Beige	398-9571
77678-03	CCL1	Knitted Band	Small	Long	Beige	398-9605
77678-01	CCL1	Knitted Band	Medium	Regular	Beige	398-9589
77678-04	CCL1	Knitted Band	Medium	Long	Beige	398-9613
77678-02	CCL1	Knitted Band	Large	Regular	Beige	398-9597
77678-05	CCL1	Knitted Band	Large	Long	Beige	398-9621
77678-06	CCL1	Silicone Band	Small	Regular	Beige	398-9639
77678-09	CCL1	Silicone Band	Small	Long	Beige	398-9662
77678-07	CCL1	Silicone Band	Medium	Regular	Beige	398-9647
77678-10	CCL1	Silicone Band	Medium	Long	Beige	398-9670
77678-08	CCL1	Silicone Band	Large	Regular	Beige	398-9654
77678-11	CCL1	Silicone Band	Large	Long	Beige	398-9688
77680-00	CCL2	Knitted Band	Small	Regular	Beige	398-9696
77680-03	CCL2	Knitted Band	Small	Long	Beige	398-9720
77680-01	CCL2	Knitted Band	Medium	Regular	Beige	398-9704
77680-04	CCL2	Knitted Band	Medium	Long	Beige	398-9738
77680-02	CCL2	Knitted Band	Large	Regular	Beige	398-9712
77680-05	CCL2	Knitted Band	Large	Long	Beige	398-9746
77680-06	CCL2	Silicone Band	Small	Regular	Beige	398-9753
77680-09	CCL2	Silicone Band	Small	Long	Beige	398-9787
77680-07	CCL2	Silicone Band	Medium	Regular	Beige	398-9761
77680-10	CCL2	Silicone Band	Medium	Long	Beige	398-9795
77680-08	CCL2	Silicone Band	Large	Regular	Beige	398-9779
77680-11	CCL2	Silicone Band	Large	Long	Beige	398-9803



## Gauntlet Ordering Information

Product Code	CCL	Option	Size	Length	Colour	PIP Code
75279-00	CCL1	With Thumb	Small	-	Beige	325-6823
75280-00	CCL1	With Thumb	Medium	-	Beige	325-6831
75281-00	CCL1	With Thumb	Large	-	Beige	325-6922
75282-00	CCL2	With Thumb	Small	-	Beige	325-6773
75283-00	CCL2	With Thumb	Medium	-	Beige	325-6781
75284-00	CCL2	With Thumb	Large	-	Beige	325-6799
75282-03	CCL2	With Thumb	Small	-	Black	363-7543
75283-03	CCL2	With Thumb	Medium	-	Black	363-7535
75284-03	CCL2	With Thumb	Large	-	Black	363-7527

## Go digital and order using JOBST Online

Bringing you an easier way to order JOBST. Order direct with payment by account or on prescription by selecting JOBST Delivered\* as the dispenser on JOBST Online.

Scan QR code to register or visit [order.jobst.com/UK](https://order.jobst.com/UK)

\*JOBST Delivered is operated by Daylong Direct, Unit 5, Heritage House, Cossall Industrial Estate, Ilkeston, Derby, DE7 5UD  
Privacy Policy: [www.daylong.co.uk](https://www.daylong.co.uk)



By completing this order form, you are confirming that you are aware of your obligation to obtain consent from the patient on the processing of their data for the production of their JOBST™ compression garment in accordance with the Data Protection Act 2018, the UK's implementation of the (EU) General Data Protection Regulation (GDPR). For more information and our Privacy Policy, visit [www.jobst.co.uk](https://www.jobst.co.uk).



JOBST®,  
an Essity brand



Essity UK Ltd  
Southfields Road, Dunstable,  
Bedfordshire, LU6 3EJ

[www.jobst.co.uk](https://www.jobst.co.uk)  
Tel: 01482 670100  
E-mail: [compression.uk@jobst.com](mailto:compression.uk@jobst.com)

📍 /JOBST UK  
📱 /JOBSTUK