



JOBST® Opaque

Rajstopy

Miejsca pomiarowe kończyny



LK - długość kończyny od podstawy do 2 palców poniżej pośladków



	I	II	III	IV	V	VI
cG	43 - 57 cm	45 - 62 cm	49 - 67 cm	53 - 72 cm	56 - 77 cm	60 - 81 cm
cD	27 - 33 cm	29 - 36 cm	32 - 39 cm	34 - 42 cm	36 - 45 cm	38 - 48 cm
cC	28 - 34 cm	30 - 37 cm	33 - 40 cm	35 - 43 cm	37 - 46 cm	39 - 49 cm
cB1	23 - 27 cm	24 - 29 cm	26 - 32 cm	29 - 35 cm	31 - 37 cm	33 - 39 cm
cB	18 - 20 cm	20 - 22 cm	22 - 24 cm	24 - 26 cm	26 - 28 cm	28 - 30 cm
Maxi option	cT > 80 cm/cH > 95 cm	cT > 90 cm/cH > 105 cm		cT > 110 cm/cH > 125 cm	cT > 120 cm/cH > 135 cm	

LK Petite: 65 - 75 cm / Regular: 75 - 85 cm