

JOBST® LOWER LIMB COMPRESSION GARMENTS



DEFY GRAVITY

The JOBST logo is located in the bottom right corner. It consists of the word "JOBST" in a white, sans-serif font, with a registered trademark symbol (®) to its upper right. The logo is set against a dark blue rectangular background.

JOBST® Medical Grade, Ready-To-Wear garments with graduated compression help to improve blood flow and reduce oedema

What type of compression does my patient need?

1 Select the required compression

CCL 1: 15-20 mmHg

- Tired aching legs
- Minor ankle, leg and foot swelling
- Minor varicose veins during pregnancy
- Long distance travel

CCL 2: 20-30 mmHg

- Tired aching legs
- Moderate to severe varicose veins during pregnancy
- Post vein surgery/post sclerotherapy
- Moderate tendency of oedema
- Superficial thrombophlebitis
- Helps prevent recurrence of ulcerations

CCL 3: 30-40 mmHg


- Severe varicose veins
- Lymphatic oedema
- Post vein surgery/post sclerotherapy
- Venous ulcers
- Manages manifestations of Post Thrombotic Syndrome (PTS), Chronic Venous Insufficiency CVI/CVD
- Aids in prevention of thrombotic syndrome

JOBST® medical legwear should not be worn and is contraindicated for the following conditions:

- Ischemia (eg. Advanced arterial disease)*
- Uncontrolled congestive heart failure
- Untreated septic phlebitis
- Phlegmasia cerulea dolens

*As part of a comprehensive assessment, a Doppler Ultrasound is used in the vascular screening process to clarify the vascular aetiology. A hand-held Doppler will assist in determining the degree, if any, of arterial disease. Using the Doppler the Ankle/Brachial Index (ABPI) can be calculated.

2 Select garment style



Some considerations when selecting the most appropriate style for your patient:

- Indication being treated/managed
- Lifestyle, mobility and age
- Ability to don and tolerate level of compression
- Skin Integrity
- Open toe or closed toe requirements
- Maternity Style only available in JOBST® Maternity Opaque

3

Measure for JOBST® medical legwear

UltraSheer, Relief, forMen, ActiveWear, Sport and Maternity Opaque

Measuring positions

cH - Hip

Circumference at largest dimension around buttocks

cH

cG - Thigh

Circumference at largest dimension

cG

cC - Calf


Circumference at largest dimension

cC

cB - Ankle

Circumference at smallest dimension

cB



Provide the following measurements

Knee High - cB & cC

Thigh High - cB, cC & cG

Waist High - cB, cG & cH

Chap - cB, cC, cG & cH

NOTE: Each leg is measured independently. Measure early in the morning when possible as swelling often increases throughout the day.

Style	Measurement		S	M	L	L Full Calf	XL	XL Full Calf
UNISEX	Ankle		18 - 21 cm	21 - 25 cm	25 - 29 cm	25 - 30 cm	29 - 33 cm	30 - 36 cm
	Calf		28 - 38 cm	30 - 42 cm	32 - 46 cm	46 - 61 cm	34 - 50 cm	46 - 61 cm
	Thigh		40 - 62 cm	46 - 70 cm	54 - 78 cm		60 - 81 cm	
	Hip		71 - 117 cm	76 - 127 cm	81 - 137 cm		102 - 166 cm	
MATERNITY	Ankle		18 - 21 cm	21 - 25 cm	25 - 29 cm			
	Thigh		40 - 62 cm	46 - 70 cm	54 - 78 cm			
	Hip		81 - 132 cm	86 - 142 cm	91 - 152 cm			

Knee high length: Petite ≤ 38cm, Regular > 38cm. Thigh high length ≥ 68.5cm.

EVERYDAY Ready-To-Wear



JOBST® Relief



Long lasting fabric



Allows air to pass



Improves blood circulation

S	S	M	M	L	L	LFC	LFC	XL	XL	XLFC	XLFC	Colour
Open Toe	Closed Toe	Open Toe	Closed Toe	Open Toe	Closed Toe	Open Toe Full Calf	Closed Toe Full Calf	Open Toe	Closed Toe	Open Toe Full Calf	Closed Toe Full Calf	
Relief Garments Knee High												
78649-00	75339-12	78649-01	75340-12	78649-02	75341-12	78649-04	75343-12	78649-03	75342-12	78649-05	75344-12	● Black
75339-10	75339-11	75340-10	75340-11	75341-10	75341-11	75343-10	75343-11	75342-10	75342-11	75344-10	75344-11	● Beige
78052-20	78052-32	78052-21	78052-33	78052-22	78052-34	78052-24	78052-36	78052-23	78052-35	78052-25	78052-37	● Black
78052-26	78052-38	78052-27	78052-39	78052-28	78052-40	78052-30	78052-42	78052-29	78052-41	78052-31	78052-43	● Black Petite
78050-04	78050-15	78051-04	78051-15	78052-04	78052-15	77975-15	77975-16	78053-04	78053-15	77976-15	77976-16	● Beige
78054-04	78054-15	78055-04	78055-15	78056-04	78056-15	77977-15	77977-16	78057-04	78057-15	77978-15	77978-16	● Beige
Relief Garments Thigh High												
78649-20	75345-02	78649-21	75346-02	78649-22	75347-02			78649-23	75348-02			● Black
78650-20	78042-07	789- < 9	78043-07	78650-22	78044-07			78650-23	78045-07			● Black
78042-08	78042-09	78043-08	78043-09	78044-08	78044-09			78045-08	78045-09			● Beige
78059-08	78059-09	78047-08	78047-09	78048-08	78048-09			78049-08	78049-09			● Beige
Relief Garments Waist High												
78476-10	78476-11	78477-10	78477-11	78478-10	78478-11			78479-10	78479-11			● Beige
78480-10		78481-10		78482-10				78483-10				● Beige
Relief Garments Chaps												
78031-05		78032-05		78033-05				78034-05				● Beige L
78031-06		78032-06		78033-06				78034-06				● Beige R
78035-05		78036-05		78037-05				78041-05				● Beige L
78035-06		78036-06		78037-06				78041-06				● Beige R

JOBST Relief knee high and thigh high garments are supplied as a pair. Only thigh high variants are finished with a dotted silicone top band.

EVERYDAY Ready-To-Wear

JOBST® UltraSheer



Sheerness for beautiful legs



Allows air to pass




For a glossy look

Length	S		M		L		XL		Colour
	Open Toe	Closed Toe	Open Toe	Closed Toe	Open Toe	Closed Toe	Open Toe	Closed Toe	
Knee High	75165-30	75165-25	75166-30	75166-25	75167-30	75167-25	75168-30	75168-25	● Natural Regular
		75408-02		75409-02		75410-02		75415-02	● Natural Petite
		75165-04		75166-04		75167-04		75168-04	● Black
		75201-19		75202-19		75203-19		75204-19	● Natural Regular
		75435-02		75436-02		75437-02		75438-02	● Natural Petite
		75201-04		75202-04		75203-04		75204-04	● Black
Thigh High		75169-22		75170-22		75171-22		75172-22	● Natural
		75169-04		75170-04		75171-04		75172-04	● Black
		75205-23		75206-23		75207-23		75208-23	● Natural
		75205-04		75206-04		75207-04		75208-04	● Black
Sensitive Thigh High		77686-12		77686-13		77686-14		77686-15	● Natural
		77686-04		77686-05		77686-06		77686-07	● Black
		77687-12		77687-13		77687-14		77687-15	● Natural
		77687-04		77687-05		77687-06		77687-07	● Black
Waist High		75173-12		75174-12		75175-12		75176-12	● Natural
		75173-04		75174-04		75175-04		75176-04	● Black
		75209-13		75210-13		75211-13		75212-13	● Natural
		75209-04		75210-04		75211-04		75212-04	● Black

JOBST UltraSheer knee high and thigh high garments are supplied as a pair.
Only JOBST UltraSheer thigh high variants are finished with a silicone top band.

JOBST® forMen



Durable, long-lasting compression



StayFresh odour control



Specially constructed cushioned sole



Anatomic toe for enhanced comfort & fit

Length	S	M	L	LFC	XL	XLFC	Colour
	Closed Toe	Closed Toe	Closed Toe	Closed Toe Full Calf	Closed Toe	Closed Toe Full Calf	
Knee High	60523-40	60523-41	60523-42	60523-44	60523-43	60523-45	● Black
	60524-40	60524-41	60524-42	60524-44	60524-43	60524-45	● Khaki
	60529-40	60529-41	60529-42	60529-44	60529-43	60529-45	● Black
	60530-40	60530-41	60530-42	60530-44	60530-43	60530-45	● Khaki
	60535-40	60535-41	60535-42	60535-44	60535-43	60535-45	● Black
	60536-40	60536-41	60536-42	60536-44	60536-43	60536-45	● Khaki
Thigh High	60423-20	60423-21	60423-22		60423-23		● Black
	60423-60	60423-61	60423-62		60423-63		● Khaki
	60425-20	60425-21	60425-22		60425-23		● Black
	60425-60	60425-61	60425-62		60425-63		● Khaki
	60427-20	60427-21	60427-22		60427-23		● Black
	60427-60	60427-61	60427-62		60427-63		● Khaki

All JOBST forMen garments are supplied as a pair.
Knee high garments are finished with SoftFit silicone top band.
Thigh high garments are finished with a micro-dotted silicone top band.

SITUATIONAL Ready-To-Wear

JOBST® ActiveWear



Extra padding provides comfort



Allows air to pass



Helps keep feet dry

Length	S	M	L	XL	Colour
	Closed Toe	Closed Toe	Closed Toe	Closed Toe	
Knee High	75147-03	75148-03	75149-03	75150-03	● Black
	75147-02	75148-02	75149-02	75150-02	○ White
	75155-03	75156-03	75157-03	75158-03	● Black
	75155-02	75156-02	75157-02	75158-02	○ White
	75411-03	75412-03	75413-03	75414-03	● Black
	75411-02	75412-02	75413-02	75414-02	○ White

JOBST ActiveWear garments are supplied as a pair without a silicone top band.

JOBST® Sport



For comfort during sport



Helps to keep feet fresh



Helps keep feet dry

Length	S	M	L	XL	Colour
	Closed Toe	Closed Toe	Closed Toe	Closed Toe	
Knee High	75289-70	75289-71	75289-72	75289-73	■ Black
	75289-80	75289-81	75289-82	75289-83	■ Royal Blue
	75289-90	75289-91	75289-92	75289-93	■ Pink

JOBST Sport garments are supplied as a pair without a silicone top band.

JOBST® Maternity Opaque



Adapts to the growing belly



Helps keep feet dry



Blend of soft yarns

Length	S	M	L	Colour
	Closed Toe	Closed Toe	Closed Toe	
Maternity Waist High	60412-20	60412-21	60412-22	● Caramel
	60412-70	60412-71	60412-72	● Caramel
	60412-50	60412-51	60412-52	● Black

SITUATIONAL Ready-To-Wear

JOBST® Activa soSoft 8-15 mmHg



Easy donning



Soft touch-
blend of soft
yarns



Helps keep
feet dry

Length	S	M	L	XL	Colour
	Closed Toe	Closed Toe	Closed Toe	Closed Toe	
Knee High	60562-00	60562-01	60562-02	60562-03	● Black
	73680-12	60562-21	60562-22	60562-23	● Sand

JOBST® Activa soSoft measurement table

Size	S	M	L	XL
Ankle Circumference	18 - 21 cm	21 - 25 cm	25 - 29 cm	29 - 33 cm
Calf Circumference	28 - 38 cm	30 - 42 cm	32 - 46 cm	34 - 50 cm
AUS Women's Shoe	4.5 - 6.5	7 - 9	9.5 - 11	11.5 - 13.5
AUS Men's Shoe	3.5 - 5.5	6 - 8	8.5 - 10	10.5 - 12

JOBST Activa soSoft garments are supplied as a pair without a silicone top band.

Shoe size alone is not an indicator for correct sizing. Measurements of calf and ankle circumference are mandatory for the proper performance of this product.

SITUATIONAL Ready-To-Wear

JOBST® Travel Sock



For a more
convenient
journey



Allows air to
pass



Natural
wearing
comfort

Length	Size 1	Size 2	Size 3	Size 4	Size 5	Colour
	Closed Toe	Closed Toe	Closed Toe	Closed Toe	Closed Toe	
Knee High	78844-07	78845-07	78846-07	78847-07	79272-07	● Black
	78844-12	78845-12	78846-12	78847-12	79272-12	● Beige

JOBST® Travel Sock measurement table

Size	1	2	3	4	5
Ankle Circumference	18 - 21 cm	20 - 23 cm	22 - 25cm	24 - 27 cm	26 - 29 cm
Calf Circumference	27 - 37 cm	30 - 40 cm	33 - 43 cm	36 - 46 cm	39 - 49 cm
AUS Women's Shoe	4.5 - 6.5	7 - 9	9.5 - 11	11.5 - 13.5	-
AUS Men's Shoe	3.5 - 5.5	6 - 8	8.5 - 10	10.5 - 12	12.5 - 14
EU Shoe	37 - 38	39 - 40	41 - 42	43 - 44	45 - 46
UK Shoe	3 - 5	5 - 7	7 - 9	9 - 11	11 - 13

JOBST Travel garments are supplied as a pair without a silicone top band.

THE RIGHT FIT. LOOK NO FURTHER.



Trusted JOBST® Partners

You can purchase JOBST Ready-To-Wear Compression garments from a range of online stockists and local pharmacies.

Scan the QR code for a list of online stockists



JOBST® How to Measure Videos

Scan the QR code to access how to measure videos for knee high, thigh high and waist high garments



JOBST®,
an Essity brand

Always read the
label and follow the
directions for use.

ESS0876_v2
03/2026

Essity Health & Medical
PO Box 337, Mount Waverley,
Victoria, Australia 3149
medical.essity.com.au
jobst.com.au

Customer Service AU
T 1300 276 633 F 1300 998 830
Customer Service NZ
T 0508 276 111 F 0508 998 830